
**FINESSE**<sup>®</sup>

*Ultra*

*Breast Biopsy Probe*

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**BIOPSY SYSTEMS**

## Instructions for Use

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.


 Refer to the FINESSE® Ultra Biopsy Driver Instructions for Use for additional important information.

 Single-use only.

### A. Device Description

The FINESSE® Ultra Breast Biopsy System is composed of a driver and a probe and allows for obtaining and storing multiple samples with a single probe insertion. The components of the system are designed to operate safely when used together for diagnostic sampling of tissue during a biopsy procedure. The driver is a self-contained, handheld, reusable electro-mechanical vacuum-assisted biopsy device that utilizes a rechargeable lithium-ion battery and comes complete with AC power accessories and carrying case. The biopsy probe is composed of a hollow outer cutting cannula, an inner needle with a notched sample chamber, and an integrated coaxial cannula. The outer cannula and integrated coaxial cannula have sharpened cutting edges and the inner needle has a sharpened trocar tip. The driver utilizes DC operated motors to create a vacuum while simultaneously moving the outer cannula backward to expose the sample chamber to draw tissue in. The outer cannula is released and moves forward over the sample chamber to cut the tissue. The tissue sample is transported via the sample chamber to the sample container. The sample chamber is returned to the original position and is available to acquire another sample. The integrated coaxial cannula may be detached after the biopsy and remain in the breast to retain a track to the biopsy site when placing a tissue marker.

**NOTE:** The integrated coaxial cannula must remain attached while acquiring tissue samples.

Catalog Number	Description	Gauge	A	B	C	D
F14105US	Ultrasound Biopsy Probe	14G	105mm	17mm	7mm	20mm
F10105US	Ultrasound Biopsy Probe	10G	105mm	21mm	9mm	20mm

### B. Indications for Use

The FINESSE® Ultra Breast Biopsy System is indicated to obtain tissue samples from the breast or axillary lymph nodes for diagnostic analysis of breast abnormalities. The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality. The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

### C. Contraindications

1. The FINESSE® Ultra Breast Biopsy System is for diagnostic use only, NOT for therapeutic use.
2. The FINESSE® Ultra Breast Biopsy System is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

### D. Warnings

1. **Patients who may have a bleeding disorder or who are receiving anticoagulant therapy may be at increased risk of complications.**
2. **As with any biopsy instrument, there is a potential for infection.**
3. **Do not reuse FINESSE® Ultra Biopsy Probes. Reusing the probe bears the risk of cross-patient contamination as biopsy probes – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the probe for an indeterminable period of time. The residue of biological material can promote the contamination of the probe with pyrogens or microorganisms which may lead to infectious complications.**
4. **Do not resterilize FINESSE® Ultra Biopsy Probes. After resterilization, the sterility of the probe is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilizing the probe increases the probability that it will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**
5. **The FINESSE® Ultra Biopsy Probe must only be used with the FINESSE® Ultra Biopsy Driver.**
6. **All breast biopsies should be performed under imaging guidance to confirm the probe position relative to the target region to be sampled and to help mitigate the occurrence of a false negative biopsy. The FINESSE® Ultra Biopsy Probe is intended for use with ultrasound imaging only.**

### E. Precautions

1. Do not use the FINESSE® Ultra Biopsy Probe without the integrated coaxial cannula. The integrated coaxial cannula may be removed after the biopsy to retain a track to the biopsy site when placing a tissue marker.
2. The FINESSE® Ultra Breast Biopsy System should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.

## F. Potential Complications

1. Potential complications are those associated with any percutaneous removal/ biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, a non-healing wound, pain and tissue adherence to the biopsy probe while removing it from the breast.
2. As per routine biopsy procedures, it may be necessary to cut tissue adhering to the biopsy probe while removing it from the breast.

## G. Equipment Required


The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- FINESSE® Ultra Biopsy Driver
- FINESSE® Ultra Biopsy Probe
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary

## H. Directions for Use

### Loading the Probe into the Driver

1. Ensure that the driver has been fully charged. The driver should be charged overnight after daily use.
2. Remove the probe from the package using aseptic technique by grasping the probe from below. Ensure that the probe remains sterile.
3. Insert the probe into the driver (Figure 1).
4. Remove the retaining pin (Figure 2).
5. Slide the probe carriage back ensuring the tab securely locks into place (Figure 3).


**NOTE: The driver does not contain an on-off switch. The driver automatically recognizes the probe type, senses when the probe has been properly inserted and performs an initialization to synchronize the driver and probe. The driver sounds three beeps and the “PRIME/PIERCE” and “SAMPLE” indicator lights display green indicating the probe is properly inserted. During the initialization the motor is audible.**

6. After initialization is complete, remove the protective sleeve from the probe.

**NOTE: Refer to the FINESSE® Ultra Biopsy Driver Instructions for Use for important information on the operation of the driver.**

### Performing a Biopsy

1. Prepare the biopsy site using appropriate aseptic techniques and adequate local anesthesia. Appropriate precautions should be followed for any percutaneous procedure (personal protective gear such as gloves and glasses/shields, etc. should be worn).
2. Make a small incision or puncture the skin with a scalpel.
3. Using ultrasound guidance, insert the probe through the incision and position the tip to the appropriate depth:
  - If the lesion is going to be pierced: Position the tip of the probe proximal to the edge of the lesion.
  - If the lesion is not going to be pierced: Position the sampling chamber to the center of the target region. (The center of the sampling chamber is approximately 17mm from the tip of the probe.)
4. If piercing is desired, prior to inserting the probe into the breast press the “PRIME/PIERCE” button to retract the probe 20-22mm. After inserting the probe press the “PRIME/PIERCE” button to advance the probe 20-22mm into the region of interest prior to taking the first tissue sample (Figure 4).

**NOTE: Image guidance should be used to confirm probe position relative to the target region to be sampled.**

**NOTE: During the initial positioning of the probe, it is recommended to advance the probe into the tissue using the “PRIME/PIERCE” option. During subsequent samplings, the physician may wish to forego this step.**

**NOTE: The probe may be primed for piercing while inside the patient if additional piercing shots are desired to access dense tissue or lesions.**

5. Begin tissue sampling by pressing the “SAMPLE” button. The sampling process proceeds automatically: A vacuum is created, the outer cannula of the probe is automatically retracted, tissue is drawn by vacuum into the sampling chamber and the outer cannula is released cutting the tissue. The tissue sample is then transported to the sample container. During this process the “SAMPLE” indicator light continuously illuminates orange and the motor will be audible. The sample container illuminates to aid in visualizing the tissue sample. After the tissue sample is deposited into the sample container the system resets to acquire another sample if desired.

**NOTE: Do not press on the sample container or remove the sample container while acquiring tissue samples.**


6. To acquire multiple samples, repeat Step 5 above.

**NOTE: Collecting multiple tissue samples may help minimize the risk of a false negative biopsy.**

7. After obtaining the final tissue sample, remove the probe from the breast and care for the incision site as appropriate.

**NOTE:** Prior to removing the probe from the breast, the integrated coaxial cannula may be detached and remain in the breast to retain a track to the biopsy site while placing a tissue marker (Figure 5).

Figure 5


8. Remove the tissue samples from the sample container by twisting the lid counter clockwise. If additional samples are required, ensure the sample container lid is in place before acquiring additional samples.

**NOTE:** The sample container may be removed completely from the probe by pulling up on the lid and may be removed before or after the probe is removed from the driver. If the probe is loaded into the driver and the sample container is removed from the probe, the “PRIME/PIERCE” and “SAMPLE” indicator lights will flash orange. After the sample container is reinserted into the probe, the “PRIME/PIERCE” and “SAMPLE” indicator lights will continuously illuminate green.

9. Remove the probe from the driver by pressing down on the locking tab, pushing the carriage completely forward and pulling the probe straight from the driver. (Figure 6).

Figure 6

**NOTE:** Once a probe has been removed from the driver following a procedure, the probe cannot be re-inserted. If additional tissue samples are required, insert a new probe into the driver.

**NOTE:** Driver will initiate mechanical reset following probe removal. Do not insert a new probe prior to completion of this reset.

10. Refer to the FINESSE® Ultra Biopsy Driver Instructions for Use for important information on cleaning and charging the driver.

**NOTE:** The probe, coaxial cannula and sample container are for single use only. After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

#### I. How Supplied

The FINESSE® Ultra Biopsy Probes are supplied sterile for single patient use only.

#### J. Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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**Bard and Finesse are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.**
**Copyright ©2011 C. R. Bard, Inc. All rights reserved.**
Copyright ©2011 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2011 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2011 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2011 C. R. Bard, Inc. Todos los derechos reservados.