TITAN III[®] Serum Protein Electrophoresis

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INTENDED PURPOSE

The Helena Serum Protein Electrophoresis procedure is intended for the separation and quantitation of serum proteins by electrophoresis on cellulose acetate.

Serum contains over one hundred individual proteins, each with a specific set of functions and subject to specific variation in concentration under different pathologic conditions. Since the introduction of moving-boundary electrophoresis by Tiselius² and the subsequent use of zone electrophoresis, serum proteins have been fractionated on the basis of their electrical charge into five classical fractions: albumin, α_1 , α_2 , beta, and gamma proteins. Each of these classical electrophoretic zones (with the exception of albumin) normally contains two or more components. Approximately fifteen serum proteins have been studied extensively because they may be measured easily.

Proteins are large molecules composed of covalently linked amino acids. Depending on electron distributions resulting from covalent or ionic bonding of structural subgroups, proteins have different electrical charges at a given pH. In the Helena Serum Protein procedure, the proteins are separated according to their respective electrical charges at pH 8.8 on a cellulose acetate plate using both the electrophoretic and electroendosmotic forces present in the system. After the proteins are separated, the plate is placed in a solution of sulfosalicylic acid and trichloroacetic acid (to immobilize the proteins) and Ponceau S (to stain the protein bands). The staining intensity is related to protein concentration. After dehydration in methanol, the plate background is then rendered transparent by treatment with a clearing solution.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. Ponceau S Stain (Cat. No. 5526)

The stain is supplied as a powder and is a mix of Ponceau S and Sulphosalicylic Acid. Dissolve the contents of the bottle in 1 litre of deionised water and mix for 30 minutes or until completely dissolved.

2. Electra[®] HR Buffer (Cat. No. 5805)

Contains a Tris-barbital Sodium Barbital buffer, pH 8.6-9.0. Dissolve one pack of dry buffer in 750ml of purified water. The buffer is ready for use when all material is dissolved and completely mixed.

3. Clear Aid (Cat No. 5005)

Contains polyethylene glycol. To prepare the clearing solution, mix 30 parts glacial acetic acid, 70 parts absolute methanol, and 4 parts Clear Aid. Stir until well mixed.

4. Titan III[®] Cellulose Acetate Plate (Cat. No. 3023 or 3024)

Contains a cellulose acetate plate. Refer to instructions for use for preparation instructions.

5. Other kit components

Each kit contains Instructions For Use.

STORAGE AND SHELF-LIFE

1. Ponceau S Stain (Cat. No. 5526)

The stain may be stored as packaged or in a tightly closed staining dish at 15...30°C. Unused stain is stable until the expiry date indicated on the label. Deterioration may be indicated by excessive evaporation or large amounts of precipitate occurring.

2. Electra® HR Buffer (Cat. No. 5805)

The packaged buffer is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C when stored tightly closed. Deterioration may be indicated by signs of dampness or discoloration or if it becomes turbid.

3. Clear Aid (Cat No. 5005)

Clear Aid is stable at 15...30°C until the expiry date indicated on the label. Deterioration may be indicated by gross contamination or discoloration. Discard if the plates appear cloudy after the clearing procedure.

4. Titan III® Cellulose Acetate Plate (Cat. No. 3023 or 3024)

The cellulose acetate plate should be stored at 15...30°C and is stable until the expiry date indicated on the label.

SAMPLE COLLECTION AND PREPARATION

Serum is the preferred specimen. The use of plasma should be avoided, as fibrinogen will appear as a distinct narrow band between the beta and gamma fractions. Cereospinal fluid may be used if concentrated approximately 100 times; urine may be used if concentrated up to 300 times, depending on original protein concentration.

Interfering factors:

1. Hemolysis may cause false elevation in the alpha₂ and beta fractions.
2. Inaccurate results may be obtained on specimens left uncovered, due to evaporation.

Storage and Stability: Fresh serum or plasma is the specimen of choice. If storage is necessary, samples may be stored covered at 2...6°C for 48 hours. Cerebrospinal fluid and urine specimens may be used after proper concentration with a concentrator.

ITEMS REQUIRED BUT NOT PROVIDED:

5% acetic acid: Add 50ml of glacial acetic acid to 950ml of purified water.

Absolute methanol, reagent grade

Cat. No. 3024 Titan® III Cellulose Acetate Plates (94 x 76 mm)

Cat. No. 3023 Titan® III Cellulose Acetate Plates (60 x 76mm)

Cat. No. 4084 Super Z Applicator

Cat. No. 4090 Super Z-12 Applicator

Cat. No. 4085 Super Z-12 Sample Well Plate (2)

Cat. No. 4096 Super Z-12 Sample Well Plate (2)

Cat. No. 4094 Super CPK Aligning Base

Cat. No. 4086 Super Z Aligning Base

Cat. No. 5114 700 Staining Set

Cat. No. 5122 1000 Staining Set

Cat. No. 5526 Ponceau S

Cat. No. 5093 Bufferizer

Cat. No. 5005 Clear Aid

- Cat. No. 5090 Zip Zone® Prep
- Cat. No. 1283 Zip Zone® Chamber
- Cat. No. 5081 Zip Zone® Chamber Wicks
- Cat. No. 5037 Blotter Pads (102 x 108 mm)
- Cat. No. 5034 Blotter Pads (76 x 102mm)
- Cat. No. 5805 Electra® HR Buffer
- Cat. No. 5053 Titan Plastic Envelopes (102 x 120 mm)
- Cat. No. 5052 Titan Plastic Envelopes (63 x 120mm)
- Cat. No. 6008 Microdispenser and Tubes
- Cat. No. 8008 Micro-Hood (220V)
- Cat. No. 5205 Electrophoresis Manual
- Cat. No. 5211 Serum Protein Report Forms
- Cat. No. 5000 Helena Marker
- Cat. No. 5109 Quality Control Chart
- Cat. No. 5015 Titan Identification Labels
- Cat. No. 5002 Glue Stick

STEP-BY-STEP PROCEDURE

1. Properly code the required number of Titan III Plates by marking on the glossy, hard (Mylar)® side with a Helena BioSciences marker. It is suggested that the identification mark be placed in one corner so that it is always aligned with sample No. 1.
2. Soak the plates for 20 minutes in diluted HR buffer. The plates should be soaked in the Bufferizer according to the instructions for use included with the Bufferizer. Alternatively, the plates may be wetted by slowly and uniformly lowering a rack of plates into the buffer such that air is not trapped in the plates. **NOTE:** The same soaking buffer may be used for soaking up to 12 plates or for approximately one week if stored tightly closed. If used for a more prolonged period residual solvents from the plate may build up in the buffer or evaporation may alter buffer concentration.
3. Pour approximately 100ml of diluted HR Buffer into each outer section of the chamber.
4. Wet two disposable wicks in the buffer and drape one over each support bridge, ensuring it makes contact with the buffer and that there are no air bubbles under the wicks.
5. Cover the chamber to prevent buffer evaporation.
6. Fill each well in the sample plate with 3 μ l of sample using the microdispenser. Expel the samples as a bead on the tip of the glass tube; then touch this bead to the well. Cover the samples with a glass slide if they are not used within 2 minutes.
7. Prime the applicator by depressing the tips into the sample wells 3 or 4 times. Apply this loading to a piece of blotter paper. DO NOT load the applicator again at this point, but proceed to the next step.
8. Remove the wetted plate from the buffer and blot once firmly between two blotters. Before placing the plate in the aligning base place a drop of water or buffer on the center of the aligning base. This prevents the plate from shifting during the sample application. Place the plate in the aligning base, cellulose acetate side up, aligning the bottom edge of the plate with the black scribe line marked "CENTER APPLICATION". The identification mark should be aligned with sample No. 1.
9. Apply the sample to the plate by depressing the applicator tips into the sample well 3 or 4 times and transfer the applicator to the aligning base. Press the button down and hold it for 5 seconds.

10. Place the plate(s), cellulose acetate side down, in the chamber. Place a weight on the plate(s) to ensure contact with the wicks. Cover the chamber securely and wait 30 seconds for the plate(s) to equilibrate.
11. Electrophorese the plate(s): 15 minutes, 180 volts. **NOTE:** Power must be applied within 5 minutes after the plate(s) has been placed in the chamber.
12. Following electrophoresis, remove the plate(s) from the chamber and place in 40-50ml of Ponceau S stain for 6 minutes. (When staining 2 or more plates, carry out the protocol vertically in a rack. **NOTE:** The stain may be reused until the plate background contains stain precipitate.
13. Destain in 3 x 2 minute washes of 5% acetic acid or until the plate background is white. The plate(s) may be dried and stored as a permanent record at this point. If a transparent background is desired (i.e. for densitometry), proceed to the next step.
14. Dehydrate by rinsing the plate in 2 absolute methanol washes of 2 minutes each wash. Allow the plate to drain for 5-10 seconds before placing in the next solution.
15. Place the plate into the clearing solution for 5-10 minutes.
16. Drain off excess solution, then place the plate, acetate side up, onto a blotter, and into a drying oven at 50...60°C for 15 minutes or until dry.

Evaluation of the Protein Bands

Scan the plates in a densitometer using a 525 nm filter and the narrow slit.

Stability of End Product

The completed, dried Serum Protein plate is stable for an indefinite period of time and may be stored in Titan Plastic Envelopes,

Calibration

The Optical Density Step Tablet (Cat. No. 1030) should be used to ensure the linearity of the instrument and a Neutral Density Densitometer Standard should be used to validate the zero adjustment and sensitivity of the instrument.

INTERPRETATION OF RESULTS^{5,7}

The fastest moving band, and normally the most prominent is the albumin band found closest to the anodic edge of the plate. The faint band next to this is alpha₁ globulin, followed by alpha₂ globulin, beta, and gamma globulins.

Calculation of the Unknown

The Helena BioSciences CliniScan™ Densitometer and other Helena BioSciences densitometers with computer accessories will automatically print the relative percent and the absolute values for each band. Alternately, the relative percent of each band can be calculated manually from the integrator scale on the scan. Refer to the Operator's Manual provided with the densitometer. The relative percent of each band is calculated by the following formula:

$$\frac{\text{No. Integration Units of the Band}}{\text{Total Integration Units}} \times 100 = \text{Relative Percent of the Band}$$

$$\text{Relative Percent of the Band} \times \text{Total Serum Protein} = \text{Absolute Value of Protein in each band.}$$

QUALITY CONTROL

Two different types of normal range quality control materials are available for use with the serum protein electrophoresis kits. Use at least one of these control sera on each plate run.

Cat. No. 5112 Electrophoresis Serum Control

Cat. No. 7024 Kentrol Serum Control-Normal

REFERENCE VALUES

The expected values for serum protein electrophoresis on cellulose acetate stained with Ponceau S were determined from a study of 51 normal subjects. These values are for illustrative purposes only. Each laboratory should establish its own range.

Protein Fraction	Range gd/L
Albumin	3.63 - 4.91
Alpha ₁	0.11 - 0.35
Alpha ₂	0.65 - 1.17
Beta	0.74 - 1.26
Gamma	0.58 - 1.74

Variations of Expected Values⁴

Studies show that values are the same for both males and nonpregnant females. (Some differences are seen in pregnant females at term and in women on oral contraceptives). Age has some effect on normal levels. Cord blood has decreased total protein, albumin, alpha₂ and beta fractions, slightly increased alpha₁, and normal or increased gamma fractions (largely of maternal origin). The gamma drops rapidly until about 3 months of age, while the other fractions have reached adult levels by this time. Adult levels of the gamma globulins are not reached until 10-16 years of age. The albumin decreases and beta globulins increase after the age of 40.

Results on normal individuals will cover age and sex-related variations and day-to-day biologic variations. Disease states in which abnormal patterns are observed include inflammatory response, rheumatic disease, liver diseases, protein-loss disorders, plasma cell dysorasis, and genetic deficiencies. Variant patterns have also been observed during pregnancy.

Further Testing Required

The serum protein electropherogram or densitometric tracing should be evaluated for abnormalities. If abnormalities are observed, appropriate follow-up studies should be initiated. These may include Immunoelectrophoresis, immunofixation, quantitation of individual component immunoglobulins, bone marrow examination, and other appropriate tests.

PERFORMANCE CHARACTERISTICS

A normal serum was run 26 consecutive times using the method described above. The following data was obtained:

Protein Fraction	Mean (Relative %)	SD (+/-)	CV (%)
Albumin	55.7	1.4	2.5
Alpha ₁	3.1	0.4	12.5
Alpha ₂	11.3	0.4	3.8
Beta	11.6	0.5	4.0
Gamma	18.1	0.6	3.5

BIBLIOGRAPHY

1. Alper, C.A., Plasma Protein Measurements as a Diagnostic Aid, N Eng J Med, 291:287-290, 1974.
2. Tiselius, A., A New Approach for Electrophoretic Analysis of Colloidal Mixtures, Trans Faraday Soc, 33:524-531, 1937.
3. Ritzmann, S.E. and Daniels, J.C., Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays in Laboratory Medicine, Harper and Row Inc., Hagerstown, 1979.
4. Ritzmann, S.E. and Daniels, J.C., Serum Protein Abnormalities: Diagnostic and Clinical Aspects, Little, Alan R. Liss Inc., 1982.
5. Killingsworth, L.M. et al, Protein Analysis, Diag Med, Jan/Feb., 47-58, 1980.
6. Henry, R.J., Cannon, D.C., and Winkelman, J.W., ed., Clinical Chemistry Principals and Techniques, 437, Harper and Row, Hagerstown, 1974.
7. Killingsworth, L.M., Plasma Protein Patterns in Health and Disease, CRC Crit Rev in Clin Lab Sci, 1-30, August, 1979.

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LITERATUR

1. Alper, C.A., Plasma Protein Measurements as a Diagnostic Aid, N Eng J Med, 291:287-290.1974.
2. Tiselius, A., A New Approach for Electrophoretic Analysis of Colloidal Mixtures, Trans Faraday Soc, 33:524-531, 1937.
3. Ritzmann, S.E. and Daniels, J.C., Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays in Laboratory Medicine, Harper and Row Inc., Hagerstown, 1979.
4. Ritzmann, S.E. and Daniels, J.C., Serum Protein Abnormalities: Diagnostic and Clinical Aspects, Little, Alan R. Liss Inc., 1982.
5. Killingsworth, L.M. et al, Protein Analysis, Diag Med, Jan/Feb., 47-58, 1980.
6. Henry, R.J., Cannon, D.C., and Winkelman, J.W., ed., Clinical Chemistry Principals and Techniques, 437, Harper and Row, Hagerstown, 1974.
7. Killingsworth, L.M., Plasma Protein Patterns in Health and Disease, CRC Crit Rev in Clin Lab Sci, 1-30, August, 1979.

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BIBLIOGRAFIA

1. Alper, C.A., Plasma Protein Measurements as a Diagnostic Aid, N Eng J Med, 291:287-290, 1974.
2. Tiselius, A., A New Approach for Electrophoretic Analysis of Colloidal Mixtures, Trans Faraday Soc, 33:524-531, 1937.
3. Ritzmann, S.E. and Daniels, J.C., Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays in Laboratory Medicine, Harper and Row Inc., Hagerstown, 1979.
4. Ritzmann, S.E. and Daniels, J.C., Serum Protein Abnormalities: Diagnostic and Clinical Aspects, Little, Alan R. Liss Inc., 1982.
5. Killingsworth, L.M. et al, Protein Analysis, Diag Med, Jan/Feb., 47-58, 1980.
6. Henry, R.J., Cannon, D.C., and Winkelman, J.W., ed., Clinical Chemistry Principles and Techniques, 437, Harper and Row, Hagerstown, 1974.
7. Killingsworth, L.M., Plasma Protein Patterns in Health and Disease, CRC Crit Rev in Clin Lab Sci, 1-30, August, 1979.

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1. Alper, C.A., Plasma Protein Measurements as a Diagnostic Aid, N Eng J Med, 291:287-290. 1974.
2. Tiselius, A., A New Approach for Electrophoretic Analysis of Colloidal Mixtures, Trans Faraday Soc, 33:524-531, 1937.
3. Ritzmann, S.E. and Daniels, J.C., Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays in Laboratory Medicine, Harper and Row Inc., Hagerstown, 1979.
4. Ritzmann, S.E. and Daniels, J.C., Serum Protein Abnormalities: Diagnostic and Clinical Aspects, Little, Alan R. Liss Inc., 1982.
5. Killingsworth, L.M. et al, Protein Analysis, Diag Med, Jan/Feb., 47-58, 1980.
6. Henry, R.J., Cannon, D.C., and Winkelman, J.W., ed., Clinical Chemistry Principles and Techniques, 437, Harper and Row, Hagerstown, 1974.
7. Killingsworth, L.M., Plasma Protein Patterns in Health and Disease, CRC Crit Rev in Clin Lab Sci, 1-30, August, 1979.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1436P 2007/06 (5)