IL Test™ Total Bilirubin
Cat. No. 0018258340 - 4 x 45 mL BILI T R1 / 8 x 19.5 mL BILI T/D R2 / 1 x 4.5 mL Sodium Nitrite
Cat. No. 0018254640 - 4 x 100 mL Total Bilirubin R1 / 4 x 74 mL T/D Bilirubin R2 / 1 x 13 mL Sodium Nitrite

Intended Use
IL Test™ Total Bilirubin is intended for the quantitative in vitro diagnostic determination of Direct Bilirubin in human serum and plasma using the ILab Chemistry Systems.
Summary

Total bilirubin measurements are used in the diagnosis and treatment of various liver diseases, and of haemolytic and metabolic disorders in adults and newborns.¹
Principle

End point analysis, modified Jendrassik-Grof assay.²

Total bilirubin + diazotized sulfanilic acid

Lithium dodecylsulfate

➔

 azobilirubin

The increase in absorbance due to the formation of azobilirubin is proportional to the total bilirubin concentration in the sample. Absorbance measurements are taken at 546 nm.
Composition

BILI T R1 containing: Lithium dodecyl sulfate 173 mmol/L.
BILI T/D R2 containing: Hydrochloric acid 180 mmol/L, and Sulfanilic acid 28.9 mmol/L.
Sodium Nitrite containing: Sodium nitrite 72.5 mmol/L.

PRECAUTION AND WARNINGS
TOTAL BILIRUBIN R1
Hazard class: Xi (Irritant)
Risk phrases: R36/37/38: Irritating to eye, respiration system and skin.
R43: May cause sensitization by skin contact.
Safety phrases: S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.
S64: If swallowed, rinse mouth with water (only if the person is conscious).

T/D BILIRUBIN R2 is not classified as hazardous
Hazard class: none
Risk phrases: none
Safety phrases: none
SODIUM NITRITE is not classified as hazardous
Hazard class: none
Risk phrases: none
Safety phrases: none
This product is For *in vitro* Diagnostic Use.

Reagent Preparation
Cat.No 0018258340
BILI T R1: ready to use. Uncap the container and place into reagent tray
BILI T/D R2: add 10 drops (0.5 mL) of Sodium Nitrite into the BILI T/D R2, bottle, cap and mix by inverting a few times, then uncap and place into reagent tray.

Cat.No 0018254640
BILI T R1: ready to use. Uncap the container and place into reagent tray
BILI T/D R2: add 35 drops (1.75 mL) of Sodium Nitrite into the BILI T/D R2, bottle, cap and mix by inverting a few times, then uncap and place into reagent tray. This reagent may be used immediately after mixing.
Note: for both Cat. No the BILI T/D R2 reagent is also used in the Direct Bilirubin assay.

Reagent storage and Stability
Unopened reagents are stable until expiration date indicated on the bottle when stored at 15-25°C.
BILI D R1: on board until expired.
BILI T/D R2: 30 days on board.

Instrument/Test Procedures
For the assay procedure refer to the relevant instrument Application Manual and to the relevant Instrument Operator Manual.

Specimen collection and Preparation
Serum, lithium heparinised plasma.
Samples are stable for 3 days at 2-8°C.^{3,4} Protect serum and plasma from exposure to direct sunlight or fluorescent light.

Calibration
Use ReferrIL™ G Calibrator (Cat. No. 0018257000). See calibrator insert sheet for lot specific concentration and calibrator traceability.
Recalibrate the instrument every 30 days or when a new lot of reagent is used. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Quality Control
Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice. Normal and abnormal controls such as SeraChem® Control Level 1 (Cat. No. 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program.
Refer to Westgard *et al.*⁵ for identification and resolution of out-of-control situations.
Calculation of Analytical Results
Refer to the Application Sheet and to the ILab Operator Manual.
To convert mg/dL to SI units (µmol/L) multiply by 17.1.¹

Reference Interval
Newborn 0 to 5 days old: 1.4 to12.0 mg/dL (23.9 - 205.2 µmol/L).¹
Adults: 0.3 to 1.2 mg/dL (5.1 - 20.5 µmol/L).¹
It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Performance Characteristics
All performance data were obtained on ILab 600/650. For performance data on other ILab systems refer to the relevant application Manual.
Limitation/Interfering Substances

No interference from lipemia up to sample absorbance of 2.0/cm at 660 nm (640 mg/dL or 7.2 mmol/L triglycerides). No interference up to 250 mg/dL (0.150 mmol/L) hemoglobin.
For a comprehensive review of interfering substances, refer to the publication by Young.⁶

Precision							
Typical precision results obtained according to NCCLS protocol. ⁷							
	Samples/Runs	Mean (mg/dL)	SD (mg/dL)	CV (%)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Within Run	3/20	0.54	0.01	2.6	3.70	0.02	0.5
Total	3/20	0.54	0.04	6.7	3.70	0.05	1.5

Method Comparison	
Comparison Method (x)	Same reagent
Comparison Instrument	ILab 900
Slope	0.955
y intercept	0.034
Range (mg/dL)	0.02 - 26.52
Mean x (mg/dL)	1.78
Mean y (mg/dL)	1.73
r	0.999
SE	0.096
n	99

Linearity		
	Linearity range (mg/dL)	Linearity range (µmol/L)
no rerun	0.1 - 31.0	2 - 530
with rerun	0.1 - 70.0	2 - 1197

Sensitivity	
milliabsorbance change per 1 mg/dL 44	
milliabsorbance change per 1 µmol/L 2.57	

Minimum Detection Limit	
0.1 mg/dL (1.71 µmol/L).	

IL Test™ Total Bilirubin
Cat. No. 0018258340 - 4 x 45 mL BILI T R1 / 8 x 19.5 mL BILI T/D R2 / 1 x 4.5 mL Sodium Nitrite
Cat. No. 0018254640 - 4 x 100 mL Total Bilirubin R1 / 4 x 74 mL T/D Bilirubin R2 / 1 x 13 mL Sodium Nitrite

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ILab 600/650 Test Parameters

Photometric Test Parameters		Serum
Test No.		4
Test Name, Test Code		T. Bilirubin, Bili T
Sample Type		Serum
Reporting Unit, Decimal Points		mg/dL (µmol/L), 1 (1)
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		End Pont, 16/22
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		546
Sampling Conditions		
Sampling 1	Sample Vol.	15
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	7
	Sample/Diluent Vol.	0/0
Sampling 3 and Sampling 4		***
Diluent Code		N/A
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
HIGH!, ABS!		Sampling 2
Reagent Volumes		
R1	Code	01041
	Rgt/Dil. Vol., Stirring	135/50, ON
	Low Vol. Warning Limit	***
	Stability (days)	0 (until expiration)
R2	Code	04052
	Rgt/Dil. Vol., Stirring	100/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	30

Bibliography / Bibliografia / Bibliografia

- Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, 3 rd edition 1999, pp. 1130, 1803.
- Pearlman and Lee, Detection of Total Bilirubin in Serum with use of Surfactant as Solubilizing Agents, Clin.Chem. 20/4 pp. 447-453 (1974)
- Henry, RJ, et al., Clinical Chemistry: Principles and Techniques , 2nd edition, Hagerstown, MD; Harper and Row, 1059, 1974.
- Tietz, NW, ed., Clinical Guide to Laboratory Tests, 3rd Edition, WB Saunders Company, Philadelphia, PA, 1995, p. 88.
- Westgard J.O., and Barry P.L., Cost Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, D.C., 1986.

- Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Wa shington, DC, 5th ed. 2000.
- Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

IL Test™ and SeraChem®, are trademarks of Instrumentation Laboratory Company.

Ranges and Evaluation Criteria		Serum
Normal Range-Male		0.3 - 1.2 (5.1 - 20.5)
Normal Range-Female		0.3 - 1.2 (5.1 - 20.5)
Normal Range-Other		**
Valid Range		0.1 - 31.0 (2 - 530)
Hemolysis/Icterus/Lipemia Limit		***
Reaction Slope		Positive
Absorbance Limit		Above, 3500
Prozone Limit		N/A
Non Linear Limit		N/A
Slope/Intercept Correction		1/0
Qualitative Report		OFF
Calibration Conditions		
Calibration		1 Point, Linear, ** Reps
Stability (days)		30
Calibrator, Concentration		Ref G, *
R-Blank Limit (mAbs)		***
Cal. Reps Range (%)		***
Min. Cal. Response (mAbs)		***
Cal. Factor Change (%)		***
M-Point Curve Fit (%)		N/A
Reagent Blank		ON
Auto R-Blank by Bottle		ON

- * Lot dependent - see package insert / Vedere l'inserito del ReferrIL G per il valore appropriato / Ver en el prospecto del ReferrIL G
- ** Operator definable data / Definibile dall'operatore / Datos definidos por el operador
- *** Optional/ Facoltativo / Opcional
- N/A Not applicable to this test / Non applicabile a questo test / No aplicable a este método

IVD

In Vitro Diagnostic Medical Device

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ATTENZIONE! Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.