IL Test™ Total Protein

Cat. No. 0018258640 - 5 x 50 mL Total Protein

Cat. No. 0018251440 - 5 x 100 mL Total Protein

Intended Use

IL Test™ Total Protein is intended for the quantitative in vitro diagnostic determination of total protein in human serum and plasma, using the ILab Chemistry Systems.

Summary

Changes in serum total protein concentration are generally caused by a change in the volume of plasma water or changes in the concentration of one or more of the specific proteins in the plasma. Total protein measurements are used in the diagnosis and treatment of a variety of liver and kidney diseases as well as metabolic and nutritional disorders.^{1,2}

Principle

End point analysis (modified biuret methodology) based on the reaction of peptide bonds with Cu⁺⁺ ions in alkaline solution to form a coloured complex.

The increase in absorbance due to the complex is proportional to the protein concentration in the sample.

Absorbance measurements are taken at 546/700 nm.

Composition

Total Protein containing: Cupric sulfate 9.0 mmol/L, Sodium Hydroxide 2.0 mol/L, Potassium iodide 9.0 mmol/L, EDTA 20 mmol/L.

PRECAUTION AND WARNINGS

TOTAL PROTEIN

Hazard class: C (Corrosive)

Risk phrases: R35: Causes severe burns.

R52/53: Harmful to aquatic organisms, may cause long term adverse effects in the aquatic environment.

Safety phrases: S24/25: Avoid contact with skin and eyes.

S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S37/39: Wear suitable gloves and eye/face protection.

S45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

S61: Avoid release to the environment. Refer to special instructions/ safety data sheet.

This product is For *in vitro* Diagnostic Use.

Reagent Preparation

TP: ready to use. Uncap the container and place into reagent tray.

Reagent storage and Stability

Unopened reagent is stable until the expiration date indicated on the bottle when stored at 15-25°C.

40 days on-board.

Instrument/Test Procedures

For the assay procedure refer to the relevant instrument Application Manual and to the relevant Instrument Operator Manual.

Specimen collection and Preparation

Serum and plasma.

Samples may be stored for 7 days at 15-25°C or 1 month at 2-8°C.⁵

Calibration

Use Referril™ G Calibrator (Cat. No. 0018257000).

See calibrator insert sheet for lot specific concentration and calibrator traceability.

Recalibrate the instrument every 21 days or when a new lot of reagent is used. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Quality Control

Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice. Normal and abnormal controls such as SeraChem® Control Level 1 (Cat. No. 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program.

Refer to Westgard et al.⁴ for identification and resolution of out-of-control situations.

Calculation of Analytical Results

Refer to the Application Sheet and to the ILab Operator Manual.

To convert mg/dL to g/L multiply by 10.²

Reference Interval

6.4 to 8.3 g/dL (64 - 83 g/L).²

It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

ENGLISH - Insert revision 10/2011

Performance Characteristics

All performance data were obtained on ILab 600/650. For performance data on other ILab systems refer to the relevant application Manual.

Limitation/Interfering Substances

No interference from lipemia up to sample absorbance of 1.0/cm at 660 nm (330 mg/dL or 3.7 mmol/L triglycerides). No interference up to 100 mg/dL (0.06 mmol/L) hemoglobin. No interference up to 20 mg/dL (340 µmol/L) bilirubin.

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Precision

Typical precision results obtained according to NCCLS protocol.⁶

	Samples/Runs	Mean (g/dL)	CV (%)	Mean (g/dL)	CV (%)
Within Run	3/20	3.8	1.1	6.9	1.0
Total	3/20	3.8	1.2	6.9	1.1

Method Comparison

Comparison Method (x)	Same reagent
Comparison Instrument	ILab 900
Slope	0.966
y intercept	0.16
Range (g/dL)	4.6 - 8.8
Mean x (g/dL)	7.3
Mean y (g/dL)	7.2
r	0.992
SE	0.12
n	98

Linearity

	Linearity range (g/dL)	Linearity range (g/L)
with rerun	1 - 13.8	10 - 138

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1. Tietz Clinical Guide to Laboratory Tests, W.B. Saunders, Philadelphia, PA, Third edition 1995, pp. 518-519.
2. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, Third edition 1999, pp. 523 and 1832.
3. Tietz Fundamentals of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, Fifth edition 2001, pp. 349.
4. Westgard J.O., and Barry P.L., Cost Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, DC, 1986.
5. Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, DC, 5th ed. 2000.

6. Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

IL Test™ and SeraChem®, are trademarks of Instrumentation Laboratory Company.

IL Test™ Total Protein

Cat. No. 0018258640 - 5 x 50 mL Total Protein

Cat. No. 0018251440 - 5 x 100 mL Total Protein

ILab 600/650 Test Parameters

Photometric Test Parameters	Serum
Test No.	13
Test Name, Test Code	T Protein, TP
Sample Type	Serum
Reporting Unit, Decimal Points	g/dL (g/L), 1 (0)
Reaction Cycle	Standard
Twin Analysis	OFF
Methodology Type, Measuring Point	End Pont, 11
Photometric Methodology	2 Wavelength
Primary/Secondary Wavelength	546/700

Sampling Conditions

Sampling 1	Sample Vol.	4
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	3
	Sample/Diluent Vol.	0/0
Sampling 3 and Sampling 4		***
Diluent Code		N/A
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
HIGH!, ABS!		Sampling 2

Reagent Volumes

R1	Code	01131
	Rgt/Dil. Vol., Stirring	200/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	40
R2	Code	N/A
	Rgt/Dil. Vol., Stirring	N/A
	Low Vol. Warning Limit	N/A
	Stability (days)	N/A

Ranges and Evaluation Criteria	Serum
Normal Range-Male	6.4 - 8.3 (64 - 83)
Normal Range-Female	6.4 - 8.3 (64 - 83)
Normal Range-Other	**
Valid Range	1.0 - 13.8 (10 - 138)
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 3500
Prozone Limit	N/A
Non Linear Limit	N/A
Slope/Intercept Correction	1/0
Qualitative Report	OFF
Calibration Conditions	
Calibration	1 Point, Linear, ** Repts
Stability (days)	21
Calibrator, Concentration	Ref G, *
R-Blank Limit (mAbs)	3500
Cal. Repts Range (%)	100
Min. Cal. Response (mAbs)	0
Cal. Factor Change (%)	100
M-Point Curve Fit (%)	0
Reagent Blank	ON
Auto R-Blank by Bottle	ON

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Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.