IL Test™ Direct Bilirubin

Cat. No 0018258440 - 4 x 45 mL BILI D R1 / 8 x 19.5 mL BILI T/D R2 / 1 x 4.5 mL Sodium nitrite

Cat. No 0018254740 - 4 x 90 mL BILI D R1 / 4 x 74 mL BILI T/D R2 / 1 x 13 mL Sodium nitrite

ENGLISH - Insert revision 02/2013

Intended Use

IL Test™ Direct Bilirubin is intended for the quantitative in vitro diagnostic determination of Direct Bilirubin in human serum and plasma using the ILab Chemistry Systems.

Summary

Direct bilirubin measurements are used in the diagnosis and treatment of various liver diseases, and metabolic disorders.¹

Principle

End point analysis, based on the Diazo reaction:
 direct bilirubin + diazotized sulfanilic acid → azobilirubin

The concentration of azobilirubin is proportional to the concentration of Direct Bilirubin in the sample. Absorbance measurements are taken at 546/660 nm.

Composition

BILI D R1 containing: Hydrochloric acid 50 mmol/L.
 BILI T/D R2 containing: Hydrochloric acid 180 mmol/L and Sulfanilic acid 28.9 mmol/L.
 Sodium Nitrite containing: Sodium nitrite 72.5 mmol/L.

PRECAUTION AND WARNINGS

BILI D R1, BILI T/D R2 and Sodium Nitrite are not classified as dangerous.
 Hazard class: none.
 Risk phrases: none.
 Safety phrases: none.

Reagent Preparation

Cat. No 0018258440
 BILI D R1: ready to use. Uncap the container and place into reagent tray
 BILI T/D R2: add 10 drops (0.5 mL) of Sodium Nitrite into the BILI T/D R2, bottle, cap and mix by inverting a few times, then uncap and place into reagent tray.

Cat. No 0018254740
 BILI D R1: ready to use. Uncap the container and place into reagent tray
 BILI T/D R2: add 35 drops (1.75 mL) of Sodium Nitrite into the BILI T/D R2, bottle, cap and mix by inverting a few times, then uncap and place into reagent tray. This reagent may be used immediately after mixing.

Note: the BILI T/D R2 reagent is also used in the Total Bilirubin assay.

Reagent storage and Stability

Unopened reagents are stable until expiration date indicated on the bottle when stored at 15-25°C.
 BILI D R1: on board until expired.
 BILI T/D R2: 30 days on board.

Instrument/Test Procedures

For the assay procedure refer to the relevant instrument Application Manual and to the relevant Instrument Operator Manual.

Specimen collection and Preparation

Serum, lithium heparinised plasma.
 Samples are stable for 3 days at 2-8°C.^{2,3} Protect serum and plasma from exposure to direct sunlight or fluorescent light.

Calibration

Use ReferrIL™ G Calibrator (Cat. No. 0018257000). See calibrator insert sheet for lot specific concentration and calibrator traceability.

Recalibrate the instrument every 30 days or when a new lot of reagent is used. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Quality Control

Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice. Normal and abnormal controls such as SeraChem® Control Level 1 (Cat. No. 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program.

Refer to Westgard *et al.*⁴ for identification and resolution of out-of-control situations.

Calculation of Analytical Results

Refer to the Application Sheet and to the ILab Operator Manual.
 To convert mg/dL to SI units (µmol/L) multiply by 17.1.¹

Reference Interval

0.00 to 0.30 mg/dL (0.0 to 5.1 µmol/L).⁷

It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Performance Characteristics

All performance data were obtained on ILab 600/650. For performance data on other ILab systems refer to the relevant application Manual.

Limitation/Interfering Substances

No interference from lipemia up to sample absorbance of 4.7/cm at 660 nm (1000 mg/dL or 11 mmol/L triglycerides). Do not use hemolyzed samples.
 For a comprehensive review of interfering substances, refer to the publication by Young.⁵

Precision

Typical precision results obtained according to NCCLS protocol.⁶

	Samples/Runs	Mean (mg/dL)	SD (mg/dL)	CV (%)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Within Run	5/10	0.44	0.00	0.9	1.33	0.01	0.5
Total	5/10	0.44	0.01	1.6	1.33	0.02	1.5

Method Comparison

Comparison Method (x)	Same reagent monochromatic	Doumas method
Comparison Instrument	ILab 600/650	Manual
Slope	1.000	1.040
y intercept	0.001	-0.060
Range (mg/dL)	0.12 - 7.80	0.09 – 7.45
Mean x (mg/dL)	2.07	2.05
Mean y (mg/dL)	2.07	2.07
r	0.999	0.999
SE	0.04	0.14
n	28	28

Linearity

	Linearity range (mg/dL)	Linearity range (µmol/L)
no rerun	0.0 -11.0	0.0 - 188
with rerun	0.0 - 23.0	0.0 - 393

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ILab 600/650 Test Parameters

Photometric Test Parameters	Serum
Test No.	5
Test Name, Test Code	D. Bilirubin, Bili D
Sample Type	Serum
Reporting Unit, Decimal Points	mg/dL (µmol/L), 2 (1)
Reaction Cycle	Standard
Twin Analysis	OFF
Methodology Type, Measuring Point	End Pont, 16/33
Photometric Methodology	2 Wavelength
Primary/Secondary Wavelength	546/660
Sampling Conditions	
<i>Sampling 1</i>	Sample Vol. 18
	Sample/Diluent Vol. 0/0
<i>Sampling 2</i>	Sample Vol. 9
	Sample/Diluent Vol. 0/0
<i>Sampling 3 and Sampling 4</i>	***
Diluent Code	N/A
Diluent Warning Limit	N/A
First Run	Sampling 1
Below/Above Normal Range	***
Panic L	***
Panic H	Sampling 2
Noise	***
Prozone	N/A
HIGH!, ABS!	Sampling 2
Reagent Volumes	
R1	Code 01051
	Rgt/Dil. Vol., Stirring 115/65, ON
	Low Vol. Warning Limit ***
	Stability (days) 0 (until expiration)
R2	Code 04052
	Rgt/Dil. Vol., Stirring 95/0, ON
	Low Vol. Warning Limit ***
	Stability (days) 30

Ranges and Evaluation Criteria

	Serum
Normal Range-Male	0.0 - 0.30 (0.0 - 5.1)
Normal Range-Female	0.0 - 0.30 (0.0 - 5.1)
Normal Range-Other	**
Valid Range	0.0 - 11.0 (0.0 - 188)
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 3500
Prozone Limit	N/A
Non Linear Limit	N/A
Slope/Intercept Correction	1/0
Qualitative Report	OFF
Calibration Conditions	
Calibration	1 Point, Linear, ** Repts
Stability (days)	30
Calibrator, Concentration	Ref G, *
R-Blank Limit (mAbs)	***
Cal. Repts Range (%)	***
Min. Cal. Response (mAbs)	***
Cal. Factor Change (%)	***
M-Point Curve Fit (%)	N/A
Reagent Blank	ON
Auto R-Blank by Bottle	ON

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- Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, 3rd edition 1999, pp. 1130, 1803.
- Henry, RJ, et al., Clinical Chemistry: Principles and Techniques, 2nd edition, Hagerstown, MD; Harper and Row, 1059, 1974.
- Tietz, NW, ed., Clinical Guide to Laboratory Tests, 3rd Edition, WB Saunders Company, Philadelphia, PA, 1995, p. 88.
- Westgard J.O., and Barry P.L., Cost Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, D.C., 1986.
- Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, DC, 5th ed. 2000.
- Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.
- Solberg, H., "IFCC Expert Panel on Theory of Reference Values, Part 5, Statistical Treatment of Collected Reference Values, Determination of Reference Limits", J.Clin.Chem.Clin.Biochem., 21:749-760, 1983

IL Test™ and SeraChem®, are trademarks of Instrumentation Laboratory Company.

In Vitro Diagnostic Medical Device
Dispositivo medico-diagnostico in vitro
Producto sanitario para diagnóstico in vitro

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ATTENZIONE! Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.