INSTRUCTIONAL INFORMATION

INTENDED USE

This product is an assayed whole blood control for evaluating the accuracy and precision of automated procedures that measure blood cell parameters.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation protocols.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument or manual methods, this whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

These whole blood reagents may contain any or all of the following: stabilized human red blood cells, a nucleated red blood cell analog, a white blood cell component consisting of human or non-human analogs and a platelet component consisting of a human or non-human analog in a preservative medium.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
- This product should not be disposed of in general waste, but should be disposed of with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
- Control products are not to be used as calibrators.

STORAGE AND STABILITY

STaK-Chex Plus Retics is stable through the expiration date when stored at 2° to 10°C before and after opening. After opening, STaK-Chex Plus Retics is stable throughout the open-vial expiration dating, as indicated on the assay sheet, when stored at 2° to 10°C.

Trending in the MCV parameter over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

- Review control product package insert and the operating procedure of the instrument.
- Check expiration date of the product on the vial. Discard outdated products.
- Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Technical Services at 800-843-0912 or online at www.streck.com.

INSTRUCTIONS FOR USE

- Remove vials of STaK-Chex Plus Retics from the refrigerator and warm to room temperature (18° to 30° C) for 15 minutes before use.
- To Mix: (Do not mix mechanically.)
 - Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
 - Mix by rapid inversion to ensure the cells are suspended.
 - Vials stored for an extended period of time may require extra mixing.
 - Gently invert the vials 8 to 10 times immediately before sampling.
- Refer to the instrument manual for the system in use for analyzing control material.
- For further assistance, contact Technical Services, Streck at 800-843-0912 or consult Technical Services online at www.streck.com.
- After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and returning to refrigeration.

QC FILE SETUP AND QC RUN PROCEDURES

MANUAL SETUP PROCEDURE Combined CBC/Diff/Reticulocyte file

- Beckman Coulter® LH 750/755, LH 780/785 and Beckman Coulter® GEN+S™**
- Select the yellow diamond icon on the Command Center to display the System Setup window.
 - Select the QA icon, select the Controls tab and click on the New Control folder icon.
 - Select the following: Source - Other; Type - CBC/Diff/Retic; Level - Low, Normal or High.
 - Select the Set Up New Lot icon to display the control setup screen.
 - Enter the 8-digit CBC/Diff/Retic lot number found in the lot number chart on the assay; this includes the leading zero. Enter the expiration date, assigned assay means and expected ranges from the assay for each level. Select the OK icon to save the data.
 - Repeat steps 3 - 5 for the remaining levels. All values will be transmitted into one control file for each level.

RUNNING CONTROLS

Beckman Coulter LH 750/755 and LH 780/785

- Note: NRBC parameters can be recovered by turning auto-print on for the controls and recording the values manually for the LH 700 series.
 - Select Run Configuration on the Command Center
 - Select Control in the QA Samples
- In the Command Center located at the bottom of the screen, select the following:
 - Process Type – Auto Analysis
 - Default Type – CDR (CBC/Diff/Retic)
 - Random Access – Select the box
- Position the vials in the cassette holder so the instrument can scan the **CBC/Diff/Retic** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Beckman Coulter GEN+S

- On the System Run Screen, select CBC/Diff/Retic as the Test Mode.
- In the Command Center, select Auto Analysis as the Process Type.
- Position the vials in the cassette holder so the instrument can scan the **CBC/Diff/Retic** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

ASSAY DISK SETUP PROCEDURE Separate CBC/Diff and Retic files

Beckman Coulter LH 750/755, LH 780/785, GEN+S, LH 500* and HmX†
If your instrument has the ability to upload assay values from a disk, please refer to your instrument operator's manual for assay disk setup.

- * LH 500 retic assay values must be entered manually
- † HmX retic assay values may be entered using the assay barcode setup procedure or entered manually.

RUNNING CONTROLS

Beckman Coulter LH 750/755, LH 780/785 and LH 500

- CBC/Diff file**
- Note: NRBC parameters can be recovered by turning auto-print on for the controls and recording the values manually for the LH 700 series.
 - Select Run Configuration on the Command Center
 - Select Control in the QA Samples
 - In the Command Center, select the following:
 - Process Type – Auto Analysis
 - Default Type – CD (CBC/Diff)
 - Random Access – Select the box
 - Position the vials in the cassette holder so the instrument can scan the CBC/Diff barcode.
 - Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Reticulocyte file

- In the Command Center, select the following:
 - Process Type – Auto Analysis
 - Default Type – R (Reticulocyte)
 - Random Access – Select the box
- Position the vials in the cassette holder so the instrument can scan the **Retic** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Beckman Coulter GEN+S

CBC/Diff file

- On the System Run Screen, select CBC/Diff as the Test Mode.
- In the Command Center, select Auto Analysis as the Process Type.
- Position the vials in the cassette holder so the instrument can scan the **CBC/Diff** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Reticulocyte file

- On the System Run Screen, select Retic Only as the Test Mode.
- In the Command Center, select Auto Analysis as the Process Type.
- Position the vials in the cassette holder so the instrument can scan the **Retic** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

ASSAY BARCODE SETUP PROCEDURE Separate CBC/Diff and Retic files

Beckman Coulter STKS™/MAXM™/HmX

Expected Range file

- From the Control Setup screen, select Expected Range and the appropriate file type and level.
- Enter the Expected Ranges for each parameter as listed on the assay sheet.

CBC/Diff Files

Enter the Low Abnormal (87) ranges into the Abnormal I file.
Enter the Normal (88) ranges into the Normal file.
Enter the High Abnormal (86) ranges into the Abnormal II file.

Retic Files

Enter the Low Abnormal (31) ranges into the Retic Level I file.
Enter the Normal (42) ranges into the Retic Level II file.
Enter the High Abnormal (53) ranges into the Retic Level III file.

- Press F10/ESC to save.

NOTE: The Streck Expected Ranges only need to be entered once as they rarely change from lot to lot.

CBC/Diff file

- From the Control Setup screen, select CBC/Diff Files.
- Enter the 6-digit lot number and expiration date using the barcode wand. The expected ranges for each level will be transmitted into QC files after the lot number is entered.
- Enter other applicable information for your laboratory.
- Enter the assigned values from the assay sheet using the barcode wand.
- Press F10 Save/ESC to save.

Reticulocyte file

- From the Control Setup screen, select Retic Files.
- Enter the 6-digit Retic lot number located near the bottom of the assay sheet. The expected ranges for each level will be transmitted into QC files after the lot number is entered.
- Enter the expiration date and other applicable information for your laboratory.
- Enter the Ret# assigned value from the assay sheet. If you wish to have the Ret# calculated, you must enter the RBC value from the assay into the Ref RBC field.
- Press F10 Save/ESC to save.

RUNNING CONTROLS

Beckman Coulter STKS/MAXM/HmX

CBC/Diff file

- Position the vials in the cassette holder so the instrument can scan the **CBC/Diff** barcode.
- Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Reticulocyte file

- Prepare a vial and change the instrument mode according to the procedure for analyzing reticulocytes.
- On the numeric keypad, press ID, then enter the 6-digit Retic lot number located near the bottom of the assay sheet. This lot number can also be found on the vial label as LOT Retic.
- Aspirate the prepared sample.

LIMITATIONS (controls assayed for WBC differential)

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. If the product is value assigned for MCV and hematocrit, assay values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Values and ranges for methods not listed on the assay and instructional information sheet must be established by the user. It is recommended that at least 10 consecutive analyses be performed on a properly calibrated instrument for each level to establish the mean.

'R' alarms may occur with control materials on automated systems. This will not affect the validity of the results. The alarms may be disregarded if the control is performing within the assay ranges.

QUALITY CONTROL PROGRAM

Streck offers STaTS®, an interlaboratory quality control program, and STaTS-Link®, which provides internet access to STaTS reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STaTS Department at 800-898-9563, or by fax at 402-333-7874.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

| | | | | |
|--|-------------------|------------------------------|------------------------|--------|
| ECREP Authorized Representative in the European Community | LOT Batch Code | Biological Risk | REF Catalog Number | Use By |
| IVD In Vitro Diagnostic Medical Device | Manufacturer | Consult Instructions For Use | Temperature Limitation | |

Glossary of symbols may contain symbols not used in the labeling of this product.

U.S. Patents 5,262,327; 5,270,208; 5,432,089; 5,672,474; 5,677,145; 5,731,205; 5,981,282; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

The brand and product names of the instruments are trademarks of their respective holders.

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