INSTRUCTIONAL INFORMATION

INTENDED USE

STaK-Chex is an assayed whole blood control for evaluating the accuracy and precision of automated and semi-automated procedures that measure blood cell parameters.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated and semi-automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, the whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

This whole blood reagent may contain any or all of the following: stabilized human red blood cells, a white blood cell component consisting of human and/or non-human analogs and a platelet component consisting of a human or non-human analog in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. The product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. The product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.

STORAGE AND STABILITY

STaK-Chex is stable through the expiration date when stored at 2° to 10°C before and after opening. After opening, STaK-Chex is stable throughout the open-vial expiration dating, as indicated on the assay sheet, when stored at 2° to 10°C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Technical Services at 800-843-0912 or online at www.streck.com.

QC FILE SETUP PROCEDURE

Coulter® GEN*SM, LH 500 and HmX

ASSAY DISK SETUP

If your instrument has the ability to upload assay values from a disk, please refer to your instrument operator's manual for assay disk setup.

Coulter STKS™/MAXM™/HmX

Expected Range file

NOTE: The Streck Expected Ranges only need to be entered once as they rarely change from lot to lot.

1. From the Control Setup screen, select Expected Range and the appropriate file type and level.
2. Enter the Expected Ranges for each parameter as listed on the assay sheet.

CBC/Diff Files

- Enter the Low Abnormal (87) ranges into the Abnormal I file.
- Enter the Normal (88) ranges into the Normal file.
- Enter the High Abnormal (86) ranges into the Abnormal II file.

3. Press F10/ESC to save.

CBC/Diff file

NOTE: The CBC/Diff File is set up with each new lot.

1. From the Control Setup screen, select CBC/Diff Files.
2. Enter the 6-digit lot number and expiration date using the barcode wand. The expected ranges for each level will be transmitted into QC files after the lot number is entered.
3. Enter other applicable information for your laboratory.
4. Enter the assigned values from the assay sheet using the barcode wand.
5. Press F10 Save/Esc to save.

INSTRUCTIONS FOR USE

1. Remove vials of control from the refrigerator and warm to room temperature (18° to 30°C) for 15 minutes before use.
2. To mix: (Do not mix mechanically)
 - a. Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
 - b. Mix by rapid inversion to ensure the cells are suspended.
 - c. Vials stored for an extended period of time may require extra mixing.
 - d. Gently invert the vials 8 to 10 times immediately before sampling.
3. Refer to the instrument manual for the system in use for analyzing control materials.
4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. If the product is value assigned for MCV and hematocrit, assay values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual lab establishes its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Values and ranges for methods not listed on the assay must be established by the user. It is recommended that at least 10 analyses be performed on a properly calibrated instrument for each level to establish the "assay" mean.

QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, and *STATS-Link*®, which provides internet access to *STATS* reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the *STATS* Department at 800-898-9563, or by fax at 402-333-7874.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

ECIREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

The brand and product names of the instruments are trademarks of their respective holders.
U.S. Patents 5,262,327; 5,270,208; 5,672,474; 5,677,145; 5,731,205; 5,981,282

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