Hollow Fiber Membrane Blood Oxygenator with Hardshell Venous Reservoir

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GB – ENGLISH - INSTRUCTIONS FOR USE

I. CONTENTS

- I. Contents
- A. Description
- B. Technical features
- C. Intended use
- D. Safety information
- E. Set up
- F. Priming and recirculation procedure
- G. Initiating bypass
- H. Operation during bypass
- I. Terminating bypass
- J. Blood recovery after bypass
- K. Use of active venous drainage with vacuum
- L. Use of cardiotomy for post-operative autotransfusion
- M. Oxygenator replacement
- N. Medical devices for use with the PRIMO₂X
- O. Return of used products
- P. Limited Warranty

A. DESCRIPTION

The PRIMO₂X is a microporous hollow fibre membrane oxygenator consisting of a gas exchange module with an integrated heat exchanger and a Hardshell Reservoir.

The PRIMO₂X is coated with PC coating. Devices coated with PC are used when a coated blood path is desired. The PC coating improves the blood compatibility of the device by reducing platelet adhesion on the coated surfaces.

The device is single use, non-toxic, non pyrogenic, supplied STERILE and packaged singly. Sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

B. TECHNICAL FEATURES

- Recommended max blood flow	8000 ml/min
- Membrane type	Microporous Polypropylene
- Membrane surface area	1.87 m ²
- Heat exchanger surface area	0.14 m ²
- Hardshell Reservoir Capacity	4300 ml
- Static priming volume (oxygenating module + heat exchanger)	250 ml
- Connections:	
Venous reservoir return	1/2" (12.7 mm)
Venous reservoir outlet	3/8" (9.53 mm)
Oxygenator venous inlet	3/8"
Oxygenating module arterial outlet	3/8"

C. INTENDED USE

The PRIMO₂X is intended for use in adult surgical procedures requiring extracorporeal gas exchange support and blood temperature control for periods of up to 6 hours.

The PRIMO₂X should be used in combination with medical devices listed in section M (Medical devices for use with the PRIMO₂X).

D. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

⚠ WARNING

WARNING indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

⚠ CAUTION

CAUTION indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

 For single use only (Do not reuse)

 Batch code (number)
(reference for product traceability)

 Use by (Expiry date)

 Manufactured by

 Date of manufacture

 Sterile - Ethylene oxide sterilised

 Non Pyrogenic

 Contains PHTHALATE

 Latex free

 Warning: Do not resterilize.

 Contents sterile only if package is not opened, damaged or broken

 Catalogue (code) number

 Attention, see instruction for use

 This way up

 Fragile; handle with care

 Quantity

 Keep away from heat

 Keep dry

The following is general safety information with the aim of advising the operator in preparing to use the device.

Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

⚠ WARNING

- The User should carefully check the device during set-up and priming for leaks. Do not use if any leak is detected.
- The device must be used in accordance with the instructions for use provided in this manual.
- For use by professionally trained personnel only
- Sorin Group Italia is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Do not expose to temperatures below 0°C (32°F) or above 60°C (140°F).
- Keep dry. Store at room temperature.
- Always give and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.
- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device contains phthalates. Considering the nature of contact with the body, the limited duration of contact and the number of treatments per patient, the amount of phthalates which might be released from the device does not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia
- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if sterile. - For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorised local representative.

⚠ CAUTION

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Inner surfaces of the system are PC coated, currently Sorin Group Italia is not aware of any contraindications to the use of this coated device

E. SET UP

1) POSITION THE HOLDER

Position the PRIMO₂X holder on the pump structure by means of the clamp at

the upper end of the arm (fig.1, ref. 1).

2) FIX THE OXYGENATOR TO THE HOLDER

⚠ WARNING

- Do not use if sterile packaging is damaged, unsealed, or has been exposed to moisture or other conditions that would compromise the sterility of the device.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.

Remove the device from the sterile packaging.

⚠ WARNING

- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc.: as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage, which may compromise the integrity and proper functioning of the device.

Attach the device to the holder.

The word "OPEN" on the water connectors locking system must be visible.

Check that the notches on the water connectors locking system are aligned.

Only then you can attach the oxygenator onto the holder.

Insert the Hansen connectors and push the device down to the holder and turn the locking lever on "CLOSED" position.

The device will be correctly positioned only when the locking lever shows "CLOSED".

3) HEATER-COOLER SET UP

Connect the water tubes to the holder by means of the female Hansen connectors (SORIN GROUP ITALIA code 09028 or equivalent).

⚠ WARNING

- The use of different connectors from those indicated may cause resistance inside the circuit and reduce the efficiency of the heat exchanger.
- Do not obstruct the hole on the heat exchanger lower cover as it is the outlet of the safety channel which helps preventing fluids crossing from one compartment to another.
- The water temperature at the heat exchanger inlet must not exceed 42°C (108 °F).
- The special positive and negative pressure relief valve is mounted on the lid of the venous reservoir. It starts to release positive pressures at 5 mmHg (0.7 kPa / 0.007 bar / 0.1 psi) and negative pressures at -80 mmHg (-10.4 kPa / -0.10 bar / -1.53 psi). The action of the valve by itself does not prevent from the over / under pressurization of the reservoir in all operating conditions.
- Do not for any reason occlude the external access hole of pressure relief valve. Remove tab inserted in pressure relief valve (Fig.2, ref.20) from the valve prior to use.
- The water pressure in the heat exchanger must not exceed 2250 mmHg (300 Kpa / 3 bar / 44 psi).

4) CHECK THE HEAT EXCHANGER

⚠ WARNING

Check the heat exchanger by circulating water inside the heat exchanger for a few minutes. There should be no leaks from the water compartment or from the safety channel hole.

5) CIRCUIT CONNECTIONS

⚠ WARNING

All connections downstream of the pump must be secured by means of ties.

VENOUS LINE: Connect a venous line of 1/2" to the connector indicated on the venous reservoir as "VENOUS RETURN" (fig. 2, ref.2).

⚠ CAUTION

If blood saturation and hematocrit monitoring is needed by using the SAT/HCT monitor, connect the 1/2" tubing of the "SAT/HCT line for PRIMO₂X" (provided with the product) to the connector indicated on the venous reservoir as "VENOUS RETURN" (Fig. 2 - ref. 2), and the 1/2" barbed tubing connector to the venous line. When making the tubing connections push the tubing at least 1/4" (7mm) past the apex of the

innermost barb of the connector and tie wrap the two connections. Then connect the SAT/HCT Monitor to the SAT/HCT fitting on the connector of SAT/HCT line for PRIMO₂X. For complete instructions, refer to the Operator's Manual for the monitor. The two luer ports of SAT/HCT line for PRIMO₂X may be used for the placement of a thermistor probe or for blood access.

When using the SAT/HCT line connect the SAT/HCT monitor to the SAT/HCT fitting on the connector.

The Venous Return connector can be rotated 360° to find the most convenient position of the venous tubing.

CARDIOTOMY SUCTION LINES: After removing the protective caps from the "filtered" inlets on the top of the venous reservoir (four 1/4" / 6.35 mm inlets, two 3/8" / 9.53 mm inlets), connect the ends of the suction tubes and rotate the turret (fig. 2, ref. 8) orienting the filtered inlets towards the suction pumps.

ARTERIAL LINE: remove the red cap on the oxygenator arterial outlet indicated as "ARTERIAL OUTLET" (fig. 2, ref. 12) and connect a 3/8" line.

PUMP LINE: the pump segment should be set up between the Venous Reservoir outlet connector (fig. 2, ref. 5) and the oxygenator venous inlet connector (fig. 2, ref. 11) taking account of the direction of rotation of the pump.

RESERVOIR VENT Remove yellow cap from the reservoir vent connector (fig. 2, ref. 7).

⚠ CAUTION

Check that the pos lock on the coronary outlet port is closed

If oxygenated blood is necessary for blood cardioplegia, remove the red pos lock and connect the 1/4" blood line of the cardioplegia circuit to the PRIMO₂X coronary outlet port using the D523C reducer (provided with the product).

The coronary outlet port has a self-sealing valve which allows connection of the D 523C reducer during extracorporeal bypass without any leakage or spillage of fluid.

⚠ WARNING

If a connection is made to the coronary outlet port during bypass, the line to be connected must be unclamped and unpressurized such that blood will flow into it upon connection.

Connect the recirculation line to one of the rear connectors of the reservoir.

6) SAMPLING SYSTEM

VERSION WITH PRE-CONNECTED SAMPLING SYSTEM

If it is preferred to take samples from a more distant position, remove the A/V sampling manifold (fig. 2, ref. 23) from its housing and attach it on the dedicated holder. The tubing lines of the sampling system allow the manifold positioning in a range of about 1 meter. Position the stopcock handles toward the access ports of the manifold (A/V shunt position).

VERSION WITH SAMPLING SYSTEM TO BE CONNECTED BY THE USER

The arterial/venous sampling system is supplied in a sterility bag. About 1 meter venous and arterial tubing lines are attached to a 3-stopcock manifold. Remove the device from the packaging.

ARTERIAL SAMPLING LINE: Remove the protective cap from the arterial sample site luer connector located on the side of the arterial outlet (fig. 2, ref. 16). Connect the red male luer of the arterial sampling line to the arterial sampling site luer.

⚠ WARNING

The arterial sampling site connector luer does not include one way valve. Verify that there is a one way valve in the sampling line to prevent accidental air introduction into the arterial line.

VENOUS SAMPLING LINE: Remove the protective cap from the venous sample site luer connector located on the venous return connector (fig.2, ref. 6). Connect the blue male luer of the venous sampling line to the venous sampling site luer.

⚠ CAUTION

Check the security of all luer connections. All accessory lines connected to the device must be connected tightly in order to prevent accidental introduction of air into the device or blood loss.

Attach the A/V sampling manifold on the dedicated holder or on the specific housing located to the top of the venous reservoir. Position the stopcock handles toward the access ports of the manifold (A/V shunt position).

7) CONNECT THE TEMPERATURE PROBES

The connection for the arterial temperature probe (fig. 2, ref. 19) is located next to the arterial outlet, while the venous probe site (fig. 2, ref. 3) is located on the venous inlet connector of the reservoir.

Use SORIN GROUP ITALIA temperature probes (code 09026) or equivalent temperature probes (ref TMPMV).

8) CLOSE THE PURGING/RECIRCULATION LINE

Close the purging/recirculation stopcock (refer to the diagram on the label) prior to priming.

9) CONNECT THE GAS LINE

Remove the green cap from the gas inlet connector indicated as "GAS INLET" (fig. 2, ref. 13) and connect the 1/4" gas line.

Ensure that the gas supply is from a suitable air/oxygen mixer such as the Sechrist, code 09046 (available from SORIN GROUP ITALIA) or a system with compatible technical features.

A capnograph connector can be found in the center of "GAS ESCAPE" (fig. 2, ref. 14) connector.

⚠ CAUTION

- The "GAS ESCAPE" system is designed to avoid any possible risk of blocking the gas outlet; such blockage could cause the immediate passage of air to the blood pathway.
- Sorin Group Italia recommends the use of a bubble trap or arterial filter on the arterial line to reduce the risk of emboli transmission to the patient.
- The user should check for occlusions in tubing during set up.

F. PRIMING AND RECIRCULATION PROCEDURE

⚠ WARNING

Do not use priming solutions containing alcohol: such solutions could compromise the proper function of the oxygenating module.

1) KEEP THE GAS FLOW OFF

2) KEEP THE OXYGENATOR PURGING/ RECIRCULATION LINE CLOSED

Check that the purging/recirculation stopcock is closed.

3) CLOSE VENOUS AND ARTERIAL LINES

Clamp the venous line. Clamp the arterial line some centimetres away from arterial outlet connector of the oxygenator.

4) CHECK THE HEAT EXCHANGER

Check for leakage from the heat exchanger, paying particular attention to possible water leaks.

5) VENOUS RESERVOIR PRIMING

Secure with ties all suction lines connected to the Venous Reservoir. Fill the Venous Reservoir with sufficient prime to ensure the intended haematocrit is obtained, taking into account:

- the static priming volume of the oxygenator is 250 ml;
- the 3/8" tube capacity is 72 ml/m;
- the 1/2" tube capacity is 127 ml/m.

Clamp the Venous Reservoir outlet.

6) CIRCUIT PRIMING

⚠ WARNING

- The pressure level inside the blood compartment of the oxygenating module shall not exceed 725 mm Hg (100 KPa / 1 bar / 14 psi).

Remove the clamp from the venous line and start with the pump at a reasonably high flow so that the pump loop tubing is primed. Reduce the flow rate to 500-1000 cc/min and prime the oxygenator.

7) OPEN VENOUS AND ARTERIAL LINES

Remove the clamp from the venous and arterial lines and increase flow up to 1500 ml/min.

8) OPEN THE PURGING/RECIRCULATION LINE

Once the steps up to and including point 7 have been carried out, increase the arterial pump speed until the flow reaches the maximum value of 8 l/min. Open the purging/recirculation stopcock briefly until all air is removed from the line in order to prime the oxygenator purging/recirculation line.

9) PURGE THE AIR CONTAINED IN THE CIRCUIT

During this phase it is necessary to tap the entire circuit tubing in order to facilitate the removal of microbubbles from the tube walls. After some minutes in which the flow is maintained at a high rate, all air will be removed.

10) PRIME OF THE SAMPLING SYSTEM

Priming of the A/V sampling system occurs automatically when the arterial, venous and central stopcock handles are positioned in such a way as to allow the prime to flow from the arterial outlet to the venous reservoir.

11) CLOSE THE PURGING/ RECIRCULATION LINE

After 3-5 minutes in which the flow is maintained at a high rate, all air will be removed and it is possible to close the purging/recirculation line stopcock.

12) CLAMP THE VENOUS AND ARTERIAL LINES

Clamp the venous and arterial lines.

⚠ CAUTION

- Do not use pulsatile flow during priming. Sudden changes in flow rate during priming can pull air across the membrane in to the blood pathway.

- SORIN GROUP ITALIA recommends the use of the pump speed control to reduce or stop the arterial flow slowly.
- Do not use the pump on/off switch until the pump speed is zero.
- Do not turn off the heater-cooler.
- Check the correct dosage of anticoagulant in the system before starting the bypass.

-The user should carry out a gross inspection for air removal.

⚠ WARNING

- If the reduction connector (D523C) and a circuit have been connected to the coronary outlet port, check the priming of the connected line.
- Do not create a negative pressure at the coronary outlet. Negative pressure in the blood compartment could cause microbubble formation.

G. INITIATING BYPASS

1) OPEN THE ARTERIAL AND VENOUS LINES

Remove first the clamp from the arterial line, then remove the clamp on the venous line. Start bypass with a blood flow rate appropriate to patient size. Monitor the blood level in the Venous Reservoir.

2) CHECK THE CORRECT OPERATION OF THE HEAT EXCHANGER

Check the temperature of the venous and arterial blood.

3) SELECTION OF THE APPROPRIATE GAS FLOW

The suggested gas/blood flow ratio in normothermia is 1:1 with a FiO₂ of 80 to 100%.

⚠ WARNING

- Always open the gas flow after the blood flow. The gas/blood flow ratio must never exceed 2:1.
- The pressure in the blood compartment must always exceed that of the gas compartment to prevent gas emboli appearing in the blood compartment.

4) BLOOD GAS MONITORING

After a few minutes of bypass operation, measure the gas content of the blood. Depending on the values found, adjust the relevant parameters as follows:

High pO₂ → Decrease FiO₂

Low pO₂ → Increase FiO₂

High pCO₂ → Increase gas flow

Low pCO₂ → Decrease gas flow

H. DURING BYPASS

1) CHECK THE VENOUS RETURN

If a higher venous return flow is necessary lower both the oxygenator and the venous reservoir with respect to the patient position.

⚠ WARNING

- The ACT (Activated Coagulation Time) must always be longer than or equal to 480 seconds in order to ensure adequate anticoagulation of extracorporeal circuit.
- If administration of anticoagulant to the patient is necessary, use the luer connector of the central stopcock on the sampling manifold.
- The minimum operating volume in the venous reservoir is 300 ml. However to ensure adequate response time in case of venous inflow obstruction, it is recommended that an adequate volume in addition to the 300 ml minimum level be maintained. Do not exceed 4300 ml in the venous reservoir.

2) ARTERIAL SAMPLING

Insert a sample syringe into the arterial sample stopcock luer. Position all stopcocks on the sample manifold to allow arterial blood to flow through the manifold. The pressure on the arterial side will allow flow. Draw the sample of blood from the arterial sample stopcock. Turn off the arterial stopcock before removing the syringe.

3) VENOUS SAMPLING

Ensure that the arterial stopcock is closed. Insert a sample syringe into the venous stopcock luer and a flush syringe into the center stopcock. Open the center stopcock and draw at least 10-15 ml of blood prior to taking the venous sample. Close the center stopcock. Return this blood through one of the filtered luer connectors positioned on the top of the reservoir. Open the venous stopcock and draw a sample of venous blood and close the stopcock before removing the syringe.

4) DRUGS INJECTION

Insert the medication syringe into the luer connector of the central stopcock. Open the central and venous stopcocks and inject into the drug into the manifold and venous sample line.

Close the central stopcock to the medication syringe and allow an arterial-venous "wash" through the stopcock manifold. Turn the stopcocks to the closed position when the drug has been delivered to the venous line.

CAUTION

Draw blood samples from the stopcocks only when the pump is running. If not, the pressure in the blood compartment would decrease causing the formation of air bubbles.

5) LOW FLOW RECIRCULATION

(Hypothermia associated with circulatory arrest).

- Reduce the gas flow to less than 500 ml/min.
- Open the recirculation line (lever of the purging/recirculation stopcock in the "RECIRCULATION" position) and clamp the Venous Reservoir inlet line (fig. 2, ref. 2).
- Reduce the flow from the arterial pump to 1500 ml/min.
- Clamp the oxygenator arterial line (fig. 2, ref. 12).
- Recirculate at a maximum flow of 1500 ml/min. throughout the patient's circulatory arrest.
- To initiate bypass after circulatory arrest, open the venous and arterial lines and slowly increase the blood flow.
- Close the recirculation line (lever in "CLOSE" position).
- Adjust gas flow.

6) CONTINUOUS AIR PURGE

The purging/recirculation stopcock in "PURGE" position features the continuous purge of blood and any incidental air from the gas exchange module during bypass.

In this condition and at full arterial blood flow, the continuous purge diverts only few mls/min of blood from the patient.

I. TERMINATING BYPASS

Bypass should only be terminated after consideration of each individual patient's condition as follows:

- Turn the gas flow off.
- Turn the heater-cooler off.
- Slowly decrease the arterial flow to zero while clamping the venous line.
- Clamp the arterial line.

CAUTION

- If extracorporeal circulation has to be subsequently reinitiated, a minimum blood flow inside the oxygenator must be maintained (maximum 1500 ml/min).
- Verify that the cardioplegia circuit connected to the coronary outlet port is clamped.

J. BLOOD RECOVERY AFTER BYPASS

- Drain as much blood as possible from the venous line into the venous reservoir, as soon as the surgeon has removed the venous cannulae from the patient.
- Slowly return the residual blood to the patient through the aortic cannulae as required by the patient's condition
- If desired, blood in the oxygenator may also be returned by adding clear prime to the venous reservoir when the blood in the reservoir as reached minimum volume. Pump the prime slowly through the oxygenator ensuring that the reservoir never is emptied.
- When the reservoir is nearly empty stop the arterial pump and clamp the arterial line.

K. USE OF ACTIVE VENOUS DRAINAGE WITH VACUUM

This method may be applied at any time during extra-corporeal circulation, provided that the steps below are followed. Using the kit code 096834, or equivalent supplied separately, and a vacuum regulation device, PRIMO₂X can be used for active venous drainage using vacuum. This technique constitutes an alternative to venous drainage by gravity and allows the use of shorter venous lines with reduced diameter and smaller-gauged cannulas.

- Open active venous drainage system kit. Maintain sterility of the system
- Connect the end with the blue cap to the vent connector of the venous reservoir, marked "VENT/VACUUM PORT" (fig. 2, ref. 7), and the end with the red cap to the vacuum line of vacuum regulating device.
- Close the clamp and the green cap on the line connected to the reservoir
- If considered necessary to interrupt or suspend this method, remove the yellow cap and open the clamp on the line.

WARNING

- Do not exceed -50 mmHg (-6,66 kPa / -0,07 bar / -0,97 psi) negative pressure in venous reservoir. Pressure greater than -50 mmHg (-6,66 kPa / -0,07 bar / -0,97 psi) increase the risk of hemolysis.
- Periodically check functioning of the vacuum regulating device and the degree of vacuum.

L. USE OF VENOUS RESERVOIR FOR POST-OPERATIVE AUTOTRANSFUSION

If desired the Venous Reservoir may be used for post-operative autotransfusion as

follows:

- Disconnect the purging/recirculation line.
- Separate the venous reservoir from the oxygenator by opening the latch and twist the reservoir clockwise while holding oxygenator stationary.
- Attach the venous reservoir on the "Post Operative Chest Drainage Holder" and use the device, by referring to its respective instructions for use:
 - D 540 AUTOTRANSFUSION CONVERSION KIT code 05053.

M. OXYGENATOR CHANGE-OUT

A spare oxygenator must always be available during bypass in the unlikely event that the oxygenator in use requires change-out. Procedures lasting longer than 6 hours or if particular situations where the safety of the patient may be compromised, (insufficient oxygenator performance, leaks, abnormal blood parameters etc.), could require change out. Follow the steps below to change out the oxygenator.

CAUTION

Use sterile technique during all replacement procedure.

- Turn the gas flow off.
- Double clamp, the venous line (5 centimetres apart).
- Turn the arterial pump off and double clamp the arterial line (5 centimeters apart) next to the arterial outlet port.
- Turn the heater-cooler off, clamp and remove the water lines.
- Disconnect the gas line, all monitoring and sampling lines.
- Cut the venous return and the arterial lines in the section between the two clamps, leaving a sufficient length of tubing to allow connection to the new oxygenator.
- Remove the oxygenator from the holder.
- Place a new oxygenator on the holder. Connect all lines (i.e. venous line to the Venous Reservoir inlet port, arterial and gas to the oxygenator, pump line to Venous Reservoir outlet and oxygenator inlet port).

WARNING

In this phase, keep the venous and arterial lines clamped.

- Open the water lines on the holder turn the heater-cooler on and check the integrity of the new heat exchanger.
- Fill the new venous reservoir with priming solution through the 1/4" or 3/8" quick prime connectors located on the lid of the venous reservoir.
- Prime the new oxygenator and remove the microbubbles, as described in the priming and circulation procedure.
- Verify all connections and secure with ties.
- Remove clamps from the venous and the arterial line, close the purging/recirculation line and start the bypass again.
Turn on the gas flow and adjust gas flow rate as required.
- The blood remaining in the replaced Venous Reservoir may be recovered by connecting its outlet port to one of the 3/8" inlet connectors of the new reservoir and drain this blood into the new reservoir.
- The blood contained in the oxygenator and heat exchanger may be poured into the new Venous Reservoir by connecting the arterial line to one of the 3/8" inlet connectors of the new reservoir.

N. MEDICAL DEVICES FOR USE WITH THE PRIMO₂X

CAUTION

The user should observe the warnings and cautions and follow instruction for Use accompanying the separate device.

- For the post-operative autotransfusion with venous reservoir use the device:
 - D 540 AUTOTRANSFUSION CONVERSION KIT code 05053.All tubing used to make the circuit connections must be of a diameter that is compatible with the dimensions of the connectors on the device (3/8", 1/4", 1/2").
 - Temperature monitoring must be carried out using SORIN GROUP ITALIA probes (code 09026) or equivalent temperature probes (ref TMPMV) compatible to YSI Series 400.
 - Use a Sechrist air/oxygen mixer (SORIN GROUP ITALIA code 09046) or a system with compatible technical features for control of the blood gas concentrations.
 - Any heating/cooling system may be used, provided that the water line connectors are Hansen 3ST (SORIN GROUP ITALIA code 09028) or equivalent.
- Currently SORIN GROUP ITALIA is not aware of any contraindications to the use of the device with occlusive or non-occlusive peristaltic pumps or with centrifugal pumps. The use of other types of pump must be agreed with SORIN GROUP ITALIA.

O. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from SORIN GROUP USA, Inc. prior to shipping .

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

The shipping address for returned goods in the US is:

SORIN GROUP USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031

P. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological haematological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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