
IL Test™ Alkaline Phosphatase
Cat. No. 0018259640 - 8 x 33 mL ALP R1 / 4 x 11 mL ALP R2
Cat. No. 0018259740 - 8 x 18 mL ALP R1 / 8 x 3 mL ALP R2

Intended Use
The IL Test™ Alkaline Phosphatase is intended for the quantitative in vitro diagnostic determination of alkaline phosphatase in human serum and plasma using the ILab Chemistry Systems.

Summary
Increases in alkaline phosphatase can be a result of liver disease, bone disease and some carcinomas.¹

Principle
Rate analysis. The method is based on a IFCC method performed at 37°C.²

Composition
ALP R1 containing: 2-amino-2-methylpropanol (AMP) 750 mmol/L, Zinc⁺⁺ 1 mmol/L, Magnesium⁺⁺ 2 mmol/L, HEDTA 2 mmol/L.
ALP R2 containing: p-nitrophenylphosphate 16 mmol/L.

PRECAUTION AND WARNINGS

ALP R1
 **Warning.** Causes skin irritation (H315). Causes serious eye irritation (H319). Wear protective gloves/protective clothing/eye protection/face protection (P280). IF ON SKIN: Wash with plenty of soap and water (P302+P352). If skin irritation occurs: Get medical advice/attention (P332+P313). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persist: Get medical advice/attention (P337+P313). Take off contaminated clothing and wash before reuse (P362).

ALP R2 is not classified as hazardous.
5.1% of the mixture consist of ingredient(s) of unknown toxicity for human health and aquatic environment.
This product is For *in vitro* Diagnostic Use.

Reagent Preparation
Cat.No 0018259640
Reagents are ready to use. Place into ALP R1 bottle the chimney provided with the kit.
Cat.No 0018259740
Reagents are ready to use. Place into ALP R1 bottle the chimney provided with the kit.

Reagent storage and Stability
Unopened reagents are stable until the expiration date indicated on the containers when stored at 2-8 °C.
0018259640 On Board: ALP R1 14 days with chimney; ALP R2 28 days.
0018259740 On Board: ALP R1 14 days with chimney; ALP R2 28 days.

Instrument/Test Procedures
For the assay procedure refer to the relevant instrument application sheet and to the instrument operator manual in use.

Specimen collection and Preparation
Serum, lithium heparinized plasma.

Freshly collected samples should be kept at room temperature and assayed as soon as possible but preferably not later than 4 hours after collection.¹

Calibration
Use ReferrIL™ E Calibrator (Cat No. 0018256300). See insert sheet for lot specific concentration and traceability.
Recalibrate every 7 days, when a new lot of reagent is used, or when adjustments are made to the instrument.

Quality Control
Use normal and abnormal controls in accordance with standard laboratory practices. Analyze at least once per day. Use at least 2 levels of IL controls such as SeraChem® Level 1 (Cat. No. 0018162412) and SeraChem® Level 2 (Cat. No. 0018162512). The target values and expected ranges can be found on their insert sheets. For identification and resolution of out-of-control situations, references such as Westgard *et al.*³ are recommended.

Calculation of Analytical Results
The reported ALP activity in a sample is calculated by the ILab as follows:

$$\text{ACTs} = \Delta\text{As} \times \text{K} \qquad \text{K} = \frac{\text{ACTc}}{\Delta\text{Ac}}$$

Where:
ACT =Alkaline Phosphatase activity in U/L
ΔA = Change in milliabsorbance per minute at 405 nm
s = sample
c = calibrator
For more details refer to the relevant Application Manual and to the ILab Operator Manual.

Activity is expressed in International Units/Liter (U/L), where International Unit is defined as: enzyme activity that converts 1 μmol of substrate in 1 minute.
To convert U/L to nKat/L multiply by 16.67.

Reference Interval
Male: 43 - 115 U/L (716 - 1917 nKat/L).²
Female: 33 - 98 U/L (550 - 1630 nKat/L).²
It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Performance Characteristics
All performances data were obtained on ILab 600/650, using calibrator ReferrIL E. For performance data on other ILab Systems refer to the relevant application Manual.

Limitation/Interfering Substances
No interference from lipemia up to sample absorbance of 3.8/cm at 660 nm (1000 mg/dL or 11.3 mmol/L triglycerides).
No interference up to 120 mg/dL (0.075 mmol/L) hemoglobin.
No interference up to 15 mg/dL (256 μmol/L) bilirubin.

For a comprehensive review of interfering substances, refer to the publication by Young.⁴

Precision
Typical precision results obtained according to NCCLS protocol.⁵

	Samples/Runs	Mean (IU/L)	CV (%)	Mean (IU/L)	CV (%)
Within Run	3/20	55	0.9	211	0.5
Total	3/20	55	1.7	211	1.2

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Bibliography / Bibliografia / Bibliografia

1. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, Third edition 1999, pp. 676-679, 1829.
2. G. Schumann, R. Klauke, F.Canalias, S. Bossert-Reuther, P. F.H. Franck, F.Javier Gella, P. J. Jørgensen, D. Kang, J.M. Lessinger, M. Panteghini and F.Cerioti. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 9:Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clin Chem Lab Med 2011;49(9):1439–1446
3. Westgard J., and Barry P., Cost- Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, D.C., 1986.
4. Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, DC, 5th ed. 2000.
5. Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

IL Test™ and SeraChem®, are trademarks of Instrumentation Laboratory Company.

IL Test™ Alakline Phosphatase
Cat. No. 0018259640 - 8 x 33 mL ALP R1 / 4 x 11 mL ALP R2
Cat. No. 0018259740 - 8 x 18 mL ALP R1 / 8 x 3 mL ALP R2

ILab 600 Test Parameters		
Photometric Test Parameters		Serum
Test No.		28
Test Name, Test Code		Alk. Phos., ALP
Sample Type		Serum
Reporting Unit, Decimal Points		U/L (nKat/L), 0 (0)
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		Rate, 20/28
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		405
Sampling Conditions		
Sampling 1	Sample Vol.	4
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	8
	Sample/Diluent Vol.	10/90
Sampling 3 and Sampling 4		
Diluent Code		Water
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
HIGH!, ABS!		Sampling 2
Sample Volume reduction		**
Reagent Volumes		
R1	Code	03281
	Rgt/Dil. Vol., Stirring	180/10, ON
	Low Vol. Warning Limit	**
	Stability (days)	14
R2	Code	03282
	Rgt/Dil. Vol., Stirring	30/10, ON
	Low Vol. Warning Limit	**
	Stability (days)	28

Ranges and Evaluation Criteria		Serum
Normal Range-Male		43 - 115 (717 - 1917)
Normal Range-Female		33 - 98 (550 - 1630)
Normal Range-Other		**
Valid Range		3 - 1000 (50 - 16670)
Hemolysis/Icterus/Lipemia Limit		***
Reaction Slope		Positive
Absorbance Limit		Above, 3500
Prozone Limit		N/A
Non Linear Limit		40
Slope/Intercept Correction		1/0
Qualitative Report		OFF
Calibration Conditions		
Calibration		1 Point Linear, 3 Reps
Stability (days)		7
Calibrator, Concentration		Ref E, *
R-Blank Limit (mAbs)		3500
Cal. Reps Range (%)		100
Min. Cal. Response (mAbs)		0
Cal. Factor Change (%)		100
M-Point Curve Fit (%)		0
Reagent Blank		ON
Auto R-Blank by Bottle		ON

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ATTENZIONE! Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.