Intended Use

For verification of GEM Premier 3500 and 3000 iQM™ cartridges prior to use with patient samples.

Summary and Principle

GEM CVP is comprised of the following four ampule configurations:

GEM CVP 1: with low pH, pO₂, Na⁺, K⁺, glucose, lactate and high pCO₂ and Ca⁺⁺ values.

GEM CVP 2: with high pH, pO₂, Na⁺, K⁺, glucose, lactate and low pCO₂ and Ca⁺⁺ values.

GEM CVP 3: with low hematocrit values.

GEM CVP 4: with normal hematocrit values.

Composition

GEM CVP 1 and 2 ampules: Aqueous buffered bicarbonate solution containing biocides, inorganic salts and organic metabolites, equilibrated with precise concentrations of carbon dioxide and oxygen. Each ampule contains 1.8 mL of solution.

GEM CVP 3 and 4 ampules: Aqueous buffered bicarbonate solution containing biocides, inorganic salts and equilibrated with carbon dioxide and oxygen. Each ampule contains 1.8 mL of solution.

Avoid contact with skin and eyes (S24/25). Do not empty into drains (S29).

Wear suitable protective clothing (S36).

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment (R52/53).

Avoid release to environment (S61).

For *in vitro* diagnostic use.

Storage and Stability

Unopened ampules are stable until the expiration date shown on the label when stored at 2-8°C, or up to 12 months at room temperature (20-25°C), providing storage does not exceed the expiration date. DO NOT FREEZE.

WARNING: GEM CVP solutions are sensitive to ambient temperature variations (gas/liquid equilibrium) and room air contamination (diffusion gradients). Ampules equilibrated at temperatures cooler than 22°C will demonstrate pO₂ values HIGHER than those stated on the insert by approximately one percent per °C. Likewise, ampules equilibrated at temperatures above 22°C will demonstrate pO₂ values LOWER than those stated by approximately one percent per °C. There is no significant temperature effect on pH and pCO₂ within the normal range of room temperature.

Directions for Use**1. EQUILIBRATE AMPULES AT 22°C FOR AT LEAST 8 HOURS PRIOR TO USE.****2. IMMEDIATELY PRIOR TO USE, HOLD THE AMPULE BY THE TOP ABOVE THE BREAK LINE AND SHAKE VIGOROUSLY FOR APPROXIMATELY 10 SECONDS.**

3. Restore liquid to bottom of vial with gentle tapping.

4. With fingers protected (e.g., tissues, gauze) carefully snap open the ampule utilizing the ampule breaker located on the GEM Premier 3500 and 3000 systems.

5. Aspirate samples immediately after opening. Any delay in measuring may cause room air contamination and result in lower pCO₂ and higher pO₂ values than stated on the insert.

6. Refer to the GEM Premier 3500 and 3000 Operator's Manual or Quick Operator's Manual for further information regarding product use.

Limitations

The values obtained from GEM CVP are sensitive to ambient temperature variations (gas/liquid equilibrium) and room air contamination (diffusion gradients). Use the solution in each ampule immediately upon opening.

Traceability

Ampules were equilibrated to 22°C and measured at 37°C. Established protocols for GEM CVP utilize NIST materials and recognized CLSI methods where available and appropriate.

Expected Values

The reference values for this product were based upon multiple determinations performed on GEM Premier 3500 and 3000 Systems.

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