BD Liquid-Based Cytology Specimen (LBC) Dilution Tubes for the BD ProbeTec™ Q^x Amplified DNA Assays

English: pages 1 – 2
Français : pages 2 – 3

Deutsch: Seiten 4 – 5
Español: páginas 5 – 6

8083954(02)
2014-08

Свържете се с местния представител на BD за инструкции. / Kontaktiraj lokalnog predstavnika BD za upute. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Нускалар ушн жергиликті BD өкілімен хабарласыңыз. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin.

INTENDED USE

The Liquid-Based Cytology Specimen (LBC) Dilution Tubes for the **BD ProbeTec** Q^x Amplified DNA Assays are designed to allow for detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* in gynecological specimens that are collected in **BD SurePath**™ Preservative Fluid or PreservCyt™ Solution and tested with the **BD ProbeTec** CT Q^x and GC Q^x Amplified DNA Assays on the **BD Viper**™ System in Extracted Mode using an aliquot that is removed prior to processing for either the **BD SurePath** or ThinPrep™ Pap test.

REAGENTS

Materials Provided: Each Liquid-Based Cytology Specimen (LBC) Dilution Tube for the **BD ProbeTec** Q^x Amplified DNA Assays contains approximately 1.7 mL of a Tris/Sodium Chloride solution and preservative.

Materials Required But Not Provided: Liquid-Based Cytology Specimen (LBC) Dilution Tube Caps for the **BD ProbeTec** Q^x Amplified DNA Assays, **BD Viper** LBC Specimen Rack, displacement pipettes, polypropylene aerosol-resistant tips capable of delivering 0.5 ± 0.05 mL, and a vortex mixer.

Storage Instructions: Store kit at 2 – 33 °C. Do not use past the expiration date.

Warnings and Precautions:

1. For *in vitro* Diagnostic Use.
2. For **BD SurePath** or PreservCyt specimens, use only the Liquid-Based Cytology Specimen (LBC) Dilution Tube for the **BD ProbeTec** Q^x Amplified DNA Assays.
3. LBC solutions contain flammable substances. Do not place specimens that have been transferred to the LBC Specimen Dilution Tubes in the **BD Viper** Lysing Rack or the Lysing Heater. Specimens transferred to the LBC Specimen Dilution Tubes should be placed in the **BD Viper** LBC Specimen Rack.
4. For testing with the **BD ProbeTec** Q^x Amplified DNA Assays on the **BD Viper** System in Extracted Mode, be sure to obtain aliquots of specimens collected in **BD SurePath** Preservative Fluid or PreservCyt Solution, prior to processing for either the **BD SurePath** or ThinPrep Pap test. Failure to do so may result in erroneous results.
5. Do not store specimens in LBC Specimen Dilution Tubes beyond the product expiration date.

SPECIMEN COLLECTION, TRANSFER, STORAGE AND TRANSPORT

LBC Specimen Collection, Storage and Transport

BD SurePath or PreservCyt specimens must be collected using either an endocervical broom or a brush/spatula combination as described in the **BD SurePath** Preservative Fluid or PreservCyt Solution product insert. LBC specimens in their original vials can be stored and/or transported for up to 30 days at 2 – 30 °C from the date of collection.

Specimen Transfer to LBC Specimen Dilution Tubes

A 0.5 mL aliquot of the **BD SurePath** or PreservCyt specimen must be transferred from the original vial to the LBC Specimen Dilution Tube prior to processing for either the **BD SurePath** or ThinPrep Pap test.

Wear gloves when handling the LBC Specimen Dilution Tube and the **BD SurePath** or PreservCyt specimen vial. If gloves come in contact with the specimen, immediately change them to prevent contamination of other specimens.

BD SurePath Specimen Transfer

NOTE: Refer to the **BD PrepStain**™ System product insert for instructions on removing an aliquot from the **BD SurePath** specimen vial prior to performing the **BD SurePath** liquid-based Pap test.

1. Label an LBC Specimen Dilution Tube with patient identification information.
2. Remove the cap from the LBC Specimen Dilution Tube.
3. Transfer 0.5 mL from the specimen vial to the LBC Specimen Dilution Tube. Avoid pipetting fluid from the bottom of the vial. Discard pipette tip. NOTE: A separate pipette tip must be used for each specimen.
4. Tighten the cap on the LBC Specimen Dilution Tube securely.
5. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the specimen and diluent are well mixed.

PreservCyt Specimen Transfer

NOTE: Refer to the ThinPrep™ 2000/3000 System Operator's Manual Addendum for instructions on removing an aliquot from the PreservCyt sample vial prior to performing the ThinPrep Pap test.

1. Label an LBC Specimen Dilution Tube with patient identification information.
2. Remove the cap from the LBC Specimen Dilution Tube.
3. Transfer 0.5 mL from the specimen vial to the LBC Specimen Dilution Tube. Avoid pipetting fluid from the bottom of the vial. Discard pipette tip. **NOTE:** A separate pipette tip must be used for each specimen.
4. Tighten the cap securely on the LBC Specimen Dilution Tube.
5. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the specimen and diluent are well mixed.

LBC Specimen Storage and Transport in the LBC Specimen Dilution Tube

Store the **BD SurePath** or PreservCyt specimens that have been transferred to the LBC Specimen Dilution Tubes at 2 – 30°C for up to 30 days. Specimens may also be stored at -20 °C for up to 90 days. Specimen storage should not exceed the expiration date of the LBC Specimen Dilution Tube.

For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Recommended time and temperature conditions for storage must be maintained during transport.

AVAILABILITY

Cat. No. Description

441444	Liquid Based Cytology Specimen (LBC) Dilution Tubes for the BD ProbeTec™ Q ^x Amplified DNA Assays, 400
441443	Liquid Based Cytology Specimen (LBC) Dilution Tube Caps for the BD ProbeTec™ Q ^x Amplified DNA Assays, 400

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD Liquid-Based Cytology Specimen (LBC) **Tubes de dilution pour les dosages BD ProbeTec Q^x** **Amplified DNA Assays**

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PreservCyt and ThinPrep are registered trademarks of Hologic, Inc.
BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD
Benex Limited