Imuno-Hemato Anti-A e Anti-B

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ENGLISH

SUMMARY
The ABO System was described in 1900 by Landsteiner and remains until today as the most important system in blood transfusion. The wrong ABO transfusion can result in death of the patient, with an intravascular hemolytic reaction, followed by biochemical and immunological changes.

The determination of blood group consists in the identification of the antigens of the ABO system existing in erythrocytes and that are genetically determined by 3 genes, respectively, the genes A, B and O.

There are two types of antibodies in ABO blood system: the natural occurrence and the immune. The natural occurrence begins to appear in serum approximately three to six months after the birth. They are the anti-A (individuals in the group B), anti-B (individuals in the group A), anti-A and anti-B (individuals in the group O) antibodies and absent in individuals of group AB. The natural antibodies represent a mixture with higher amount of immunoglobulins of class M (IgM) than the immunoglobulins class G (IgG). The immune occur through prior alloimmunizations, which can be through heteroimmunization through substances of animal or bacterial origin, or through alloimmunization, through pregnancy or ABO incompatible transfusion. These antibodies are usually referred to as hemolysins, and most of them are IgG.

Blood typing of an individual should always be performed through two tests: the direct and the reverse. In the direct test, the antisera is tested against 5% suspension of red blood cells of the individual, whose objective is to identify the antigens present on the surface of the red blood cell. Individuals from the group A will have antigens A in their erythrocytes, from group B will have antigen B, from AB group will have antigens A and B, and from group O does not have either antigens A or B on its surface. In the reverse test, the objective is to identify the free antibodies in the serum, testing the sera of individuals against a 5% suspension of known red blood cells, one of red blood cell of type A, negative Rh, and another of type B, Rh negative. Individuals in group A will have anti-B antibodies, in group B will have anti-A antibodies, in AB group will not have antibodies, and in the group O will have anti-A and anti-B circulating antibodies.

The sera of newborns does not yet contain the natural antibodies anti-A and anti-B, and can also have the presence of maternal anti-A and anti-B antibodies and, therefore, the reverse test should not be performed in this age group, because it will lead to incorrect results. The use of direct and reverse test leads to correct results and prevent complications, such as the presence of subgroups A and B that can produce weak or negative reactions in direct test.

The anti-serum **Imuno-HEMATO Anti-A and Anti-B from WAMA** are monoclonal IgM class antibodies, prepared from the supernatant of in vitro cell line culture of hybridized mouse secreting immunoglobulins. The anti-serum offers specificity (antigenicity), potency (titer) and avidity (quickness to establish the agglutination) equal to or greater than the international reference standards, ensuring optimum quality in the identification of A and B antigens on the blood typing tests.

PRINCIPLE OF THE METHOD
The procedure of the assay is based on the principle of agglutination, where red blood cells of individuals possessing A and/or B antigen on their surface will agglutinate in the presence of anti-A and/or anti-B antiserum, indicating a positive result. The test is considered negative when no agglutination appears.

KIT PRESENTATION
Imuno-Hemato Anti-A
[REF**]** **1054010-A – 10ml (anti-A)**
1. Anti-A antiserum (1 x 10 mL)
2. Instructions for use

[REF**]** **1054050-A – 50ml (anti-A)**
1. Anti-A antiserum (5 x 10 mL)
2. Instructions for use

[REF**]** **1054100-A – 100ml (anti-A)**
1. Anti-A antiserum (10 x 10 mL)
2. Instructions for use

Imuno-Hemato Anti-B
[REF**]** **1054010-B - 10ml (Anti-B)**
1. Anti-B antiserum (1 x 10 mL)
2. Instructions for use

[REF**]** **1054050-B - 50ml (Anti-B)**
1. Anti-B antiserum (5 x 10 mL)
2. Instructions for use

[REF**]** **1054100-B - 100ml (Anti-B)**
1. Anti-B antiserum (10 x 10 mL)
2. Instructions for use

REAGENT PREPARATION AND STABILITY
• **ANTISERUM ANTI-A:** ready to use. Contains monoclonal IgM antibodies, anti-A antigen of the ABO System. Reagent containing blue dye. Stable at 2 – 8 °C until expiration date printed on the bottle. Improper Storage affects the effectiveness of the reagent. Contains NaN3 < 0.095% as a preservative. DO NOT FREEZE. For in vitro diagnostic use only.

• **ANTISERUM ANTI-B:** ready to use. Contains monoclonal IgM antibodies, anti-B antigen of the ABO System. Reagent containing yellow dye. Stable at 2 – 8 °C until expiration date printed on the bottle. Improper Storage affects the effectiveness of the reagent. Contains NaN3 < 0.095% as a preservative. DO NOT FREEZE. For in vitro diagnostic use only.

Note: The kit has the same performance after first use, and is stable until the expiration date described on the label if it is kept in the indicated temperature (2-8 °C).

SPECIMENS
Use fresh blood sample. Avoid hemolyzed samples.
If it is not possible to perform the tests after collection, the sample must be kept refrigerated at 2 – 8 °C. In accordance with the anticoagulant used, the samples should be tested within the following time limits:

- Heparine or Oxalate Sodium : 2 days
- EDTA or Sodium Citrate : 14 days
- ACD or CPD : 28 days

The long term storage of the sample may cause deterioration of the antigen and result in incorrect classification.

MATERIAL REQUIRED BUT NOT PROVIDED
• Glass slide Technique : glass slide, 50 µl pipettor and tips or Pasteur pipette, mixing rod.
• Test Tube Technique: test tubes (10 x 75 mm or 12 x 75 mm), 50µl micropipette and tips, centrifuge, isotonic 0.85 or 0.9 % saline solution.
• Container for material disposal.

PROCEDURE
Regardless of the technique used, let the reagent and the sample reach room temperature before starting the tests.

I. GLASS SLIDE TECHNIQUE
1. On a clean glass slide, drip 1 drop (50 µl) of Anti-A or Anti-B antiserum (1) and 1 drop (50 µl) of whole blood to be tested.
2. Mix evenly the reagent and blood cells using a mixing rod, over an area of approximately 2.5 cm in diameter.
3. With gentle rotation movements and tilt of the glass slide check for the formation of agglutination within 1 minute (in the presence of antigen A, the reaction starts in seconds). Nonspecific Reaction can occur due to drying of the reaction mixture or if the blade is heated.

II. TUBE TECHNIQUE
1. Prepare a suspension in 5% of washed erythrocytes to be tested in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution, pipette 50 µl of washed erythrocytes + 1 ml of saline into a 12 x 75 mm test tube. Other volumes can be used).
2. In a labeled 10 x 75 mm or 12 x 75 mm test tube, pipette 50 µl of 5% erythrocytes suspension.
3. Pipette in the same tube 50 µl of the Anti-A or Anti-B antiserum.
4. Mix by gentle agitation and incubate at room temperature (15 -30 °C) for 1 to 15 minutes.
5. Centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.

6. Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

III. MICROPLATE TECHNIQUE

- Prepare a suspension of 2-3% of washed erythrocytes to be tested in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution and, in a 12 x 75 mm test tube, , pipette 25 µl of washed erythrocytes + 1 ml of saline- other volumes can be used).
- In a microplate with U-bottom, pipette 50 µl of suspension at 2.5% in each cavity of the microtitre plate, according to the number of samples to be tested and the distribution plan on board.
- Pipette in each cavity 50 µl of the anti-A or anti-B antiserum (1).
- Mix by gentle agitation, preferably with a microplate shaker, and incubate at room temperature (15- 30° C) for 15 minutes.
- Centrifuge the microplate 140 r.c.f. for 1 minute.
- Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

IMPORTANT:

→ Although the slide glass technique is used by many laboratories, it is suggested to use the tube or microplate technique, with direct and reverse test, because such procedure avoids some interpretation errors when only the search of antigen is used.

INTERPRETATION OF THE RESULTS

POSITIVE REACTION: visible agglutination of erythrocytes, indicating the presence of A or B antigen or both in the tested sample.

NEGATIVE REACTION: Absence of agglutination of erythrocytes, indicating absence of antigen A and/or B in the tested sample.

DISCREPANT REACTION: if the results obtained in the reverse test do not correlate with the direct test, it necessary to perform other investigations.

INTERPRETATION OF REACTION FOR BLOOD TYPING				
DIRECT TEST - REAGENT		REVERSE TEST - SUSPENSION		BLOOD TYPES
Anti -A	Anti -B	A1	B	
+	-	-	+	A
-	+	+	-	B
+	+	-	-	AB
-	-	+	+	O

TEST SPECIFICITY:

The performance of the ***Imuno-HEMATO Anti-B and Anti-B*** was evaluated in 1125 samples collected with the recommended anticoagulants. The evaluation showed 100% specificity for each reagent when compared with the results of a Kit used on the market.

PRECISION OF THE TEST:

a. Intra-Assay Precision

4 samples were selected from quality control (1 for each blood type) to perform the repeatability test. 10 repetitions of the test was done by the same operator on the same day. According to the tests carried out, it was verified that the kit ***Imuno-Hemato Anti-A and Anti-B*** Family repeated satisfactorily the results with the samples tested under the same conditions.

b. Inter-Assay Precision

4 samples were selected from quality control (1 for each blood type) to perform the test of reproducibility; 20 tests were performed with the same samples in different days and with different technicians.

According to the tests carried out, it was verified that the kit ***Imuno-Hemato Anti-A and Anti-B*** Family repeated satisfactorily the results with the samples to be tested even under different conditions.

LIMITATIONS OF USE

- False-positive and negative results can occur due to contamination of the test materials, inadequate concentration of suspensions of cells, inadequate incubation time or temperature, or errors of procedure.
- The centrifugation times mentioned are recommendations. A appropriate centrifugation time must be determined by each laboratory, and will be the one that produces the strongest reaction of the antibody with the antigen in positive reactions, and will also to allow for easy resuspension of RBC buttons.
- Weaker reactions will occur more frequently with stored blood than with fresh blood.
- ABO antigens are not fully developed at birth, therefore, weak reactions can occur when using red blood cells from the umbilical cord or of neonates. Cord Samples contaminated with Wharton's jelly can give false-positive results.
- Blood Samples of subgroups A and B may let to false-negative results or weak reactions.
- It is recommended to you use control with red blood cells with knowingly present antigen A and negative RBCs, daily, before starting the tests to confirm the reactivity and specificity of the reagent.
- Reverse test should always be performed to search for circulating antibodies in individuals older than 6 months of age, in order to confirm that the results of the direct test correlates with the reverse. If they do not correlate, other investigations should be carried out.
- Red blood cells coated with alloantibodies or auto-antibodies with specificity equal or similar to the reagent, which is, cells that are positive in the direct antiglobulin test (Coombs test) can give weak responses. In extreme cases, false-negative results may occur.

PRECAUTIONS AND WARNINGS

- The reagents should be kept between 2 - 8°C.
- The expiry date refers to the last day of the month indicated on the labels of the vial and of the kit box.
- Avoid exposing the reagents to high temperatures and direct sunlight.
- Do not freeze the reagent, as this may cause irreversible damage.
- The **Imuno-HEMATO Anti-A and Anti-B** from **WAMA**, though derived from animal source, cannot assume to be free of infectious agents. Therefore, it is recommended to treat all reagents as potentially infectious, as well as having the care in disposing of these materials.
- Allow the reagents to reach the room temperature (20 -25°C) before starting the test.

- The reagents must be capped immediately after use.
- Do not use the reagent after the expiration date.
- The procedures are only for manual testing. When using automatic or semi-automatic instruments, procedure contained in the operator's manual provided by the instrument manufacturer must be followed.
- Discard the material following local regulations.
- Use the Good Laboratory Practices (GLPs) in storage, handling and disposal of materials.

WARRANTY

WAMA Diagnóstica guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. **WAMA** and its distributors are not liable for failures on the performance of the kits under these conditions.

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- Levine, P.: The discovery of Rh haemolytic disease. Vox Sanguinis 47: 187-190, 1984.
- Jones, J., Scott, M.M. and Voak, D.: Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine, 5: 171-184, 1995.
- Avent, N.D. and Reid, M.E.: The Rh blood group system: a review. Blood, 95: 375-387, 2000.
- Zallen, D. T., Christie, D. A. and Tansey, E. M.: The Rhesus Factor and Disease Prevention. Witness Seminar held by the Wellcome Trust Centre for the History of Medicine at UCL, London, June, 2003.
- Scott, M. L.: The complexities of the Rh system. Vox Sanguinis, 87: 58-62, 2004.
- Westhoff, C. M.: The Rh blood group system in review: A new face for the next decade. Transfusion, 44: 1663-1673, 2004.
- Moulds, M.K.: Review: monoclonal reagents and detection of unusual or rare phenotypes or antibodies. Immunohematology, 22 (2): 52-63, 2006.
- Van Kim, C., Colin, Y. and Cartron, J.: Rh proteins: key structural and functional components of the red cell membrane. Blood Reviews, 20: 93-110, 2006.
- Westhoff, C.M.: The structure and function of the Rh antigen complex. Seminars in Hematology, 44: 42-50, 2007.

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