IL Test™ Glucose (Oxidase)
Cat. No. 0018259140 - 10 x 30 mL Glucose OX
Cat. No. 0018250840 - 5 x 95 mL Glucose OX

Intended Use
IL Test™ Glucose (Oxidase) is intended for the quantitative in vitro diagnostic determination of glucose in human serum, plasma, urine or cerebrospinal fluid (CSF) using the ILab Chemistry System.

Summary
Glucose measurements are used in the diagnosis and treatment of disorders of carbohydrate metabolism such as diabetes mellitus, hypoglycemia and hyperglycemia. Refer to the publication “Effects of Diseases on Clinical Laboratory Tests”¹ for a summary of causes of increased and decreased glucose.

Principle
End point analysis. Trinder methodology: Glucose oxidase (GOD) / Peroxidase (POD)^{2,3}.
$$\text{Glucose oxidase} \\ \text{β-D-Glucose} + \text{O}_2 + \text{H}_2\text{O} \longrightarrow \text{Gluconic acid} + \text{H}_2\text{O}_2$$
$$\text{Peroxidase} \\ 2\text{H}_2\text{O}_2 + \text{phenol} + 4\text{-aminoantipyrine} \longrightarrow \text{red quinoneimine} + 4\text{H}_2\text{O}$$

The increase in absorbance generated by the red dye is proportional to the glucose concentration in the sample.

Composition
Glucose Ox containing: Glucose oxidase >12500 U/L, Peroxidase >2500 U/L, Phenol 7.5 mmol/L, 4-aminoantipyrine 1.25 mmol/L, Phosphate buffer and less than 0.1% of sodium azide.

PRECAUTION AND WARNING
Glucose Ox is not classified as dangerous.
Hazard class: none.
Risk phrases: none.
Safety phrases: none.
This product is For *in vitro* Diagnostic Use

Reagent Preparation
Cat. No 0018259140: On ILab 600/650 reagent is ready-to-use. Uncap and place Glucose OX bottle into instrument reagent tray.
On ILab Aries: using a volumetric pipette add 15 mL of deionized or distilled water to one bottle of Glucose OX. Invert to mix then place into ILab Aries reagent tray.
Cat. No 0018250840: On ILab 600/650 reagent is ready-to-use. Uncap and place Glucose OX bottle into instrument reagent tray.
On ILab Aries: using a volumetric pipette mix in empty bottle 30mL of Glucose OX with 15 mL of deionized or distilled water. Invert to mix then place into ILab Aries reagent tray.

Reagent storage and Stability
Unopened reagents are stable until expiration date indicated on the bottle when stored at 2-8°C.
ILab 600/650 : 35 days on board.
ILab Aries: 50 days on board.

Instrument/Test Procedures
For the assay procedure refer to the relevant instrument Application Manual and to the relevant Instrument Operator Manual.

IL Test™ Glucose (Oxidase)
Cat. No 0018259140 - 10 x 30 mL Glucose OX
Cat. No 0018250840 - 5 x 95 mL Glucose OX

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Specimen collection and Preparation
Serum, lithium heparinized plasma, urine, CSF.
Serum: Specimens should be separated from cells as soon as possible after collection to avoid loss due to glycolysis. Glucose is stable in serum for up to 3 days at 2-8°C if kept free of bacterial contamination.⁴
Urine: 24 hours urine collections must be made in amber tinted containers and remain refrigerated (2-8 °C) during collection period. Prior to collection add 5 mL Glacial Acetic Acid to container. Final pH of 4-5 is needed to inhibit bacterial activity. Stable for 3 days at 2-8 °C if kept free of bacterial contamination.⁴
CSF: analyze immediately to avoid falsely low results.⁴ Bring samples to room temperature and mix well before analysis.
Note: Analyze urine and CSF samples without dilution.

Calibration
Use ReferrIL™ G Calibrator (Cat. No. 0018257000). See calibrator insert for lot specific concentration.
Urine
Since for either the urine or the CSF calibration the ILab600/650 uses the K-Acceptance method, the results are calculated from the serum/plasma calibration. A separate calibration is not necessary. For other ILab systems refer to the relevant Application Manual.
On ILab 650 and ILab Aries recalibrate every 14 days, when a new lot of reagent is used, or when adjustments are made to the instrument. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Quality Control
Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice. Normal and abnormal controls such as SeraChem® Control Level 1 (Cat. No 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program. Refer to Westgard et al.⁵ for identification and resolution of out-of-control situations.

Calculation of Analytical Results
Refer to the relevant Application Sheet and to the relevant ILab Operator Manual.
To convert mg/dL to SI units (mmol/L) multiply by 0.055.⁶

Reference Interval
Serum: 74 - 106 mg/dL (4.1 - 5.9 mmol/L).⁴
Urine: less than 0.5 g/24h (2.8 mmol/24h).⁴
CSF: 40 - 70 mg/dL (2.2 - 3.9 mmol/L).⁴
IL recommends that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Performance Characteristics
All performance data were obtained on ILab 600/650. For performance data on other ILab systems refer to the relevant application Manual.

Limitation/Interfering Substances
No interference from lipemia up to sample absorbance of 4.9/cm at 660 nm (1100 mg/dL or 12 mmol/L triglycerides). No interference up to 250 mg/dL (0.150 mmol/L) hemoglobin. No interference up to 5 mg/dL (85 µmol/L) bilirubin.
NOTE: urine and CSF measurements by this method are subject to interferences by uric acid and ascorbic acid. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*⁷

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ENGLISH - Insert revision 09/2010					
Precision					
Typical precision results obtained according to NCCLS protocol. ⁸					
	Samples/Runs	Mean (mg/dL)	CV (%)	Mean (mg/dL)	CV (%)
Serum					
Within Run	3/20	73	1.3	248	0.8
Total	3/20	73	1.5	248	0.9
Urine					
Within Run	3/20	56	1.4	201	1.3
Total	3/20	56	2.6	201	1.9
CSF					
Within Run	3/10	46	1.5	90	0.8
Total	3/10	46	1.5	90	1.0

Method Comparison					
	Serum	Urine	CSF	Serum	Serum
Comparison Method (x)	Same Reagent		Hexokinase PN0018250740*	Hexokinase*	
Comparison Instrument	ILab 900		ILab 650	Competitor Instrument*	
Slope	0.944	0.949	0.932	0.9979	1.024
y intercept	7.9	4.3	-0.1	-4.13	-0.1542
Range (mg/dL)	51 - 393	0 - 801	40 - 226	30 - 358	51 - 259
Mean x (mg/dL)	107	172	94	146	118
Mean y (mg/dL)	109	168	87	141	121
r	0.997	0.997	0.999	0.999	0.999
SE	4.6	11.5	0.9	4.68	2.40
n	137	80	20	48	53

* the IL Glucose Oxidase reagent correlates well against both IL Glucose Hexokinase and a Glucose Hexokinase available commercial method.

Linearity		
	Serum	
	Linearity range (mg/dL)	Linearity range (mmol/L)
no rerun	2 – 512	0.1 – 28.2
with rerun	2 – 3000	0.1 – 165.0
	Urine, CSF	
	Linearity range (mg/dL)	Linearity range (mmol/L)
no rerun	2 - 512	0.1 – 28.2
with rerun	2 - 5120	0.1 – 281.6
Sensitivity		
milliabsorbance change per 1 mg/dL 0.8		
milliabsorbance change per 1 mmol/L 14.5		
Minimum Detection Limit		
2 mg/dL (0.1 mmol/L).		

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Calculation of Analytical Results

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Limitation/Interfering Substances

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ILab™ 600/650 Test Parameters		
Photometric Test Parameters	Serum	Urine/CSF
Test No.	12	54
Test Name, Test Code	Glucose, GLUOX	Glucose UR/CSF, GLUOX U/CSF
Sample Type	Serum	Urine/other
Reporting Unit, Decimal Points	mg/dL (mmol/L), 0 (1)	mg/dL (mmol/L), 0 (1)
Reaction Cycle	Standard	Standard
Twin Analysis	OFF	OFF
Methodology Type, Measuring Point	End point, 1/3	End point, 1/3
Photometric Methodology	2 Wavelength	2 Wavelength
Primary/Secondary Wavelength	510/600	510/600
Sampling Conditions		
Sampling 1	Sample Vol.	3
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	2
	Sample/Diluent Vol.	0/0
Sampling 3	Sample Vol.	6
	Sample/Diluent Vol.	10/90
Sampling 4		
Diluent Code	***	***
	Water	Water
Diluent Warning Limit	N/A	N/A
First Run	Sampling 1	Sampling 1
Below/Above Normal Range	***	***
Panic L	***	***
Panic H	Sampling 2	Sampling 2
Noise	***	***
Prozone	N/A	N/A
HIGH!, ABS!	Sampling 3	Sampling 3

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1. Friedman R., Young D.S.: Effects of Disease on Clinical Laboratory Tests, AACC Press, Washington, D.C. 1989

2. Komberg A., Horecker B.L.: Methods in Enzymology, Colowickand Kaplan, eds., Vol. 1, Academic Press, New York 323, 1955.

3. Bathelmal W., Czock R., Enzymatic detriminations of glucose in the blood, cerebrospinal fluid and urine, Klin.Wochennschr., 40:585,1962.

4. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, 3 rd edition 1999, pp. 776-82.

5. Westgard J., Barry P., Cost Effective Quality Control: Managing the Quality and Productivity of Analitycal Processes, AACC Press, Washington, D.C. 1986

6. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, 3 rd edition 1999, pp. 1815.

7. Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press Washington, D.C., 1990.

8. Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

IL Test™ and SeraChem®, are trademarks of Instrumentation Laboratory Company.

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Instrumentation Laboratory SpA - V.le Monza 338 - 20128 Milano (Italy)
ATTENZIONE! Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.