# Imuno-Hemato

## Anti-D (Rho)

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Monoclonal antibodies IgG/IgM (Blend) reagent for qualitative determination of D (Rho) Antigens, Rh factor, in human erythrocytes in slide, tube and microplate by agglutination technique.

Monoclonal antibodies IgG/IgM (Blend) reagent for qualitative determination of D (Rho) Antigens, Rh factor, in human erythrocytes in slide, tube and microplate by agglutination technique.

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#### ENGLISH

#### SUMMARY

The Rh system was discovered by Wiener, Landsteiner and Levine in 1940 when they were studying of sensitization in rabbits with red blood cells of Rhesus monkeys.

After the antigens of the ABO system (A and B), the D antigen (Rh) is the most important antigen in blood typing, once it is related to the development of the perinatal hemolytic disease, hemolytic transfusion reactions, and autoimmune hemolytic anemias. The phenotype D negative occurs in 15% of Caucasians and 9 to 10% of the blacks.

The Rh is a complex system and currently comprises of f4f4 antigens, where the main ones are the D, C/c, E/e. The D antigen is the most immunogenic, although all others may induce the production of antibodies that are related with the majority of the problems associated with alloimmunization by transfusion and pregnancy.

Some individuals may present negative D when directly tested with anti-D antisera, but when the test continues to the human antiglobulin phase (Coombs Serum), variable agglutination occurs. These individuals were originally called D<sup>+</sup> positive. Currently, the term D<sup>+</sup> was replaced by "weak D", and they should be considered as individuals D (Rh) positive.

Some individuals may have partial D antigen, characterized by the absence of some of its epitopes. This category of partial D (Rh) antigens is differentiated one from others by the presence or absence of one or more epitopes (DII, DIII and DIV, DV, DVI, DVII, DFR, DBT, RoHan).

The anti-serum Imuno-HEMATO Anti-D (Rh) from WAMA mixture of human monoclonal antibodies, prepared from the supernatant of in vitro culture of heterohibridoma cells, where a strain of cells secrete IgG antibodies and another strain secretes IgM antibodies, which react specifically with the corresponding antigen, providing assurance in the determination of the presence or absence of the D antigen on human red blood cells. The antibody reacts negatively with erythrocytes that have partial D antigen, category DVI, with the IgM antibodies and positively with the IgG antibodies in indirect Coombs tests.

#### PRINCIPLE OF THE METHOD

The procedure of the assay is based on the principle of agglutination, where red blood cells of individuals possessing D antigen on their surface will agglutinate in the presence of anti-D antiserum, indicating a positive result. The test is considered negative when no agglutination appears. In such cases, the test should continue to search for antigen "weak D".

#### KIT PRESENTATION

REF 1055010-D - 10ml

- 1. Anti-D antiserum (1 x 10 ml)

2. Instructions for use

REF 1055050-D - 50ml

- 1. Anti-D antiserum (5 x 10 ml)

2. Instructions for use

#### REAGENT PREPARATION AND STABILITY

• **ANTISERUM ANTI-D:** ready to use. Containing a blend of IgM and IgG anti-antigen D monoclonal antibodies.

Stable at 2 – 8 °C until expiration date printed on the bottle. Improper Storage affects the effectiveness of the reagent. Contains sodium azide < 0.095% as preservative. Also contains sodium chloride, bovine albumin and extenuating macromolecular architecture. Supplied in dropper bottles with a volume of 50µl per drop. DO NOT FREEZE. For in vitro diagnostic use only.

Note: The kit has the same performance after first use, and is stable until the expiration date described on the label if it is kept in the indicated temperature (2 – 8 °C).

#### SPECIMENS

&lt;p

D antigen in the tested sample.

**NEGATIVE REACTION:** visible agglutination of erythrocytes, indicating the presence of D antigen in the tested sample. In this case, the search for "weak" D antigen must be performed through the indirect Coombs test.

#### INDIRECT COOMBS TEST TO SEARCH FOR "WEAK" D

If you are performing the technique in tube, continue the search for the "weak" D in the same tube; otherwise follow the procedure:  
1. Prepare a suspension in 5% of washed erythrocytes to be tested in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution, pipette 50µl of washed erythrocytes + 1ml of saline into a 12 x 75 mm test tube. Other volumes can be used).  
2. In a labeled 10 x 75 mm or 12 x 75 mm test tube, pipette 50µl of 5% erythrocytes suspension.  
3. Pipette in the same tube 50µl (1 drop) of the anti-D antiserum (1).  
4. Mix by gentle agitation and incubate at 37°C for 15-30 minutes.  
5. Wash the cells 3 times in 0.85% or 0.9% saline: fill the tube with saline, resuspending the cells, mix well through gentle agitation, centrifuge at 2,000 RPM (800 to 1,000 g) for 1 minute and discard the supernatant. Repeat the washing 3 times.  
6. After discarding the supernatant in the last wash, pipette 100µl (2 drops) of the anti-human polyclonal Coombs serum from WAMA, mix and centrifuge at 2,000 RPM (800 to 1,000 g) for 1 minute.  
7. Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

**INTERPRETATION:** visible agglutination of erythrocytes indicates the presence of "weak" D antigen in the sample tested, and must be reported as D(Rho) positive ("weak" D). Absence of agglutination of erythrocytes indicates absence of D antigen in the tested sample and, therefore, D(Rho) negative.

**CAUTION:** to perform the test searching for Rh antigen and the Indirect Coombs for "weak" D, an Rh control must always be used. We recommend the use of the Rh serum control from WAMA. The purpose for using serum control is to identify false-positive reactions through the presence of autoimmune antibodies, cold agglutinins or formation of royleas. The false positive reaction detected in the human antiglobulin phase (Coombs Serum) generally indicates that the red blood cells are positive for in the Direct Coombs test. For the correct interpretation of the test with Rh Control, see the reagent package insert.

#### TEST SPECIFICITY:

The performance of the **Imuno-HEMATO Anti-D (Rho)** was evaluated in 1125 samples collected with the recommended anticoagulants. The evaluation showed 100% specificity for each reagent when compared with the results of a Kit used on the market.

#### PRECISION OF THE TEST:

##### a. Intra-Assay Precision

2 samples were selected from quality control (1 positive and 1 negative) to perform the repeatability test. 20 repetitions of the test was done by the same operator on the same day.

According to the tests carried out, it was verified that the kit **Imuno-Hemato Anti-D (Rho)** repeated satisfactorily the results with the samples to be tested under the same conditions.

##### b. Inter-Assay Precision

20 samples were selected from quality control (16 negatives and 4 negatives) to perform the test of reproducibility: 3 tests were performed with the same samples in different days and technicians and with two different lots.

According to the tests carried out, it was verified that the kit **Imuno-Hemato Anti-D (Rho)** repeated satisfactorily the results with the samples to be tested even under different conditions.

#### LIMITATIONS OF USE

False-positive and negative results can occur due to contamination of the test materials, inadequate concentration of suspensions of cells, inadequate incubation time or temperature, or errors of procedure.

The centrifugation times mentioned are recommendations. A appropriate centrifugation time must be determined by each laboratory, and will be the one that produces the strongest reaction of the antibody with the antigen in positive reactions, and will also to allow for easy resuspension of RBC buttons.

Weaker reactions will occur more frequently with stored blood than with fresh blood.

• Positive and negative control should be used when carrying out the tests.

• In the glass slide method, "weak" D antigens and certain types of partial D may not be recognized.

Improper Storage affects the effectiveness of the reagent.

• It's recommended to always use an Rh control to avoid false positive results.

• Red blood cells coated with alloantibodies or auto-antibodies, whose cells are positive in the direct antiglobulin test (Coombs test) can give weak responses. In extreme cases, false-negative results may occur.

The described procedures are for use with the Coombs sera of **WAMA**. As a principle, Coombs sera from other manufacturers may be used, but the Laboratory must follow the procedures of these manufacturers. Also the laboratory must follow approved validation procedures to demonstrate the compatibility of these products.

#### PRECAUTIONS AND WARNINGS

- 1. The reagents should be kept between 2 - 8°C.
- 2. The expiry date refers to the last day of the month indicated on the labels of the vial and of the kit box.
- 3. Avoid exposing the reagents to high temperatures and direct sunlight.
- 4. Do not freeze the reagent, as this may cause irreversible damage.
- 5. The **Imuno-HEMATO Anti-D (Rho)** from **WAMA**, though derived from animal source, cannot assume to be free of infectious agents. Therefore, it is recommended to treat all reagents as potentially infectious, as well as having the care in disposing of these materials.
- 6. Allow the reagents to reach the room temperature (20 - 25°C) before starting the test.
- 7. The reagents must be capped immediately after use.
- 8. Do not use the reagent after the expiration date.
- 9. The procedures are only for manual testing. When using automatic or semi-automatic instruments, procedure contained in the operator's manual provided by the instrument manufacturer must be followed.
- 10. Discard the material following local regulations.
- 11. Use the Good Laboratory Practices (GLPs) in storage, handling and disposal of materials.

#### WARRANTY

**WAMA Diagnóstica** guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. **WAMA** and its distributors are not liable for failures on the performance of the kits under these conditions.

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