Imuno-Hemato Soro de Coombs

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ENGLISH

SUMMARY

Coombs, Mourant and Race, a few years after the discovery of anti-Rh antibodies(1945) , described the antiglobulin test, which uses human antiglobulin serum, prepared in rabbits, to detect "incomplete" Rh antibodies Today the test is used to detect in vivo sensitization of erythrocytes, helping in the diagnosis of newborn hemolytic disease and autoimmune hemolytic anemia, making it possibly the most widely used test in immunohaematology.

Generally the antibodies involved in transfusion reactions are of two types, complete and incomplete. The complete antibodies agglutinate red blood in saline solution and the incomplete sensitize red cells without agglutination. Usually IgM antibodies, and IgG1 and IgG3 types of IgG antibodies fix the complement. Thus, the cellular lysis, in vivo, is mediated through the complement system.

The anti-human globulin serum (Coombs reagent) is used in two types of tests: direct and indirect. In the direct, the anti-human globulin is used for detecting antibodies adsorbed to red blood cells in vivo, and, therefore, antibodies bound to red blood cells of the patient, which is used in the study of hemolytic disease of the newborn in negative Rh mothers, in the research of autoimmune hemolytic anemias, in drug induced hemolysis, and in post-transfusion hemolytic reactions. In the indirect test, the anti-human globulin is used for detecting antibodies adsorbed onto the red blood cells in vitro, and, therefore, free antibodies in serum or plasma, which is used in the research and identification of irregular antibodies and compatibility test.

The ***Imuno-HEMATO SORO DE COOMBS*** from **WAMA** is a serum of broad spectrum, human antigama and non gamma globulin, characterized by a mixture of antibodies against immunoglobulins and complement factors, prepared by immunization of rabbits with purified human IgG (polyclonal anti-IgG) and by in vitro culture of mouse hybridoma (cell line Bric-8) secreting anti-C3 monoclonal antibody, IgM class (monoclonal anti-C3d). The reagent is green in color.

PRINCIPLE OF THE METHOD

The procedure of the test is based on the principle of agglutination, where washed red blood cells from individual with antibodies or complement factors adsorbed to their surface (Coombs) or antibody-free (Indirect Coombs) will agglutinate in the presence of anti-human globulin Serum poly specific (Coombs), indicating a positive result. The test is considered negative when no agglutination appears.

KIT PRESENTATION

[REF**]** **1058010-D - 10ml**

- Coombs Serum (1 x 10 ml)
- Instructions for use

[REF**]** **1058050-D - 50ml**

- Coombs Serum (5 x 10 ml)
- Instructions for use

[REF]	1058100-D - 100ml
[REF]	1058100-D - 100ml

REAGENT PREPARATION AND STABILITY

• **COOMBS SERUM (1):** ready for use. Stable at 2 – 8 °C until expiration date printed on the bottle. Improper Storage affects the effectiveness of the reagent. Contains sodium azide < 0.095% (w/v) as a preservative, sodium chloride and bovine serum albumin. DO NOT FREEZE. For in vitro diagnostic use only.

Note: The kit has the same performance after first use, and is stable until the expiration date described on the label if it is kept in the indicated temperature (2-8 °C).

SPECIMENS

Use fresh blood sample. Avoid hemolyzed samples.

If it is not possible to perform the tests after collection, the sample must be kept refrigerated at 2 – 8 °C. In accordance with the anticoagulant used, the samples should be tested within the following time limits:

- Heparin or Sodium Oxalate
- EDTA or Sodium Citrate
- ACD or CPD

MATERIAL REQUIRED BUT NOT PROVIDED

- Test Tubes (10 x 75 mm or 12 x 75 mm)
- Micropipette of 50 µl, 100 µl and 1000 ul and tips.
- Centrifuge
- Physiological saline Solution 0.85 or 0.9%
- Water Bath
- Container for material disposal.

PROCEDURE

Regardless of the technique used, let the reagent and the sample reach room temperature before starting the tests.

I. ANTI-HUMAN GLOBULIN TEST (DIRECT COOMBS)

- Prepare a suspension in 5% of washed erythrocytes to be tested in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution, pipette 50 µl of washed erythrocytes + 1 ml of saline into a 12 x 75 mm test tube. Other volumes can be used).
- In a labeled 10 x 75 mm or 12 x 75 mm test tube, pipette 50 µl of 5% erythrocytes suspension.
- Complete the tube with physiological saline solution and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Decant the supernatant by inverting the tube, resuspend the erythrocyte pellet and complete again with physiological saline solution. Repeat this procedure 3 times.
- After the last wash, decant the supernatant completely, inverting the tube and drying the edge with filter paper.
- Pipette 100 µl (2 drops) of Coombs Sera (1).
- Mix by gentle agitation and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

II. INDIRECT ANTI-HUMAN GLOBULIN TEST (INDIRECT COOMBS)

- Prepare a suspension in 5% of washed erythrocytes compatible with the patient serum (usually from type O, positive Rh), which should be washed in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution, pipette 50 µl of washed erythrocytes + 1 ml of saline into a 12 x 75 mm test tube. Other volumes can be used). There is commercial RBCs to perform the test.
- In a labeled 10 x 75 mm or 12 x 75 mm test tube, pipette 50 µl of 5% erythrocytes suspension.
- Pipette in the same tube 100 µl of the serum to be tested.
- Mix by gentle agitation and incubate at 37 ° C for 15-30 minutes.
- Remove the tube from the Water-Bath, complete with physiological saline solution and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Decant the supernatant by inverting the tube, resuspend the erythrocyte pellet and complete again with physiological saline solution. Repeat this procedure 3 times.
- After the last wash, decant the supernatant completely, inverting the tube and drying the edge with filter paper.
- Pipette 100 µl (2 drops) of Coombs Sera (1).
- Mix by gentle agitation and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

INTERPRETATION OF RESULTS (DIRECT AND INDIRECT COOMBS)

POSITIVE REACTION: Visible agglutination of erythrocytes.

NEGATIVE REACTION: Absence of agglutination of erythrocytes.

III.SURVEY OF IRREGULAR ANTIBODIES IN INDIRECT ANTI-HUMAN GLOBULIN (INDIRECT COOMBS)

Its purpose is to detect antibodies of blood groups capable of causing hemolytic reactions. It is carried out in 4 phases: the 1st, 2nd, and 3rd stages are more efficient in the detection of IgM, while the 4th phase detects IgG.

1ST PHASE: Saline solution at Room Temperature

- Prepare a suspension in 5% of washed erythrocytes compatible with the patient serum (usually from type O, positive Rh), which should be washed in physiological saline solution (wash the cells 3 times in 0.85% or 0.9% saline solution, pipette 50 µl of washed erythrocytes + 1 ml of saline into a 12 x 75 mm test tube. Other volumes can be used). There is commercial RBCs to perform the test.
- In a labeled 10 x 75 mm or 12 x 75 mm test tube, pipette 50 µl of 5% erythrocytes suspension.
- Pipette in the same tube 100 µl of the serum to be tested.
- Mix by gentle agitation and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently re suspend the button of erythrocytes and examine for the presence of agglutination and/or hemolysis.

2ND PHASE: Protein Environment at Room Temperature

- Pipette in the same tube referred to above 100 µl (2 drops) of Bovine Albumin at 22% (1).
- Mix by gentle agitation and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently re suspend the button of erythrocytes and examine for the presence of agglutination and/or hemolysis.

3RD PHASE: Protein Environment at 37°C

- Incubate the tube at 37°C for 15-30 minutes.
- Centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently re suspend the button of erythrocytes and examine for the presence of agglutination and/or hemolysis.

4TH PHASE: Human Antiglobulin (Indirect Coombs)

- Complete the tube with physiological saline solution and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Decant the supernatant by inverting the tube, resuspend the erythrocyte pellet and complete again with physiological saline solution. Repeat this procedure 3 times.
- After the last wash, decant the supernatant completely, inverting the tube and drying the edge with filter paper.
- Pipette 100 µl (2 drops) of **SORO DE COOMBS** from **WAMA**.
- Mix by gentle agitation and centrifuge at 1,000 RPM (± 100 to 125 g) for 1 minute.
- Gently resuspend the button of erythrocytes and examine for the presence of agglutination.

TEST INTERPRETATION:

POSITIVE TEST: presence of agglutination and/or hemolysis in one or more phases of the test, characterizing presence of irregular antibodies.

NEGATIVE TEST: absence of agglutination and/or hemolysis in all phases of the test, characterizing absence of irregular antibodies.

TEST SPECIFICITY:

The performance of the ***Imuno-HEMATO Soro de Coombs*** was evaluated in 86 samples collected with the recommended anticoagulants. The evaluation showed 100% specificity for each reagent when compared with the results of a Kit used on the market.

Coefficient of Variation (inter / intra assay test) / Application Limit:

a. Intra-Assay Precision

2 samples were selected from quality control (2 negative) to perform the repeatability test. 20 repetitions of the test was done by the same operator on the same day.

According to the tests carried out, it was verified that the kit ***Imuno-Hemato Soro de Coombs*** repeated satisfactorily the results with the samples to be tested under the same conditions.

b. Inter-Assay Precision

10 semples were selected from quality control (10 negatives) to perform the test of reproducibility; 3 tests were performed with the same samples in different days and technicians and with same lot.

According to the tests carried out, it was verified that the kit ***Imuno-Hemato Soro de Coombs*** repeated satisfactorily the results with the samples to be tested even under different conditions.

LIMITATIONS OF USE

- False Results can occur due to contamination of the test materials, inadequate concentration of suspensions of cells, inadequate incubation time or temperature, or errors of procedure.
- The centrifugation times mentioned are recommendations. A appropriate centrifugation time must be determined by each laboratory, and will be the one that produces the strongest reaction of the antibody with the antigen in positive reactions, and will also to allow for easy resuspension of RBC buttons.
- Weaker reactions will occur more frequently with stored blood than with fresh blood.
- It is recommended to use fresh RBCs and sera, preferably not chilled.
- The inadequate washing of erythrocytes can cause false negative through the presence of free serum globulins, neutralising the Coombs serum.
- It is recommended to always use a positive and negative control for the quality control of the technique.

PRECAUTIONS AND WARNINGS

- The reagents should be kept between 2 - 8 °C.
- The expiry date refers to the last day of the month indicated on the labels of the vial and of the kit box.
- Avoid exposing the reagents to high temperatures and direct sunlight.

- Do not freeze the reagent, as this may cause irreversible damage.
- The ***Imuno-HEMATO Soro de Coombs*** from **WAMA**, though derived from animal source, cannot assume to be free of infectious agents. Therefore, it is recommended to treat all reagents as potentially infectious, as well as having the care in disposing of these materials.
- Because there is sodium azide as a preservative in the reagents, the disposal of reagents must be followed by large volumes of water to prevent azide buildup in pipes, because this can react with lead or copper, forming highly explosive salts. Furthermore, the azide is toxic when ingested.
- Allow the reagents to reach the room temperature (20 -25°C) before starting the test.
- The reagents must be capped immediately after use.
- Do not use the reagent after the expiration date.
- The procedures are only for manual testing. When using automatic or semi-automatic instruments, procedure contained in the operator's manual provided by the instrument manufacturer must be followed.
- The procedures are only for manual testing. When using automatic or semi-automatic instruments, procedure contained in the operator's manual provided by the instrument manufacturer must be followed.
- Discard the material following local regulations.
- Use the Good Laboratory Practices (GLPs) in storage, handling and disposal of materials.

WARRANTY

WAMA Diagnóstica guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. **WAMA** and its distributors are not liable for failures on the performance of the kits under these conditions.

WARRANTY

WAMA Diagnóstica guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. **WAMA** and its distributors are not liable for failures on the performance of the kits under these conditions.

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WARRANTY

WARRANTY

- Levine, P.: The discovery of Rh haemolytic disease. Vox Sanguinis 47: 187-190, 1984.
- Jones, J., Scott, M.M. and Voak, D.: Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine, 5: 171-184, 1995.
- Avent, N.D. and Reid, M.E.: The Rh blood group system: a review. Blood, 95: 375-387, 2000.
- Zallen, D. T., Christie, D. A. and Tansey, E. M.: The Rhesus Factor and Disease Prevention. Witness Seminar held by the Wellcome Trust Centre for the History of Medicine at UCL, London, June, 2003.
- Scott, M. L.: The complexities of the Rh system. Vox Sanguinis, 87: 58-62, 2004.
- Westhoff, C. M.: The Rh blood group system in review: A new face for the next decade. Transfusion, 44: 1663-1673, 2004.
- Moulds, M.K.: Review: monoclonal reagents and detection of unusual or rare phenotypes or antibodies. Immunohematology, 22 (2): 52-63, 2006.
- Van Kim, C., Colin, Y. and Cartron, J.: Rh proteins: key structural and functional components of the red cell membrane. Blood Reviews, 20: 93-110, 2006.
- Westhoff, C.M.: The structure and function of the Rh antigen complex. Seminars in Hematology, 44: 42-50, 2007.

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