***IL Test™ Uric Acid - 0018480100***

**Intended Use**


IL Test™ Uric Acid is intended for the quantitative *in vitro* diagnostic determination of uric acid in human serum, using the ILab 300 Plus Chemistry System.

**Summary**

Uric acid measurements are used in the diagnosis and treatment of gout and impaired renal function.<sup>1</sup>

**Principle**

End point analysis (Bichromatic). Uricase / Peroxidase (POD)<sup>1</sup>.

TIBA: 3-hydroxy- 2,4,6- triiodobenzoic acid

The increase in absorbance generated by the red dye is proportional to the uric acid concentration in the sample. Absorbance measurements are taken at a primary wavelength of 510 nm and a blanking wavelength of 620 nm.

**Composition**

UA Dry: 6 vials of lyophilized powder containing: uricase, peroxidase, 3-hydroxy- 2,4,6- triiodobenzoic acid, 4-aminoantipyrine, buffer.  
Uric Acid: 6 bottles 33 mL containing: purified water

**PRECAUTIONS AND WARNINGS**

- Irritating to eyes, respiratory system and skin (R36/37/38).
- Do not breathe dust (S22).
- Avoid contact with skin and eyes (S24/25).
- In case of contact with eyes rinse immediately with plenty of water and seek medical advice (S26).
- Do not empty into drains (S29).
- Wear suitable protective clothing (S36).

This product is For *in vitro* Diagnostic Use.

**Reagent Preparation**

Using the plastic adapter enclosed in the kit connect one vial of UA Dry with one bottle of Uric Acid: reconstitute the powder with the diluent, invert several times, do not shake and do not allow to foam, refill the Uric Acid bottle with the solution. Allow the working reagent to stand for 30 minutes and mix gently before use. Place the bottle into the ILab 300 Plus reagent rack in the refrigerated area (positions 14 to 33). After daily use cap the bottle and store at 2-8°C.

**Reagent Storage and Stability**

Unopened reagents are stable until the expiration date indicated on the containers when stored at 2-8°C. When used as described in Reagent Preparation the reagent is stable for 10 days.

**Instrument/Test Procedures**

For the assay procedure refer to the instrument application sheet enclosed in this kit and to the instrument Operator Manual.

**Specimen collection and preparation**

Serum. Samples may be stored for 3-5 days at 2-8°C or 6 months at -20°C.<sup>2</sup>

**Calibration**

Use ReferrL™ G Calibrator (Cat. No. 0018257000). See calibrator insert sheet for lot specific concentration. Recalibrate the instrument every 15 days or when a new lot of reagent is used. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

**Quality Control**

Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice. Normal and abnormal controls such as SeraChem® Control Level 1 (Cat. No. 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program. Refer to Westgard *et al.*<sup>3</sup> for identification and resolution of out-of-control situations.

**Calculation of Analytical Results**

Refer to the Application Sheet and to the ILab 300 Plus Operator Manual. To convert mg/dL to SI units (μmol/L) multiply by 59.5.<sup>2</sup>

**Limitations/Interfering Substances**

No interference from lipemia up to sample absorbance of 5.9/cm at 660 nm (1000 mg/dL or 11 mmol/L triglycerides). No interference up to 250 mg/dL (0.150 mmol/L) hemoglobin. No interference up to 25 mg/dL (425 μmol/L) bilirubin. For a comprehensive review of interfering substances, refer to the publication by Young.<sup>4</sup>

**Reference Interval**

Serum  
Male: 3.5 to 7.2 mg/dL (208 – 430 μmol/L).<sup>2</sup>  
Female: 2.6 to 6.0 mg/dL (155 – 360 μmol/L).<sup>2</sup>  
It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

**ENGLISH - Insert revision 02/2004**

**Performance Characteristics**

**Precision**

Typical precision results obtained according to NCCLS protocol <sup>5</sup>.

	Samples/Runs	Mean (mg/dL)	CV (%)	Mean (mg/dL)	CV (%)
Within Run	5/10	6.0	1.1	10.5	1.2
Total	5/10	6.0	4.7	10.5	3.4

**Method Comparison**

Comparison Method (x)	Same reagent
Comparison Instrument	ILab 600
Slope	1.038
Intercept	0.04
Range (mg/dL)	3.1 – 9.4
Mean x (mg/dL)	5.4
Mean y (mg/dL)	5.6
r	0.996
SE	0.13
n	56

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**Bibliography / Bibliografía / Bibliografia / Βιβλιογραφία**

1. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, Third edition 1999, pp. 1245-1250.
2. Tietz, NW, ed., Clinical Guide to Laboratory Tests, 3<sup>rd</sup> Edition, WB Saunders Company, Philadelphia, PA, 1995, p. 624.
3. Westgard J., and Barry P., Cost- Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, D.C., 1986.
4. Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, D.C., 1990.
5. Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

Description: URIC ACID	Reference Range							
	LOW VALUES				HIGH VALUES			
Unit: mg/dL (µmol/L)	Male: 0.2 (12)	0.2 (12)	0.2 (12)	3.5 (208)	7.2 (430)	20 (1190)	20 (1190)	20 (1190)
Decimals: 1 (0)	Female: 0.2 (12)	0.2 (12)	0.2 (12)	2.6 (155)	6.0 (360)	20 (1190)	20 (1190)	20 (1190)
LIS Code: **	Children: 0.2 (12)	0.2 (12)	0.2 (12)	**	**	20 (1190)	20 (1190)	20 (1190)
Unit Factor: 1.0		Low Alert	Very Low	Low	Normal Values	High	Very High	High Alert
Slope: 1.00	Rerun: No	No					No	No
Intercept: 0.00								

		Parameter					
Reaction Type: End Point		Predilut.→	S.+R.1→	Reag.2→	Reag.3→	Incubation→	Read
Direction: Up	Times (sec):		0	0	0	458	0
E.P. Limit (abs): 0.0200	Dil./Rgt. Code:		URIC	0			) = Kinetic
Depl Limit (abs): N/A	Lot Number:					Filter 1 (nm): 510	
First Limit (abs): N/A	Ratio/Vol. (ul):	1/1	300	0	0	Filter 2 (nm): 620	
Linear Factor: N/A	Rinse (ul):		0	0	0	Bicr. Factor: 1.00	
Fit: N/A	Sample (ul):		20				

Lin Limit. Low: 0.2 (12) RBL Min (abs): 0.0000  
High: 20 (1190) Max (abs): 1.5000  
Rerun when Over: 1/2

Calculation Model: Standard RBL Stability (days): 1  
Factor: N/A Calibration Stab. (days): 15  
Sample Blank: No Dinamic Controls (min): None

\*\* Operator definable data / Definibile dall'operatore /  
Datos definidos por el operador / Δεδομένα Χειριστού  
N/A Not applicable to this test / Non applicabile a questo test /  
No aplicable a este método / ΔΕΝ ΕΦΑΡΜΟΖΕΤΑΙ ΣΕ ΑΥΤΟ ΤΟ ΤΕΣΤ

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Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.