
Combur⁹ Test

REF 04510038 191

REF 04510046 040

cobas[®]

Σ 50

Σ 100

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Use fresh urine that has not been centrifuged. The urine specimen should not stand for more than 2 hours before testing. In case of longer standing, mix before use.

Use only clean, well-rinsed vessels to collect urine.

Do not add preservatives to the urine.

Materials provided

- REF 04510038191, package with 50 test strips
- REF 04510046040, package with 100 test strips

Materials required (but not provided)

- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document.

1. Use fresh urine that has not been centrifuged. Thoroughly mix the urine sample. The sample should be at room temperature when the test is performed and should not have been standing for more than 2 hours.

2. Take a test strip out of the container. Close the container again with the original desiccant stopper immediately after removal of the strip. This is important as otherwise the test areas may become discolored due to environmental influences such as moisture or nitrite gases in the air. Incorrect results may be obtained.

3. Briefly (about 1 second) dip the test strip into the urine making sure that all test areas are moistened.

4. When withdrawing the test strip, wipe the edge against the rim of the vessel to remove excess urine.

5. After 60 seconds (60-120 seconds for the leukocyte test area), compare the reaction colors of the test areas with the colors on the label and assign always the value of the nearest color block.

▪ Compare the 9th (blood) test area with both color scales as separate color scales are given for erythrocytes and hemoglobin.

Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

Quality control

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

For quality control, use commercially available urine controls, or other suitable control material.

Important note for reporting results (for professional users).

According to the regulations from the German Medical Association for quality assurance of medical laboratory analyses dated 11/23/2007, the decision to classify a laboratory test result to either part B1 or B2 depends on the way the test results are expressed in the report (scale level).

The specification in the report defines whether a determination is quantitative or qualitative and therefore which legal requirements for quality assurance (B1 for quantitative or B2 for qualitative) need to be followed. Examples of qualitative characteristics include titer levels, concentrations/color ranges (+ to +++) or a defined range of values. A characteristic of a quantitative value is when the value has a corresponding measured unit value.

Limitations - interference

Leukocytes: Formaldehyde (stabilizer) and medication with imipenem, meropenem and clavulanic acid may cause false-positive reactions. If the urine specimen has a pronounced intrinsic color (for example due to the presence of bilirubin or nitrofurantoin), the reaction color may be intensified due to an additive effect. Urinary protein excretions in excess of 500 mg/dL and urinary glucose excretions in excess of 1 g/dL may diminish the intensity of the reaction color, as can cephalaxin and gentamicin if administered in high daily doses, or boric acid if used as a preservative.

Nitrite: Prolonged urinary retention in the bladder (4-8 hours) is essential in order to obtain an accurate result. Administration of antibiotics or chemical drugs should be discontinued 3 days before the test. Large amounts of ascorbic acid decrease the sensitivity of the test. **Attention:** Nitrogen oxides present in the atmosphere may have an influence on the stability of the nitrite test pad.

Protein: False-positive readings may be found after infusion of polyvinylpyrrolidone (blood substitute), or when the urine collection vessel contains chlorhexidine or traces of disinfectants possessing quaternary ammonium groups.

Glucose: The effect of ascorbic acid has been largely eliminated so that at glucose concentrations of 100 mg/dL and above even high ascorbic acid concentrations are not likely to give false-negative results.

Ketone bodies: Phenylketones and phthalain compounds produce red colors on the test patch; they are, however, quite distinct from the violet colors produced by ketone bodies and can lead to false-positive results. Captopril, mesna (2-Mercaptoethanesulfonic acid sodium salt) and other substances containing sulphydryl groups may produce false-positive results.

Urobilinogen: Nitrite concentrations above 5 mg/dL or formaldehyde (stabilizer) above 200 mg/dL may cause a decrease in the color reaction.

Bilirubin: Large amounts of ascorbic acid lower the sensitivity of the test.

Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these two parameters.

Drugs that turn red in an acid environment (e.g. phenazopyridine) may produce false-positive readings or reddish colorations on the test patches for nitrite, protein, urobilinogen and bilirubin.

Blood: The values appearing on the printout refer to intact erythrocytes. At concentrations of about 5-50 Ery/ μ L, significant hemolysis (such as may occur on prolonged standing of the urine) leads to values which are higher than the corresponding concentrations given for intact erythrocytes. Ascorbic acid has virtually no effect on the test. In women the test for blood may be falsified from 3 days before to 3 days after a period. It is therefore advisable not to perform the test during this time. After physical activity, e.g. strenuous jogging, raised values for erythrocytes and protein may occur without being signs of disease.

False-positive readings, particularly for erythrocytes, glucose and protein can result from residues of strongly oxidizing disinfectants in the specimen collection vessel.

Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is therefore advisable to repeat the test after discontinuing a particular drug.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

By visual reading

Measuring range

Leukocytes: NEG – ~ 500 LEU/ μ L (3+), **Nitrite:** NEG – POS (1+), **pH:** 5-9, **Protein:** NEG - 500 mg/dL (5 g/L; 3+), **Glucose:** NORM – 1000 mg/dL (55 mmol/L; 4+), **Ketone bodies:** NEG – 150 mg/dL (15 mmol/L; 3+), **Urobilinogen:** NORM – 12 mg/dL (200 μ mol/L; 4+), **Bilirubin:** NEG – ~ 6 mg/dL (100 μ mol/L; 3+), **Blood:** NEG – ~ 250 ERY/ μ L (4+).

Lower limits of measurement

Lower detection limit

Leukocytes: 10-25 LEU/ μ L, **Nitrite:** 0.05 mg/dL (11 μ mol/L), **Protein:** 6 mg albumin/dL, **Glucose:** 40 mg/dL (2.2 mmol/L), **Ketone bodies:** For aceto-acetic acid 5 mg/dL (0.5 mmol/L), **Urobilinogen:** 0.4 mg/dL (7 μ mol/L), **Bilirubin:** 0.5 mg/dL (9 μ mol/L), **Blood:** intact erythrocytes: 5 ERY/ μ L, Hemoglobin or hemolysed erythrocytes: corresponding to 10 ERY/ μ L.

Accuracy

Leukocytes: ≥ 90 % compared with counting chamber, **Nitrite:** ≥ 90 % for 10⁷ gram-positive organisms compared with Griess test, **pH:** ≥ 95 % compared with pH-meter, **Protein:** 90 % compared with radial immuno-diffusion, **Glucose:** ≥ 90 % compared with hexokinase method, **Ketone bodies:** ≥ 85 % compared with photometric enzymatic determination of acetate, **Urobilinogen:** ≥ 95 % compared with Watson & Henry method, **Bilirubin:** ≥ 85 % compared with total bilirubin determination by Jendrassik's method (direct bilirubin), **Blood:** ≥ 90 % compared with counting chamber.

Expected values

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

For further information, please refer to the appropriate operator's manual for the instrument concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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