IL Test™ CK
Cat. No. 0018257340 - 6 x 24 mL CK R1 / 3 x 12 mL CK R2
Cat. No. 0018482100 - 8 x 10 mL CK R1 / 4 x 5 mL CK R2

Intended Use

IL Test™ CK is intended for the quantitative in vitro diagnostic determination of Creatine Kinase in human serum, using the ILab Chemistry Systems.

Summary

Creatine kinase activity measurements are used in the diagnosis and treatment of myocardial infarction and of skeletal muscle disease such as muscular dystrophy, especially in Duchenne type.¹

Principle

Rate analysis according to IFCC recommendations:²

Creatine phosphate + ADP

Creatine Kinase

Creatine + ATP

ATP + Glucose

Hexokinase

glucose-6-phosphate + ADP

glucose-6-phosphate + NADP⁺

G-6-PDH

gluconate-6-phosphate + NADPH + H⁺

The rate of increase in absorbance, due to the NADPH, is directly proportional to the Creatine Kinase activity in the sample.

Absorbance measurements are taken at 340 nm.

Composition

CK R1 containing: Imidazole buffer pH 6.5 60 mmol/L, Glucose 27 mmol/L, N-acetyl cysteine 27 mmol/L, Hexokinase > 5 KU/L, NADP 2.7 mmol/L, Magnesium acetate 14 mmol/L, EDTA-Na₂ 2 mmol/L, and less than 0.1% sodium azide.

CK R2 containing: Imidazole buffer 160 mmol/L, Creatine phosphate 160 mmol/L, ADP 11 mmol/L, AMP 28 mmol/L, Glucose-6-phosphate dehydrogenase > 14 KU/L, Diadenosine Pentaphosphate 55 µmol/L, EDTA-Na₂ 2 mmol/L, and less than 0.1% sodium azide.

Contains sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedure.

PRECAUTIONS AND WARNINGS

CK R1 and CK R2 are not classified as dangerous.

Hazard class: none

Risk phrases: none

Safety phrases: none


This product is For *in vitro* Diagnostic Use.

Reagent Preparation

The reagents are ready to use.

Cat. N°0018482100: for ease of handling uncap the bottles and screw onto them the appropriate device (Cat. 0014557200). Place the bottles, connected to the devices, in the ILab 600/650/Taurus/Aries reagent tray.

See the drawing below.

Reagent Storage and Stability

Unopened reagents are stable until the expiration date indicated on the containers when stored at 2-8°C.

On board stability: 15 days.

Instrument/Test Procedures

For the assay procedure refer to the relevant instrument Application Manual and to the relevant Instrument Operator Manual.

Specimen collection and preparation

Serum.

Samples should be analyzed immediately or refrigerated as soon as possible after collection. Protect from light. Specimens are stable for 1 day at 2-8 °C or 1 month at -20°C.³

Calibration

No calibration is required for this test, the enzyme activity is measured by the rate of increase in NADPH concentration (ΔA/min), using a molar absorption coefficient at 340 nm of 6.3x10³ per 1 cm pathlength.

Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Quality Control

Each laboratory should establish its own mean and standard deviation and adopt a quality control program to monitor laboratory testing. Controls should be analyzed in accordance with good laboratory practice.

Normal and abngmal controls such as SeraChem® Control Level 1 (Cat. No 0018162412) and SeraChem® Control Level 2 (Cat. No. 0018162512) are recommended for a complete quality control program.

Refer to Westgard *et al.*⁴ for identification and resolution of out-of-control situations.

Calculation of Analytical Results

The ILab calculates the creatine kinase activity as follows:

Enzymatic activity (U/L) = ΔA/min x Factor

Where:

ΔA/min: Change in milliabsorbance per minute at 340 nm

Factor: 4.127

The factor is obtained by taking into account: total reaction volume, sample volume, NADPH molar absorption coefficient at 340 nm of 6.3x10³ per 1 cm pathlength. The absorbance measurement in the ILab is automatically converted to 1 cm pathlength by the instrument.

For more details refer to the relevant Application Sheet and Operator Manual.

Activity is expressed in International Units/Liter (U/L), where International Unit is defined as: enzyme activity that converts 1 µmol of substrate in 1 minute.

To convert U/L to nKat/L multiply by 16.67.

Reference Interval

Serum

Male: 38 - 171 U/L (633 - 2850 nKat/L).^{1,2}

Female: 26 - 145 U/L (433 - 2420 nKat/L).^{1,2}

It is recommended that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Performance Characteristics

All performance data were obtained on ILab 600/650. For performance data on other ILab systems refer to the relevant application Manual.

Limitations/Interfering Substances

No interference from lipemia up to sample absorbance of 10/cm at 660 nm (2000 mg/dL or 22.6 mmol/L triglycerides).

No interference up to 300 mg/dL (0.18 mmol/L) hemoglobin.

No interference up to 14 mg/dL (240 µmol/L) bilirubin.

For a comprehensive review of interfering substances, refer to the publication by Young.⁵

Precision

Typical precision results obtained according to NCCLS protocol.⁶

	Samples/Runs	Mean (U/L)	CV (%)	Mean (U/L)	CV (%)
Within Run	5/10	129	1.1	654	0.5
Total	5/10	129	3.5	654	1.8

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- Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, PA, Third edition 1999, pp. 657-664.
- IFCC Primary Reference Procedures for the Measurements of Catalytic Activity Concentration of Enzymes at 37°C. Part 2: Reference Procedure for the Measurements of Catalytic Concentration of Creatine Kinase. Clin Chem Lab Med 2002, 40(6):635-642, 2002.
- Tietz, NW, ed., Clinical Guide to Laboratory Tests, 3rd Edition, WB Saunders Company, Philadelphia, PA, 1995, p. 180.
- Westgard J., and Barry P., Cost- Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, Washington, D.C., 1986.
- Young D.S., Effect of Drugs on Clinical Chemistry Laboratory Tests, AACC Press , Washington, D.C., Fifth edition, 2000, Volume 1, pp. 3/234-3/239.
- Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

IL Test™, SeraChem®, ReferrILTM are trademarks of Instrumentation Laboratory Company.

ILab 600/650 Test Parameters		
Photometric Test Parameters		Serum
Test No.		23
Test Name, Test Code		CK, CK
Sample Type		Serum
Reporting Unit, Decimal Points		U/L (nKat/L), 0 (0)
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		Rate, 20/32
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		340
Sampling Conditions		
Sampling 1	Sample Vol.	8
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	8
	Sample/Diluent Vol.	10/90
Sampling 3 and Sampling 4		***
Diluent Code		Water
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
HIGH!, ABS!		Sampling 2
Sample Volume Reduction		**
Reagent Volumes		
R1	Code	01231
	Rgt/Dil. Vol., Stirring	160/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	15
R2	Code	01232
	Rgt/Dil. Vol., Stirring	40/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	15

Ranges and Evaluation Criteria	Serum
Normal Range-Male	38 - 171 (633 - 2850)
Normal Range-Female	26 - 145 (433 - 2420)
Normal Range-Other	**
Valid Range	1 - 1000 (17 - 16670)
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 2500
Prozone Limit	N/A
Non Linear Limit	40
Slope/Intercept Correction	1/0
Qualitative Report	OFF
Calibration Conditions	
Calibration	Factor: 4.127 (68.8), 3 Reps
Stability (days)	N/A
Calibrator, Concentration	N/A
R-Blank Limit (mAbs)	***
Cal. Reps Range (%)	N/A
Min. Cal. Response (mAbs)	N/A
Cal. Factor Change (%)	N/A
M-Point Curve Fit (%)	N/A
Reagent Blank	ON
Auto R-Blank by Bottle	ON

* See ReferrIL G package insert for appropriate value / Vedere l'inserto del ReferrIL G per il valore appropriato/ Ver en el prospecto del ReferrIL G

** Operator definable data / Definibile dall'operatore / Datos definidos por el operador

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ATTENZIONE! Materiale di proprietà della IL S.p.A. che non deve essere assolutamente modificato o danneggiato.