BAUSCH+LOMB

EN

BAUSCH+LOMB

FR

BAUSCH+LOMB

DE

BAUSCH+LOMB

IT

BAUSCH+LOMB

ES

BAUSCH+LOMB

PT

BAUSCH+LOMB

DA

BAUSCH+LOMB

NL

BAUSCH+LOMB

NO

BAUSCH+LOMB

IS

BAUSCH+LOMB

SV

BAUSCH+LOMB

EL

INSTRUCTION FOR USE PHACO PACKS

⚠ These instructions for use are not intended to substitute for the necessity of reading and understanding the owner’s manual on the Bausch & Lomb Microsurgical System. The owner’s manual provided with the Bausch & Lomb Microsurgical System includes in-depth materials intended to familiarize the surgical team with the proper operation of the equipment.

INDICATIONS FOR USE

The Phaco Packs are intended for use with the Bausch & Lomb Stellaris® System and Stellaris PC systems for the phacoemulsification of an opacified crystalline lens during anterior segment surgery.

SET-UP INSTRUCTIONS

• Proper aseptic delivery and usage of components during setup and procedure is required to assure patient safety.

- ⚠ Balanced Salt Solution wetting of the sleeve prior to assembly is advised to assist assembling the infusion sleeve to the needle.

• When appropriate, prime the system per the owner’s manual to remove any air bubbles assuring proper operation of equipment and accessories. Verify proper tubing connection as indicated below.

Tableting Components	Color Scheme on Tubing Components
Irrigation Lines	Green or Green Color Stripe
Aspiration Lines	Clear

OPERATIONAL INSTRUCTIONS

Operate the system in accordance to the owner’s manual. Bausch & Lomb assumes no responsibility for complications which may arise as a result of improper operation or use of components contained within this pack or available separately.

DISASSEMBLY INSTRUCTIONS

Dispose of contaminated single use products after use in accordance with applicable regulations. Follow cleaning and sterilization instructions for nonincluded limited reuse products prior to use.

⚠ CAUTION: ANTERIOR SEGMENT SURGERY ONLY

• This product should only be used under the guidance of a trained and registered physician.

• Do not re-sterilize and do not reuse.

• Only use this product with Bausch & Lomb products and Bausch & Lomb Microsurgical Systems. Use of this product with non Bausch & Lomb products may affect system performance and create hazards.

• Make certain that the handpiece and any limited reuse accessories are sterilized as specified before use. Allow the handpiece to cool for 20 minutes after sterilization and before use. The electrical connection must be completely dry before connecting to system.

• Make certain to thread the tip onto the handpiece in a clock-wise direction and secure tightly with needle wrench. Improper needle connection could result in damage to equipment.

• A loose needle may lead to improper calibration and could cause shedding of metal fragments into the eye resulting in serious permanent patient injury.

• Do not use any item in this pack if received in an unsatisfactory manner. Customer inquiries or concerns should be promptly directed to Consumer Affairs of Bausch & Lomb or to local sales representative if more convenient.

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© Bausch & Lomb Incorporated.