UNISTRIP™
GENERIC
Blood Glucose Test Strips

Test Strip / 24850, 25520

Intended Use

The UniStrip1™ Test Strips are only for use with the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini® and OneTouch® UltraSmart® meters purchased before October 2012, set at calibration code 49, for measuring glucose (sugar) in whole capillary blood. The UniStrip1™ is meant for self-testing of blood glucose as an aid to monitor the effectiveness of diabetes control.

The UniStrip1™ Test Strip can only be used on the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini® and OneTouch® UltraSmart® meters purchased before October 2012.

They are for single patient use only and should not be shared.

They are used to quantitatively measure glucose in fresh capillary whole blood samples taken from the finger, palm, or forearm.

Testing is done outside the body (in vitro diagnostic use).

They are indicated for use by people with diabetes in their home as an aid to monitor the effectiveness of diabetes control.

Not intended for the diagnosis of or screening for diabetes mellitus and is not intended for use on neonates.

UniStrip1™ Test Strips allow alternate site testing (AST) from the fingertip, palm and/or the forearm. (Use of palm AST is not to be done with OneTouch® Ultra® meter). Alternative site testing should only be done during steady-state times (when glucose is not changing rapidly).

Operating Conditions

UniStrip1™ Operating Temperature	43°-111°F or 6°-44° C
UniStrip1™ Storage Temperature	39° - 86° F or 4° - 30°C (Do Not Freeze)
Relative Humidity	10-90%
UniStrip1™ Strip Vial Open	
Altitude	(10,000 feet) 3,275 meters

Taking care of your test strips:

Storage:

- Store the test strip vial in a cool dry place. Do Not Freeze.
- Keep away from direct sunlight, heat, cold or moisture (water).
- Failure to follow storage guidelines may result in wrong readings.
- Store your test strips in their original vial only.

Handling Your Test Strips:

- Do not take a test strip out of the test strip vial until you are ready to use it.
- Close the test strip vial tightly as soon are you remove a test strip. Use each test strip right away after you remove it from the vial.
- Do Not use test strips that are damaged, wet, or bent. Do not bend, cut or alter the strip in any way.
- You may touch the test strip on its surface with clean dry hands.

Expiration:

- Write the discard date (date opened plus 3 months) on the vial label when you first open it.
- Do Not use test strips beyond the expiration (printed on vial label) or discard date, whichever comes first. Your blood sugar level readings may be wrong if you use expired test strips.

Disposal

- Do Not return the used test strip to the vial.
- Used test strips may be considered biohazardous waste in your area. Be sure to follow your healthcare professional's recommendations for proper disposal.

Accuracy:

- Make sure your meter and test strips are about the same temperature before you test.
- You may not have correct readings if your strips are left in the heat, cold or get wet.
- Test strips are for single use only. Never reuse a test strip that had blood or control solution applied to it.

⚠WARNING:

- Keep the test strip vial away from children. Test strips are a choking hazard.**
- Do not swallow test strips.**
- The test strip vial may contain drying agents that are harmful if inhaled or swallowed and may cause skin or eye irritation.**
- Do not ingest or swallow any items.**

Alternate Site Testing (AST)

- Use AST only two hours or more after a meal, two hours or more after taking insulin, or two hours or more after exercise.
- Do not use AST if you are aware that your glucose level is not as stable as usual.
- Do not use AST if you think you have hypoglycemia (low blood sugar).
- Do not use AST if you think you have hyper glycemia (high blood sugar).
- Do not use AST when you think your blood glucose may be rising or falling quickly.
- Do not use AST if your AST results do not match the way you feel.
- AST value should not be used in insulin does calculations.
- AST values should never be used to calibrate continuous glucose monitors (CGMS).
- Do not use palm alternate site for testing with OneTouch® Ultra® meter.

Blood Glucose Test Procedure

For instructions on performing a blood test (including blood sample collection), refer to the User Guide that came with your system.

⚠CAUTION: Some meters require coding. For these meters, matching the code on the meter and the code on the test strip vial is vital to obtain correct results. Each time you test, check to make sure the code numbers match. For instructions on coding the meter, refer to the User Guide that came with your system.

Test Results

Blood glucose management requires the help of a doctor or nurse. Together you can set your own range of expected blood glucose values, arrange your testing times, and discuss the meaning of your blood glucose results.

Low Glucose Values

If your test result is below 20 mg/dL (1.1 mmol/L), a warning message will appear showing a low glucose level. This may be a sign of severe low blood glucose. Treat this condition right away, as directed by your doctor or nurse. Although this message could be due to a test error, it is safer to treat first, and then do another test.

High Glucose Values

If your test result is above the maximum value, repeat your test with a new strip to confirm this reading. If the high results should appear again, call your doctor or nurse right away.

***Reference values:** Please consult your physician to establish your target ranges. The American Diabetes Association (ADA) suggests the following targets for most nonpregnant adults without diabetes. More or less stringent goals may be appropriate for each individual.*

Glycemic control Time of Day ADA Suggested Reading

Preprandial plasma glucose 70–100 mg/dl (before a meal)

Postprandial plasma glucose <140 mg/dl (after a meal)

Source: American Diabetes Association. Standards of medical care in diabetes - 2012. 2012;35 (supp.1 Diabetes Care):S11-563.

If You Get Unexpected Results

If your blood glucose result is below 70 mg/dL (3.9 mmol/L), showing low blood glucose, or above 180 mg/dL (10.0 mmol/L), showing high blood glucose, you should contact and follow your doctor or nurse's treatment advice. If you keep on getting unexpected results, check your system with UniStrip™ control solution. If you are feeling symptoms that are not consistent with your blood glucose test results AND you have followed all instructions described in your User Guide, call your health care professional. Never ignore symptoms. Never make major changes to your diabetes control program without speaking to your healthcare professional.

Checking the System

A control solution test is performed to check that the meter and test strips are working together properly and that you are doing the test correctly. For instructions on how and when to check the system by doing a control solution test, refer to the User Guide that came with your system.

Use only UniStrip™ Control Solution and follow the instructions on the package insert.

Limitations of Procedure

UniStrip1™ Test Strips give correct results when the following limitations are observed:

- DO Not use for screening purposes.
- Do Not use for the diagnosis of diabetes or for testing of neonates.
- Test strips are for single use only. Do Not reuse. Do not share test strips.
- The test strips are specific to D-glucose and do not react to other sugars, which may be present in blood.
- Use only fresh capillary whole blood. Do Not use serum or plasma.
- Hematocrit is the percentage of red blood cells in the blood. Extremes in hematocrit may affect test results. Hematocrit levels below 30% may cause falsely high readings. Hematocrit levels over 55% may cause falsely low readings. If you do not know your hematocrit level, consult your doctor or nurse.
- UniStrip1™ Test Strips may be used at altitudes up to 10,000 feet (3048 meters) without an effect on test results. *Healthcare professionals- please note these extra limitations of procedure:*
- Do Not use other anticoagulants or preservatives.
- Interferences: Acetaminophen, salicylates, uric acid, ascorbic acid (vitamin C), and other reducing substances (when occurring in normal blood or normal therapeutic concentrations) do not affect results by more than 10%. However, more than normally high concentrations in blood may cause inaccurately high results.
- Patients going through oxygen therapy may yield falsely low results.
- Test results may be falsely low if the patient is severely dehydrated, in shock, or in a hyperosmolar state (with or without ketosis). These test strips should not be used with meters to test critically ill patients.
- Lipemic samples: Cholesterol levels up to 700 mg/dL (18.1 mmol/L) and triglycerides up to 3000 mg/dL (33.9 mmol/L) do not affect the results. Grossly lipemic patient samples have not been tested. They are not advised for testing with any OneTouch® Ultra® Meters.

Test Principle

The OneTouch® Ultra® Meters are plasma-calibrated to allow easy comparison of results with laboratory methods. Glucose in the blood sample mixes with special chemicals on the UniStrip1™ Test Strip and a small electrical current is produced. This current is measured by the OneTouch® Ultra® Meters and displayed as your blood glucose result. The strength of the current changes with the amount of glucose in the blood sample.

Reagent Composition

Each UniStrip1™ Test Strip contains: Glucose oxidase (Aspergillus niger) 20IU ; Potassium ferricyanide 0.12mg ; Non-reactive ingredients 1.8mg

Performance Characteristics

The performance of UniStrip1™ Test Strips has been tested both in laboratory and in clinical tests.

Measurement Range:

The UniStrip1™ with the OneTouch® Ultra® Family of Meters measures from 20 to 600 mg/dL (1.1-33.3mmol/L).


Introduction to Accuracy Study

To assess the accuracy of the UniStrip1™ Test Strips, a clinical study was performed using two each of the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and three different lots of UniStrip1™ Test Strips. Different glucose concentrations were included in the study, ranging from 47 mg/dL to 456 mg/dL. Readings were taken with the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and the YSI 2300 blood glucose analyzer. Each OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meter performed one test at each glucose concentration with each lot of test strips.

Accuracy Study Results

OneTouch® Ultra® Blood Glucose Meter: Finger capillary whole blood					
For glucose concentration<75 mg/dL					
Within±5 mg/dL (Within±0.28 mmol/L)	Within±10 mg/dL (Within±0.56 mmol/L)	Within±15 mg/dL (Within±0.83 mmol/L)			
11/18	(61%)	17/18	(94%)	18/18	(100%)


For glucose concentration ≥ 75 mg/dL							
Within±5%	Within±10%	Within±15%	Within±20%				
38/82	(46%)	71/82	(87%)	79/82	(96%)	82/82	(100%)

Regression Analysis of OneTouch® Ultra® Blood Glucose Meter

OneTouch® Ultra®2 Blood Glucose Meter: Finger capillary whole blood					
For glucose concentration<75 mg/dL					
Within±5 mg/dL (Within±0.28 mmol/L)	Within±10 mg/dL (Within±0.56 mmol/L)	Within±15 mg/dL (Within±0.83 mmol/L)			
11/18	(61%)	17/18	(94%)	18/18	(100%)


For glucose concentration ≥ 75 mg/dL							
Within±5%	Within±10%	Within±15%	Within±20%				
40/82	(49%)	70/82	(85%)	80/82	(98%)	82/82	(100%)

Regression Analysis of OneTouch® Ultra®2 Blood Glucose Meter

OneTouch® UltraMini® Blood Glucose Meter: Finger capillary whole blood					
For glucose concentration<75 mg/dL					
Within±5 mg/dL (Within±0.28 mmol/L)	Within±10 mg/dL (Within±0.56 mmol/L)	Within±15 mg/dL (Within±0.83 mmol/L)			
11/18	(61%)	17/18	(94%)	18/18	(100%)


For glucose concentration ≥ 75 mg/dL							
Within±5%	Within±10%	Within±15%	Within±20%				
42/82	(51%)	66/82	(80%)	80/82	(98%)	81/82	(99%)

Regression Analysis of OneTouch® UltraMini® Blood Glucose Meter

OneTouch® UltraSmart® Blood Glucose Meter: Finger capillary whole blood					
For glucose concentration<75 mg/dL					
Within±5 mg/dL (Within±0.28 mmol/L)	Within±10 mg/dL (Within±0.56 mmol/L)	Within±15 mg/dL (Within±0.83 mmol/L)			
10/18	(56%)	17/18	(94%)	18/18	(100%)

For glucose concentration ≥ 75 mg/dL							
Within±5%	Within±10%	Within±15%	Within±20%				
42/82	(51%)	71/82	(87%)	81/82	(99%)	82/82	(100%)


Regression Analysis of OneTouch® UltraSmart® Blood Glucose Meter


Introduction to Linearity Regression Study

To assess the linearity regression of the UniStrip1™ Test Strips, a clinical study was performed using five each of the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and three different lots of UniStrip1™ Test Strips. Different glucose concentrations were included in the study, ranging from 15 mg/dL to 620 mg/dL. Readings were taken with the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and the YSI 2300 blood glucose analyzer. Each OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meter performed one test at each glucose concentration with each lot of test strips.


Linearity Regression Analysis Results					
	Slope	Y-Intercept	R ²	r	Audit
OneTouch® Ultra® Blood Glucose Meter:	0.9709	2.0425	0.9979	0.9989	pass
OneTouch® Ultra®2 Blood Glucose Meter:	0.9704	1.9647	0.9974	0.9986	pass
OneTouch® UltraMini® Blood Glucose Meter:	0.9699	2.7701	0.9972	0.9985	pass
OneTouch® UltraSmart® Blood Glucose Meter:	0.9796	0.5354	0.9979	0.9989	pass


Linearity regression analysis of OneTouch® Ultra® Blood Glucose Meter

Linearity regression analysis of OneTouch® UltraMini® Blood Glucose Meter

Linearity regression analysis of OneTouch® Ultra®2 Blood Glucose Meter

Linearity regression analysis of OneTouch® UltraSmart® Blood Glucose Meter

Introduction to Precision Study

To assess the precision of the UniStrip1™ Test Strips, a clinical study was performed using five each of the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and three different lots of UniStrip1™ Test Strips. This study is designed to determine repeatability (Within-Run) and intermediate precision (Between-Run). The Within-Run tests were performed in one day at five blood glucose concentrations ranging from 48 mg/dL to 319 mg/dL. The Between-Run tests were performed over ten days at three blood glucose concentrations ranging from 46 mg/dL to 298 mg/dL. Readings were taken with the OneTouch® Ultra®, OneTouch® Ultra®2, OneTouch® UltraMini®, and OneTouch® UltraSmart® meters, and the YSI 2300 blood glucose analyzer. The FDA requires that the SD must be < 5 mg/dL at glucose concentrations < 75 mg/dL, and the CV must be < 5% for glucose concentrations > 75 mg/dL. Additionally, the FDA requires that >95% of results fall within Zone A of the Error Grid Analysis, and <5% of results fall within Zone B.

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