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This procedure gives the weight calibration method by substitution to obtain both real and conventional mass values

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This procedure applies for the mass and weights calibration with nominal value from 1 mg to 200 kg. For this calibration a weighted device is used with the substitution method application and the use of a same nominal value mass.

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Balance: Instrument measures the body apparent mass from the net force on it, when it is in the earth gravity field and there are other perturbations. The net force is the composition of all series of forces that act on the body; between them are the followings:

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Magnetic forces arise from the interaction between the magnetic field around the balance and the weight. These are no quantifiable forces.

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Mass calibration: Set of operations that settle down, under specified conditions, the relationship among the mass value of an unknown mass weight, object of calibration called “sample”, and the mass value of a known mass weight mass called “standard”

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Weight accuracy class: Weights class satisfy determined metrological requisites having the deviations within specified limits

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Density: The body’s density is the ratio between its mass and its volume

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Mass error: Difference between the mass value and its nominal value

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Mass: Is a physical quantity determines a body’s amount of matter

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Reference mass: Mass o mass set used as reference standard to determine the mass of one or other masses (test ones) or weighted instruments.

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Weight: Materialized mass measurement with physical and metrological requirements: form, dimensions, material, surface quality, nominal value and maximum permissible error. [2]

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Additional weight: low mass nominal value weight with respect the test mass that it is used to reduce the differences indicated between the reference mass and the test mass. It is possible to put on additional weights both with the reference mass as well as with the test mass.

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Mass conventional value: For a body at 20 °C, its conventional mass or the mass conventional value resulting from weighing in the air is the reference weight mass with a density of 8 000 kg·m⁻³ "in equilibrium" with it, into an air density of 1,2 kg·m⁻³. [3]

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k : factor de cobertura / factor of coverage

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The substitution method is followed for the mass calibration. This method compares two masses with equal nominal value. One mass is the reference (A) and the other one is the test (B), in order to get the unknown value from the reference value. The steps to follow are different in according to the weighing device used. In every one step is intended to put on the balance's plate in a correlative way the masses A and B so the balance drift error is cancel along cycle. Henceforth, when the reference mass or the test mass is called, it is possible that a weight is called, ever that the mass is standardized.

- For the mono-plate balances the sequence to follow is the double substitution known as A/B/B/A, some times, the cycle can include a witness weight, so the cycle is: A/B/B/T/A. The witness weight is only used to high accuracy calibrations.*
- For the automatic mass comparators, the sequence to follow can be the sequence (a) or as well as $A_1/B_1/A_2/B_2/\dots A_n/B_n/A_{n+1}$*
- For the equal arms balance, a mass with sensitivity m_s is added at the cycle end without remove the reference mass and the sequence to follow is A/B/B/A/A+ m_s . In a similar way, a witness weight can be introduced: A/B/B/T/A/A+ m_s . The witness weight is only used to high accuracy calibrations.*

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
Figura 1

Figure 1

The mass determination's physical principle is to compare the gravity force (F_1) exerted by a test mass (object of calibration) on the balance's load receptor plate with the force exerted on the same sense by a known value mass (reference) with the same nominal value; For this comparison, other forces take part: $\sum \vec{F}_i$ from the electric, magnetic or other any source interactions with the mass. Generally, these interactions should be lower in several orders than the other ones so they can be neglected in

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the every day calibrations, otherwise with the maximum accuracy calibrations or if, it is certain that their order is not neglected.

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The most significant influence force and usually is not neglected is the air's buoyancy (F_2), (this force depends on the weight volume and on the air density).

The equilibrium is achieved as the forces on the body are: $\vec{F}_1$ (weight), $\vec{F}_2$ (air buoyancy) and $\vec{F}_3$ normal force from the balance's plate, as shown the figure 1.

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The equilibrium equation is gived by:

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Balance's reading in balance's units I These values are taken as mass values as the balance's sensitivity is not determined

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By substitution for the forces in equation (4b):

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Equation (10) is the basic weighing equation.

$$(m_p - \rho_{ap} V_p)\vec{g} = (m_r - \rho_{ar} V_r)\vec{g} + \Delta\vec{F}_{pr} \quad (11)$$

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For the comparison time: ρ_{ap} and ρ_{ar} are equal to ρ_a

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Where, $\Delta\vec{F}_{pr}$ is an additional small force needed to compensate both gravity forces exerted on the compared masses. This force depends on the reading difference $\Delta I = I_p - I_r$ between the balance and the gravity value. The equation is: $\Delta\vec{F}_{pr} = f \cdot \Delta I \cdot \vec{g}$. The proportionality factor f , is the inverse of the balance's sensitivity coefficient and is a factor to convert the readings values into mass values.

$$\Delta\vec{F}_{pr} = f \cdot \Delta I \cdot \vec{g} = \frac{m_s - \rho_{as} V_s}{\Delta I_s} \cdot \Delta I \cdot \vec{g} \quad (12)$$

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Where, m_s is the sensitivity weight's mass value. V_s is the sensitivity weight's volume. ΔI_s is the balance reading's variance as the sensitivity weight is put on it.

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By assuming that the masses under comparison are nearly equal to the nominal mass value the sensitivity factor would be the same for both masses, generally f equals 1 (For maximum accuracy calibrations, this is continuously checked).

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By comparing the forces act on m_r y m_p , the equation would be.

This force depends on the balance's reading difference ΔI ($\Delta I = I_p - I_r$)

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The right hand of the equations (13) and (14) is the mass difference between the test mass and the standard weight (reference)

The weights calibration procedure by substitution method allows to calculate the test mass by direct comparison with a standard weight, by means of an established weighted sequence and the repetition of a number of n cycles.

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Several weighing devices will be used, they are described in the Mass Laboratory's inventory, according to their measurement range and accuracy and the nominal value of the weights to be calibrated and the required uncertainty. It is needed to verify that they are perfectly levelness and that there are no vibrations, avoiding external influences that produce interactions disturb the reading.

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The mass standards must be properly calibrated and their updated certificate of calibration

Mass standards set described in the Mass Laboratory's inventory with accuracy class higher than the masses to be calibrated 's one and with the suitable nominal values to make a direct comparison. The standard's drift must never be higher than the standard's uncertainty.

The standard's metrological requirements and physical characteristics must be the ones that are described in the OIML Recommendation R 111-1 [1].

5.1.2 Instrumentos para control de condiciones ambientales / *Instruments for the environmental conditions control*

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*Climatic instruments that are described in the Mass Laboratory's inventory are used to measure the temperature, the pressure and the air humidity (or the dew point temperature) in order to determine the air density. See technical instruction: Calculation of the air density (**LM01-MA-105**).*

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The comparison equipment (mass comparators or balances) are selected according to the mass set to be calibrated. The comparison equipment must be connected to the electric energy power supply at least 24 hours before the calibration. Before the start of the calibration process, verify that there is no dust on the plate, so must do a cleaning by blowing air with the air pear or cleaning with a paintbrush.

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*Before the calibration, the weight set or the weight must be checked in order to determine if there are impurities, if so of easy removing, it will do a surface cleaning. See technical instructive **LM01-MA-IO2**. If the presence of particles is detected during the weight changes, it will blow with the air pear.*

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During the comparison the reference weight will be let on a green plate and the test weight on a red plate, this for the non automatic comparison. It does not apply for the mass comparator because the comparison process is automatic.

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To record the test weight, the standards and the environmental conditions data. These data are directly recorded within the mass comparators' software.

The environmental conditions are recorded to determine the air density, between each measurement cycle in the case of high accuracy weight and only at the beginning and at the end of the calibration in the case of lower accuracy weights.

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Operations sequence

- Environmental conditions data taking, at beginning and at the end for each measurement cycle or at beginning and the end of the calibration depending on the accuracy class*
- Weighing cycle realization according the sequence of the calibration method (items 4.2 and 5.4)*
- n -weighing cycles repetition*
- Real mass value or conventional mass value estimation of the mass under calibration and estimation of their associated uncertainties.*

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*The registers **LM01-MA-F01** or **LM01-MA-F04** must be used to the calibration in the case of data taking by hand. Once upon the environmental conditions and/or the start time of the calibration are recorded, at least four (4) pre-weighing are made by using the test weight, these data are not recorded. The weights must be let in a softly way without hit the plate.*

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It is needed to wait to the reading stabilization when the (reference) standard weight is put on the balance's plate. The reading I_{r1} is recorded, to withdraw the weight and then the test weight is put on, wait for reading stabilization, the data I_{p1} is recorded; to withdraw and put on again the test weight and the reading I_{p2} is recorded. To close the cycle, again put the standard weight and the reading I_{r2} is recorded. For higher accuracy weights, the cycle ends with the witness weight reading I_T . It must no adjust the balance during a cycle. This adjustment can be done between cycles.

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For the calibration of weights Class E_2 repeat five-times the cycle, four-times for the calibration of weights Class F_1 , three-times for the calibration of weights Class F_2 and two ones for the calibration of weights Class M_1 and lower ones. This indicates the minimum number of cycles, which no excludes that a greater number of cycles can eventually be applied.

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For the case of the equal arms balances, a counter-weight with the same nominal value of the weight to be calibrated and the standard weight is put on the left plate. The standard weight is put on the right plate, wait for the reading stabilization, the reading L_{N1} is recorded, the standard weight is withdrawn, and then the test weight is put, wait for the reading stabilization, the reading L_{P1} is recorded; to withdraw and put on again the test weight and the reading L_{P2} is recorded; the test weight is withdraw, again put the standard weight and the reading I_{r2} is recorded. To close the cycle, without blocking the balance, the sensitivity weight (with its nominal value greater than the balance's division of scale) is put and the reading I_{r2+ms} is recorded.

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It must be checked that the environmental conditions are within the tolerance indicated in the technical manual and that the standard deviations are not greater than population's standard deviation.

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For each cycle, it has a difference between the test weight and the standard weigh, I_{p-r} , which one is determined by the following:

From the balance's readings values for every weight are firstly obtained the mean readings.

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Find out the difference between the test and standard weights readings mean values

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According to the test weight's accuracy class and as the process is n -times repeated then the mean difference between the test and standard weights readings is obtained from the arithmetic mean of the n cycles differences.

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As the ABA method is used, the mean value is obtained by the following equation:

$$\bar{I}_{p-r} = \frac{\sum_i^n \left(I_{p_i} - \frac{I_{r_i} + I_{r_{i+1}}}{2} \right)}{n} \quad (18)$$

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The calibration determines the real mass value or the conventional mass value. The calibration is performed by the substitution method directly determining the test conventional mass value or rather its real mass value and from this the conventional mass value.

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By applying the equation (14) to the calibration cycles set, the test real mass value is obtained, since the reference real mass value, its volume and the test volume are known, as well as the air density at the calibration time as shown the equation (19). Taking into account the reference's drift, so for the calibration time, the mass value should include in addition to the nominal value error, the drift from the last calibration. Also, the mean reading difference value must be corrected by the weighing instrument's eccentricity δe , equation (20). As commonly, these values: standard's drift and eccentricity are taken zero as corrections, but non zero as uncertainty components

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The same test weight's mathematical treatment is applied to the witness standard when it is included in the calibration process. So, replace m_p by m_t and the reading difference by $\bar{I}_{t-r}$. For control issues must be observed that the historical difference between the witness standard and the reference standard is maintained.

In terms of mass values, a standard's error denotes to the nominal value deviation, equation (21)

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The Document OIML D – 28: 2004 defines the conventional mass in the following way:

$$m_c \left(1 - \frac{\rho_0}{\rho_{ref}} \right) = m \left(1 - \frac{\rho_0}{\rho} \right) \quad (22)$$

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In the case of this Laboratory, it must be careful with this model application because the air density differs more than 10% of the reference value 1,2 kg·m⁻³, therefore, the proposed correction in the document OIML R111-1:2004 must be done, see numeral 5.4.1.3, where the air density measured in the Mass Laboratory is included.

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To perform a conventional mass calibration for a test weight (whose volume is V_p), must be determined the difference between the test weight conventional mass m_{cp} and the reference weight conventional mass m_{cr} . The reference weight volume is V_r .

From equations (14) and (22), the conventional mass difference is:

$$\Delta m_c = f \cdot (l_p - l_r) + m_{cr} \left(\frac{1}{\rho_r} - \frac{1}{\rho_p} \right) \cdot (\rho_a - \rho_p) \quad (23)$$

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Taking into account, the effects from sensitivity, eccentricity, magnetic susceptibility and standard's drift, results the mathematical model for the test weight conventional mass determination.

The factors are no taken as effective correction factors, but uncertainty components.

$$m_{cp} = \bar{l}_{p-r} + \underbrace{\delta l_{ba} + \delta_e + \delta_{ma}}_{\delta m_{ba}} + \underbrace{(m_{cr} + \delta m_{cr})}_{\substack{\text{masa convencion} \\ \text{del patron} \\ \text{corregida}}} \cdot (1+C) \quad (25a)$$

$$m_{cp} = m_{cr} + \bar{l}_{p-r} + m_{cr} \cdot C + \delta m_r + \delta m_{ba} \quad (25b)$$

$$m_{cp} = m_{cr} + \bar{l}_{p-r} + \delta m_{em} + \delta m_r + \delta m_{ba} \quad (25c)$$

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The buoyancy correction is $\delta m_{em} = (m_{cr} + \delta m_{cr}) \cdot C \cong m_{cr} \cdot C$; however, he mass drift value will be taken for the uncertainty calculation.

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The values of $U(m_{cp})$ are taken from the OIML R 111-1:2004 permissible maximum errors tables. This is valid if the uncertainty value assigned to the conventional mass be for $k=2$.

Rearrange the equation (24) terms, the conventional mass value as a function of the volumes is given in a good approximation by:

$$m_{cp} = m_{cr} + \bar{I}_{p-r} + (\rho_a - \rho_0)(V_p - V_r) + \delta m_r + \delta m_{ba} \quad (27a)$$

$$m_{cp} = m_{cr} + \bar{I}_{p-r} + \delta m_{em} + \delta m_r + \delta m_{ba} \quad (27b)$$

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This mathematical model (25b) or (27b) is in the spreadsheets, weights calibration balances comparators, weights calibration discs comparators and weights calibration in balances.

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The results are reported indicating the nominal mass value plus its error, this one can be related to real mass or conventional mass, and the expanded uncertainty.

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The calibration result is expressed in the certificate of calibration. The calibrated weight' characteristics as manufacturer, labels, form, material, density, environmental conditions during the calibration and the reference weights' SI traceability data are recorded, too.

The results table shows the real mass and the conventional mass displayed as nominal mass and error, the expanded uncertainty with a coverage factor of $k = 2$ for a level of confidence of 95,45%. (Assuming the Central Limit Theorem, see GUM numeral 6.3 [7]).

It will have into account the following factors to classify the weights as OIML weights: form, surface finish, labeled, error, uncertainty, and the error and uncertainty combination as described OIML R 111-1, numeral 5.3.1. [2] (Annex 2 LM01-MA-P01)

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For the uncertainty of measurement calculation, the following documents are the references: Guide EA-4/02 (it corresponds DKD-3) [6], Recommendation OIML R 111 - 1: 2004 and the Colombian Technical Guide GT 51: 1997 – ICONTEC (corresponds GUM [7]).

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From equation (25), the following uncertainties sources are considered:

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- ✓ Reference standard
- ✓ Reference standard's time variation
- ✓ Comparison mean (balances or mass comparator)
- ✓ Comparison process
- ✓ Buoyancy correction

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Applying the uncertainties propagation law to equation (25) and assuming that there is not a significant correlation between the uncertainty components:

$$u^2(m_{cp}) = u^2(m_{cr}) + u^2(\delta m_r) + u^2(\delta m_{em}) + u^2(\bar{l}_{p-r}) + u^2(\delta m_{ba}) \quad (32)$$

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Each equation (32) 's uncertainty source is analyzed.

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7.1.1.1 The reference standard's mass uncertainty, $u(m_r)$, is taken from its certificate of calibration, $U(m_r)$, with the appropriate coverage factor k .

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In case of carries out the calibration with some standards which are correlated, (e.g. if they have been calibrated by the same laboratory using the same equipment or the same standards):

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7.1.1.2 The standard's drift uncertainty, $u(\delta m_r)$, is calculated as follows:

From the standard's successive calibrations, a standard's mass variation Δm_r is calculated. This is the maximum mass difference between two successive calibrations. The drift uncertainty is:

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If there is no a historic record of the standard's successive calibrations, the maximum drift uncertainty is given by:

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From the proposed equation in OIML R111-1:2004 [2] y [8]:

$$u(\delta m_{em}) = \sqrt{\left[m_{cr} \frac{(\rho_p - \rho_r)}{\rho_p \rho_r} u(\rho_a) \right]^2 + [m_{cr} (\rho_a - \rho_0)]^2 \frac{u^2(\rho_p)}{\rho_p^4} + m_{cr} (\rho_a - \rho_0) [(\rho_a - \rho_0) - 2(\rho_{a1} - \rho_0)] \frac{u^2(\rho_r)}{\rho_r^4}} \quad (37)$$

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*The air density uncertainty is calculate from the technical instructive **LM01-MA-I05 - Air density calculation**, taking into account the uncertainty of climatic instruments used in the measurements.*

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The values of test weights ' density and their uncertainty are getting from the certificates of calibration for density or volume; otherwise from the manufacturer's specifications or from the tabulated density values for the weight material (Annex 3 LM01-MA-P01). If there is a values interval where it is possible to found the density value, $[\rho_{p \max}; \rho_{p \min}]$, the uncertainty is estimated from this interval as follows:

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The volume expansion uncertainty due the temperature variation is neglected.

- From the certificates of calibration, reference weights' values and their uncertainties of density, for the corresponding coverage factor.

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7.1.1.4 Difference of reading $u(\bar{I}_{p-r})$ uncertainty. Uncertainty component type A.

The uncertainty component obtained from the estimation of the population variance s_p^2 for an instrument under statistical control and for a large number of essays is considered of the type A. This parameter characterizes the dispersion better than the standard deviation from a limited number of observations, $s(\bar{I}_{p-r}) = s(\Delta m)$. In this case:

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Where, n is the number of cycles as it is described in the numeral 5.4. If s_p is unknown then it will be taken the value of $s(\bar{I}_{p-r})$ instead.

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7.1.1.5. Uncertainty component due to the comparison instrument (balance or mass comparator)

Applying the propagation of uncertainties law to equation (16) and considering that the variables I_p and I_r are correlated since they are taken by the same instrument, the uncertainty component will be:

$$u^2(\delta m_{ba}) = u_B^2(I_p) + u_B^2(I_r) - 2u_B(I_p)u_B(I_r) \cdot r((I_p);(I_r)) \quad (43)$$

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The uncertainties $u_B(I_p)$ and $u_B(I_r)$ have mainly three components, one component due to the instrument's resolution, $u(d)$, other component due to eccentricity effects and other component due to the magnetic effects. The eccentricity effect are considered small significant since the mass comparators have automatic centered systems or the eccentricity effect is no applicable; otherwise, this component must be taken into account. In a similar way, the magnetic effect is no considered because standardized weights are used and they satisfy the magnetic requisites. (It may be considered that its effect is included within the standard deviation.) The covariance between $u_B(I_p)$

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and $u_B(I_r)$ duplicates the instrument's resolution uncertainty component and nulls the other uncertainty components for the instrument's calibration. Therefore:

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To calculate $u(d)$, a rectangular distribution with amplitude d is taken.

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For the real mass case, the drift components of the standard, $u(\delta m_r)$, of the balance or mass comparator, $u(\delta m_{ba})$ and of the reading difference $u(\bar{I}_{p-r})$ do not change, meanwhile for the standard must be take the real mass uncertainty $u(m_r)$ instead of the conventional mass uncertainty $u(m_{cr})$.

For the buoyancy, it is taken the equation:

$$u(\delta m_{em}) = \sqrt{(V_p - V_r)^2 u^2(\rho_a) + \rho_a^2 (u^2(V_p) - u^2(V_r))} \quad (46)$$

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The test volume uncertainty component $u(V_p)$ is taken from the certificate of calibration, $U(V_p)$, with the appropriate coverage factor k .

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The standard volume uncertainty component $u(V_r)$ is taken from the certificate of calibration, $U(V_r)$, with the appropriate coverage factor k

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For any case, the uncertainty associated to the volume expansion by temperature variation is neglected.

7.2 Resumen de componentes de incertidumbre / Uncertainty components summary

The following mathematical model is applied once the uncertainties sources are identified. The data are replaced in the model to get the combined uncertainty.

$$u^2(m_{cp}) = u^2(m_{cr}) + u^2(\delta m_r) + u^2(\delta m_{em}) + \frac{s_p^2}{n} + \frac{d^2}{6} \quad (49)$$

COMPONENTE component x_i	ESTIMADO estimated x_i	INCERTIDUMBRE ESTANDAR Standard uncertainty $u(x_i)$	DISTRIBUCIÓN DE PROBABILIDAD Probability distribution	COEFICIENTE DE SENCIBILIDAD Sensitivity coefficient	CONTRIBUCIÓN contribution $u(y)$	CONTRIBUCIÓN CUADRÁTICA Square contribution $u(y)$
PATRON Standard (m_{cr})			normal	1		
VARIACION DEL PATRON Standard's drift (δm_r)			rectangular	1		
EMPUJE (em)			normal	1		
COMPARACION Comparison (S_p)			normal	1		
MEDIO DE COMPARACION Comparison mean (d)			rectangular	1		

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The nominal value uncertainty is null, $u(m_n)$, because it is a constant. value.

Expanded uncertainty calculation

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If the uncertainty value for the standard's real mass is no significantly different from the uncertainty value for the conventional mass, it is considered that the test's mass uncertainty value and the real mass no significantly change and hence the same value is expressed, because the other terms do not produce significant changes. [9]

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* ΔT = Initial difference between weight temperature and laboratory temperature.

Anexo 2 / Annex 2

Table 1. Maximum permissible errors for weights ($\pm \delta w$ in mg)

Nominal value*	Class E ₁	Class E ₂	Class F ₁	Class F ₂	Class M ₁	Class M _{1,2}	Class M ₂	Class M _{2,3}	Class M ₃
5 000 kg			25 000	80 000	250 000	500 000	800 000	1 600 000	2 500 000
2 000 kg			10 000	30 000	100 000	200 000	300 000	600 000	1 000 000
1 000 kg		1 600	5 000	16 000	50 000	100 000	160 000	300 000	500 000
500 kg		800	2 500	8 000	25 000	50 000	80 000	160 000	250 000
200 kg		300	1 000	3 000	10 000	20 000	30 000	60 000	100 000
100 kg		160	500	1 600	5 000	10 000	16 000	30 000	50 000
50 kg	25	80	250	800	2 500	5 000	8 000	16 000	25 000
20 kg	10	30	100	300	1 000		3 000		10 000
10 kg	5.0	16	50	160	500		1 600		5 000
5 kg	2.5	8.0	25	80	250		800		2 500
2 kg	1.0	3.0	10	30	100		300		1 000
1 kg	0.5	1.6	5.0	16	50		160		500
500 g	0.25	0.8	2.5	8.0	25		80		250
200 g	0.10	0.3	1.0	3.0	10		30		100
100 g	0.05	0.16	0.5	1.6	5.0		16		50
50 g	0.03	0.10	0.3	1.0	3.0		10		30
20 g	0.025	0.08	0.25	0.8	2.5		8.0		25
10 g	0.020	0.06	0.20	0.6	2.0		6.0		20
5 g	0.016	0.05	0.16	0.5	1.6		5.0		16
2 g	0.012	0.04	0.12	0.4	1.2		4.0		12
1 g	0.010	0.03	0.10	0.3	1.0		3.0		10
500 mg	0.008	0.025	0.08	0.25	0.8		2.5		
200 mg	0.006	0.020	0.06	0.20	0.6		2.0		
100 mg	0.005	0.016	0.05	0.16	0.5		1.6		
50 mg	0.004	0.012	0.04	0.12	0.4				
20 mg	0.003	0.010	0.03	0.10	0.3				
10 mg	0.003	0.008	0.025	0.08	0.25				
5 mg	0.003	0.006	0.020	0.06	0.20				
2 mg	0.003	0.006	0.020	0.06	0.20				
1 mg	0.003	0.006	0.020	0.06	0.20				

* The nominal weight values in Table 1 specify the smallest and largest weight permitted in any class of R 111 and the maximum permissible errors and denominations shall not be extrapolated to higher or lower values. For example, the smallest nominal value for a weight in class M₃ is 100 mg while the largest is 5 000 kg. A 50 mg weight would not be accepted as an R 111 class M₃ weight and instead should meet class M₁ maximum permissible errors and other requirements (e.g. shape or markings) for that class of weight. Otherwise the weight cannot be described as complying with R 111.

Anexo 3 / Annex 3

Table 5 Minimum and maximum limits for density ($\rho_{\min}$, $\rho_{\max}$)

Nominal value	$\rho_{\min}, \rho_{\max}$ (10^3 kg m^{-3})							
	Class of weight (for class M_3 , no value is specified)							
	E_1	E_2	F_1	F_2	M_1	M_{1-2}	M_2	M_{2-3}
$\geq 100 \text{ g}$	7.934 – 8.067	7.81 – 8.21	7.39 – 8.73	6.4 – 10.7	≥ 4.4	> 3.0	≥ 2.3	≥ 1.5
50 g	7.92 – 8.08	7.74 – 8.28	7.27 – 8.89	6.0 – 12.0	≥ 4.0			
20 g	7.84 – 8.17	7.50 – 8.57	6.6 – 10.1	4.8 – 24.0	≥ 2.6			
10 g	7.74 – 8.28	7.27 – 8.89	6.0 – 12.0	≥ 4.0	≥ 2.0			
5 g	7.62 – 8.42	6.9 – 9.6	5.3 – 16.0	≥ 3.0				
2 g	7.27 – 8.89	6.0 – 12.0	≥ 4.0	≥ 2.0				
1 g	6.9 – 9.6	5.3 – 16.0	≥ 3.0					
500 mg	6.3 – 10.9	≥ 4.4	≥ 2.2					
200 mg	5.3 – 16.0	≥ 3.0						
100 mg	≥ 4.4							
50 mg	≥ 3.4							
20 mg	≥ 2.3							

Table B.7 Method F2 - List of alloys most commonly used for weights

Alloy/material	Assumed density	Uncertainty ($k = 2$)
Platinum	21 400 kg m^{-3}	$\pm 150 \text{ kg m}^{-3}$
Nickel silver	8 600 kg m^{-3}	$\pm 170 \text{ kg m}^{-3}$
Brass	8 400 kg m^{-3}	$\pm 170 \text{ kg m}^{-3}$
Stainless steel	7 950 kg m^{-3}	$\pm 140 \text{ kg m}^{-3}$
Carbon steel	7 700 kg m^{-3}	$\pm 200 \text{ kg m}^{-3}$
Iron	7 800 kg m^{-3}	$\pm 200 \text{ kg m}^{-3}$
Cast iron (white)	7 700 kg m^{-3}	$\pm 400 \text{ kg m}^{-3}$
Cast iron (grey)	7 100 kg m^{-3}	$\pm 600 \text{ kg m}^{-3}$
Aluminum	2 700 kg m^{-3}	$\pm 130 \text{ kg m}^{-3}$