Selectra Slitter

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Target Group

This manual is intended for cardiologists, electrophysiologists and surgeons who are familiar with, and have practical experience with, cardiac resynchronization therapy and implantation of leads in the coronary venous system.

Description

The Selectra slitter tool is designed for use in conjunction with the BIOTRONIK Selectra CS lead introducer system.

It is used to remove Selectra guiding catheters from CS leads with an outer diameter between 4.6 and 5.8 F and can cut up to three Selectra guiding catheters.

BIOTRONIK's Selectra CS lead introducer system is a combination of CS guiding catheters and implant accessories used to facilitate access to the coronary venous system for suitable leads and catheters.

Intended Use, Indications

Selectra accessories are used in conjunction with the CS lead introducer system to facilitate lead implantation in the left side of the heart via the coronary sinus.

Contraindications

Use of the Selectra CS lead introducer system is contraindicated for the following:

- Patients with an existing or possible occlusion of the coronary vessels or unsuitable anatomy of the coronary veins
- Patients with active systemic infection

No Reuse

The Selectra slitter tool is intended for single-use only.

It must not be resterilized.

Reuse of implant accessories can result in infections, embolisms and damage to the product.

Refer to Other Manuals

Also refer to the manuals for the CS guiding catheters and the accessories kit for the Selectra CS lead introducer system, as well as the manual for the lead used.

Warnings and Precautions

- To inhibit blood clotting, BIOTRONIK recommends intravenously administering 50-100 IU heparin per kg body weight if the patient is not otherwise anticoagulated.
- Only use water soluble contrast media.

Potential Negative Side Effects

Possible complications are as follows:

- Allergic reactions to contrast media
- Hematoma
- Bleeding
- Infection
- Embolism
- Pneumothorax
- Cardiac tamponade
- Myocardial damage
- Venous or cardiac perforation
- Local tissue response, fibrotic tissue formation
- Damage to the heart valves
- Vascular occlusion
- Chronic nerve damage

Storage

The product should only be stored at temperatures between 5°C and 55°C (41°F to 131°F).

Packaging and Sterility

The Selectra slitter tool is packaged in two sterile bags so that the exterior of the inside bag is also sterile. Ethylene oxide is used for sterilization purposes.

The two sterile bags are packaged in a container that bears a quality control seal and a product information label. The label lists the model type, technical data, use by date, and information about sterility and storage of the package and contents.

- 1 To ensure sterility, the container should be checked for damage prior to opening.
- 2 If it is suspected that the sterile packaging has been opened or damaged, return the product to BIOTRONIK.

Scope of Delivery

- 1 Slitter tool (sterile)
- 1 Technical manual (non-sterile)

Opening the Package

In a non-sterile area:

- Remove the outer sterile bag from the box and open it.

In a sterile area:

- Remove and open the inner sterile bag.

Suggested Slitter Tool Preparation

Prior to use, verify the accessory's functionality and ensure that its size and shape are appropriate for the procedure to be performed.

Disposal

Used implant accessories must be disposed of as medical waste in an environmentally friendly and proper manner.

Accessories do not contain material that requires further controls for disposal.

Handling the Slitter Tool

Note: The CS guiding catheters that pertain to BIOTRONIK's Selectra CS lead introducer system must be removed with a slitter tool.

- 1 Take the slitter tool in one hand, and the lead with the catheter that is to be removed in the other - see Figure 1.

Note: The slitter tool is suitable for both left-handed and right-handed use.


Fig. 1: Holding the slitter tool, clamp not open

- Legend:
- 1 Grip plate of the slitter tool
 - 2 Blade
 - 3 Clamp for the lead

- 2 To open the clamp, press and hold the release key - see Figure 2.
- 3 With the other hand, insert the lead into the clamp - see Figure 2.


Fig. 2: Insert the lead into the open clamp

- 4 Disengage the release key - the lead is held in the clamp.
- 5 Position the handle of the guiding catheter on the blade of the slitter tool - see Figure 3.
- 6 Support the hand that is holding the slitter tool (if possible), holding the catheter handle level with the slitter tool. Pull the catheter against the slitter tool, slitting it open.


Fig. 3: Hold the slitter tool level with the catheter handle.
Pull the catheter parallel to the edge of the slitter tool against the blade.

If the lead lies in a telescopic system consisting of an inner and an outer catheter, the inner catheter is removed first, followed by the outer catheter, using the slitter tool. The slitter tool is suitable for slitting up to three catheters.

Note: When using the slitter tool, make sure that the guiding catheter moves straight along the slitter tool (Fig. 3).

If the slitting procedure requires an unusual amount of effort, or if other signs of problems occur while slitting the catheters, resume slitting or reclamp the lead or replace the slitter tool with a new one and then continue slitting.

- 7 If present, remove the introducer sheath appropriately.

Disclaimer

BIOTRONIK's Selectra CS lead introducer system (including all individual accessories available separately), used in conjunction with devices for implantation in the coronary sinus, has been qualified, manufactured, and tested according to proven and accepted standards and procedures. Nonetheless, the physician must be aware that the accessories may easily be damaged through improper handling or use. Except as set forth in BIOTRONIK's limited warranty, BIOTRONIK makes no expressed or implied warranties for its accessories.

Legend for the Label

Manufacturing date

Use by

Temperature limit

BIOTRONIK order number

Lot number

Sterilized with ethylene oxide

Do not resterilize

Do not reuse

Non-sterile

Consult instructions for use

Contents

Do not use if packaging is damaged

CE mark

Slitter tool (Selectra)

Clamp of slitter tool (Selectra) for leads with a diameter of between 4.6 and 5.8 F

Caution:
Federal (U.S.A.) law restricts this product to sale by, or on the order of, a physician.

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All rights reserved.
BIOTRONIK, Inc.