Contains/Contient/Contiene/Enthält/ Contiene:

Aliphatic Urethane Methacrylate

Warning/May cause an allergic skin reaction./Wear protective gloves, protective clothing, eye protection and face protection./If skin irritation or rash occurs: Get medical attention./Take off contaminated clothing and wash before reuse./Dispose of contents/container in accordance with local and national regulations.

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© 1986, 2014 DENTSPLY International Inc. All rights reserved. Printed in U.S.A.

90256 Rev. 10 (12/14)

DENTSPLY

Trubyte®

# Triad® Gel

Visible Light Cure Material  
Matériau photopolymérisable  
Materiale fotopolimerizzabile a  
luce visibile  
Lichthärtendes Material  
Material fotopolimerizable

DIRECTIONS FOR USE  
CONDITIONS D'EMPLOI  
ISTRUZIONI PER L'USO  
VERARBEITUNGSSANLEITUNG  
INSTRUCCIONES DE EMPLEO

**Indications for Use:**

Triad Gel Material is a visible light cure material indicated for the fabrication of orthodontic and other intra-oral appliances including, but not limited to, retainers, functional appliances, and occlusal splints/night guards. Triad Gel Material is indicated where materials with high flow characteristics are desired, such as when encapsulating and retaining orthodontic wires, mesh, etc.

Triad Gel Material is a no-mix urethane dimethacrylate-based material that contains no methyl methacrylate monomer.

**Contraindications:**

Triad Gel Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).

**Warnings:**

Triad Gel Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

**Precautions:**

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
3. All cure times cited are appropriate for Triad 2000 and Triad II curing units. Consult your Technical Manual or call Technical Support for Triad I cure times.
4. Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.
5. Unused Triad Materials should be polymerized prior to disposal.



**Adverse Reactions:**

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

**STEP-BY-STEP INSTRUCTIONS****I. BASIC FUNCTIONAL APPLIANCE:**

1. Prepare wires, retainers, retentive mesh, etc. in the normal manner. Allow approximately 0.5 mm (0.02") space between the metal retention and the cast surface (Figure 1A).


Thoroughly clean metal portions of appliances by brushing with methyl methacrylate monomer liquid or ethyl alcohol. Allow to air dry while the cast is being prepared.

**Fig. 1B**

Cross-section view of Fig. 1A.

**Fig. 1A**


2. Prepare the cast by applying a thin layer of Triad Model Release Agent (MRA) to only those areas that will come in contact with Triad Gel. Wipe off excess MRA.
3. Position the wires, etc. on the cast. Secure with sticky wax on buccal and facial areas not covered with MRA (Figure 1B).
4. Extrude Triad Gel under and around the wires, then flow around teeth, into interproximal and other areas.

**Fig. 2A**

- a. Cure the appliance on the cast for four (4) minutes in the Triad VLC Unit.
- b. Remove the appliance from the cast.

**Fig. 2B**

- c. Paint all Triad material surfaces with Triad Air Barrier Coating (ABC).
- d. Cure the appliance off the cast, tissue side up, for three (3) minutes in the Triad VLC Unit.
5. Remove the appliance from the curing unit and scrub thoroughly with warm water and a brush to remove all traces of Triad ABC.
6. Finish in the normal manner.

**Fig. 2C**

## **II. PALATAL EXPANDERS:**

1. Prepare the cast as described under "Basic Functional Appliance."
2. With expander positioned, apply Triad Gel to one-half of the palate.
3. Cure 30 seconds in the Triad VLC Unit. Apply Triad ABC to the exposed Triad Gel Material surface, and cure one (1) additional minute on the cast.
4. Ensure that Triad ABC is dry before proceeding.
5. Apply Triad MRA to the dry Triad ABC coated surface, especially the palate edge what will contact the second half of the expander.
6. Apply Triad Gel Material to the rest of the palate to form the second half of the expander appliance.
7. Complete the appliance as previously described in "Basic Functional Appliance."

**B. Occlusal Splints** can be made easily on the cast with Triad Gel following conventional design procedures.

## **C. Patient Identification Procedures:**

1. For new appliances:

- a. Apply a thin layer of Triad Gel to the prepared cast surface in the desired location and cure for 30 seconds.
- b. Moisten "onion skin" paper containing the patient's typed identification data with heat-cured methyl methacrylate monomer (MMA), and place the paper on the thin, partially-cured Triad Gel Material surface.

- c. Continue the Triad Gel Material application as previously described in "Basic Functional Appliance."

2. For existing appliances:

- a. Create a shallow well in a non-tissue bearing surface of the appliance.
- b. Moisten "onion skin" paper containing the patient's identification data using heat-cured MMA. Position the paper in the prepared area.

- c. Cover with a layer of Triad Gel Material and cure for 30 seconds.
- d. Cover with Triad ABC

and cure an additional 3 minutes.

- e. After curing, scrub gently with warm water and a brush to remove all traces of Triad ABC.
- f. If needed, trim excess material and polish in a normal manner.

## **D. Additions, Modifications, and Repairs**

Additions of Triad Gel Material can be made at any time during construction before Triad ABC has been applied. After application of the Triad ABC, Triad Gel Material can be added to the cured surface using the following procedure:

1. Remove the Triad ABC by scrubbing thoroughly with a brush and warm water.
2. Freshen the Triad material surface by lightly grinding the area where the addition is to be made.
3. Apply a thin layer of Triad VLC Bonding Agent to the prepared surface, and allow to bench set for 2 minutes.
4. Cure the bonding agent in the Triad VLC Unit for 2 minutes.
5. Add Triad Gel Material and process as previously described in "Basic Functional Appliance."

\*Heat-cured methyl methacrylate monomer (Lucitone Liquid) can be substituted for Triad VLC Bonding Agent. Refer to the repair instructions in the Triad VLC System Technique Manual and follow these procedures through completion.

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