Medtronic

5866-46

Sizing sleeve

Manchon de calibrage

Manguito ajustable

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Description/intended use

The 5866-46 sizing sleeve is designed to allow IS-1¹ unipolar or bipolar lead connectors to fit into Medtronic® pulse generators with unipolar connectors (5 mm) or bipolar connectors (5 mm bifurcated).

All four of the sizing sleeves are made of silicone rubber.

Caution: The 5866-46 sizing sleeves are not intended for use with Medtronic pulse generators that have a universal (5/6 mm) unipolar connector. They are too short for the universal (5/6 mm) connector cavity. Any attempt to remove these sizing sleeves from the universal (5/6 mm) connector cavity may require using an instrument which could damage the connector block.

The Model 5866-46 is intended for single use only.

¹ IS-1 refers to an International Connector Standard (ISO 5841-3) whereby pulse generators and leads so designated are assured of a basic mechanical fit.

Contraindications

If the reliability of the implanted lead is in doubt, use of the 5866-46 sizing sleeve is contraindicated.

Contents of package

2 short sizing sleeves (clear)

2 long sizing sleeves (white opaque)

Documentation

Instructions for use

Procedure to fit IS-1 lead(s) into unipolar (5 mm) or bipolar (5 mm bifurcated) Medtronic pulse generator connectors.

1. If the sizing sleeves are to be used with implanted lead(s), disconnect the lead from the pulse generator using suitable wrenches. Examine the lead for signs of connector pin corrosion, pitting, or any other visible deterioration and test the leads for adequate pacing and sensing thresholds.

Caution: If the lead(s) are electrically or mechanically unreliable, do not use them.

2. Insert the short sizing sleeve (clear) into the top pulse generator connector port. Insert the long sizing sleeve (white opaque) into the bottom connector port (see Figure 1). If necessary, use sterile water as an insertion lubricant.

Caution: Do not use medical adhesive to secure the sizing sleeve; this could make it impossible to remove the sizing sleeves.


Figure 1. Insert the short and long sizing sleeves into the pulse generator connector.

3. Gently insert the male connector end of the lead(s) into the pulse generator connectors (see Figure 2). If necessary, use sterile water as an insertion lubricant.

The IS-1 connector pins may not be readily visible through the connector block. Care should be taken to assure full insertion of the lead.

Caution: Inserting a bipolar lead (IS-1 BI) into pulse generator connector port containing a sizing sleeve will insulate the ring electrode and cause the lead to function only in a unipolar mode.


Figure 2. Insert the lead(s) male connector into the pulse generator.

4. Secure the pulse generator lead assembly by tightening all of the setscrews as instructed in the pulse generator technical manual.

Warnings

An implanted lead forms a direct low-resistance current path to the myocardium. Care must be taken during surgical procedures to insulate the terminal pins of an implanted lead from leakage currents that may arise from nearby line-powered equipment.

Since the sizing sleeve is a functional extension of the pacing lead, it is recommended that the pacing lead

manual be reviewed for all appropriate warnings, complications, and precautions.

Precautions

Necessary Hospital Equipment

Defibrillating equipment should be kept nearby for immediate use during the implantation procedure and whenever arrhythmias are possible or intentionally induced.

Handling the Sizing Sleeves

Lead insulators attract small particles such as lint and dust; therefore, protect the sizing sleeves from materials shedding these substances.

Ensuring Package Integrity

Medtronic has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized.

Carefully inspect the package before opening.

It is not recommended that you use the product after its expiration date.

Potential complications

The potential complications related to the use of the 5866-46 sizing sleeves include, but are not limited to the following:

- intermittent or continuous loss of pacing and/or sensing,
- erosion,
- infection, inflammation, irritation or body rejection, and
- skeletal muscle and nerve stimulation.

Special notice

Medtronic Sizing Sleeves are used with leads, which are implanted in the extremely hostile environment of the human body. Sizing sleeves may fail to function for a variety of causes, including, but not limited to: medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, leads or sizing sleeves may be easily damaged before, during, or after insertion by improper

handling or other intervening acts. Consequently, for leads used with sizing sleeves, no representation or warranty is made that failure or cessation of function will not occur or that the body will not react adversely to the implantation of leads or that medical complications (including perforation of the heart) will not follow.

Disclaimer of warranty

The following Disclaimer of Warranty Applies to United States Customers Only:

ALTHOUGH THE MEDTRONIC MODEL 5866-46 SIZING SLEEVE HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following Disclaimer of Warranty Applies to Customers Outside the United States:

Although the Medtronic Model 5866-46 Sizing Sleeve, hereafter referred to as "Product" has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of this **DISCLAIMER OF WARRANTY**.

Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **DISCLAIMER OF WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with

applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© Medtronic, Inc. 2010
The following are trademarks of Medtronic:
Medtronic B.V.
Medtronic Australasia Pty Ltd