Seraclone® Anti-k (KEL2)

For the detection of red blood cell characteristic k

Tradeform

REF	808 126	2 ml
PRES	Preservative: 0.1 % NaN ₃	
ACT	Activity / Min. titer	

IVD Diagnostic reagent for in vitro use only

To be used by trained laboratory personnel only. The **ACT** is printed on the label. Symbols are used according to standard DIN EN 980. Symbols not contained in standard DIN EN 980 are described.

Intended use

Seraclone® Anti-k (KEL2) is used for the determination of the k (Cellano) antigen of red blood cells using the tube test.

Summary Kell-System

Since the description of the antigen K in 1946 by Coombs *et al* and its allele k in 1949 by Levine *et al*, the Kell blood group system has been shown to be increasingly complex and over 20 antigens are now known to be associated with the system. These are probably controlled from a series of closely linked loci so that Kell antigens, like CDE in the Rh system, are inherited as haplotype.

The antigens of the Kell blood group system are of further interest in that they tend to occur either very frequently or relatively infrequently and show considerable ethnic variation (e.g. k (KEL2): Caucasians 99.8%, Blacks 100%)² and K (KEL1): Caucasians 9%, Blacks 2 %, Asians rare, Arabs as high as 25%)². Kell system antibodies are capable of causing hemolytic transfusion reactions and hemolytic disease of the newborn.

Test principle

The test principle is a hemagglutination test to detect the respective red cell antigens. The blood type k is defined by presence or absence of the k characteristic on red cells. The antigen k reacts specifically with the corresponding antibody in Seraclone® Anti-k (KEL2) by direct agglutination. The reagent is used in tube testing.

Reagent

As reactive component Seraclone® Anti-k (KEL2) contains a murine monoclonal antibody of the immunoglobulin class IgM which is derived from hybridoma cell lines created by fusing mouse antibody producing B lymphocytes with mouse myeloma cells and demonstrate consistent specificity and reproducibility characteristic for monoclonal antibodies.

Antibodies are diluted in citrate buffer pH 5.2 to 5.3 with bovine albumin and macromolecular potentiators.

Seraclone® Anti-k (KEL2) clone Lk1

Materials required but not provided

- Pipettes (volume 50 µl)
- Phosphate Buffered Saline (PBS) pH 7.3 ± 0.2 or isotonic saline
- Glass tubes
- Laboratory centrifuge

Specimen collection

Fresh samples of clotted, EDTA or citrate anticoagulated whole blood collected following general blood sampling guidelines are acceptable. The specimen should be tested as soon as possible after collection. If testing is delayed, specimens should be stored at 2 to 8°C.

Blood specimens exhibiting gross hemolysis or contamination should not be used. Clotted samples or those collected in EDTA may be tested within 14 days from collection. Donor blood stored in citrate anticoagulant may be tested until the expiration date of the donor unit.

Test procedure

1. Prepare a 3 to 5% suspension of red blood cells to be tested in saline or PBS.
2. Place one drop (approximately 50 µl) reagent into an appropriately labelled tube.
3. Add one drop (approximately 50 µl) of red blood cell suspension into the tube and mix.
4. Incubate for 5 to 10 minutes at room temperature.
5. Centrifuge for :
 - a. 20 seconds at 800 to 1000 x g or
 - b. at a time and speed appropriate for the centrifuge calibration.
6. Gently dislodge the red blood cell button and observe for macroscopic agglutination.

Interpretation of results

Evaluation of the reaction strength should be carried out according to the Technical Manual¹

Reaction strength	Agglutination
4+	One single agglutinate
3+	Several large agglutinates.
2+	Medium-size agglutinates, clear background.
1+	Small agglutinates, turbid background
+/-	Barely visible agglutination, turbid background
-	No agglutination

Stability of the reaction

Following centrifugation, all tube tests should be read immediately and results interpreted without delay. Time delays may cause a dissociation of the antigen-antibody complexes resulting to false negative or more often weak positive reactions.

Quality control

The reactivity of all blood typing reagents should be confirmed on each day of use.

To confirm the reactivity or specificity of Bio-Rad monoclonal Seraclone® Anti-k (KEL2) Blood Grouping Reagent should be tested with antigen-positive (preferably from heterozygous individuals) and antigen-negative red blood cells. The reagent is acceptable for use if it reacts only with antigen-positive red blood cells.

Performance characteristics and limitations of the method

- Red cells from individuals of the Kell phenotype Kk, Kp (a+b+) show a substantially weakened expression of the k antigen.
- The Kell and Cellano antigenicity can be weakened in case of antigen variants which occur rarely, for example the McLeod syndrome. The McLeod syndrome is associated with the McLeod phenotype, whereby the Kell antigens are slightly expressed and Km (KEL20) and Kx are missing.¹ Excessive agitation may disrupt weak agglutination and produce false negative results.
- In case of ambiguous results it is recommended to wash red blood cells at least 2 times.
- Stored red blood cells may exhibit weaker reactions.
- Some conditions that may cause false results are:
 - Contamination of sample or reagent,
 - Improper storage or preparation of red blood cells
 - Autoantibodies
 - Antibodies to antibiotics
 - Cold antibodies
 - Positive Direct Anti-Globulin Test
 - Strongly lipemic or icteric samples

In case of questionable results of unknown origin our Bio-Rad Service will assist you (phone: +49-6103-3130-611).

Shelf life

Storage temperature of Seraclone Anti-k (KEL2) is 2...8°C. After opening the reagent the product can be stored until its expiry date under proper storage conditions. The expiry date is printed on the label.

Warning and precautions

- Do not use damaged vials
- Turbidity or other visible changes may indicate bacterial contamination. Do not use if turbid.
- All handling techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user and needs to be validated by the user.
- Used tests must be discarded as hazardous material. Waste information is provided in the safety data sheet.
- Because of biotechnological manufacturing of monoclonal antibodies the risk of contamination with infectious agents can be considered as insignificant. Nevertheless, all test reagents must be regarded as potential infectious. Appropriate safety precautions are recommended.
- The packaging of this product contains natural rubber latex which may cause allergic reactions.
- The product contains sodium azide (NaN₃), which may react with lead or copper plumbing to form explosive azides. If discarded in the sink, flush with large amounts of water to prevent the build-up of explosive metal azides.

References

- (1) Technical Manual, 17th edition, Section 1, American Association of Blood Banks
- (2) The Blood Group Antigen Facts Book, 3rd edition, Marion Reid, Christine Lomas-Francis and Martin L. Olsson

Key: Underline = Addition or significant change ◀= Deletion of text

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