DENTSPLY
TRUBYTE

Triad[®]

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DENTSPLY
TRUBYTE

Manufactured by:
DENTSPLY International Inc.
570 West College Ave.
York, PA 17405-0872
1-800-786-0085

www.trubyte.dentsply.com or www.dentsply.com

DeguDent GmbH
Postfach 1364 • 63403 Hanau
GERMANY
Telefon +49/6181/5950

Indications for Use:

Triad Reline Material is indicated for relining complete and partial removable dentures and for denture repairs.

Triad Reline Material is an urethane dimethacrylate-based acrylic that contains no methyl methacrylate monomer.

Contraindications:

1. Triad Reline Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).
2. Triad Reline Material is contraindicated for patients with infected, inflamed or lacerated oral tissues. These conditions should be alleviated before treatment.

Warnings:

1. Triad Reline Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

Precautions:

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation,

masks and vacuum systems should be used.

3. All cure times cited are appropriate for Triad 2000 and Triad II curing units. Consult your Technique Manual or call Technical Support for Triad I cure times.
4. Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.
5. For best results, remove the denture 1 ½ - 2 hours prior to the relining procedure to allow sufficient tissue relaxation.
6. To avoid compromising occlusal relationships, prepare maxillary and mandibular relines sequentially.
7. If patient has inadequate salivary flow, tissue lubrication may be needed.
8. For patients with excessive undercuts, the relining procedure could lock the denture in place. To avoid this problem, reduce the intraoral cure time with a hand-held light to 5-second bursts in the undercut areas.
9. Triad Reline sheet materials will cause some tissue distortion during the impression phase. Allow sufficient time for tissue to relax prior to curing the relining.
10. Do not use Triad materials that crack or craze when handled.
11. Thermoplastic denture bases may be distorted in the Triad Unit.

12. Unused Triad Materials should be polymerized prior to disposal.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.
3. Some patients may experience an unpleasant taste from uncured Triad Materials.

STEP BY STEP INSTRUCTIONS

Direct Reline Procedures

The maxillary and mandibular dentures are prepared identically but each reline is accomplished as a separate procedure in order to maintain proper occlusal relationships.

1. A. With a bur remove excessive undercuts and reduce the height of the peripheral border approximately 2 mm (0.08"), leaving slight labial-facial sloping. With a #8 round bur, establish a bead line at the termination of the labial-facial slope to produce a definitive finish line.
- B. With the bur, relieve the entire tissue-bearing surface approximately 2 mm (0.08"). (Consider the thickness of the original denture base as you trim the tissue surface). This space will accommodate the Triad Reline Material and minimize any alteration of the vertical dimension of occlusion.

NOTE: If you wish to restore lost vertical dimension, only "freshen" the entire internal surface with the bur for a clean bond. This will allow the Triad Reline Material to add vertical dimension and achieve an opening effect.

- C. Paint the entire tissue-bearing surface of the denture to be relined with Triad VLC Bonding Agent. *Benchset for 2 minutes and cure for 2 minutes in the Triad VLC Light Curing Unit* before applying Triad Reline Material.

2. Have the patient rinse thoroughly with an astringent mouth wash and use gauze "sponges" to remove any surface debris from the mucosa. A 3% H₂O₂ (hydrogen peroxide) rinse may be helpful.

NOTE: To increase the "flowability" of the Triad Reline Material, before opening the light proof pouch, submerge the material in 110°F (43°C) water for a few minutes. This will result in a softer, more flowable material.

3. Remove the Triad Reline Material from the light-safe packet and carefully adapt the material to the "primed" denture base with finger pressure and/or instrumentation. Use a small moist gauze sponge to aid this adapting procedure, as necessary.
4. Excess material can be removed readily with a sharp blade. Save a small amount of this excess material in its original light safe pouch to use for possible additions to the reline.

5. Place the denture in the patient's mouth, seat with moderate pressure, and manipulate it to establish denture border limits and muscle attachments.
6. Guide the patient into a proper centric position.
7. When you are satisfied with the occlusal relationships, instruct the patient to hold teeth together with moderate pressure for 3 minutes. Then instruct the patient to manipulate lips and tongue to further refine the impression.
8. After the initial seating, remove the denture and rinse in tepid water. Remove excess material and/or apply small amounts of reline material to any defects. Mould and reinsert. Permit the patient to hold as previously described for 3 to 5 additional minutes.
9. Before curing, allow several minutes for the tissues to relax with the denture in place. Instruct the patient to keep the teeth lightly together and/or talk with your assistant. This relaxation time will vary depending on the character of the tissue.
10. Using a high intensity hand-held visible light source, such as the DENTSPLY Caulk Max Light Unit, set the reline material with 10-second bursts along the entire flange and palatal seal areas. Spot cure the reline in approximately 12 positions.
11. Continue with the hand-held curing light onto the palatal surface, tracing the lingual aspects of the teeth approximately 1/2 inch below the tooth line. Five exposure positions should be sufficient.
12. Finally, cure the mid-palatal area in a similar manner, as required.
13. Carefully remove the denture and examine for areas of distortion. If acceptable, immediately place it in the *Triad VLC Light Curing Unit* and cure for 30 seconds to set the reline material.
14. After this preliminary cure, remove the denture from the curing unit and thoroughly rinse under running water to remove any saliva residue. Thoroughly dry the denture.
15. Apply Triad Air Barrier Coating (ABC) to all exposed surfaces of the Reline Material.
16. Place the relined denture in the Triad Curing Unit, tissue side up, and cure for an additional 10 minutes. Thin dentures and partials with metal frameworks should be cured using several shorter cure cycles. Cool in water between cycles. Thermoplastic denture bases may be distorted in the Triad Unit.
17. After curing, rinse the denture with warm water and gently brush to remove all traces of Triad ABC.
18. Trim, polish, and adjust occlusal relation in the standard manner.
19. Generally, a posterior palatal seal (post dam) is not necessary. However, if additional retention is desired, slightly

reduce the posterior palatal seal area by grinding. Paint with Triad VLC Bonding Agent, and add a small amount of the uncured Triad Reline Material to the prepared surface. Mould with finger pressure, re-insert and repeat steps 9 through 17.

NOTE: 24 to 48 hours after insertion some patients may exhibit a slight "plaque-like" deposit on the tissue bearing surface of the denture. Such areas may be gently removed with a denture brush and pumice. This occasionally results if tissues were not adequately cleaned before the impression was made with the Triad Reline Material.

Laboratory Reline Procedures

1. Remove undercuts and excess material from tissue side of the denture. Make the impression and pour the cast. (Plaster is recommended as a cast material, however, dental stone may also be used.)
2. Index the cast and lubricate the base with petrolatum. Mount the cast in the reline jig with the teeth and cast indexed in stone or plaster.
3. Carefully separate the denture from the cast and grind approximately 1 mm (0.04") of denture base material from all tissue surfaces of the denture. The tissue surface should be clean and free of undercuts. To permit excess material to escape more easily, drill relief holes in the palatal and/or flange areas and shorten the periphery one

or two millimeters. Grind buccal-labial flange areas to slightly roughen the surfaces where new reline material will be added.

4. Paint the cast with Triad Model Release Agent (MRA) to aid separation of the Triad reline, when cured. Before applying the Triad Reline Material, apply Triad VLC Bonding Agent to the entire prepared denture surface to be relined. *Benchset for 2 minutes and cure for 2 minutes in the Triad VLC Light Curing Unit.* Alternative Bonding Method: Before relining prime the surface of the reline area by painting it with at least one coat of Lucitone® Liquid rather than Triad VLC Bonding Agent.
5. Adapt Triad Reline Material to the cast or adapt the material into the prepared denture. (If the alternative bonding method was used earlier, be sure to redampen the reline surface of the old denture with Lucitone Liquid to enhance the bond between the old and new materials.)
6. With the cast and dentures positioned in the jig, very slowly close the jig to allow the Triad Reline Material to flow and escape from peripheral areas and/or palatal or flange vent holes. Bring the cast and denture together to the original tooth/cast relationship.

NOTE: To increase the "flowability" of the Triad Reline Material, before opening the light proof pouch, submerge the material in 110°F (43°C) water for a few

minutes. This will result in a softer, more flowable material.

7. Remove excess Triad Reline Material at the periphery and adapt the remaining material to the proper contour. Coat all exposed reline material surfaces with Triad Air Barrier Coating (ABC).
8. Remove the cast from the reline jig with the denture in place. *Cure the relined denture on the*

cast for four (4) minutes in the Triad VLC Curing Unit.

9. Carefully remove the denture from the cast. Paint the tissue surface with Triad Air Barrier Coating (ABC) and *cure the relined denture, tissue side up, for an additional eight (8) minutes.* Rinse with warm water and a soft brush to remove all traces of Triad ABC. Finish and polish in the usual manner.

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