Medtronic

RF MARINR® UNIPOLAR

Steerable electrode catheter for intracardiac mapping and ablation

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! USA

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

CE
0123

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners. Medtronic, RF Marinr

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

Figure 7. To increase the deflected curve, pull the Curve Radius Control. /

Figure 7. Pour augmenter la courbure, tirer le curseur du rayon de courbure. / **Abbildung 7.** Zum Erhöhen des Krümmungsradius den Krümmungsradiusregler zurückziehen. / **Figura 7.** Para aumentar la deflexión de la curva, tire del control de radio de la curva. / **Afbeelding 7.** Om de buigingscurve te vergroten, trekt het besturingsmechanisme voor de curvestraal naar u toe. / **Figura 7.** Per aumentare il raggio di curvatura, retrarre il dispositivo di comando del raggio di curvatura. / **Figur 7.** För att öka böjningens kurva, dra i spetsens kurvradieglage. / **Εικ. 7.** Για αύξηση της καμπύλης εκτροπής, τραβήξτε το χειριστήριο ελέγχου ακτίνας της καμπύλης. / **Figur 7.** Træk vinkelradiusstyret tilbage for at øge defleksionsvinklen. / **Figura 7.** Para aumentar a curva deflectida, puxe o controlo do raio da curva



Figure 8. To apply lateral deflection, rotate the Lateral Deflection Control. /

Figure 8. Pour appliquer une déflexion latérale, faire tourner le curseur de déflexion latérale. / **Abbildung 8.** Zur seitlichen Krümmung den Lateralkrümmungsregler drehen. / **Figura 8.** Para aplicar deflexión lateral, gire el control de deflexión lateral. / **Afbeelding 8.** Om laterale buiging toe te passen, draait u het besturingsmechanisme voor laterale buiging. / **Figura 8.** Per deflettere lateralmente la punta, ruotare il dispositivo di deflessione laterale. / **Figur 8.** Vinkla spetsen lateralt genom att vrida reglaget för lateral vinkling. / **Εικ. 8.** Για εφαρμογή πλευρικής εκτροπής, περιστρέψτε το χειριστήριο ελέγχου πλευρικής εκτροπής. / **Figur 8.** Anvend lateral defleksion ved at dreje det laterale defleksionsstyr. / **Figura 8.** Para aplicar uma deflexão lateral, rode o controlo de deflexão lateral

Description

The Medtronic RF Marivr Unipolar steerable ablation catheter is a flexible, radiopaque catheter constructed of extruded polymer over stainless steel braid. The RF Marivr Unipolar catheter is designed for intracardiac radio frequency (RF) ablation via the tip electrode and separate dispersive electrode, when connected to a Medtronic RF power generator. The RF Marivr Unipolar catheter may also be used for intracardiac recording or stimulation. A unipolar mapping catheter with a reference electrode on the proximal shaft requires less filtering than an external reference electrode, such as Wilson's Central Terminal. Reduced filtering can permit the detection of smaller amplitude signals. With unipolar mapping, signals may be localized from a single point source rather than two sources as in bipolar mapping.

*Not included on all models

Contents of package

The RF Marivr Unipolar catheter is supplied sterile. The package contains the following items:

- 1 RF Marivr Unipolar catheter
- Product documentation

Intended use

The RF Marivr Unipolar catheter is intended for use with the Medtronic RF power generator to deliver RF energy for intracardiac ablation of accessory atrioventricular (AV) conduction pathways associated with tachycardia for the treatment of AV nodal re-entrant tachycardia and for creation of complete AV block in patients with a difficult to control ventricular response to an atrial arrhythmia.

Contraindications

Contraindications for the RF Marivr Unipolar include the following conditions:

- Active sepsis
- Known sensitivity to heparin
- Blood clotting abnormalities

- Venous filtering device (Greenfield Filter)
- Use of catheter in an obstructed or damaged vessel

The catheter is not recommended for use in patients who cannot undergo standard anticoagulation protocol for a left-sided cardiac procedure, or who have had a recent coagulopathy or embolic event.

The transeptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.

The retrograde transaortic approach is contraindicated in patients with aortic valve replacement.

Warnings and precautions

General

Related product literature – Do not attempt to operate the Medtronic ablation system or connect the ablation catheter to a Medtronic RF power generator prior to completely reading and understanding the Medtronic ablation system technical manual and the ablation catheter instructions for use.

System compatibility – Use the catheter with only a Medtronic RF power generator and accessories. The safety and use with other RF power generators or accessories has not been tested. Use only Medtronic cables.

Qualified users – The catheter should be used only by or under the supervision of physicians trained in ablation procedures using this catheter and the Medtronic RF ablation generator.

Required use environment – Cardiac ablation procedures should be performed only in a fully equipped electrophysiology laboratory.

Ablation therapy hazards

Serious adverse events – A number of serious adverse events have been documented for catheter ablation procedures, including pulmonary embolism; myocardial infarction; cerebrovascular accident; cardiac damage, perforation, and tamponade; perforation of the vasculature; and death. See the “Adverse events” section for additional potential adverse events.

Left-sided ablation procedures – Patients undergoing left-sided ablation procedures should be closely monitored during the post-ablation period for clinical manifestations of infarction.

Distal pair electrode spacing of >2 mm – Catheters with distal pair electrode spacing greater than 2 mm should not be used in the ablation of septal accessory pathways or in the treatment of AV nodal re-entrant tachycardia because of the potential for creating inadvertent complete AV block.

Catheter manipulation and placement – Provide adequate fluoroscopic visualization during catheter manipulation and placement. During a transaortic approach, avoid placement of the ablation catheter within the coronary vasculature. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. Catheter placement and RF power application within a coronary artery have been associated with myocardial infarction and death.

X-ray and fluoroscopic exposure – Due to the x-ray beam intensity and the duration of the fluoroscopic imaging during ablation procedures, patients and laboratory staff may be subjected to acute radiation injury and increased risk for somatic and genetic effects. The long-term effects of protracted fluoroscopy have not been established.

- Minimize x-ray exposure.
- Carefully consider the use of the device in pregnant women and prepubescent children.

AV conduction – Closely monitor AV conduction during RF energy delivery in patients undergoing AV node modification or septal accessory pathway ablation. These patients may be at risk for complete atrioventricular (AV) block. Immediately terminate energy delivery if partial or complete AV block is noted.

Leakage current – Use only isolated amplifiers, pacing equipment, and ECG equipment (IEC 601-1 Type CF equipment, or equivalent) or patient injury or death may occur. Leakage current from any connected device to the patient must not exceed 10 μ A under any circumstances.

Catheter removal – See your Medtronic ablation system technical manual for information concerning catheter removal following generator shutdown.

Long-term risk – The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.

Emergency use – In the event of an emergency, do not use the RF Mariner Unipolar catheter to sustain the patient's life. Patient death or injury may occur.

Flammable materials – Do not allow flammable material in the area where RF ablation procedures are performed. The risk of igniting flammable gases, flammable agents used for cleaning or disinfecting, or other materials is inherent in the application of RF energy.

Placement of monitoring electrodes – Place all physiological monitoring electrodes as far away as possible from the patient return electrodes and their leads to avoid RF interference, which affects the ability to interpret patient electrograms (EGMs).

Catheter storage and handling

Storage conditions – Store the catheter in normal operating room temperatures and humidity levels and in a manner that protects the integrity of the package and the sterile barrier. Keep dry. The catheter has the following storage temperature limits: -30 °C (-22 °F) to 60 °C (140 °F), up to 80% relative humidity.

Sterile package inspection – Inspect the sterile packaging and catheter prior to use. If the sterile packaging or the catheter exhibits damage, do not use the catheter. Contact your local Medtronic representative.

For single use only – This device is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize this device for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Sterilization – Medtronic has sterilized the package contents using irradiation before shipment. This device is for single use only and is not intended to be resterilized.

Disposal of catheter – Dispose of the catheter according to hospital biohazard requirements. If returning the catheter, contact your local Medtronic representative.

Catheter handling and care –

- Do not excessively bend or kink the catheter. Excessive bending or kinking may damage internal electrode wires and/or distal tip shaping capabilities.
- Do not allow moisture onto the connectors on the catheter, Medtronic RF power generator, or cables. If the connectors get wet, the system may not function correctly.
- Do not wipe the catheter with organic solvents, such as alcohol, to maintain optimal patient safety and catheter electrode integrity.

Procedure-related equipment use

Power limit – Medtronic recommends limiting the power output to no more than 50 Watts on the RF power generator.

Low battery – Do not begin an ablation procedure if the LOW BATTERY display is illuminated on the Medtronic RF power generator. Replace the battery.

Handling RF catheters during therapy – Do not touch the ablation electrode of the RF catheter and the dispersive electrode at the same time, especially while operating the Medtronic RF power generator. Operator injury may occur.

Use of RF energy near implanted devices – Implantable devices, such as pacemakers and implantable cardioverter-defibrillators (ICDs), may be inhibited or otherwise affected by RF energy. Refer to the RF generator operator manual and the appropriate implantable device technical manual for additional information.

Be aware of the following when using RF energy near implanted devices:

- Keep external sources of pacing and defibrillation available during ablation.
- Deactivate ICDs during the delivery of RF energy.
- Exercise extreme caution when delivering ablation energy in close proximity to atrial or ventricular implanted leads.
- Perform complete implantable device testing before and after ablation.

Magnetic resonance imaging (MRI) – Catheter materials are not compatible with magnetic resonance imaging (MRI).

DC energy – Do not deliver DC energy through the RF power generator or the catheter. Neither the RF power generator nor the catheter were designed to deliver DC energy. No testing has been performed on the RF power generator or the catheter regarding the delivery of DC energy.

Contact with other electrodes – Verify that the catheter tip electrode is not in contact with electrodes of other catheters in the heart, or heating of the other electrodes may occur during RF energy delivery.

Sudden rise in impedance – Discontinue energy delivery if a sudden rise in impedance is observed during the ablation procedure. Examine the catheter ablation electrodes, and remove coagulum if present.

Catheter performance – Bench testing of Medtronic RF catheters demonstrated that they could withstand 25 RF energy deliveries without any reduction in performance.

Clinical summary

Medtronic performed clinical studies using RF ablation systems. Enrollment in the study was 683 patients. Please refer to the Atakr Ablation System Technical Manual for a description of the clinical background. The following summarizes the complications and adverse events related to the clinical study.

There were 18 deaths reported among the 683 patients enrolled in the study. None of the 18 deaths was judged by the investigators as definitely related to the use of the Atakr RFCA System, but two deaths were judged as possibly device related and a third as possibly procedure related.

Device related complications were reported in 30 of the 683 patients (4%) including unintentional AV block in 13 patients (1.9%), pericardial effusion in 8 patients (1.2%), and cardiac tamponade in 3 patients (0.4%).

Permanent pacing was required in 1.6% (2/128) of septal accessory pathway, 3.1% (1/32) of left posterior pathways, and 1.7% (4/238) of AV nodal modification patients who experienced inadvertent partial or complete AV block during the study.

Use the following instructions to view, download, print, or order the Atakr Technical Manual from the Medtronic website:

1. Navigate your web browser to <http://www.medtronic.com/manuals>.
2. Select United States for the location.
3. In the search field on the left side of the screen, type "Atakr".
4. Click [Search].

Adverse events

Potential adverse events associated with cardiac catheter ablation procedures include, but are not limited to, the following conditions:

| | |
|--------------------------------------|--|
| Angina | Pericardial effusion |
| Arrhythmias | Pericarditis |
| Atrial lead dislodgement | Phlebotrombosis |
| Bacterial endocarditis | Pleural effusion |
| Bowel obstruction | Pneumonia |
| Brachial plexus injury | Polymorphic ventricular tachycardia |
| Bradycardia | Proarrhythmia |
| Cardiac tamponade | Pulmonary embolism |
| Cerebrovascular accident | Puncture site hematoma |
| Coronary artery spasm | Respiratory depression |
| Coronary sinus perforation | Temperature elevation |
| Death | Thrombophlebitis |
| Deep vein thrombosis | Thrombus/embolic event |
| Femoral artery laceration | Transient ischemic attack |
| High CPK | Thromboembolism |
| Hypotension | Unintended AV or bundle branch block requiring new pacemaker |
| Infected IV site | Valvular insufficiency |
| Infection | Vasovagal reaction |
| Inferior vena cava clot | Ventricular fibrillation |
| Myocardial infarction | Ventricular flutter |
| Nerve or blood vessel injury | Ventricular tachycardia |
| Non-clinical ventricular tachycardia | |

Instructions for use

Read the Medtronic ablation system technical manual and Sequencr switch box instructions for use before using or connecting the catheter.

Caution: Medtronic recommends limiting the power output to no more than 50 Watts on the RF power generator.

General

1. Inspect the catheter package before opening. The contents of this package have been sterilized before shipment. Do not use the contents if the package is open or damaged.
 2. Use aseptic technique to remove the catheter from the package and place it on a sterile working area (Figure 1).
 3. Inspect the catheter carefully for electrode integrity and overall condition. Do not use the catheter if the electrodes or tip are loose, misshapen, or otherwise visibly damaged.
 4. Create vascular access using sterile techniques. The catheter may be used from femoral, brachial, subclavian, or jugular access sites.
 5. Connect the catheter to the catheter cable.
 - a. Match connector colors for proper cable to catheter connection.
 - b. Connect the cable (Figure 2): Align the double arrows on the cable's plastic connector with the external button on the catheter handle connector.
 - c. Press the connectors together. Do not force the connection.
- Note:** To disconnect the catheter cable from the catheter (Figure 3), pull back on the grip ring to release the lock before removing the connector.
6. Advance the catheter into the desired area of the heart using fluoroscopic and ECG guidance.

Handle Control Operation

1. Pull back on the Tip Deflection Control (Figure 4) to deflect the tip of the catheter. The tip can be deflected up to 180°. To straighten the tip, push the Tip Deflection Control forward (Figure 5).
2. Rotate the Tip Deflection Control clockwise for greater friction and counter-clockwise for less friction (Figure 6) to control the tip position. To maintain tip position, rotate the Tip Deflection Control clockwise to the locked position.
3. Pull back on the Curve Radius Control (Figure 7) to increase the deflected curve radius. To decrease the deflected curve radius, push the Curve Radius Control forward.
4. Rotate the Curve Radius Control clockwise for greater friction and counter-clockwise for less friction (Figure 6) to adjust the force needed to maintain the curve radius.
5. Turn the Lateral Deflection Control clockwise or counter-clockwise to rotate the tip laterally up to 45° in either direction (Figure 8).
6. Before withdrawing the catheter, re-align marks on the Lateral Deflection Control, push the Tip Deflection Control forward and verify by fluoroscopy that the tip is in a neutral position.

Specifications

| | |
|-----------------------------|--|
| Catheter size | 2.3 mm (0.09 in; 7 French) |
| Usable length | 110 cm (43.0 in.) |
| Number of electrodes | 4 unipolar |
| Electrode spacing | 2 mm, 15 mm, 250 mm |
| Electrode size | 4.0 mm (0.16 in.) |
| Curve reach | |
| | Model: 075802 40 - 60 mm (1.57 in. - 2.36 in.) |
| | Model: 075812 45 - 70 mm (1.77 in. - 2.75 in.) |

Medtronic limited warranty

For complete warranty information, see the accompanying warranty document.

Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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© Medtronic, Inc. 2012
Medtronic B.V.
Medtronic Australasia Pty Ltd