Braided Sheath Cutter

STERILE EO

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Pressure Products
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760
www.pressure-products.com

Merit Medical Systems, Inc.
65 Great Valley Parkway
Malvern, PA 19355 - USA
Tel +1-801-253-1600
USA Customer Service:
Tel +1-800-356-3748

EC REP

Merit Medical Ireland, Ltd
Parkmore Business Park West
Galway, Ireland

European Customer Service:
Tel +31 43 358 82 22

CE
0344

Pressure Products
ACCESSING THE HEART OF INNOVATION


The Braided Sheath Cutter is covered by USA Patent 6,497,681.
SafeSheath and CSG are registered trademarks of Pressure Products Medical Supplies, Inc.










Technical Manual

en Cutter

da Skærer
de Cutter
el κόπτη
es Cortador
fi Leikkuri
fr Couteau
it Dispositivo di taglio
ja カッタ
nl Cutter
no Kniv
pt Cortador
sv Skärkniv

en	Directions for Use	2
da	Brugsanvisning	2
de	Gebrauchsanleitung	2
el	Οδηγίες Χρήσης	3
es	Instrucciones de uso	4
fi	Käyttöohje	4
fr	Mode d'emploi	5
it	Istruzioni per l'uso	5
ja	使用説明書	6
nl	Gebruiksaanwijzing	6
no	Bruksanvisning	7
pt	Instruções de Utilização	7
sv	Bruksanvisning	8

en-Directions for Use

To aid in removal of a braided sliceable introducer after it has been used to place a pacemaker lead or catheter.

Precautions

- Do not alter this device in any way.
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- As specified in the following chart, use an introducer sheath and pacing lead within the parameters listed in the chart below.

Largest Introducer Size	Largest Lead Size
4.0mm/12F	2.0mm/6F

Suggested Directions for Use

The Cutter is to be used after pacing lead or catheter is in the desired location and the sheath introducer is ready to be removed. It is suggested a lead stylet or stiffening wire be retained in the pacing lead during sheath removal.

If cutting a SafeSheath CSG Braided Series Sheath:

1. Open package and place contents on sterile field.
2. It is suggested that the braided sheath be withdrawn to the furthest extent possible prior to splitting the sheath valve and cutting the sheath.

NOTE: The lead tip and lie of the lead along the RA floor should be carefully observed fluoroscopically as the braided sheath is withdrawn and cut, to avoid inadvertent lead dislodgment.

3. Manually split the handles of the sheath valve. One half of the valve will remain attached to the braided sheath. Dispose of the other valve section. (see Figure 1)
4. Position the distal 'C' section of the Cutter under the sheath tubing, distal to the sheath hub. (see Figure 2)
5. Place the blade hook into the 'V' notch of the tubing at the proximal end of the sheath. (see Figure 3)
6. **Thread the pacing lead over the notch at the back end of the Cutter and compress the lead with your thumb,** while holding the Cutter in a 'pistol grip' fashion. (see Figure 4)
7. While holding the Cutter and the lead with one hand, slowly withdraw the braided sheath over the Cutter blade with the other hand. As the sheath is withdrawn it is cut away from the lead in the process. (see Figure 4)

If cutting a braided Lateral Vein Introducer (LVI) within a CSG sheath:

1. Open package and place contents on sterile field.
2. With a stabilizing stylet or stiffening wire in place, and while holding the lead in place, under fluoroscopic control, withdraw the LVI into the distal end of the SafeSheath CSG. (see Figure 5)

NOTE: The lead tip should be carefully observed fluoroscopically as the LVI sheath is withdrawn and cut, to avoid inadvertent lead dislodgment. The outer CSG sheath should be stabilized by an assistant grasping the CSG distal to the CSG hub as the LVI is withdrawn and cut to prevent inadvertent dislodgment of the CSG outer sheath from the mid CS.

3. Manually split the handles of the LVI sheath valve. One half of the valve will remain attached to the braided LVI. Dispose of the other valve section. (see Figure 6)
4. Position the distal 'C' section of the Cutter under the LVI sheath tubing distal to the sheath hub. (see Figure 7)
5. Place the blade hook into the 'V' notch of the LVI tubing at the proximal end of the LVI Guide. (see Figure 3)
6. While maintaining the LVI Guide, pacing lead and CSG Sheath in position, hold the Cutter in 'pistol grip' fashion and cut the LVI by advancing the Cutter toward

the CSG sheath until the distal 'C' portion of the Cutter engages the CSG hub. (see Figure 8)

NOTE: Make sure the pacing lead, LVI guide and CSG sheath are not advanced or withdrawn with the Cutter during this step of the procedure.

7. **After the Cutter has engaged the CSG hub, thread the pacing lead over the notch at the back end of the Cutter and compress the lead with your thumb.** (see Figure 9)
8. While holding the Cutter engaged with the CSG hub and lead with one hand, slowly withdraw the LVI sheath over the Cutter blade with the other hand. As the LVI is withdrawn it is cut away from the lead in the process. (see Figure 9)
- NOTE:** The outer CSG sheath should continue to be stabilized distal to the CSG hub as the LVI is withdrawn and cut to prevent inadvertent dislodgment of the CSG outer sheath from the mid CS.
9. After the LVI is cut and removed the CSG is removed as described above under the heading "If cutting a SafeSheath CSG Braided Series Sheath," or in the appropriate DFU.

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