<b>TrueLok™</b>		<b>TrueLok™</b>	
EXTERNAL FIXATION SYSTEM		EXTERNT FIXERINGSSYSTEM	
3		86	
<b>TrueLok™</b>	EN	<b>TrueLok™</b>	SV
SISTEMA DI FISSAZIONE ESTERNA		ΣΥΣΤΗΜΑ ΕΞΩΤΕΡΙΚΗΣ ΣΤΑΘΕΡΟ ΟΙΗΣΗΣ	
11		94	
<b>TrueLok™</b>	IT	<b>TrueLok™</b>	EL
SYSTÈME DE FIXATION EXTERNE		創外固定システム	
19		102	
<b>TrueLok™</b>	FR	<b>TrueLok™</b>	JP
EXTERNES FIXATIONSSYSTEM		外固定系統	
27		110	
<b>TrueLok™</b>	DE	<b>TrueLok™</b>	CN
SISTEMA DE FIJACIÓN EXTERNA		ZEVNÍ FIXAČNÍ SYSTÉM	
36		117	
<b>TrueLok™</b>	ES	<b>TrueLok™</b>	CS
EKSTERNT FIKSATIONSSYSTEM		ZEWNĘTRZNY SYSTEM MOCOWANIA	
46		125	
<b>TrueLok™</b>	DA	<b>TrueLok™</b>	PL
ULKOINEN KIINNITYSJÄRJESTELMÄ		SISTEM ZUNANJEGA FIKSATORJA	
54		133	
<b>TrueLok™</b>	FI	<b>TrueLok™</b>	SL
SYSTEM FOR EKSTERN FIKSERING		외고정 시스템	
62		141	
<b>TrueLok™</b>	NO		KR
EXTERN FIXATIESYSTEEM			
70			
<b>TrueLok™</b>	NL		
SISTEMA DE FIXAÇÃO EXTERNA			
78			

ORTHOFIX Srl

Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy  
Telephone 0039-0456719000 - Fax 0039-0456719380

0123

PQ TLK I 06/12

# TrueLok™

## External fixation system RX Only

### Important information for the surgeon performing the operation Read instructions before use

#### General information

The Orthofix TrueLok™ is a modular circular external fixation system consisting of rings, partial rings, posts, plates, hinges, distraction elements, connection elements, half-pins, wires, and fixation bolts. The rings and partial rings are manufactured from aluminum. The linear distractors are manufactured from aluminum, stainless steel, and nylon. The angular distractors are manufactured from aluminum, stainless steel, and Techron PPS. All other elements are manufactured from stainless steel. The system is used with 1.8 mm K-wires and bone screws manufactured by Orthofix Srl. This external fixation system is a metal bone fixation device for stabilisation and correction during the normal healing process. The purpose of the fixation devices is to keep the bone segments aligned during the healing process; they are not designed to support weight bearing loads. The device may only be applied by a physician with adequate surgical knowledge of the system and its relative mechanical limitations. The selection of the correct frame and implant size for each patient is fundamental for the success of the procedure. The metal screws and wires are subject to repeated stress during use and their resistance will be limited by the size and shape of the bones being treated. Unless care is taken in patient selection, proper placement of the frame hinges, wires and screws, and careful postoperative management to minimise the stresses on the frame, metal fatigue may occur, with consequent breakage or bending before the healing process is complete, which may result in further injury or the need to remove the device prematurely. While using the TrueLok™ system for bone fixation, the surgeon should keep the following in mind:

- Patient weight: an obese or overweight patient places loads on the device that could cause breakage or bending.
- Patient occupation or activity: the risk of bending or breakage of an internal or external bone fixation device during post-operative rehabilitation may increase if the patient carries out activities that involve lifting or heavy muscular activity, as these movements subject the device to forces that could cause it to break.
- Mental condition of the patient: the risk of breakage of a fixation device is greater in older patients, or patients with mental deficiency, alcoholics or drug addicts or patients who, for other reasons, may ignore the necessary restrictions and precautions to be observed while using the device.
- Patient sensitivity to certain material: if it is suspected that the patient is sensitive to the used material, suitable tests must be carried out before selecting or implanting the material.
- Debilitated patients: debilitated persons who have difficulty in using weight support devices can be at risk during post-operative rehabilitation.

## **INDICATIONS**

The TrueLok™ system is intended for limb lengthening by metaphyseal or epiphyseal distractions, fixation of open and closed fractures, treatment of non-union or pseudoarthrosis of long bones and correction of bony or soft tissue defects or deformities.

## **CONTRAINDICATIONS**

The Orthofix TrueLOK system is not designed or sold for any use except as indicated. Use of the system is contraindicated in the following situation:

- Patients with mental or physiological conditions who are unwilling or incapable of following postoperative care instructions.

## **WARNINGS**

Each implantable device is disposable. After use, the system must be discarded as it could be damaged and cause problems if reused.

The fixation device must never be subjected to bending, cutting or scratching as these could reduce the fixation frame's resistance to stress and increase the risk of bending or breakage.

Specialized instrumentation is required to apply the fixation devices.

Instruments subjected to prolonged use or excessive force are susceptible to fracture and should be examined prior to surgery.

The patient must be informed that a second minor surgery for the removal of the fixation system is required.

## **WITH REGARD TO USING SCREWS AND WIRES**

It is essential that special care be taken to ensure that the screws do not penetrate joints or growth cartilage in pediatric patients. The length of both the entire screw and the threaded section must be selected in accordance with the dimensions of the bone and soft tissues.

Screws should always and only be inserted using a driver or T-wrench, regardless of whether or not pre-drilling has been performed. It is important to apply moderate force to penetrate the first cortex. Diaphyseal bone screws should always be inserted in the centre of the bone axis, in order to avoid weakening. In all cases, the surgeon should be mindful of the amount of torque required to insert the screw. Should this appear more difficult than usual, it is recommended the screw be removed and cleaned and the bone drilled once more using the drill bit, even if this has already been used.

Excessive penetration of the second cortex must be avoided because this may cause soft tissue damage. Bone screws must never be inserted so far that the smooth shaft penetrates the first cortex, in order to avoid the risk of bone damage.

The diameter of the screw must be selected in accordance with the diameter of the bone: for a bone diameter greater than 20 mm, 6.0 mm, 6.0 - 5.0 or 6.0 - 5.6 mm should be used; for a bone diameter of between 15 and 20 mm, 5.0 mm screws should be used and for a bone diameter between 10 mm and 15 mm, 4.0 mm bone screws are recommended.

Self-drilling screws with a thread diameter of 5.0 mm or more must never be inserted with a power tool, and should be inserted manually or with the use of a hand drill. Self-drilling screws with a smaller thread diameter may be inserted either manually or with the use of a low-speed electric drill.

## **ORTHOFIX SCREWS**

Standard Orthofix screws have a conical thread that is tapered from 5.0 to 6.0 mm; XCaliber screw threads taper from 5.6 mm to 6.0 mm. The use of blunt tipped bone screws requires pre-drilling using suitable drill bits and drill guides. The matching grooves on the screws and drill bits help the surgeon choose the correct drill bit. Blunt drill bits can cause heat damage to the bone and should never be used.

XCaliber bone screws are usually cut to length, either before insertion or after they have been inserted, the clamp applied and the clamp locking screws firmly tightened. They should never be cut after insertion and prior to application of the clamp, since some of the cutting force might be transferred to the bone.

XCaliber bone screws are self-drilling. In cancellous bone, direct insertion using a hand drill or T-wrench is recommended; in diaphyseal bone it is recommended that pre-drilling be performed. If the bone is strong, use the 4.8 mm drill bit and drill guide; where the bone is of poor quality or in the metaphyseal region, where the cortex is thin, use a 3.2 mm drill bit.

## **TrueLok™ SCREWS (Half Pins)**

**TrueLok™ HALF PINS SHOULD BE USED ONLY WITH THE TrueLok™ SYSTEM**

TrueLok™ screws are not tapered and have a thread that has the same diameter throughout.

TrueLok™ screws are available in 4, 5, and 6mm diameters, each having a total length of 180mm and a thread length of 55mm. If the screw has to be inserted in areas wider than 50mm, XCaliber bone screws should be used. TrueLok™ screws are self-drilling, self-tapping, and have a standard quick connect end for ease of use in the operating room.

The following steps outline the recommended method of screw insertion and fixation:

1. A half pin fixation bolt is inserted into an appropriate hole of the ring.  
The half pin fixation bolt acts as a guide for half pin insertion
2. A K-wire is passed through the hole of the fixation bolt along the intended direction of the half pin insertion.  
Make a stab incision in the skin at this level
3. Use scissors or hemostats to make a track through the soft tissues to the bone with blunt dissection only.
4. The desired half pin is inserted through the fixation bolt and soft tissue track, and drilled into the first bony cortex. Insertion through the second cortex should be done manually using the dedicated Orthofix T-wrench.
5. The half pin fixation bolt is then secured firmly to the external support with a nut.

If TrueLok™ screws are cut to length, screw covers should be used to protect the cut end. TrueLok™ screws which have been cut to length could be removed manually with the Universal Chuck with T-Handle (17955) or with a power drill.

Care should be taken of the soft tissues during screw or wire insertion. Never insert a K-wire more than once, since the tip can be blunted, and excessive heating of the bone might occur.

Use suitable Orthofix tools to correctly insert bone screws and K-wires.

Wherever a K-wire or a guide wire is used to guide a cannulated reamer, a drill bit or a screw into position:

- a) Always use NEW K-wire or wire guides.
- b) The wire should always be checked prior to insertion in order to exclude the presence of scratches or bends.
- c) The surgeon must check the tip of the wire as frequently as possible during the introduction of any tool or device onto a wire, so as to avoid the wire being inadvertently inserted more than is necessary.
- d) Each time a tool or device is handled, the surgeon must check that no bone or other residues have accumulated on the wire or in the tool or device, which might force the wire forward.

## K-WIRE INSERTION

When used with a ring, regardless of whether operating on a hybrid or circular system:

- a) K-wires should be introduced into the area of greatest soft tissue risk.
- b) They should be inserted by pushing through the soft tissue and drilling through the bone; they should never be drilled through soft tissue.
- c) Wires should be inserted with awareness of the safe corridors to avoid damaging vital structures.
- d) Used wire should always be disposed of if removed (since the tip can be blunted, and excessive heating of the bone might occur).
- e) To avoid causing injury the ends of wires should be protected with special covers or bent at the ends as soon as they are tensioned.

To tension K-wires, the handles of the wire tensioner must be fully open and the arms fully inserted over the wire down to the wire locking nut, ensuring that at least 6 cm of wire protrudes. When wires are mounted on a circular ring, they must be tightened to a minimal tension of 1200 N. When wires with a central olive are used, to stabilize a fragment, the tension must be reduced to 800/1000N. Grooved washers may be used in numerous positions on a ring or distant from it. The tension possible varies with the position, up to a maximum of 1000 Newtons. When it is attached directly to a ring, care must be taken to maintain the level of tension, to avoid the wire becoming kinked and damaged.

All equipment should be examined carefully prior to use in order to ensure it is working properly. If it is believed that a component or tool is defective, damaged or suspicious, DO NOT USE IT. Hybrid Fixation frames for use in the progressive correction of deformities must be pre-assembled and checked prior to application in order to ensure they provide the correction required and that their hinges are at the correct level.

The fixation must be applied at such a distance from the skin as to allow post-surgical swelling and cleansing, bearing in mind that the stability of the fixation depends on the distance between it and the bone. Additional equipment might be required for the application and removal of the fixation devices, such as wire-cutters, hammers and electric drills.

Do not use Magnetic Resonance imaging in segments where a fixator has been applied. Check the status of the screws and the fixation at regular intervals.

Ensure the screw and wire sites are kept meticulously clean.

The patient must be instructed regarding the use and maintenance of the fixation and care of the screw sites.

The patient must be instructed to report any adverse or unexpected effects to their family physician.

The fracture or bone gap must be checked periodically during treatment, making any necessary adjustments to the fixation. An excessive or persistent gap can delay consolidation.

In patients undergoing callus distraction (usually 1 mm per day, i.e. 1/4 of a turn of the distractor-compressor every 6 hours) the regenerate bone must be checked regularly and monitored radiologically.

Removal of the device: the surgeon should decide when it is time to remove the fixation device. Before final removal, temporary removal of a section of the frame is advised to check the strength of a healed fracture, osteotomy or lengthened segment.

## POTENTIAL ADVERSE EFFECTS

Loosening, bending, rupture or fracture of the fixation device(s) or loss of fixation or migration which may result in nerve, soft tissue, or organ damage, including perforation through the skin or other hemorrhaging. Loss of anatomical position with non-union or malunions with rotation or angulation. Corrosion with localized tissue reaction or pain. Infection, local or systemic. Pain, discomfort or abnormal sensations of the nervous system

resulting from the presence of the device. Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction. Bone loss or reduced bone density due to a reduction in the tension applied to the bone.

1. Nerve or vascular damage following the insertion of wires or screws.
2. Deep or surface infections of the bone screw site, osteomyelitis, septic arthritis, including chronic drainage of the bone screw sites following removal of the device.
3. Oedema or possible compartmental syndrome.
4. Joint constriction, subluxation, dislocation or loss of motor movement.
5. Premature bone callus consolidation during distraction.
6. Possible tension affecting the soft tissues and/or the fixation during manipulation of the callus (e.g. corrections of deformities and/or elongation).
7. Lack of satisfactory bone regeneration, the development of mal-union or nonunion.
8. Fracture of regenerated bone, or at the bone screw holes, following removal of the device.
9. Loosening or breakage of the bone screws.
10. Bone damage due to erroneous bone screw selection.
11. Bone deformities or talipes equinus.
12. The persistence or recurrence of the initial condition subject to treatment.
13. New surgery to replace a component or the entire fixation frame.
14. Abnormal growth cartilage development in skeletally immature patients.
15. Foreign body reactions due to bone screws or components of the fixation frame.
16. Tissue necrosis secondary to the insertion of bone screws.
17. Pressure on the skin caused by external components when the free space is insufficient.
18. Limb dysmetria.
19. Excessive surgical bleeding.
20. Intrinsic risks associated with anesthesia.
21. Unmanageable pain.
22. Secondary bone sequestration due to rapid perforation of the cortex with accumulation of heat and bone necrosis.
23. Vascular disorders, including thrombophlebitis, pulmonary embolism, wound hematoma, avascular necrosis.

**Warning: This device is not approved for fixing or attachment by means of screws to posterior elements (pedicles) of the cervical, thoracic or lumbar spinal column.**

### **Important**

Not all surgical cases are resolved with a positive result. Additional complications may arise at any time due to improper use, for medical reasons or as a result of failure of the device, requiring surgery in order to remove or replace the external fixing device.

Presurgical and surgical procedures, including knowledge of the surgical techniques and the correct selection and positioning of the external fixation devices, are important aspects for the successful application of Orthofix external fixation devices by the surgeon.

Proper selection of the patient and their ability to comply with the physician's instructions and follow the prescribed treatment regimen will have a significant influence on the results. It is important to conduct patient screening and select the optimal therapy, while bearing in mind the requirements and/or limitations in terms of

physical and/or mental activity. Should a candidate for surgery have any contraindications or be predisposed to any contraindications, DO NOT USE any Orthofix external fixation device.

## Materials

The Orthofix TrueLok™ External Fixation System is comprised of components made from stainless steel, aluminium alloy and plastic. The components coming into contact with the patient's body are the percutaneous screws (bone screws), the K-wires, the drill bits, the guides used during insertion of the screws, the trocars and bone depth measurement devices; said components are made from surgical device grade stainless steel. Certain of the bone screws used in the Orthofix External Fixation are supplied with a thin coating of hydroxyapatite (HA), applied as a spray, over the threaded portion of the shaft.

## STERILE and NON STERILE products

Orthofix supply certain external fixation devices as STERILE and others as NON STERILE. It is recommended the product label be checked in order to determine the sterility of each device.

### Sterile

The devices or kits supplied as STERILE carry a label indicating said status.

**The contents of the package are STERILE unless the package is opened or damaged. Do not use the device if the package has been opened or damaged.**

### Non Sterile

Unless indicated otherwise, the components of Orthofix external fixation devices are provided as NON STERILE. Orthofix recommends that all non sterile components be properly cleaned and sterilized, following the specific cleaning and sterilization procedures recommended. The integrity and performance of the product is only guaranteed in the case the package is undamaged.

## Cleaning and Maintenance (Please note, see instructions PQ ISP)

Prior to use, any NON STERILE product should be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water. After cleaning, the device and/or system is duly rinsed in sterile distilled water and dried using a clean non-woven fabric cloth.

Prior to sterilization, it is recommended all components be inspected, since any damage to the surface of the metal can reduce the fatigue strength and the robustness and lead to corrosion. Should the components be damaged in any way, they must be replaced with new ones.

It is recommended the fixation device be assembled in order to ensure that all the components are present.

**ALL DEVICES LABELED "SINGLE-USE ONLY" MUST NOT BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR THE SAFETY AND EFFICACY OF SINGLE USE COMPONENTS WITH THE FIRST PATIENT.**

The institution or the user will be responsible for any subsequent use of said devices.

### Sterilization

The recommended sterilization cycle is:

Method	Cycle	Temperature	Exposure time
Steam	Pre-vacuum	132°-135°C [270°- 275°F]	Minimum 10 minutes

Sterility cannot be guaranteed if the sterilization tray is overloaded.

Do not overload the sterilization tray or include additional systems or instruments.

**Orthofix is only responsible for the safety and efficacy in the case of first use of the external fixation devices by the patient.**

**The responsibility for any subsequent use will be entirely that of the institution or the attending physician.**

#### **PLEASE NOTE**

In the case where the normal sensitivity of the limb is altered, abating the usual proprioceptive feedback, the fixation system may be subjected to greater loads than usual. Under similar circumstances, the patient will have to be informed of the risk of excessive stress on the fixation system, and the physician should be especially attentive regarding problems associated with overloading which might lead to loosening, bending or breakage of the components. In such situations, it is recommended the rigidity of the fixation system be increased compared to what is usually necessary.

#### **Implant removal**

The external fixation devices are designed to stabilize the fracture site during the rehabilitation and/or limb lengthening process. Once the limb rehabilitation and/or lengthening is complete, the devices must be removed. Removal may also be recommended in other cases, for example if there are signs of corrosion, reaction or localized pain; implant breakage, curvature, fracture or loosening; or bone loss.

#### **Ordering information**

Contact your local Orthofix sales representative.

**Attention:** US federal law restricts the sale of this device exclusively to or on the order of a physician or a hospital.

TrueLok™ is a registered trademark of Orthofix Inc.

**Note:** the physician is fully responsible for the correct selection of the patient, proper training, experience in selecting and applying TrueLok™, as well as for the selection of the subsequent post-operative procedures. For complaints or additional information regarding the product and its use, contact Orthofix Srl.

## RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

### IMPLANTABLE DEVICE\*

The "SINGLE USE" implantable device\* of Orthofix is identified through symbol  reported on the product label.

After the removal from the patient, the implantable device\* has to be dismantled.

The re-use of implantable device\* introduces contamination risks for users and patients.

The re-use of implantable device\* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\*): Implantable device

Any device which is intended:

Any device intended to be totally / partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is also considered an implantable device.

### NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol  reported on the label or are indicated in the "Instructions For Use" supplied with the products.

The re-use of "SINGLE USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

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All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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