DENTSPLY
TRUBYTE

Triad[®]

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DENTSPLY
TRUBYTE

Manufactured by:
DENTSPLY International Inc.
570 West College Ave.
York, PA 17405-0872
1-800-786-0085

www.trubyte.dentsply.com or www.dentsply.com

CE REP

DeguDent GmbH
Postfach 1364 • 63403 Hanau
GERMANY
Telefon +49/6181/5950

Indications for Use:

TRIAD Denture Base Material is a visible light cured (VLC) resin designed for use in the fabrication of full and partial dentures, baseplates, orthodontic appliances and appliance repairs. This VLC material is indicated as an alternative to traditional heat-cured and auto polymerizing resins.

Triad Denture Base is a urethane dimethacrylate-based pre-mixed dough that contains no methyl methacrylate monomer.

Contraindications:

Triad Denture Base Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).

Warnings:

Triad Denture Base Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

Precautions:

1. Users with special skin problems, cuts or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation,

masks and vacuum systems should be used.

3. All cure times cited are appropriate for Triad 2000 and Triad II VLC curing units. Consult your Technique Manual or call Technical Support for Triad I cure times.
4. Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.
5. Unused Triad Materials should be polymerized prior to disposal.
6. Do not use Triad Materials that crack or craze when handled.
7. Thermoplastic denture bases may be distorted in the Triad Unit.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

STEP-BY-STEP INSTRUCTIONS

Functional Orthodontic Appliances

1. Apply a coating of Triad Model Release Agent (MRA) to the cast.
2. Position and secure wires on the cast.

3. Adapt Triad Denture Base Material to the cast with finger pressure and the modeling tool. (Adapt a small amount of Denture Base Material under and/or around the wires before applying a sheet or rope of the material).
4. Trim to the desired shape with a sharp blade.
5. Apply Triad Air Barrier Coating (ABC) to all exposed surfaces of the VLC material.
6. Place the appliance in the Triad VLC Unit for a four (4) minute preliminary cure.
7. Following the preliminary cure, paint the tissue side of the appliance with Triad Air Barrier Coating and cure, tissue side up, for six (6) minutes.
8. Gently scrub the appliance with soap, warm water, and a soft brush to remove all traces of Triad Air Barrier Coating.
9. Trim and polish as required.

Denture Repairs

NOTE: This method is appropriate for both VLC dentures and acrylic dentures.

1. Prepare a cast and coat it with Triad Model Release Agent (MRA).
2. Prepare the fracture area by opening the break to approximately 1/4 inch. Bevel the edge of the opening with greater width on the outer side of the denture. Prime roughened surfaces of the repair areas with Triad VLC Bonding

Agent or Lucitone® Liquid. *If Bonding Agent is used, bench set for two (2) minutes, then cure in the Triad VLC Unit for two (2) minutes.*

3. If teeth are being replaced, clean wax from surfaces. Lightly grind areas that were in contact with wax. Paint bonding surfaces of teeth and the denture base with Triad VLC Bonding Agent. *Bench set two (2) minutes, then cure in the Triad Curing Unit for two (2) minutes.*
4. Apply Triad Denture Base Material to the appropriate areas. Coat all surfaces with Triad Air Barrier Coating (ABC) and cure for four (4) minutes. Carefully remove the denture from the cast. Paint the tissue surface with Triad Air Barrier Coating and cure in the Triad VLC Unit for six (6) minutes, tissue side up.
5. Scrub softly with warm water and a brush to remove all traces of Triad ABC.
6. Finish and polish.

NOTE: VLC dentures may also be repaired with Lucitone 199 Repair Material or DENTSPLY Repair Material. Follow instructions furnished with the repair materials.

Record Bases

Use Triad Denture Base Material to construct Record Bases using the following technique which is similar to making impression trays.

1. Utilizing the master cast, block out undesirable undercuts —

refer to custom tray technique — however, do not add any wax for spacing.

2. Paint a thin coating of Triad Model Release Agent (MRA) on the cast and relief wax to aid in the removal of the record base from the cast.
3. Carefully adapt Triad Denture Base Material to the cast. Do not overthin material on ridge portion or hamular notch areas of the cast.
4. *Process in the Triad VLC Unit for two (2) minutes to set the Triad Material.* Separate the record base from the cast and remove all remaining wax from the record base. Apply Triad Air Barrier Coating (ABC) to the entire surface of the record base and cure in the Triad VLC Unit, tissue side up, for an additional six (6) minutes.
5. Carbide Cutters are recommended for adjusting and finishing borders to desired outlines. Polish as usual.

Complete Denture Procedure Using Individual Artificial Teeth And The Triad System

NOTE: The use of mechanical retention is strongly recommended for all denture teeth used in implant prostheses with Triad Denture Base.

The following technique assumes that a set-up and wax-up have been accomplished.

1. Form three notches in the land area of the cast, one in the anterior midline area and one on each side in the posterior,

distal to hamular notch areas. Wax-up may be lightly sealed to the cast. Wax out undercuts in the areas of the locating notches.

2. Apply Triad Model Release Agent (MRA) to incisal/occlusal areas of the teeth, to the notches formed in the land area of the cast, and to wax areas which may be covered with Triad VLC material.
3. Adapt a rope of Triad Denture Base to adequately cover incisal/occlusal surfaces of the teeth, providing definite support and positive seating for the teeth. The Denture Base Material should not be adapted gingivally beyond the greatest convexity of the tooth. Extend each end of the rope to engage each notch in the land area of the cast. Also attach the Denture Base Material in the anterior area, extending from the notch in the land area to the material adapted over the anterior teeth.
4. Paint the VLC material matrix with Triad Air Barrier Coating (ABC) and cure in the Triad VLC Unit for five (5) minutes.

The amount of heat generated during this five minute curing cycle will soften the wax to permit separation of the matrix and teeth from the wax-up and cast.

5. Thoroughly eliminate any traces of wax remaining on the teeth or the matrix with clean boiling water and detergent. Remove any remaining wax from the cast or Triad integral baseplate, if used. The cast, teeth and matrix (and integral baseplate) should

- have a final rinse with clear boiling water.
- (Teeth should remain in the matrix if sufficient coverage and support of the buccal/labial/lingual surfaces was provided during fabrication of the matrix.)
6. Incorporate a posterior palatal seal on the cast at this time.
 7. Apply a coating of Triad Model Release Agent (MRA) to the cast.
 8. During adaptation of the Triad Denture Base Material, maintain as uniform a thickness as possible over the entire cast. Provide sufficient space between teeth and cast for the denture base by grinding ridge lap areas. *Cure the Triad baseplate at least two (2) minutes in the TRIAD VLC Unit.*
 9. Paint Triad VLC Bonding Agent on the ridge lap areas of the teeth and also about 2 mm onto the facial and lingual collars and interproximal areas. Bonding areas on the baseplate may also be painted with Triad VLC Bonding Agent. *Bench set for two (2) minutes, then cure in the Triad VLC Unit for two (2) minutes. (Adjust height of the curing table so lamps are properly focused on teeth to achieve maximum cure.)*
 10. Apply Triad Denture Base Material to form a union between the teeth and the baseplate. Relate the teeth to the baseplate by repositioning the matrix. Excess Denture
- Base Material should be removed before curing.
11. *Tack teeth to the denture base with a two (2) minute cure cycle.*
 12. Following the cure cycle, remove the matrix.
 13. Paint a band of Triad VLC Bonding Agent on the junction line between the teeth and the Denture Base Material. *Allow to bench set two (2) minutes, then cure for two (2) minutes.*
 14. Contour the lingual, then buccal surfaces with Triad Denture Base Material. Careful contouring and trimming at this stage will save time later during finishing.
 15. When contoured, paint the denture with Triad Air Barrier Coating (ABC) and *cure five (5) minutes in the Triad VLC Unit.*
- Remove the cast from the denture and paint the tissue side of the appliance with Triad ABC. *Cure in the Triad VLC Unit for five (5) minutes.*
16. When cured, wash dentures to remove all traces of Air Barrier Coating.
- NOTE:** Recommended total cumulative curing time is a minimum of ten (10) minutes. Half of the curing should be done tissue side up, half completed tissue side down.
17. Finish and polish by normal methods.

Removable Partial Dentures

1. Apply a coating of Triad Model Release Agent (MRA) to denture base areas of the cast.
 2. Apply Triad Denture Base Material to the edentulous areas. Trim to the general outline with a sharp blade.
 3. Seat framework firmly, embedding in the uncured Triad Denture Base Material. Remove excess material.
 4. *Cure the denture base material in the Triad VLC Unit for two (2) minutes.*
 5. Adjust the teeth to fit the denture base. Trim ridge lap areas of the teeth to provide as much space as possible for the base material between the teeth and the retentive mesh.
 6. Clean the tooth surfaces to be bonded to the denture material. Apply Triad VLC Bonding Agent to the tooth surfaces. *Bench set two (2) minutes and cure two (2) minutes in the Triad VLC Unit.*
7. Adapt Triad Denture Base Material to each tooth. Position teeth and remove excess Denture Base Material. Tack teeth to the Triad Denture Base Material using a hand-held light.
 8. Paint a band of Triad VLC Bonding Agent between the teeth and the Denture Base Material. *Bench set two (2) minutes and cure in the Triad VLC Unit for two (2) minutes.*
 9. Complete the contouring with Triad Denture Base Material.
 10. Paint the denture material with Air Barrier Coating (ABC) and *cure for two (2) minutes in the Triad VLC Unit.*
 11. Remove the cast from the partial denture. Paint tissue surfaces with Triad ABC and *cure in the Triad VLC Unit, tissue side up, for six (6) minutes.*
 12. Wash, trim, and polish.

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