Selectra Catheters

Guiding Catheters for the Selectra Telescopic CS Lead Introducer System	Technical Manual • en
Katétrý pro teleskopický systém Selectra k zavádění elektrod do koronárního sinu	Technická příručka • cs
Kateter til CS-elektrode-teleskop-indføringsystemet Selectra	Brugermanual • da
Katheter für das CS-Elektroden-Teleskop-Einführungssystem Selectra	Gebrauchsanweisung • de
Catéter guía para el sistema telescópico de introducción de electrodos por el SC Selectra	Manual técnico • es
Ohjainkatetri cs-johdon teleskooppiselle sisäänvientijärjestelmälle Selectra	Käyttöohje • fi
Gaines pour le système d'implantation télescopique de sonde SC Selectra	Manuel technique • fr
Catetere guida per il sistema di introduzione per seno coronarico Selectra	Manuale tecnico di istruzione • it
Katheter voor het CS-elektroden-telescoop-inbrengsysteem Selectra	Gebbruikershandleiding • nl
Cewnik do teleskopowego systemu wprowadzania elektrod CS Selectra	Instrukcja obsługi • pl
Cateter para o sistema telescópico de introdução de eletrodos para seio coronário Selectra	Manual técnico • pt
Kateter för CS-elektrod-teleskopinföringsystemet Selectra	Bruksanvisning • sv

en	• English	2
cs	• Český	7
da	• Dansk	12
de	• Deutsch	17
es	• Español	22
fi	• Suomi	27
fr	• Français	32
it	• Italiano	37
nl	• Nederlands	42
pl	• Polski	47
pt	• Português	52
sv	• Svenska	58

Target group

This manual is intended for cardiologists, electrophysiologists and surgeons who are familiar with and have practical experience with cardiac resynchronization therapy and implantation of leads in the coronary venous system.

Description

BIOTRONIK's Selectra CS lead introducer system is a combination of guiding catheters and implantation accessories used to facilitate access to the coronary veinous system for suitable leads and catheters.

The Selectra CS lead introducer system consists of several individually available guiding catheters with various different curve shapes and the Selectra accessory kit.

The catheters are available as inner and outer catheters which jointly form a telescope system.

Please refer to the Selectra accessory kit manual to get an overview of the different catheter forms and designs that are available.

Using the Selectra system has the following advantageous options:

- Quick probing of the ostium of the coronary sinus
- Easy access to the coronary sinus
- Introduction of contrast media for angiography
- Easy and fast positioning and repositioning of the leads

Observe other manuals

This manual is intended to supplement the Selectra accessory kit manual.

Please refer to the accessory kit manual for all instructions pertaining to the safe handling of the guiding catheter as part of using the Selectra lead introducer system.

Moreover, the manual of the lead to be implanted is to be observed.

Intended use

In conjunction with the Selectra accessory kit, Selectra guiding catheters are used to facilitate lead implantation in the left side of the heart via the coronary sinus.

Product line for the Selectra CS lead introducer system

The Selectra CS lead introducer system consists of the accessory kit described in this technical manual and a selection of compatible guiding catheters with different curve shapes and working lengths. Table 1 shows our product line.

Figure	Item	Order number
	Outer guiding catheters for Selectra	
	Selectra Amplatz 6.0-45	375 519
	Selectra Amplatz 6.0-55	375 520
	Selectra BIO 2-45	375 523
	Selectra BIO 2-55	375 524
	Selectra Extended Hook-45	375 527
	Selectra Extended Hook-55	375 528


Figure	Item	Order number
	Selectra Hook-45	375 529
	Selectra Hook-55	375 530
	Selectra MPEP-45	375 531
	Selectra MPEP-55	375 532
	Selectra Multipurpose Hook-45	375 533
	Selectra Multipurpose Hook-55	375 534
	Selectra Right-45	375 535
	Selectra Right-55	375 536
	Selectra Straight-45	375 537
	Selectra Straight-55	375 521
Inner guiding catheters for Selectra		
	Selectra IC-50-65	375 545
	Selectra IC-50-75	375 546
	Selectra IC-90-65	375 547
	Selectra IC-90-75	375 548
	Selectra Accessory Kit	375 518
More accessories		
	Introducer Set HLS-1009 SafeSheath 9 F	343 594
	Steerable EP catheter ViaCath NG 4/S/5mm (not available in the U.S.)	351 197
	Venogram balloon catheter Corodyn P1 (not available in the U.S.)	336 074

Table 1: Available components and accessories for the Selectra 7 F CS lead introducer system

Compatible leads and accessories

The individually available Selectra guiding catheters with a nominal inner diameter of 7 French (outer catheter) and 5 French (inner catheter) can only be used in combination with the accessory kit of the Selectra CS lead introducer system.

The outer guiding catheters and the implantation accessory are designed for implanting suitable leads with polyurethane coatings or coatings with comparable anti-friction properties and a diameter of less than 6.6 F.

The inner catheters can be used to implant leads with polyurethane coatings or coatings with comparable anti-friction properties and a diameter of less than 5 F.

Not for multiple usage

This implant accessory is intended for single use only!

Do not resterilize!

Reuse of implant accessories can result in infections, embolisms and damage to the device.

Contraindications

Use of the Selectra CS lead introducer system is contraindicated for the following:

- Patients with an existing or possible occlusion of the coronary vessels or unsuitable anatomy of the coronary veins
- Patients with active systemic infection

Potential negative side effects

Possible complications are as follows:

- Allergic reactions to contrast media
- Hematoma
- Bleeding
- Infection
- Embolism
- Pneumothorax
- Cardiac tamponade
- Myocardial damage
- Venous or cardiac perforation
- Local tissue response, fibrotic tissue formation
- Damage to the heart valves
- Vascular occlusion
- Chronic nerve damage

Storage

The product should only be stored at temperatures between 5 and 55°C (41°F to 131°F).

Packaging and sterility

The Selectra guiding catheter is delivered in a blister with a lid. This blister is placed in two sterile bags so that the inside bag and the blister are also sterile on the outside. Ethylene oxide is used for sterilization purposes.

The sterile bags with the blister are packaged in a box that bears a quality control seal and a product information label. The label lists the model type, technical data, use by date, and information about sterility and storage of the package and content.

- 1 To ensure sterility, the container should be checked for damages prior to opening.
- 2 If it is suspected that the sterile packaging has been opened or damaged, return the product to BIOTRONIK.

Scope of delivery

- 1 guiding catheter (sterile)
- 1 dilator (sterile) - for the outer catheter only
- 1 technical manual (non-sterile)

Opening the packaging

In a non-sterile area:

- Remove the sterile bag from the box and open it.

In a sterile area:

- Remove and open the inner sterile bag.
- Remove and open the blister.

Use and handling of the catheters

Use of the Selectra guiding catheter in combination with other implantation accessories and a suitable CS lead is described in the manuals for the Selectra accessory kit and the lead.

Caution: Rinse thoroughly before use.**- Avoid excessive sliding friction on the interior wall of the catheter**

- Before use, rinse with sufficient physiological NaCl solution to activate the hydrophilic inner coating of the catheter.

In addition, this prevents air from being trapped.

Disclaimer

BIOTRONIK's Selectra CS lead introducer system (including all available accessories), used in conjunction with devices for implantation in the coronary sinus, has been qualified, manufactured, and tested according to proven and accepted standards and procedures. Nonetheless, the physician must be aware that the accessories may easily be damaged through improper handling or use. Except as set forth in BIOTRONIK's limited warranty, BIOTRONIK makes no expressed or implied warranties for its accessories.

Technical data

Common features of the Selectra catheters are as follows:




- Color of the shaft:
Different shades of blue for segments with distinct rigidity.
- Ergonomic handle with Luer lock connector side port and with integrated hemostatic silicone valve in main port.
- Handle with integrated valve can be slit open using a slit tool.
- Hydrophilic inner coating for improved gliding properties.

Outer guiding catheters	Models "...-45"	Models "...-55"
Working length	45 cm	55 cm
Inner diameter	2.44 mm (7.3 F)	
Inner material	Hydrophilic coating	
Outer diameter	2.91 mm (8.7 F)	
Material	PEBAX	
X-ray contrast	Barium sulfate in the entire catheter area	
Radiopaque marker at the distal end	Black marker zone (width: 20 mm)	
Soft tip length	1.0 mm (blue)	

Inner guiding catheters	Models "...-65"	Models "...-75"
Working length	65 cm	75 cm
Inner diameter	1.83 mm (5.4 F)	
Inner material	Hydrophilic coating	
Outer diameter	2.31 mm (6.9 F)	
Material	PEBAX	
X-ray contrast	Barium sulfate in the entire catheter area	
Radiopaque marker on the distal end	Black marker zone (width: 20 mm)	
Soft tip length	1.0 mm (blue)	

Dilator	Models "...-45"	Models "...-55"
Working length	53 cm	63 cm
Material	PE	
Outer diameter	2.3 mm (6.9 F)	

Legend for the label

	Manufacturing date of the device		Contents
	Use by		Do not use if package is damaged!
	Storage temperature		European approval mark
	BIOTRONIK order number		Inner diameter
	Lot number		Outer diameter
	Sterilized with ethylene oxide		Total length
	Do not resterilize!		Working length
	Single use only. Do not reuse!		Dilator
	Non-sterile		Guiding catheter (Selectra), various curve shapes (example)
	Consult the instructions for use		Caution: Federal (U.S.A.) law restricts this product to sale by, or on the order of, a physician.

CE mark

This product conforms with EC Directive 90/385/EEC relating to active implantable medical devices. It is therefore designated with the CE mark. The product can be used in all European Union countries as well as in countries that recognize the above-mentioned directive.

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All rights reserved.
BIOTRONIK, Inc.