BioMonitor

Implantable Cardiac Monitor with Home Monitoring
Implantovatelný monitor srdečního rytmu s funkcí Home Monitoring
Implanterbar hjertemonitor med Home Monitoring
Implantierbarer Herzmonitor mit Home Monitoring
Monitor cardiaco implantable con Home Monitoring
Home Monitoring -järjestelmällä varustettu asennettava sydänmonitori
Moniteur cardiaque implantable avec Téléc@rdiologie
Monitor cardiaco impiantabile con Home Monitoring
Implanteerbare cardiale monitor met Home Monitoring
Wszczepialny monitor rytmu serca z Home Monitoring
Monitor cardíaco implantável com Home Monitoring
Implanterbar hjärtmonitor med Home Monitoring

Technical Manual • en
Technická příručka • cs
Brugermanual • da
Gebrauchsanweisung • de
Manual técnico • es
Käyttöohje • fi
Manuel technique • fr
Manuale tecnico di istruzione • it
Gebruikershandleiding • nl
Instrukcja obsługi • pl
Manual técnico • pt
Bruksanvisning • sv

BIOTRONIK

excellence for life

en	• English	2
cs	• Český	11
da	• Dansk	20
de	• Deutsch	29
es	• Español	39
fi	• Suomi	49
fr	• Français	58
it	• Italiano	68
nl	• Nederlands	78
pl	• Polski	88
pt	• Português	98
sv	• Svenska	108

Product Description

Intended Medical Use

Intended use

BioMonitor is the name of an implantable cardiac monitor for the monitoring of heart rhythm.

Its primary purpose is to provide early detection and diagnostics of symptoms of arrhythmias, such as atrial fibrillation and the causes of syncope, which can be manifested clinically.

Note: BioMonitor does not have a pacing function.

Form of diagnosis

The heart rhythm is continuously monitored; the possible detection types are atrial fibrillation, high ventricular rate, asystole, or bradycardia. Depending on the preset parameters, subcutaneous ECGs and other data are recorded.

The patient can also initiate a recording of a subcutaneous ECG.

BIOTRONIK Home Monitoring® enables physicians to perform a multi-year diagnosis management at any time.

Required expertise

In addition to having basic medical knowledge, the user must be thoroughly familiar with the operation of the implantable cardiac monitor and the specified conditions for its use. Only qualified medical specialists having the special knowledge are permitted to implant the BioMonitor and make diagnoses.

Guidelines of cardiological societies

It is recommended that the indications published by the German Cardiac Society (Deutsche Gesellschaft für Kardiologie, Herz- und Kreislau fforschung) and the ESC (European Society of Cardiology) are observed. This also applies to the guidelines published by the Heart Rhythm Society (HRS), the American College of Cardiology (ACC), the American Heart Association (AHA), and other national cardiology associations.

Indications

- Clinical symptoms or increased risk of cardiac arrhythmias
- Temporary symptoms that may indicate cardiac arrhythmia

Contraindications

There are no known contraindications.

However, the particular patient's state of health determines whether a subcutaneous device will be tolerated long-term.

System Overview

Device family

BioMonitor is a cardiac monitor. This device is not available in all countries.

Parts

The system consists of the following parts:


- Device with electrodes integrated in the header
- Programmer and M-50 permanent magnet
- Current device program

Cardiac monitor

The device's housing is made of biocompatible titanium, welded from the outside and therefore hermetically sealed. It is coated in silicone. The ellipsoid shape facilitates ingrowth into the pectoral muscle area.

The header contains 2 subcutaneous electrodes. An area on the bottom of the housing that is recessed from the silicone sheath serves as the third pole. The device labeling provides information about the BioMonitor.

The cardiac monitor's sensing principle is based on three vectors:

The resulting signals are filtered and summed to yield a meaningful signal:

Programmer

Implantation and follow-up are performed with a portable BIOTRONIK programmer. The programmer is used during implantation to transfer the current device program to the device. In addition, the programmer is used to set parameter combinations, as well as for interrogation and saving of data from the device. Subcutaneous ECG, markers, and functions are displayed simultaneously on the color display.

Note: It is not permitted to use the device's ECG display for diagnostic purposes because it does not meet all requirements for diagnostic ECG devices (IEC 60601-2-25).

Telemetry

Telemetric communication between the device and the programmer can be carried out by applying a programming head (PGH).

BIOTRONIK Home Monitoring®

The BIOTRONIK cardiac monitor provides a complete diagnosis management system:

- With Home Monitoring, diagnostic information as well as technical data of the device are automatically and wirelessly sent to a stationary or mobile transmitter via an antenna in the device header. The data are encrypted and sent from the transmitter to the BIOTRONIK Service Center via the cellular phone network.
- The received data are deciphered and evaluated. Each physician can set the criteria for evaluation to be used for each patient and can configure the time of notification via E-mail, SMS or fax.
- A clear overview of the results of this analysis is displayed for the attending physicians on the protected Internet platform Home Monitoring Service Center (HMSC).
- Data transmission from the device is performed with a daily device message.
- Device messages which indicate special events in the heart or in the device are forwarded at the set time.
- A test message can be initiated at any time using the programmer to immediately check the Home Monitoring function.

Technical manuals

The following technical manuals provide information about usage of the device systems:

- Technical manual for the device
- Technical manual for the HMSC

- Technical manuals for the programmer
- Technical manual for device programs as online help on the user interface and as a PDF file in the Manual Library at www.BIOTRONIK.com

BioMonitor order numbers

Not all device types are available in all countries:

Device	Order number
BioMonitor	394119

Scope of delivery

The storage package includes the following:

- Sterile container with device
- Serial number label
- Patient ID card
- Technical manual for the device

The sterile container includes the following:

- Device

Diagnostic Functions

General overview

- Automatic functions make it easy and fast to implant, configure, and check the BioMonitor.

Detection and data storage

- The signals are recorded and stored once a detection type is set.
- Multiple detection types can be set simultaneously.
- Individual episodes of up to 40 s can be stored. Recordings triggered by the patient of up to 7.5 min can be stored.
- The device can record a total of 35.8 min of episodes in subcutaneous ECGs.
- When performing follow-ups using the programmer, the subcutaneous ECG is indicated with markers after applying the programming head.

Home Monitoring functions

Important medical information include, among others, the following:

- Ongoing atrial and ventricular arrhythmia
- Current statistics
- High definition subcutaneous ECGs and sending of these recordings with device messages

General Safety Instructions

Operating Conditions

Care during shipping and storage

- Devices are not to be stored or transported close to magnets or sources of electromagnetic interference.
- Observe the use by date (see the battery data).

Temperature

Extremely low and high temperatures affect the service time of the battery in the device.

- The following temperatures are permitted for transport, storage, and use:
-10°C to 45°C.

Sterile delivery

The device is delivered after sterilization using gas. Sterility is guaranteed only if the blister and quality control seal have not been damaged.

Sterile container

The device is packaged in 2 separately sealed blisters. The inner blister is also sterile on the outside so that it can be transferred in a sterile state during implantation.

Single use only

The device is only intended for one-time use.

- Do not use the device if the package is damaged.
- The device must not be resterilized and reused.

Possible Complications

General information on medical complications

Complications for patients and device systems generally recognized among practitioners also apply to BIOTRONIK devices.

- Complications may include, for example, foreign body rejection phenomena, local tissue reactions, migration of the device, or infections. Primary sources of complication information include current scientific and technological knowledge.
- It is impossible to guarantee the efficacy of implantable cardiac monitors.

Possible technical failures

Technical failure of a device system cannot be entirely ruled out. Possible causes can include the following:

- Device component failures
- Battery depletion

Electromagnetic interference (EMI)

Any device can be sensitive to interference, for example, when external signals are sensed as intrinsic rhythm.

- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimal.
- Under unfavorable conditions, especially during diagnostic procedures, sources of interference may induce such a high level of energy into the device that the data recording can be influenced or the device may be damaged.

Possible Risks

Contraindicated procedures

Due to possible damage to the patient or the device, the use of the following procedures is contraindicated:

- Therapeutic ultrasound: damage to the patient via excess warming of body tissue near the device
- Transcutaneous electrical nerve stimulation
- Hyperbaric oxygen therapy
- Applied pressures higher than normal pressure

Potentially risky therapeutic and diagnostic procedures

If energy from an external source is introduced to the body for diagnostic or therapeutic purposes, then the device can be subjected to interference and its function can be impaired.

Harmful effects may occur, for example, during electrocautery, HF ablation, HF surgery, or lithotripsy.

If potentially risky procedures cannot be avoided, the following should be observed at all times:

- Electrically insulate the patient.
- Do not introduce energy near the device.
- Monitor the patient during and after every intervention.

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. The sensing properties may change as a result.

- Place adhesive electrodes anterior-posterior or perpendicular to the axis formed by the device to the heart at least 10 cm away from the device and from implanted leads.

Radiation therapy

The use of radiation therapy is contraindicated due to possible damage to the device and the resulting impaired functional safety. If this type of therapy is to be used anyway, prior risk/benefit analysis is absolutely necessary. The complexity of influencing factors such as different sources of radiation, a variety of devices and therapy conditions makes it impossible to issue directives that guarantee radiation therapy without an impact on the device. The EN 45502 standard pertaining to active implantable medical devices requires the following information in combination with therapeutic ionic radiation:

- Adhere to instructions for potentially risky therapeutic and diagnostic procedures.
- Shield device against radiation.
- After applying radiation, double-check the device system to make sure it is functioning properly.

Note: Please contact BIOTRONIK with questions on the risk/benefit analysis.

Magnetic resonance imaging

Magnetic resonance imaging is contraindicated due to the high frequency fields and the associated magnetic flux density: damage or destruction of the device by strong magnetic interaction and damage to the patient by excessive warming of the body tissue in the area surrounding the device.

Implantation

Implantation Procedure

Having parts ready

The following parts that correspond to the requirements of EC Directive 90/385/EEC are required:

- BIOTRONIK device
- BIOTRONIK programmer and approved cables
- External multi-channel ECG device
- Keep spare parts for all sterile components.

Keeping an external defibrillator ready

In order to be able to respond to unforeseeable emergencies or possible technical failures of the device:

- Keep an external defibrillator and paddles or patch electrodes ready.

Unpacking the device

WARNING

Inadequate function due to defective device

If an unpacked device is dropped on a hard surface during handling, electronic parts could be damaged.

- Use a replacement device.
- Send the damaged device to BIOTRONIK.

- Peel the sealing paper off of the outer blister at the marked position in the direction indicated by the arrow. The inner plastic container must not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments.
- Take hold of the inner blister by the gripping tab and take it out of the outer blister.
- Peel the sealing paper off of the sterile inner blister at the marked position in the direction indicated by the arrow.

Implantation site

- Depending on the patient's anatomy, the BIOTRONIK cardiac monitor can be implanted subpectorally or subcutaneously on the left side.

Preventing leakage currents

Leakage currents between the tools and the device must be prevented during implantation.

- Electrically insulate the patient.

Implanting

1	Form the device pocket.
2	Insert the device.
3	Position the device optimally.
4	Guide the fixation suture through the 2 openings in the header and fixate the device in the prepared device pocket.
5	Check the signals.
6	Close the device pocket.

Position the device optimally

While the device is being positioned, the signals are recorded by the programming head. It must therefore be placed in the surgical field in sterile condition.

- Use the programming head with a suitable sterile cover.

⚠ CAUTION

Poor signal quality

Unreliable signals may impair detection and even lead to false detections.

- Try different positions for the implanted device until the best signals are displayed on the programmer.

Applying the programming head

The programming head features a diagram of the device. This is used to assist in positioning the head to ensure proper telemetry.

- Make sure the PGH is positioned correctly.

Establishing telemetry contact

- When the programming head is applied, time remains for device interrogation. All detection parameters are disabled during this time.

Activate diagnostics

- Load the suitable device program in the programmer.

Follow-up

Follow-up intervals

Follow-ups must be performed at regular, agreed intervals.

- Following the ingrowth phase, approximately 3 months after implantation, the first follow-up should be carried out by the physician using the programmer (in-office follow-up).
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

Follow-up with BIOTRONIK Home Monitoring®

Monitoring using the Home Monitoring function does not serve to replace regular in-office appointments with the physician required for other medical reasons. Follow-up supported by Home Monitoring can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted if symptoms worsen or if new symptoms arise despite use of the Home Monitoring function.
- Device messages are transmitted regularly.
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device system are sufficient. If not, an in-office follow-up has to be carried out.

Possible early detection due to information gained via Home Monitoring may necessitate an additional in-office follow-up. For example, the data may provide early indication of a foreseeable end of service time (ERI). Furthermore, the data could provide indications of previously unrecognized arrhythmias.

Follow-up with the programmer

Use the following procedure for in-office follow-up:

1	Record and evaluate the ECG
2	Interrogate the device.
3	Check the sensing function.
4	Evaluate the status and automatically measured follow-up data.
5	Possibly evaluate statistics and subcutaneous ECGs.
6	If necessary, customize the program functions and parameters.
7	Transmit the program permanently to the device.
8	Print and document follow-up data (print report).
9	Finish the follow-up for this patient.

Patient Information

Patient ID card

A patient ID card is included in delivery.

- Provide the patient with the patient ID.
- Request that patients contact the physician in case of uncertainties.

Prohibitive signs

Places with prohibitive signs must be avoided.

- Draw the patient's attention to prohibitive signs.

Possible sources of interference

Electromagnetic interference should be avoided in daily activities. Sources of interference should not be brought into close proximity with the device.

- Draw the patient's attention to special household appliances, safety locks, anti-theft devices, strong electromagnetic fields, cell phones, and transmitters among other things.
- Request patients to do the following:
 - Use cell phones on the opposite side of their body from the device.
 - Keep the cell phone at least 15 cm away from the device both during use and when storing.

Magnet application by patients

Patients need to know the following:

- Which magnet is used?
BIOTRONIK M-50 permanent magnet
- When may the magnet be used?
In case of symptoms such as severe dizziness, indisposition, or palpitations, as well as after fainting.
- How long is the magnet placed on the BioMonitor?
1 to 2 s.
- What happens when the magnet is applied?
The subcutaneous ECG of the last 7.5 min is stored and the last 40 s are transmitted at the daily transmission time via Home Monitoring.
- What has to happen after magnet application?
The patient has to contact the physician for a follow-up.

Replacement Indications

Possible battery levels

- BOS: Beginning of Service: > 70% charge
- ERI: Elective Replacement Indication
- EOS: End of Service

Elective Replacement Indication (ERI)

The device can no longer monitor the heart rhythm when in ERI charging status. ERI can be detected by Home Monitoring.

EOS replacement indication

The end of service has been reached.

Explantation and Device Replacement

Explantation

- Interrogate the device status.
- Remove the device using state-of-the-art technology.
- Explants are biologically contaminated and must be disposed of safely due to risk of infection.

Device replacement

Basic principles:

- The device must not be resterilized and reused.

Cremation

Devices should not be cremated.

- Explant the device before the cremation of a deceased patient.

Disposal

BIOTRONIK takes back used products for the purpose of environmentally safe disposal.

- Clean the explant with an at least 1% sodium hypochlorite solution.
- Rinse off with water.
- Send the cleaned explant to BIOTRONIK.

Parameters

Arrhythmia Detection Parameters

Diagnostics: Atrial fibrillation

The following can be set:

Parameter	Range of values	Standard	Factory
Atrial fibrillation (AF)	ON; OFF	ON	OFF
AF sensitivity	Low; medium; high	Medium	—

Additional AF expert parameters can be set:

Parameter	Range of values	Standard	Factory
RR variability limit	6.25; 12.5; 18.75%	12.5%	—
Analysis window onset/end	8/16; 16/24; 24/32 detection/end	8/16	—
Onset intervals	5 ... [2] ... 23 minimum number of consecutive RR cycles with an RR interval outside of the variability limit	5	—
Resolution intervals	1 ... [2] ... 7 maximum number of consecutive RR cycles with an RR interval outside of the variability limit	1	—
Confirmation time	1 ... [1] ... 5 min	2 min	—

Diagnostics: High ventricular rate

The following can be set:

Parameter	Range of values	Standard	Factory
High ventricular rate	ON; OFF	ON	OFF
HVR limit	150 ... [10] ... 200 bpm	180 bpm	—
HVR counter	4; 8; 12; 16	8	—

Diagnostics: Bradycardia

The following can be set:

Parameter	Range of values	Standard	Factory
Brady zone limit	OFF; 30 ... [5] ... 80 bpm	40 bpm	OFF
Brady duration	5 ... [5] ... 30 s	10 s	—
Brady rate decrease	OFF; 10 ... [10] ... 50%	30%	—
Brady sensitivity	Low; medium; high	Medium	—

Additional bradycardia expert parameters can be set:

Parameter	Range of values	Standard	Factory
Baseline intervals	32; 48; 64	48	—
Rate-drop intervals	4; 8; 16	16	—

Diagnostics: Asystole duration

The following can be set:

Parameter	Range of values	Standard	Factory
Asystole duration	OFF; 2 ... [1] ... 10 s	3 s	OFF

Diagnostics: Patient trigger

The following can be set:

Parameter	Range of values	Standard	Factory
Patient trigger	ON; OFF	ON	OFF

Global parameters

The following settings apply to all detection types:

Parameter	Range of values	Standard	Factory
Start resting period	00:00 ... [60] ... 23:00 hh:mm	2:00 hh:mm	—
Duration of resting period	00:30 ... [30] ... 12:00 hh:mm	4:00 hh:mm	—

Home Monitoring Parameter Settings

Home Monitoring: HM episode trigger

Home Monitoring can be set for all detection types:

Parameter	Range of values	Standard	Factory
Home Monitoring	ON; OFF	ON	OFF
HM episode trigger: AF; HVR; bradycardia; asystole; patient trigger	ON; OFF	ON	OFF

Note: Depending on the set detection type, the oldest, newest, and the longest episode is stored. With patient trigger, the oldest and the 2 most recent episodes are stored.

Global parameters

The following settings apply to all detection types:

Parameter	Range of values	Standard	Factory
Time of transmission	STD; 00:00 ... (30) ... 23:30 hh:mm	STD	
Periodic subcutaneous ECG	OFF; 1; 2; 30; 60; 90; 120; 180 days	30 days	OFF

Technical Data

Mechanical Characteristics

Measurements for the housing

Device	W x H x D [mm]	Volume [cm ³]	Mass [g]
BioMonitor	53.3 x 42.7 x 7.1	12.5	26

Note: D = housing without header

X-ray identification

BIO VP

Materials in contact with body tissue

- Housing: titanium
- Leads: titanium, fractally coated
- Housing sheath: silicone
- Silicone plug (screw cover): Silicone
- Header: epoxy resin

Electrical Characteristics

Components and input values

Electrical characteristics determined at 37°C, 500 Ω:

Circuit	Hybrid electronics with VLSI-CMOS chip
Input impedance	> 50 kΩ

Housing shape

The BioMonitor device housing has a flattened, ellipsoid shape.

Telemetry

Telemetry data for Home Monitoring:

Nominal carrier frequency	Maximum power of transmission
403.62 MHz	< 25 µW -16 dBm

Industry Canada

Telemetry data

- This device must not interfere with stations operating in the 400.150 - 406.000 MHz band in the meteorological aids, meteorological-satellite, and earth exploration-satellite services and must accept any interference received, including interference that may cause undesired operation.
- This device will be registered with Industry Canada under the following number: IC: 4708A-PRIMUS
- The code IC in front of the certification/ registration number only indicates that the technical requirements for Industry Canada are met.

Battery Data

Battery characteristics

The following data is provided by the manufacturers:

Manufacturer	GREATBATCH, INC. Clarence, NY 14031 USA	LITRONIK GmbH 01796 Pirna Germany
Battery type	GB 2596	LiS 3150M
System	Ag/SVO/CFx QMR®	LiMnO ₂
Device type	BioMonitor	
Battery voltage at BOS	3.0 V	3.1 V
Open-circuit voltage	3.0 V	3.1 V
Nominal capacity	1.3 Ah	1.2 Ah
Remaining capacity at ERI	0.25 Ah	0.15 Ah

Power consumption

- Battery current: 9 µA
- Average additional battery current in transmission mode: 4 µA

Average service time

Average service times are precalculated using the battery manufacturer's technical specifications and the setting of different detection parameters.

For BioMonitor: 4 years

Shortening of the service time after long storage period

Depending on the storage period, the service time from the beginning of service BOS to the replacement time ERI decreases as follows:

- After 1 year: by 8 months
- After 1.5 years: by 12 months

Legend for the Label

The label icons symbolize the following:

	Manufacturing date		Use by
	Storage temperature	REF	Order number
SN	Serial number	PID	Product identification number
CE	CE mark		
	Contents		Follow the instructions for use!

STERILE EO	Sterilized with ethylene oxide		
	Do not resterilize		Single use only. Do not re-use!
	Do not use if packaging is damaged		Non-sterile

	Transmitter with non-ionizing radiation at designated frequency
---	---

	Implantable cardiac monitor
Example	

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All rights reserved. Specifications subject
trademark of BIOTRONIK SE & Co. KG