
Modèle présenté : Table UNIVERSIS 6090 équipée d'une extension tête, dossier et une plaque jambes.
Description type: UNIVERSIS 6090 table equipped with head and back extension, and leg plate.
Modelo presentado: Mesa UNIVERSIS 6090 equipada con una extensión para la cabeza, el respaldo y las piernas.

UNIVERSIS 6090

Tables d'opérations
conventionnelles motorisées
Conventional motorized operating tables
Mesas convencionales motorizadas

Notice technique
Technical manual
Manual Técnico

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GUARANTEE

FHSurgical warrants the equipment against all defects of materials or manufacture for 12 months after delivery.

*The warranty covers only equipment supplied by **FHSurgical**.*

It is limited to replacement of parts acknowledged to be defective. The labour, travel, living and insurance expenses of the crew replacing parts acknowledged to be defective are incumbent on the customer.

*Title to replaced parts reverts to **FHSurgical** and they must be returned to it within one month of replacement; otherwise **FHSurgical** reserves the right to bill the value of the replaced part.*

Warranty repairs and replacements do not extend the initial warranty period.

This guarantee does not cover the normal wear and tear of parts, nor does it apply to consumables (in particular bulbs, fuses, neon tubes, UV tubes, semi-conductors, valves and diodes, rubber seals, mattress, cushions, etc.).

In order to benefit from the terms of this guarantee, the Purchaser must inform us of any defect by letter or by telex followed by letter of confirmation within three days of discovering the defect.

*The Buyer must submit to **FHSurgical** all evidence that the defect is incumbent on **FHSurgical** and afford **FHSurgical** every facility to investigate the same.*

The customer may buy spare parts simultaneously with the equipment.

***FHSurgical** may use the Buyer's spare parts subject to replacement thereof.*

*On request, **FHSurgical** will supply the principal electrical diagrams and lists of parts to duly accredited technicians.*

*The warranty determines if the supplies are altered or repaired by the Buyer or other parties without **FHSurgical**'s prior written consent.*

No damages are recoverable under this warranty, which is personal to the Buyer and determines forthwith if it disposes of the equipment.

The warranty does not cover:

- incidents caused by fortuitous events or force majeure as described in article X of the general conditions of sales,

- repair or replacement required by normal wear of the equipment,
- damage or accidents due to the customer's act such as negligence, lack of supervision, faulty connection and use/procedures, disregard of the operating or maintenance instructions or improper use due inter alia to overload of any kind.

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SYMBOLS USED

1 - DESCRIPTION / CHARACTERISTICS

- 1-1 MECHANICAL AND ELECTRICAL PARTS
- 1-2 TECHNICAL CHARACTERISTICS OF THE CHARGER
- 1-3 INFORMATION CONCERNING LEAD-ACID BATTERIES
- 1-4 HYDRAULIC PARTS

2 - MAINTENANCE

- 2-1 MAINTENANCE SCHEDULE
- 2-2 SPARE PARTS
- 2-3 REPAIR TOOLS
- 2-4 RECOMMENDED SPARE PARTS

3 - TROUBLESHOOTING PROCEDURES

- 1) REMOVAL OF THE COVERS
- 2) REMOVAL OF THE BELLOW
- 3) PURGING OF THE HYDRAULIC SYSTEM
- 4) FIXING THE TABLE RELATIVE TO THE GROUND
- 5) REMOVAL OF THE HYDRAULIC UNIT
- 6) REPLACEMENT OF PRESSURE SWITCH
- 7) REMOVAL OF HYDRAULIC DISTRIBUTOR ASSEMBLY
- 8) REMOVAL OF LATERAL TILT CYLINDER
- 9) REMOVAL OF TRENDLENBURG / REVERSE TRENDLENBURG CYLINDER
- 10) REMOVAL OF UP/DOWN CYLINDER
- 11) REMOVAL OF BACK REST CYLINDER
- 12) REMOVAL OF FRAME
- 13) COLUMN PLAY ADJUSTMENT
- 14) LOCKING DEVICE
- 15) REPLACEMENT OF THE SEALED LEAD BATTERIES
- 16) SERVICE OF THE CONTROL UNIT
- 17) REMOVAL OF RECEPTION BOARD
- 18) RESET ADJUSTMENTS
- 19) REPLACEMENT OF TRENDLENBURG / REVERSE TRENDLENBURG POTENTIOMETER
- 20) REPLACEMENT OF LATERAL TILT POTENTIOMETER
- 21) REPLACEMENT OF LEFT BACK REST POTENTIOMETER

4 - PARTS LIST AND DRAWINGS

- DRAWING 1.A BASE
- DRAWING 1.B BASE
- DRAWING 2.A DRUM
- DRAWING 2.B DRUM
- DRAWING 3.A SEAT
- DRAWING 3.B SEAT
- DRAWING 4. POTENTIOMETERS
- DRAWING 5. BACKREST
- DRAWING 6. COLUMN
- DRAWING 7. ELECTRICAL S/A
- DRAWING 8. LOCKING DEVICE

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
- 1 - Base
- 2 - Column
- 3 - Back plate
- 4 - Seat plate
- 5 - Access-holder rails
- 6 - Cassette rail
- 7 - Hand control unit
- 8 - Battery charger
- 9 - Protection fuse (12 A)
- 10 - Hand control unit bracket
- 11 - Extension blocking system
- 12 - Controls connection
- 13 - Ground securing shoes
- 14 - Cushion self-attachment strip
- 15 - Connection panel

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1 - DESCRIPTION/ CHARACTERISTICS

Description

The UNIVERSIS 6090 is a mobile, motorized electro-hydraulic operating table, self-powered from batteries.

This conventional table with multi-section universal top comprises (FIG. 1):

- a base (1),
- a column (2) with sealed bellows,
- a back plate (3) and a seat plate (4) on which a variety of extensions can be mounted. The table is equipped with 25 x 10-section accessory-holding rails (5) and cassette rails (6) designed for cassettes with formats of up to 36 x 43 cm,
- a wire-connected control unit (7),
- a battery charger (8),
- a protection fuse (9).

On option, the table can be controlled by a foot pedal for up/down and Trendelenburg/Reverse Trendelenburg movements.

To avoid dangers linked to electrostatic charges, we recommend that you only use FHSurgical antistatic cushions, specifically designed for this product and on which the patient must be in direct contact. The operating table must be set on an antistatic floor or connected to an equalization device (equipotential connector).

Motorized movements (FIG. 2)

The UNIVERSIS 6090 has the following motorized controls:

- up/down (maximum height 1090 mm, minimum height 740 mm),
- Trendelenburg/Reverse Trendelenburg (30°),
- right and left tilts (20°),
- back rest swivel + 80°/- 40°
- reset to zero (tilts, back, low position) nephrectomy/reverse nephrectomy positioning movements,
- table secured to ground using control unit.

NOTE:

- The UNIVERSIS 6090 motorization functions are "smart" functions. In particular, when resetting or controlling nephrectomy positions, the patient's orientation on the table and Trendelenburg angles are taken into account. Movements are controlled progressively with Trendelenburg angle continuously compensated.

- When using high frequency surgical devices, defibrillators or cardiac monitors-defibrillators, manufacturer's instructions must be strictly applied.

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- The table is designed to operate safely for patients **weighting up to 150 kg** for all the operating configurations described in this manual.

Higher-power tables are available for interventions on heavier patients provided the stability of the product has been verified.

For example, in configuration 1 (head extension + short backrest + Tab6090 + leg extension) described below, the table can receive patients of up to 300kg.


The table must be secured in place (onto table shoes) to perform any surgical procedure

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Electrical power supply -Connections

- Charger input: 100/120-230 V* 50 Hz/60 Hz.
- Protection by 16 A fuse.
- Table designed for self-sufficient operation by built-in 24 V battery system.
Battery charger with European power cable, length 5.80 m.

Insulation

The UNIVERSIS 6090 table is equipped with an earthing/grounding connector on the connection panel for prevention of electrical shock.

* The charger supply voltage is not automatically selected. The supply voltage is checked and modified on the selector in the charger. This operation should be performed by a technician.

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Battery recharge :

Special external charger
CHA 6090(Europe)

Observing the technical specifications defined in these manuals implies that the floor in the area of use for the operating table satisfies the minimum characteristics. The stability of the table has been proven to comply with these characteristics.

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1-2 Technical characteristics of the charger

- Charge graph type ----- IUI
- Rated input voltage ----- 110 / 220 V
- Max. input voltage ----- 130 / 260 V (not permanent)
- Min. input voltage ----- 90 / 180 V
- Frequency ----- 50 / 60 Hz
- Charging voltage for sealed lead battery ----- 28V ± 1%
- First output charge current ----- 1.8A average ± 5%
- End of charge current (for 3H) ----- 0.2A average ± 10%
- Charge current if U < 20 V ----- 7A average ± 10%
- Temperature compensation ----- - 2mV / ° C
- Operating temperature ----- - 10° / + 45° C
- Storage temperature ----- - 20° / + 70° C
- Relative humidity ----- 90 %
- Dimensions ----- 220 x170 x155mm
- Weight ----- 3.5 kg
- 3-wire mains input cable with standardized connector
- Battery output connector 5.8 m cable with connector

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The IUI type charge graph is perfectly suited for batteries used in cycles. It can totally release the battery capacity using a program and time-limited end of charge current.

At the end of the charging, the use of a floating voltage allows the charger to be left connected to the battery without risking damage.

If the battery voltage does not rise to its reference value, an electronic safety device cuts off the charger in 12 hours following start up.

During operation, if the battery voltage falls below 20 V, the charger supplies a 7 A current to sustain the normal use of the table connected to it.

The current regains its rated value of 1.8 A automatically once the voltage has risen to a value > 24V.

1-3 Information concerning lead-acid batteries

The lead-acid batteries mounted on tables TAB 6090 are sealed batteries type YUASA NP17-12 (12V/17 Ah).

The batteries operate in all positions except when used permanently in the reverse position. These «NP» batteries are reversal-type batteries.

No maintenance is required.

Depending on the extent of the discharge, more than 1000 cycles (charge/discharge) are possible with this type of battery.

*The NP type batteries have a service life of 5 years when used strictly in accordance with **FHSurgical** specifications.*

At a storage temperature of 20°C, the self-discharge rate is around 3% per month.

The table below indicates the normal storage time versus ambient temperature.

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Recommendations for use

Operating temperature range: + 5 °C to + 35 °C.

Clean the batteries with a moist rag.

Never leave any splashes of oil or organic solvents such as gas-oil or thinner on the batteries.

Never clean the batteries with a rag moistened with any of these products.

Touching the conducting parts can cause electrical shocks. Always wear rubber gloves before performing any inspection or maintenance on the batteries.

Fitting the operating table with a combination of batteries of different capacities, ages or manufacturers can result in damage to the batteries and to the operating table.

Do not interrupt use of an operating table for an extended period of time without first recharging the batteries. This will avoid reducing their service life.

At the end of the battery service life, the batteries should be recovered for recycling.

Recycling of batteries

In accordance with article 9 of decree No.97-1328 dated December 30, 1997, "users of batteries and accumulators other than households are required to collect, or have collected, value or have re-valued, eliminate or have eliminated, all used batteries and accumulators". FHSurgical can take care of collecting the used batteries when our technician comes in to replace these components.

An operating table must be re-charged every evening.
If the table is not used, it is best to leave the charger connected.

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a) Safety points

Check the following :

- Upper and lower jack attachments.
- Shaft securing.
- Accessory rail securing screws.

b) Settings

- Movement speed.
- Reset.
- Column play
- Leg supports*.

c) Checks

Hydraulic parts

- Fluid level.
- No leaks.
- Draining, if necessary.
- Condition of the anti-static pads and of the fixation on the ground.
- Hydraulic pressure check (110 bar).

Mechanical parts

- Locking devices for fixed parts.
- Cleanliness of the jack attachments.
- Cleanliness and lubrication of the column bearings.
- Engagement of the leg and thigh supports*.
- Wear of the leg support pinions*.
- Operation of the leg support safety ratchets*.
- Castor condition.

Electrical parts

- Testing of the various movements.
- Earth conductivity (0.2 Ohms max.).
- Charger operation.
- Battery condition.
- Battery voltage: 11.8V min.
- Voltage difference between 2 batteries: < 0.15 V

d) External appearance

- Cover tightness.
- Condition of the keyboard and of the icons on the control unit.
- Condition of the extension cable.
- Cleaning using a damp sponge and soap.

* Depending on the option.
Refer to the technical instructions of the option.

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NOTE: in this procedure, the terms left and right are determined by the left or right side of the patient lying on his back (over the "T" section of the operating table base) "Normal use" configuration "One".

1) REMOVAL OF THE COVERS

FIG.3

- Place the table in the low position, if possible.
- Remove the 8 screws located on the lower part of the bellows.

Telescopic covers

a) Removing the upper half-cover (A)

- Separate the 2 upper half-covers using a small screwdriver and then disassemble them.

b) Removing the middle covers (B)

- Separate the 2 middle covers using a small screwdriver and then disassemble them.

c) Removing the lower covers (C)

- Remove the 4 screws located on the lower half-covers.
- Disassemble the half-covers.

Base cover

- Lift the table cover from **F1** first and then from **F2**. The cover is secured on the table base with four metallic clips.

2) REMOVAL OF THE BELLOWS

- Remove the 8 screws.

NOTE : the bellows must be removed prior to any work on the following parts:

- tilting, Reverse Trendelenburg/ Trendelenburg jacks,
- shaft potentiometers,
- charger control units and foot pedal plugs.

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* couple de serrage 10 Nm.

3) PURGING OF THE HYDRAULIC SYSTEM

- Remove the covers following the procedure described in procedure 1.
- Set the table in the fixed position relative to the ground.

For all jacks:

- Drain is done automatically.
- For that, do the movement regarding exchanging fixing jack several times.

Checking the oil level

- Immobilize the table then place table horizontal, in its highest position, then on the tank:

Drawing 1.B

- Remove the screw **(55)** printed "oil".
- Visually gauge the oil level (middle rib on plastic tank)
- Fill if necessary
- Replace the screw "Oil"

4) FIXING THE TABLE RELATIVE TO THE GROUND

a) Replacing a pad

NOTE : any operation on the whole fixing jack requires the replacement of 2 bellows attachment collars, seals and hollow bolts of the hydraulic connections. The pad must be secured on the jack rod using weak threadlocking product.

Position operating table on these pads.
Place, under the base, closed from the pads, two 8 cm height wood squared chocks.

b) Replacing the fixing jack (Drawing 1.B)

- Remove the covers following the **procedure**.
- Loosen the hydraulic connection **(30)**.
- Remove the 3 securing screws **(110)** of the jack unit
- Replace the faulty parts.
- Put the parts back into place following the reverse order.

NOTA : The 3 securing screws (110) are tightened to 8 N.m and stick with low threaded glue.

- Lock the connection screw * **(30)** into position.
- Perform an automatic drain cycle by alternately pressing the pushbuttons used to lock and unlock the table on the floor.

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c) Emergency release of the fixing to the ground

- Remove the covers following the procedure N° 1.

Drawing 1.B

When the table is immobile and the ground fixing system is no longer operating, it is possible to release the hydraulic pressure of the jacks manually.

- Remove the plug "Oil" (55) from the tank.
- Insert the flexible hose of the OPM 050* tool in the tank opening
- Unscrew the plug from the connector "pressure plug" (15).
- Screw the connector of the tool onto the connector "pressure plug"
- Check the fitting of the pads
- Disassemble the tool
- Reassemble in reverse order.

5) REMOVAL OF THE HYDRAULIC UNIT

Drawing 1.A

- Remove :
 - the fuse,
 - the covers following the procedure N° 1,
 - the pig iron.
- Cut the 2 white secure collars,
- Disconnect the electrical power supply units of the the hydraulic pump(wires + &-) from the relay (slot1&A1),
- Remove the 2 hydraulic connections,
- Remove the 3 THc pump supportscrews (170),
- Remove the relay holder and the electrovalve black square connector,
- Replace the hydraulic unit.
- Re-assemble in reverse order and **check the oil level as instructed in procedure No. 3 «Checking the oil level».**

NOTE : before fitting a new unit, check that there is a ventilation hole on the "oil" cap.

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6) REPLACEMENT OF PRESSURE SWITCH

Drawing 1.B

Function: Indicates that hydraulic pressure used to lock table on floor has been reached (100 ± 5 bars). This is displayed by a green LED on the control unit.

Disassembly:

- Remove the main fuse (see user manual).
- Raise the cover on the base as instructed in procedure No. 1.
- Slide the cover along the electrical cable.
- Disconnect the cable.
- Remove the pressure switch using a 24mm ring spanner.
- Replace the assembly.
- Re-assemble in reverse order.

7) REMOVAL OF HYDRAULIC DISTRIBUTOR ASSEMBLY

Drawing 2.A

- Remove the covers as instructed in procedure No. 1
- Secure the table in the low position (OPM 036 *).
- Set the back rest in the low position.
- Secure the tilt movements using straps OPM 037 *.
- Disconnect the 8 connectors on the electric distributors.
- Disconnect the 10 hydraulic couplings.
- Remove the 2 black covers and remove the 2 nuts located above the distributor block assembly.
- Remove the 2 screws **(190)** located under the distributor block assembly, then remove the distributor blocks.
- Replace the distributor assembly.
- Re-assemble by replacing all of the seals and hollow screws on the hydraulic couplings.

NOTE: Hydraulic couplings must be torque-tightened to 10 N.m (OPM 040 *).

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8) REMOVAL OF LATERAL TILT CYLINDER

Drawing 2.B

- Remove the covers as instructed in procedure No. 1.
- Lighten the top as much as possible by removing its various components.
- If possible, set the table to the left tilt position and lock the seat in this position using straps OPM 037*.
- Unscrew the upper and lower couplings and fit hydraulic covers on the pipes and couplings.
- Remove nut (170) and remove ball joint from pin (55).
- Remove retaining ring (180).
- Remove pin (30) using a pin drift.
- Recover the two spacers (40).
- Replace the faulty parts.
- Re-assemble in reverse order.

NOTE: Locknut (170) should be replaced each time cylinders are disassembled. Drain and check the oil level in the hydraulic unit.

* Repair tools, see para. 2.3.

9) REMOVAL OF TRENDELBURG / REVERSE TRENDELBURG CYLINDER

Drawing 2.A

- Remove the covers and bellows as instructed in procedure Nos. 1 and 2.
- Lighten the top as much as possible by removing its various components.
- If possible, place the table in Trendelenburg or reverse Trendelenburg position and lock the seat in this position using straps OPM 037*.
- Remove the upper and lower couplings and fit hydraulic covers on the pipes and couplings of the lower part of the cylinder.
- Release the locknut on the cylinder rod (100).
- Remove the retaining ring (290).
- Remove pin (80) using a pin drift.
- Recover the two spacers (90).
- Cut the mounting clamp.
- Uncouple the cylinder from its mounting point by unscrewing the cylinder assembly.
- Replace the cylinder.
- Re-assemble in reverse order.

* Repair tools, see para. 2.3.

NOTE: Perform an automatic drain operation and check the oil level on the hydraulic unit.

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10) REMOVAL OF UP/DOWN CYLINDER

NOTE: Replace hollow screws and seals on hydraulic couplings after each disassembly procedure.

- Remove covers as instructed in procedure No. 1.
- Secure the drum with wedge OPM 036*.
- Remove the fuse.

* Repair tools, see para. 2.3.

FIG. 6

- Loosen the roller holders.
- Secure the guide plates (5) with 2 screws Chc M5-16 (4).
- Loosen the nuts (1) and the adjustment screws (2) cylinder side, then re-block the two screws Chc M5-16.

Drawing 2.B

NOTE: The screws are only used for disassembly of the column; the screws are removed after reassembly and adjustment of the lateral plate.

- Raise the cover off the socket (see procedure No. 1), then turn it by 90°.
- Disconnect the orange connector in the base.
- Remove the gray connector.
- Disconnect the line A hose from the pump, then the double-coupling.
- Disconnect the wires.
- Remove the 2 screws securing the reeving (screws located in base).
- Disconnect the 2 hoses supplying the up/down cylinder on the distributor.
- Remove the plate/frame assembly.
- Remove the frame cover.
- On the distributor, disassemble the pipes (open the clamps required for disassembly) and replace the hydraulic couplings by 2 yellow covers.
- Remove screw Chc M12 (205).
- Lift the frame/drum assembly.
- Remove the wedge.
- Remove the 4 nuts from the column (160).
- Remove the column (20).
- Loosen the locknut (207) on the cylinder rod.
- Unscrew and remove the cylinder.
- Cut the clamp.
- Fit the new cylinder with hoses.
- Fit two new single couplings with seals and torque-tighten to 10 N.m.
- Fit the new cylinder.
- Reinstall the column.
- Install the nuts and torque-tighten to 35 N.m.
- Re-position shim OPM 036.
- Reinstall the frame/drum assembly.
- Re-position screw Chc M12 without blocking it.
- Perform the lateral adjustment (see procedure No. 13).
- Perform the longitudinal adjustment (see procedure No. 13).
- Torque-tighten screw Chc M12 to 27 N.m.

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- Re-connect the hydraulic couplings.

NOTE: The hydraulic couplings on pipes By3 and Ay4 are torque-tightened to 10 N.m. (see OPM 040, repair tools).

- Secure the gray connector.
- Connect the orange connector.
- Install the fuse.
- Set the table to the high position.
- Install the reeving guide.
- Remove wedge OPM 036.
- Perform 5 up/down movements to drain the hydraulic pipes.
- Check the oil level (see procedure No. 3).

Drawing 3.A

- Install the cover on the frame (50).
- Install the top-frame assembly.
- Install the covers.

11) REMOVAL OF BACK REST CYLINDER

Drawing 3.A

- Set the back rest to the up position at + 80°.
- Remove the top-frame assembly.
- Remove the covers (40) (drawing 3.B).
- Remove the 2 screws (310 and 330) then the accessory bar (10).
- Disconnect the cylinder hydraulic supply.
- Fit the hydraulic cover supplied with the new cylinder.
- Remove the stainless steel flange.
- Loosen the screws (70 and 90) (drawing 3.B).
- Remove the pin (60).
- Disengage the cylinder (120) and push it toward the back rest to take it out of its guide.
- Remove the rollers (110) (drawing 3.B).
- Return the cylinder to its initial position and take it out through the top (leg side).
- Remove the back rest link rods (20) (drawing 3.B).
- Replace the cylinder.
- Re-assemble in reverse order.

NOTE: Replace hollow screws and seals on hydraulic couplings after each disassembly procedure. Perform automatic drain operation and check oil level on hydraulic unit

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12) REMOVAL OF FRAME

- Remove the fuse.
- Remove the covers and bellows as instructed in procedure Nos.1 and 2.
- Disconnect the battery supply connector (located in base).

Drawing 3.A

- Remove the top-frame assembly.
- Remove the frame cover (50).
- Remove the top supports (40).
- Disconnect the hydraulic couplings from the back rest cylinders (install covers).
- Disconnect all the electrical connections connecting the frame to the drum.
- Remove the lateral tilt potentiometer (see procedure No. 20).

Drawing 2.A

- Remove the fastener (10) by removing the 2 screws (160).
- Remove the screw (165) by loosening screw (166).
- Remove the sleeve (167).
- Remove the retaining rings (209) and remove the pin from the Trendelenburg/reverse Trendelenburg cylinder.
- Tilt over the frame to the Trendelenburg position to disengage the cylinder from the tilt shaft (330).
- Remove the shaft (330).

Drawing 2.B

- Remove the nut (170) and take the lateral cylinder off pin (55).
- Lift and replace the frame.
- Re-assemble in reverse order.

Removal of headset

- Remove the frame as instructed in the procedure above.
- Remove the Trendelenburg/reverse Trendelenburg potentiometer (see procedure No. 19) and pin (240).
- Remove the screws and counter-screws (150) and (155).
- Remove the pins (240) located on either side using, if necessary, extractor OPM 017.
- Replace the headset.
- Re-assemble in reverse order.

NOTE: Drain the unit and check the oil level in the hydraulic unit. Re-calibrate each time potentiometer is replaced or disassembled.

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13) COLUMNPLAY ADJUSTMENT

- Remove the covers following the procedure 1.

Longitudinal play (Fig. 4 & 5) page 28

- Unfasten the nuts (2).
- Loosen the screws (3).
- Re-tighten the screws (1) progressively until a torque of 4 N.m is obtained (tighten the screws alternately).
- Fasten the screws (3).
- Fasten the lock-nuts (2).

NOTE : these adjustments are to be carried out on the upper and lower roller holders.

Side play (Fig. 6)

(The adjustment screws are located on either side of the column. **The adjustment is carried out on the right hand side (jack), the left hand side being used for reference.**)

- Check that the column is centred on the ramp.
- Check that the column is vertical using a spirit level.
- Block the nuts (1).
- Bring the screws (2) into contact.
- Fasten the screws (2) by 1/4 turns then refasten the lock-nuts (1).

14) LOCKING DEVICE

NOTE : before carrying out any disassembly operation make sure that you possess a locking device repair set as certain parts cannot be recovered after disassembly.

Drawing 8

Disassembly of handle (7)

- Fully loosen the clamp.
- Remove the screw (6).
- Separate the handle from pin (3).

Removal of the locking device

- Remove handle.

For the frame (Drawings 3A and 3B)

- Remove the cover (40) then the fixing screws (190, 210, 220 et 230).

For the rails

Drawings 5 and 8

- Screw in the (3) fully until it just out past the level of the rails.
- Remove the 8 fixing screws.
- Remove the locking device assembly.
- Replace the faulty parts.
- Put the parts back into place in the reverse order.

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15) REPLACEMENT OF THE SEALED LEAD BATTERIES

- Remove the fuse.
- Using OPM 037, raise the assembly formed by the stainless steel cover (**30, drawing 1.A**), and the sliding cover (**50, Drawing 6**).

Drawing 1.B

- Remove the 2 screws (**100**) and the cover (**60**).
- Replace the faulty sealed lead batteries.

16) SERVICE OF THE CONTROL UNIT

Replacement of the key board

For tables with serial No. between AR 0003 and AR 0099:

- Open the unit by removing the 6 screws located on the lower part of the unit.
- Remove the 2 screws located under the cover of the battery housing.
- Separate the two parts of the unit.
- Disconnect the power supply cable.
- Remove the 8 nuts holding the keyboard.
- Replace the faulty part.
- Put the unit back together in the reverse order.

Replacement of the extendable power cable

- Open the unit following the procedure described above.
- Disconnect the pin of the extendable cable.
- Replace the extendable cable.
- Put the unit back together in the reverse order.

Replacement of lower remote control board

- Open the housing as instructed in the procedure above.
- Disconnect the pin on the flexible cable.
- Remove the 6 screws securing the board.
- Replace the board.
- Configure the board switches if necessary (see diagram on page opposite).
- Re-assemble in reverse order.

For tables with serial No. starting with AR 0100

Replacing the extendable cable

- Remove the 6 screws.
- Remove the keypad.
- Disconnect the extendable cable.
- Reassemble in reverse order.

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17) REMOVAL OF RECEPTION BOARD

- Remove the upper half-cover, leg side.
- Remove the fuse.

Caution! 24 V present at terminal strip. This voltage can be cut off by disconnecting the orange connector located at the bottom of the reeving assembly.

Drawing 2.A

- Uncouple connectors (P6.1; P6.2; P6.3; P7.1; P7.2) located on the sides of the electronic unit.
- Remove the 3 screws (170).
- Replace the board.
- Fit the 3 screws (170) with washers (250).
- If necessary, configure the board switches (see diagram on page 32).
- Connect connectors (P6.1; P6.2; P6.3; P7.1; P7.2) and check that the cables have not been damaged when handled.
- Adjust the RESET on the new board (see procedure No. 18).
- Install the half-cover.

Caution! No replacement of the board under warranty will be accepted if the WARRANTY sticker has been removed or deteriorated.

18) RESET ADJUSTMENTS

Preliminary steps

NOTE: The left and right designations are defined with respect to the left and right of the patient lying on his back (patient back on table back rest). Table in "Normal utilization" configuration on control unit. To obtain a correct reset, the adjustments must be performed with all the potentiometers connected to the board.

- Remove the covers as instructed in procedure No. 1.
- Using the control unit, switch on the table then perform a "Lock on floor" operation. Once the table is bearing on its shoes, the bearing surface is considered as the horizontal plane reference.
- Check the reference using a spirit level.
- With the table in the maximum up position and "secured to floor", use a spirit level and the control unit to position the seat horizontally with respect to the floor then switch OFF the table.

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Trendelenburg/reverse Trendelenburg adjustment (FIG. 7)

- Set the table ON/OFF switch to ON and switch off LED 6 corresponding to the «Trendelenburg/reverse Trendelenburg» movement on the board by turning potentiometer 2 using a small screwdriver. The new RESET position is now recorded.

Adjustment of lateral tilts (FIG. 7)

- Set the table ON/OFF switch to ON and switch off LED 5 corresponding to the «lateral tilt» movement on the board by turning potentiometer 1 using a small screwdriver. The new RESET position is now recorded.

Adjustment of back rest left arm (FIG. 7)

- Set the table ON/OFF switch to ON and switch off diode LED 7 corresponding to the «back rest left arm movement» by turning potentiometer 3 using a small screwdriver. The new RESET position is now recorded.

19) REPLACEMENT OF TRENDLENBURG/REVERSE TRENDLENBURG POTENTIOMETER

NOTE: Before starting the intervention, be sure you have a Trendelenburg / reverse Trendelenburg potentiometer repair kit in your possession (see para. 2.2).

Removal of potentiometer

Drawing 3.A

- Using a spirit level, set the «Trendelenburg / reverse Trendelenburg movement» to the level position.
- Remove the frame cover.

FIG. 8

- Using a hex head wrench, remove screw (4) then unlock screw (3) to separate the potentiometer shaft (8) from the drum shaft.
- Remove the 3 screws (5) and the potentiometer.
- Disconnect the 3 plugs from the terminal strip on the board (30, 31, 32 – connector P7.1).
- Open the clamps required to remove the potentiometer assembly.

Assembly and adjustment of potentiometer

Figs. 8 and 11

- Align the two marks on the potentiometer shaft and body.

FIG. 8

- Carefully insert the potentiometer.

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Secure the potentiometer base with 3 screws (5).

- *Block the potentiometer shaft (8) with screws (3) and (4).*
- *Connect the potentiometer to the terminal strip on the board (connector P7.1):*
Red or brown wire → terminal 32
Braid or green wire → terminal 30
Black or white wire → terminal 31
- *Disconnect connector P7.1. from the terminal strip.*
- *Using an ohmmeter, check that the potentiometer mid-point is properly adjusted (see FIG. 11).*
- *Fit the clamps required to secure the assembly.*

If mid-point is misadjusted:

- *Remove the counter-screw (9) and loosen the screw (10).*
- *Turn the shaft/potentiometer assembly (11) by screws (5) to obtain the mid-point.*
- *Block screw (10) and counter-screw (9).*
- *Reconnect P7.1 to the terminal strips.*
- *Fit the clamps required to secure the assembly.*
- *Perform the RESET adjustment.*

20) REPLACEMENT OF LATERAL TILT POTENTIOMETER

NOTE: Before starting the intervention, be sure you have a lateral tilt potentiometer repair kit in your possession (see para. 2.2).

Disassembly of potentiometer FIG. 9

- *Using a spirit level, set the "lateral tilt" movement to the level position*
- *Remove the protective cover (10) by removing the two screws (9).*
- *Separate the potentiometer shaft (8) with respect to the seat by removing screw (4) and unlocking screw (3).*
- *Remove the potentiometer by removing the 3 screws (5).*
- *Disconnect the 3 plugs from the terminal strip on the board (33, 34, 35 Connector P7.1).*
- *Open the clamps required to remove the assembly.*

Assembly and adjustment of potentiometer FIG. 9 and FIG. 11.

- *Align the two marks on the potentiometer shaft and body.*

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FIG. 9

- *Carefully insert the potentiometer.*
- *Secure the base of the potentiometer using 3 screws (5).*
- *Block the potentiometer shaft (5) using screws (3) and (4).*
- *Connect the potentiometer to the terminal strip of the board (connector P7.1):*
Red or brown wire → terminal 35
Braid or green wire → terminal 33
Black or white wire → terminal 34
- *Disconnect the connector P7.1 from the terminal strip.*
- *Using an ohmmeter, check that the potentiometer mid-point is properly adjusted (see FIG. 11).*
- *Fit the clamps required to secure the assembly.*

If mid-point is misadjusted:

Drawing 2.A and FIG. 9

- *Remove the fastener (10) by removing the two screws (160).*
- *Remove the counter-screw (11) and loosen screw (12).*
- *Turn the shaft/potentiometer assembly (13) by screws (5) to obtain the mid-point.*
- *Block the screw (12) and counter-screw (11).*
- *Reconnect P7.1 to the terminal strips.*
- *Fit clamps required to secure the assembly.*
- *Perform the RESET adjustment.*

**21) REPLACEMENT OF LEFT
BACK REST POTENTIOMETER**

NOTE: Before starting the intervention, be sure you have the left back rest potentiometer repair kit in your possession (see para 2.2).

**Removal of potentiometer
Drawing 3.B**

- *Using a spirit level, set all the table movements level (Trendelenburg / reverse Trendelenburg, lateral tilt and back rest).*
- *Switch off the table.*
- *Remove the top from the frame.*
- *Remove the cover (40).*

FIG. 10

- *Remove the screw (12) immobilizing pin (11) on the arm.*
- *Remove the 3 screws (5) and remove the flange (6).*
- *Remove the 3 screws (3).*
- *Remove the flange (2) by removing 3 screws (1).*
- *Remove the pin (11).*
- *Remove the screw (10) from the pin then loosen screw (9).*
- *Remove the potentiometer from the shaft (11).*
- *Disconnect the 3 end-pieces No. 27, 28 and 29 of the potentiometer from the left arm of connector P7.1.*
- *Open the clamps required to remove the assembly.*
- *Replace the potentiometer.*

FIG. 9

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FIG. 10

- Carefully insert the potentiometer in the shaft **(11)**.
- Tighten screw **(9)** then reinstall and lock screw **(10)**.
- Insert shaft **(11)** in housing.
- Secure back rest pin of arm by installing and tightening screw **(12)**.
- Block 3 screws **(3)**.
- Connect potentiometer to terminal strip on board (connector P7.1):
Red or brown wire → terminal 27
Braid or green wire → terminal 28
Black or white wire → terminal 29
- Disconnect connector P7.1 from terminal strip.
- Using ohmmeter, check that potentiometer mid-point is properly adjusted (see FIG. 11).
- Fit clamps required to secure the assembly.

If mid-point is misadjusted:

- Loosen 3 screws **(3)**.
- Turn potentiometer to obtain mid-point.
- Block 3 screws **(3)**.
- Perform RESET adjustment.

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