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- ◆ Never place on the pan an object heavier than the maximum capacity of the balance, the sensor could be damaged.
- ◆ Neither submerges the balance nor spill liquids on it.
- ◆ If you use batteries, remove them when you are not using the balance for a long period of time.
- ◆ If any liquid comes into contact with the electric parts of the balance please do immediately disconnect it from the net and send it to the technical service, as soon as possible, for its checking and adjustment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the balance and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

NOTE: According to legislation currently in force, regarding "Non-automatic weighing instruments" in which balances are included, by means of writ dating from 22nd December 1994 (BOE 3/1/95), Nahita balances series 5032 must not be used for:

- *Commercial transactions.*
- *Calculation of taxes, tariffs, duties, remunerations, indemnities and other similar canons.*
- *Judicial surveys.*
- *Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories.*
- *Determination of the price or total amount in retail sales and in the preparation of pre-packaged products.*

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita digital balance, series 5033. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita digital balance series 5032, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the balance avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the balance.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

8. TROUBLESHOOTING

TROUBLE	CAUSE	SOLUTION
Digits are not displayed	Check the balance receives electrical current. The adapter can be damaged. The fuse is blown.	Tighten the power plugs. Change the current adapter. Replace the fuse.
The value on display does not stabilizes	The balance is in a place with vibrations. The door of the balance is not completely closed. There is an odd object under the pan causing destabilization.	Move the balance to a place without vibrations. Close all the doors of the balance.. Check that the balance is clean inside and that there are not strange objects.
The displayed value is wrong	Balance is not calibrated. Balance has not been tared before weighing. Balance is not well leveled.	Calibrate the balance. Tare the balance before each weighing. Level the balance until the bubble situated at the back side is centred.

9. MAINTENANCE

To get the best results and a higher duration of the pH-meter it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Prevent the balance from sudden movements and knocks, as well as from direct sunlight or air flows. The balance is a precision instrument, you must handle it carefully.
- ◆ The balance is supplied with an adapter. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never unplug the adapter by pulling the wire, do it from the base.
- ◆ Place the balance in a location with adequate ventilation to prevent heat accumulation. Never use the balance in a wedged location as for example a shelf.
- ◆ Never use sharp objects as pens, etc... to press the buttons of the control panel; only use your fingers.

7.3.3. AUTO PRINT C

Data are transmitted continuously when plugging the balance without giving any additional order to neither the balance nor the computer or printer.

PRINT MODE

Data will only be transmitted if PRINT key (1.6) is pressed.

COMMAND MODE

Data will be transmitted from the balance to an external device when they receive an order. Balance can be controlled from the external device by the following orders in the same way as by the frontal panel of the balance.

ORDERS (Commands)

<O><CR><LF> ON/OFF order (4F 0D 0A); the same as ON/OFF function from the frontal panel.

<P><CR><LF> Printing order (50 0D 0A); the same as PRINT function from the frontal panel.



<T><CR><LF> Tare order (54 0D 0A); the same as TARE function from the frontal panel.

<M><CR><LF> Mode order (4D 0D 0A); the same as MODE function from the frontal panel.


<C><CR><LF> Calibration order (43 0D 0A) the same as CAL function from the frontal panel.

Connection with external services:

1. With computer

2. With print

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish.....	2-17
English	18-32
French	33-48

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	19
2. DESCRIPTION.....	20
3. TECHNICAL SPECIFICATIONS.....	21
4. INSTALLATION / SETTING UP.....	21
5. CALIBRATION.....	23
6. FUNCTIONS OF THE BALANCE	25
7. SOFTWARE.....	28
8. TROUBLESHOOTING	31
9. MAINTENANCE AND CLEANING	31

1. USES OF THE INSTRUMENT


Nahita precision electronic balances are the result of many years of research, design and development in the field of analytical balances. They are basic instruments in every laboratory in weighing processes requiring the highest accuracy. Nahita balances, developed with the most advanced technology, are reliable and durable instruments with a solid, clean and ergonomic design suitable for educational, industrial or researching usages.

2. DESCRIPTION

Figure 1: front panel

- 1.1 Digital display
- 1.2 ON/OFF
- 1.3 CAL key: calibration
- 1.4 MODE key: select the weighing mode
- 1.5 TARE key
- 1.6 PRINT key
- 1.7 STAB: Stability indicator
- 1.8 %: percentage function indicator
- 1.9 COUNT: piece count function indicator
- 1.10 CT (carats) unit indicator
- 1.11 Oz (ounces) unit indicator
- 1.12 Kg (Kilograms) unit indicator
- 1.13 g (grams) unit indicator

code 55032300, 55032320

code 55032250



Figure 2

- 2.2 Level
- 2.1 Socket
- 2.3 RS 232 output

DATA RECORDING

- ◆ All the above features (*table 3*) can be modifying directly from the screen.
- ◆ Press TARE and ON/OFF keys simultaneously with the balance turned off.
- ◆ Release the key and 888888888 will be displayed.
- ◆ Press “PRINT” key (-).
- ◆ “C1-0” (CX-Y) will be displayed.
- ◆ Press “MODE” key to change the X digit between values from 1 to 7 (*see table 3*).
- ◆ By pressing “PRINT” key change the Y value. This digit can take different values depending on X digit value (*see table 3*).
- ◆ Press ON/OFF key to save the new parameters in the memory.
- ◆ The balance is ready to work with the new parameters, that will be recorded when you turn on/off the balance.

7.2 RS-232C OUTPUT. SPECIFICATIONS

- ◆ Transmission way: double, asynchronous and bi-directional
- ◆ Data format: digital transmission speed: 1200, 2400, 4800 and 9600.
- ◆ Informative bit: 8
- ◆ Bit of parity: none
- ◆ Bit of stop: 1
- ◆ Code: ASCII
- ◆ Data output: Auto-print, “PRINT” key, data emission
- ◆ Data output format: ST, +001.000xxg <CR><LF><LF>
US, +001.000xxg <CR><LF><LF>

Note:

- ST = Stable (53H, 54H) (Stable)
- US = Unstable (55H, 53H) (Unstable)
- X = Space (20H) (Space)
- <CR> = Carriage return (0DH)
- <LF> = Line Feed <0AH>

7.3 DATA OUTPUT

AUTO-PRINT MODE

There are three different ways in Auto Print mode:

7.3.1. AUTO PRINT A

Note: the pan must be clean and screen must show Zero after each weighing in AUTO PRINT A mode.

Data transmission starts when balance begins to stabilize (stability indicator appears (1.7)).

7.3.2. UTO PRINT B

Data transmission happens whenever balance stabilizes.

7. SOFTWARE**7.1 PARAMETERS**

Nahita balances series 5032 are supplied with a Software to adapt their features, as best as possible, to users' necessities:

TABLE 3

FUNCTION	PARAMETERS (CX-Y)	NAME	MEANING
C1	0~3	Ajustment type Refresh screen	
C2	0~3	Adjustment type Data exit RS-232	
C3	C3-0* C3-1	Zero and tare	ON OFF
C4	C4-0* C4-1 C4-2 C4-3	Digital transmission speed	2400 1200 4800 9600
C5	C5-0* C5-1 C5-2 C5-3	Data exit	Zero set Balance set Print Continuous exit
C6	C6-0* C6-1	Back illumination	OFF ON
C7	C7-0* C7-1 C7-2 C7-3	Stability band's amplitude allowed	1 2 3 4
C8	C8-0 C8-1 C8-2 C8-3*	Average time	1 second 2 seconds 3 seconds 4 seconds

*Factory settings.

3. TECHNICAL SPECIFICATIONS

<i>Code</i>	55032250	55032300	55032320
<i>Capacity</i>	500 g	1000 g	2000 g
<i>Readability</i>	0.01 g		
<i>Reproducibility</i>	± 0.01 g		
<i>Linearity</i>	± 0.02 g		
<i>Pan diameter</i>	150 mm		
<i>External dimensions</i>	195x310x85mm		
<i>Weight</i>	5 kg		
<i>Alimentation</i>	220 V 50/60 Hz		

Among its most outstanding characteristics are the following ones:

- ◆ Stable and rapid weighing
- ◆ Easy working, simple and light line
- ◆ Subtractive tare
- ◆ Network connection by means of adapter (included)
- ◆ Automatic calibration
- ◆ Weighing (Ct, Oz, g and Kg), and pieces count and percentage functions
- ◆ Below weighing function (animal weighing, hydrostatic weighing, etc)
- ◆ Bubble level and screwed legs
- ◆ RS 232 interface for connecting a computer or a printer

4. INSTALLATION*Preliminary inspection*

- ◆ Unwrap the balance, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or administrator immediately, so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied

Please check that all the accessories are enclosed with the equipment:

- Pan
- Cleaning cloth
- Electric current adapter
- Instructions manual
- Installation

Before using the balance, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

- ◆ Please put the balance on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the balance near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Placing the AC adapter
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz \pm 10% provided with earth wire and to the other end to the balance connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz \pm 10%.

- ◆ Place the stainless steel weighing pan.
- ◆ Check that the level is in the right position (centred), if not, centre the bubble with the knobs situated at the frontal legs of the balance.
- ◆ The balance is ready to use.

Setting up

NOTE: Before using the balance, it must be calibrated according to "Calibration" point of this manual. Recalibration must be performed whenever observing important variations or mistakes in weighings.

- ◆ Make sure that there is nothing on the pan and switch the balance on (1.2). The display will light up (1.1) and the balance will make and autotest that finishes when 0.00g is visualized on display.
- ◆ Select the desired function: weighing, percentage or pieces count.
- ◆ By pressing MODE key (1.4) select the measuring unit of weight: grams (g), kilograms (Kg), ounces (oz) or carats (ct). The following table shows the conversion of these units into grams.

displayed.

Note: To set a new reference, press CAL key and proceed as previously described (from point 4).

6.3 PERCENTAGE MODE

- ◆ In percentage mode the balance perform measurements based on a reference value.
- ◆ By pressing "MODE" key, select percentage mode (%).
- ◆ "0.00" will be displayed.
- ◆ Press "CAL" key.
- ◆ By default, "--100--" will be displayed; this value indicates the weight percentage corresponding to the reference sample.


NOTE: To modify this value, press MODE key (to increase) or PRINT key (to decrease). You can choose between 50, 100, 150, 200, 250, 300, 350 and 400 %.

- ◆ Place on the pan the sample whose weight is going to be considered as the 100 % (or the selected value).
- ◆ Press "TARE" key, and after a few seconds, the corresponding percentage will be displayed.
- ◆ Removed the sample from the pan.
- ◆ "0.00 %" will be displayed.
- ◆ Place the sample to be compared on the pan. As a result of the comparison between the reference sample and the problem sample, the deviation in percentage will be displayed.

Note: To set a new reference, press CAL and proceed as previously described (from point 4).

weightings; this technique is widely used for weighing magnetic material and its use can be extended to industry.

- ◆ For below weighing, proceed as follows:
- ◆ Remove the 2 screws that hold the rhombus-shaped lid placed at the bottom of the balance (see fig.4)
- ◆ Remove the lid to see the hook to hold weighing device (not included with the balance)
- ◆ Put the thread of the weighing device through the hook
- ◆ Place the balance on a table with a hole to let the thread hanging through freely
- ◆ Hang the weighing device at the other end of the thread (see fig.5)
- ◆ Proceed in the same way as for normal weighing

6.2 PIECE COUNTING MODE

In counting mode the balance calculates the number of pieces placed on the pan. Proceed as follows:

- ◆ By pressing "MODE" key, select Counting mode (1.9).
- ◆ "0" will be displayed.
- ◆ Press "CAL" key (1.3).
- ◆ By default, "10" will be displayed; this number indicates the number of pieces that should be placed on the pan and that will be taken as a reference for a later count.

NOTE: To modify this value, press MODE key (to increase) or PRINT key (to decrease). You can choose between 5, 10, 15, 20, 25, 30, 35 and 40 pieces.

- ◆ Place 10 pieces of those to be counted, or the selected number of pieces, on the pan.
- ◆ Press TARE key and after a few seconds, the corresponding number of pieces will be displayed.
- ◆ Remove the pieces from the pan.
- ◆ "0" will be displayed.
- ◆ Place on the pan the amount of pieces to be counted and the number of pieces will be

UNITS CONVERSION

Abreviation	Name	Conversion
Kg	Kilogram	1000 g
g	Gram	1 g
mg	Miligram	0.001 g
oz	Ounce	28.3495231 g
ct	Carat	0.1999694 g

- ◆ The balance is ready to weigh.

NOTE: To avoid possible oscillations in weighing, it is important to avoid air flow and vibrations on the table in which it is situated the balance while weighing.

- ◆ Switch the balance off by pressing On/Off key once operation has finished.

5. CALIBRATION

Whenever observing important variations in weighing or not correct weightings, balance should be calibrated. Some factors such as variation of gravity according to latitude, change of working place as well as an inappropriate use of the equipment lead to the necessity of recalibrate the balance.

Important!! Do not move the balance during calibration

Consult the weights you need to calibrate your balance based on its model (not included with the equipment).

CALIBRATION WEIGHTS:

Code Model	Max. capacity	Calibration weight
55032250	500 g	500 g
55032300	1000 g	1000 g
55032320	2000 g	2000 g

Calibration weights should be previously calibrated.

Choose the calibration weight you need according to the model of the balance to be calibrated; you need a weight with the same value as maximum capacity of the balance (see table 2). Proceed as following:

- ◆ Turn the balance on by pressing On/Off key (1.2); the pan should be empty.
- ◆ The balance will perform an autotest and finally stabilizes on 0.00g

Note: If 0.00g does not appear on display, press TARE key (1.5) until it appears and the stabilization indicator (1.7) lights up.

- ◆ Press "CAL" key (1.3) to enter calibration mode.
- ◆ "CAL 0" will be displayed.
- ◆ Press "TARE" key (1.5).
- ◆ "CAL F" will be displayed.
- ◆ Place the calibration weight on the pan.
- ◆ Press "TARE" key (1.5) again.
- ◆ "CAL END" will be displayed.
- ◆ The maximum capacity of the balance will appear on display, it is also the weight of the standard calibration weight used.
- ◆ Removed the calibration weight from the pan.
- ◆ 0.00g will be displayed, indicating that balance calibration has been well performed.

It is necessary to check the results of calibration:

Place the calibration weight on the pan when 0.00g is displayed, and check that the displayed value coincides with the weight value with a maximum deviation of $\pm 0.01\text{g}$; otherwise, balance is not correctly calibrated, and calibration should be performed again.

Note: If "CAL no" is displayed during calibration, it means that calibration weights are not the correct ones because their weight is up $\pm 10\%$ on the maximum capacity of the balance.

Change the weights for ones below these limits and do the calibration again. If "CAL no" message remains in spite of using the appropriate weights, please ask Nahita authorised technical service.

5.1 BALANCE ADJUSTMENT

If you need more accuracy and precision in measures and you also know the uncertainty of calibration weights, you can adjust manually the weighing values of the balance. Taking as an example that the calibration weight value is 200.5 g, calibration should be performed as follows:

- ◆ Turn the balance on by pressing ON/OFF key (1.2) with the pan empty.
- ◆ Press "CAL" key (1.3) when 0.00g is displayed.
- ◆ "CAL 0" will be displayed.
- ◆ Press "CAL" key (1.3) again.
- ◆ The maximum capacity of the balance will be displayed (200.0 g in the example).
- ◆ Adjust this value until reaching the exact value of the calibration weight to be used, 200.5 g in the example. For this purpose, press "MODE" or "PRINT" key either to increase or decrease the value initially displayed.

- ◆ Once the calibration weight value has been adjusted, press "TARE" (1.5) key.
- ◆ "CAL 0" will be displayed.
- ◆ Press "TARE" key (1.5) again.
- ◆ "CAL F" will be displayed.
- ◆ Place the calibration weight on the pan.
- ◆ Press "TARE" key (1.5).
- ◆ "CAL END" will be displayed.
- ◆ The adjusted value (200.5 g) will be displayed.
- ◆ Remove the calibration weight from the pan and 0.00g will be displayed; the balance is correctly calibrated and ready to use.

6. FUNCTIONS OF THE BALANCE

6.1 WEIGHING

SIMPLE WEIGHING

- ◆ Take into account the temperature of the balance; wait at least one hour for a total stabilization of the temperature if balance has been moved.
- ◆ Calibrate the balance if necessary, according to "Calibration" point.
- ◆ Turn the balance on by pressing ON/OFF key (1.2) with the pan empty.
- ◆ When 0.00g and the stabilization indicator (1.7) are displayed, the balance is ready to weight.
- ◆ When tare is necessary because a recipient is going to be used in weighing, place the recipient empty on the pan and press "TARE" key (1.5). Wait until 0.00g is stabilized; the balance is ready to weight.

Note:

1. *The value of subtractive tare is kept in balance's memory until TARE key (1.5) is pressed once again.*
2. *If you need to weight more than one substance in the same recipient and you require having their separated weights, add the first substance and, after its weight has been displayed, press "TARE" key (1.5). 0.00g will be displayed, so the next substance can be weight independently.*

If you need to weight a substance contained in a recipient, place the full recipient on the pan. When its weight stabilizes, press "TARE" key (1.5) and 0.00g will be displayed. With the aid of a spatula or similar, transfer the substance to be weight to another recipient. When transferring finishes the substance's weight (with negative sign) will be displayed.

BELOW WEIGHING

Nahita precision balance series 5032 is provided with a built-in hook to perform below