## Breezy 300

Gebruiksaanweisung

Directions for use

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Manuale per l'uso

Manual do utilizador

Manual de instrucciones

 Breezy®

**INNOVATE.™**  
**BUILD FOR THE FUTURE**

# BREEZY 300

BREEZY 305

BREEZY 312

BREEZY 315

BREEZY 310

BREEZY 321

BREEZY 341

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6

Fig. 7

Fig. 8

Fig. 9

|                 |    |
|-----------------|----|
| DEUTSCH .....   | 5  |
| ENGLISH .....   | 14 |
| FRANÇAIS .....  | 23 |
| DUTCH .....     | 32 |
| ITALIANO.....   | 41 |
| PORTUGUÊS ..... | 50 |
| ESPAÑOL.....    | 59 |

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# INTRODUCTION

Thank you for choosing a SUNRISE MEDICAL product. In view of the broad range of components and accessories available for all SUNRISE MEDICAL wheelchairs it is in your interest to consult a health professional or authorised SUNRISE MEDICAL distributor. These wheelchairs are designed for persons with a physical handicap, to be used only for transport. These wheelchairs are suitable for use both indoors and outdoors (Class B product).

The manufacturer reserves the right to modify, without prior notice, weights, dimensions, or any other technical features described in this manual. All figures, measures, and capacities included in this manual are approximate and are not specifications

Models: 305, 312, 315, 310, 321, 341

Sunrise Medical LTD.  
Sunrise Business Park  
High Street, Wollaston  
West Midlands DY8 4PS  
England  
Phone: 1384 44 66 88  
Fax: 1384 44 66 99

**IMPORTANT:**  
DO NOT USE YOUR WHEELCHAIR  
UNTIL THIS MANUAL HAS BEEN  
READ AND UNDERSTOOD.

Sunrise Medical is ISO 9001 certified, which ensures quality at all stages of the development, and production of this wheelchair.

As the manufacturer, SUNRISE MEDICAL, declares that  
the lightweight wheelchairs conform to the 93/42/  
EEC / 2007/47/EEC guideline.

## INDEX

|                        |    |
|------------------------|----|
| A.CHARACTERISTICS..... | 17 |
| B.SAFETY .....         | 18 |
| C. HANDLING.....       | 25 |
| D. ASSEMBLY .....      | 25 |
| E. ADJUSTMENTS .....   | 27 |

|                                  |    |
|----------------------------------|----|
| F.CONTROLS AND MAINTENANCE ..... | 27 |
| G.SPECIFICATIONS.....            | 28 |
| H. MEASURES .....                | 28 |
| I. GUARANTEE .....               | 88 |

## A. CHARACTERISTICS

1. Handgrip
2. Back upholstery
3. Seat upholstery
4. Armrest upholstery
5. Push rim
6. Propulsion wheel
7. Wheel lock
8. Footrest
9. Footplate
10. Fork
11. Castor

## Use

Lightweight wheelchairs are exclusively for a user who is unable to walk or has limited mobility, for their own personal use in- and outdoor.

The maximum weight limit (includes both the user and any weight of accessories fitted to the wheelchair) is marked on the serial number label, which is affixed to the crossbar or stabiliser bar below the seat.

Warranty can only be taken on if the product is used under the specified conditions and for the intended purposes.

The intended lifetime of the wheelchair is 5 years. Please DO NOT use or fit any 3rd party components to the wheelchair unless they are officially approved by Sunrise Medical.

## Area of application

The variety of fitting variants as well as the modular design mean that it can be used by those who cannot walk or have limited mobility e.g. because of:

- Paralysis
- Loss of extremity (leg amputation)
- Extremity defect/deformity
- Joint contractures/joint injuries
- Illnesses such as heart and circulation deficiencies, disturbance of equilibrium or cachexia as well as for elderly people who still have strength in the upper body.

When considering provision, please also note the body size, weight, physical and psychological constitution, the age of the person, living conditions and environment.

## General safety notes and driving restrictions

The engineering and construction of this wheelchair has been designed to provide maximum safety. International safety standards currently in force have either been fulfilled or exceeded in parts. Nevertheless, users may put themselves at risk by improperly using their wheelchairs. For your own safety, the following rules must be strictly observed.

Unprofessional or erroneous changes or adjustments increase the risk of accident. As a wheelchair user, you are also part of the daily traffic on streets and pavements, just like anyone else. We would like to remind you that you are therefore also subject to any and all traffic laws.

Be careful during your first ride in this wheelchair. Get to know your wheelchair.

Before each use, the following should be checked:

- Quick-release axles on the rear wheels
- Velcro on seats and backrests
- Tyres, tyre pressure and wheel locks.

Before changing any of the adjustments of this wheelchair, it is important to read the corresponding section of the user's manual.

It is possible that potholes or uneven ground could cause this wheelchair to tip over, especially when riding uphill or downhill. When riding over a step or up an incline frontally, the body should be leaning forward.

**Danger!**

NEVER exceed the maximum load of 125 kg for driver plus any items carried on the wheelchair. If you exceed the maximum load, this can lead to damage to the chair, or you may fall or tip over, lose control and may lead to serious injury of the user and other people.

**Danger!**

When it is dark, please wear light clothing or clothing with reflectors, so that you can be seen more easily. Make sure that the reflectors on the side and back of the wheelchair are clearly visible. We would also recommend that you fit an active light.

**Danger!**

To avoid falls and dangerous situations, you should first practice using your new wheelchair on level ground with good visibility.

**Danger!**

When getting on or off the wheelchair, do not use the footboards. These should be flipped up beforehand and swung to the outside as far as possible. Always position yourself as close as possible to the place where you wish to sit.

**Danger!**

Only use your wheelchair properly. For example, avoid travelling against an obstacle without braking (step, kerb edge) or dropping down gaps.

**Danger!**

The wheel locks are not intended to brake your wheelchair. They are only there to ensure that your wheelchair does not begin rolling unintentionally. When you stop on uneven ground, you should always use the wheel locks to prevent such rolling. Always apply both wheel locks; otherwise, your wheelchair could tip over.

**Danger!**

Explore the effects of changing the centre of gravity on the behaviour of the wheelchair, for example on inclines, slopes, all gradients or when overcoming obstacles. Do this with the secure aid of a helper.

**Danger!**

With extreme settings (e.g. rear wheels in the most forward position) and less than perfect posture, the wheelchair may tip over even on a level surface.

**Danger!**

Lean your upper body further forward when going up slopes and steps.

**Danger!**

Lean your upper body further back when going down slopes and steps. Never try to climb and descend a slope diagonally.

**Danger!**

Avoid using an escalator which may lead to serious injury in the event of a fall.

**Danger!**

Do not use the wheelchair on slopes  $> 10^\circ$ . The Dynamic safe slope is dependant on the chair configuration, the users abilities and the style of riding. As the users abilities and style of riding cannot be pre-determined then the max safe slope cannot be determined. Therefore this must be determined by the user with the assistance of an attendant to prevent tipping. It is strongly recommended that inexperienced users have Ant-tips fitted.

**Danger!**

It is possible that potholes or uneven ground could cause this wheelchair to tip over, especially when riding uphill or downhill.

**Danger!**

Do not use your wheelchair on muddy or icy ground. Do not use your wheelchair where pedestrians are not allowed.

**Danger!**

To avoid hand injuries do not grab in between the spokes or between the rear wheel and wheel lock when driving the wheelchair.

**Danger!**

In particular when using lightweight metal handrims, fingers will easily become hot when braking from a high speed or on long inclines.

**Danger!**

Only attempt stairs with the help of an attendant. There is equipment available to help you, e.g. climbing ramps or lifts, please use them. If there is no such equipment available, then the wheelchair must be tipped and pushed, never carried, over the steps (2 helpers). We recommend that users over 100 kg in weight do not use this stairway manoeuvre!

In general, any anti-tip tubes fitted must be set beforehand, so that they cannot touch the steps, as otherwise this could lead to a serious tumble. Afterwards the anti-tip tubes must be set back to their correct position.

**Danger!**

Make sure that the attendant only holds the wheelchair using securely mounted parts (e.g. not on the footrests or the sideguards).

**Danger!**

When using the lifting ramp make sure that the anti-tip tubes fitted are positioned outside the danger area.

**Danger!**

Secure your wheelchair on uneven ground or when transferring (e.g. into a car) by using the brakes.

**Danger!**

If and whenever possible, during a journey in a specially fitted vehicle for disabled people, vehicle occupants should use the seats in the vehicle and the appropriate restraint system. This is the only way to ensure that occupants will have the maximum protection if there is an accident. When using safety elements offered by SUNRISE MEDICAL and using a specially designed safety system, lightweight wheelchairs can be used as a seat when being transported in a specially fitted vehicle. (See the Chapter on "Transportation").

**Danger!**

Depending on the diameter and setting of the castors, as well as the centre of gravity setting of the wheelchair, the castors may begin to wobble at high speeds. This can lead to the castors being blocked and the wheelchair may tip over. Therefore, please make sure that the castors are adjusted correctly (see the Chapter "Castors"). In particular, do not travel on an incline without brakes, travel at a reduced speed. We recommend that novice users use anti-tip tubes.

**Danger!**

Anti-tip tubes should prevent the chair tipping over backwards unintentionally. Under no circumstances should they take the place of transit wheels, and be used to transport a person in a wheelchair with the rear wheels removed.

**Danger!**

When reaching for objects (which are in front of, to the side or behind the wheelchair) make sure that you do not lean too far out of the wheelchair, as if you change the centre of gravity there is a risk of tipping or rolling over. The hanging of additional load (back pack or similar items) onto your chair backposts can affect the rearward stability of your chair, especially when used in combination with recliner backrests. This can cause the chair to tip backwards causing injury.

**Danger!**

For thigh amputees you must use anti-tip tubes.

**Danger!**

Before setting off, check that your tyre pressure is correct. For rear wheels it should be at least 3.5 bar (350 kPa). The max. pressure is indicated on the tyre. The knee-lever brakes will only work if there is sufficient tyre pressure and if the correct setting has been made (see the Chapter on "Brakes").

**Danger!**

If the seat and back sling are damaged, you must replace them immediately.

**Danger!**

Be careful with fire, in particular with burning cigarettes. Seat and back slings can be set alight.

**Danger!**

If the wheelchair is subject to direct sunlight for a long period of time, then parts of the wheelchair (e.g. frame, legrests, brakes and sideguard) may become hot (>41°C).

**Danger!**

Always make sure that the quick-release axles on the rear wheels are set properly and lock in. If the button on the quick-release axle is not pressed in, the rear wheel cannot be removed.

**Note !**

The effect of the knee-lever brake as well as the general driving characteristics are dependant on tyre pressure. The wheelchair is significantly lighter and easier to manoeuvre when the rear wheels are pumped up correctly and both wheels have the same pressure.

Note !

Make sure that your tyres have sufficient tread! Please note that you are subject to any and all traffic laws when driving in public traffic.

Note !

Always be careful with your fingers when working or adjusting the wheelchair!

The products shown and described in this manual may not be exactly the same in every detail as your own model. However, all instructions are still entirely relevant, irrespective of detail differences.

The manufacturer reserves the right to alter without notice any weights, measurements or other technical data shown in this manual. All figures, measurements and capacities shown in this manual are approximate and do not constitute specifications.

## TRANSPORTATION IN VEHICLES

This product is suitable to use as a seat in vehicles as per certificate ISO/FDIS 7176-19:2008 (E) and using the restraint system Q strain Q-5001-T2.

A wheelchair secured in a vehicle will not provide the equivalent level of safety and security of a vehicle seating system. It is always recommended that the user transfers to the vehicle seating. It is recognised that this is not always practical for the user to be transferred and in these circumstances where the user must be transported whilst in the wheelchair the following advice should be followed.

1. Confirm that the vehicle is suitably equipped to transport a passenger in a wheelchair, and ensure the method of access/egress is suitable for your wheelchair type. The vehicle should have the floor strength to take the combined weight of the user, the wheel chair and accessories.
2. Sufficient space should be available around the wheelchair to enable clear access to attach, tighten and release the wheelchair and occupant tie down restraints and safety belts.
3. The occupied wheelchair must be located in a forward facing position and secured by the wheelchair tie down and occupant restraint straps (WTORS tie downs) meeting the requirements of ISO 10542 or SAE J2249, in accordance with the WTORS manufacturers' instructions.
4. The wheelchair's use in other positions within a vehicle has not been tested e.g. transportation in a side facing position must not be carried out under any circumstances. (Fig. 10)

Warning There is a risk of serious injury or death if this advice is ignored.

5. The wheelchair should be secured by a Tie Down Restraint system, conforming to ISO 10542 or SAE J2249 with non-adjustable front straps and adjustable rear straps, which typically use Karabiner clips/S hooks and tongue and buckle fittings. These restraints generally comprise of 4 individual straps that are attached to each corner of the wheelchair.
6. The tie-down restraints should be fitted to the main frame of the wheelchair as indicated in the diagram on the following page, and not to any attachments or accessories, e.g. not around the spokes of wheels, brakes or footrests.
7. The tie-down restraints should be attached as close as possible at an angle of 45 degrees and tightened securely in accordance with the manufacturer's instructions.
8. Alterations or substitutions must not be made to the wheelchair securement points or to structural and frame or components without consulting the manufacturer. Failure to do so will invalidate the ability of a Sunrise Medical wheelchair to be transported within a vehicle.
9. Both pelvic and upper torso restraint belts must be used to restrain the occupant to reduce the possibility of head and chest impacts with the vehicle components and serious risk of injury to the user and other vehicle occupants. (Fig. 11) The upper torso restraint belt should be mounted to the vehicle "B" pillar - failure to do so will increase the risk of serious abdominal injuries to the user.

Fig.10



Fig. 11


10. A head restraint suitable for transportation (see label on headrest) must be fitted and suitably positioned at all times during transportation.
11. Postural supports (lap straps, lap belts) should not be used or relied on for occupant restraint in a moving vehicle unless they are labelled as meeting the requirements specified in ISO 7176-19:2001 or SAE J2249.
12. The safety of the user during transportation depends upon the diligence of the person securing the tie-down restraints and they should have received appropriate instructions and/or training in their use.
13. Wherever possible remove and stow safely away from the wheelchair all auxiliary equipment, for example:  
Crutches  
Loose cushions  
Tray Tables
14. Articulating/elevating leg rest should not be used in the elevated position when the wheelchair and user are being transported and the wheelchair is restrained using Wheelchair Transport and Occupant Restraints.
15. Reclining backrests should be returned to an upright position.
16. The manual brakes must be firmly applied.

## Occupant Restraint Instruction

1. The pelvic restraint belt must be worn low across the front of the pelvis so that the angle of the pelvic belt is within the preferred zone of 30 to 75 degrees to the horizontal.

A steeper (greater) angle within the preferred zone is desirable i.e. closer to, but never exceeding 75 degrees.  
(Fig 12)

Fig.12

2. The upper torso restraint belt must fit over the shoulder and across the chest as illustrated Fig 13 and Fig 14. Restraint belts must be adjusted as tightly as possible consistent with user comfort. Restraint belt webbing must not be twisted when in use. The upper torso restraint belt must fit over the shoulder and across the shoulder as illustrated in Fig 13 and Fig 14.
3. Restraints should be mounted to the vehicle "B" pillar and should not be held away from the body by wheelchair components such as armrest or wheels.

Fig.13



Fig. 14

4. The attachment points to the chair are the inner front side frame just above the castor and the rear side frame. The straps are fitted around the side frames at the intersection of the horizontal and vertical frame tubes. (See Figs 15-20)
5. The tie down symbol (Fig 16) on the wheelchair frame indicates the position of the wheelchair restraint straps. The straps are then tensioned after the front straps have been fitted to secure the wheelchair.


## Transportability – positioning of wheelchair tie down restraints on wheelchair

Fig. 15

Fig. 16 Tie down label

The wheelchair secured with front and rear wheelchair tie down restraints. Positioning of the restraint straps is shown in more detail below.

Fig. 17

Fig. 18

Position of the front wheelchair tie down restraint and the tie down label.

Fig. 19

Fig. 20

Position of the rear wheelchair tie down restraint and the tie down label on the wheelchair.

## C. HANDLING

1. Opening: Rest both hands on each of the seat upholstery support rails and press down until the seat and back upholsteries become taut. Lower the footplates into the using position. As a safety measure, the footplates will always be in an upright position before sitting or getting out of the chair in a forward movement. Warning: Do not place your fingers between the seat upholstery and the seat tubes when opening or folding the wheelchair

2. Folding: First flip up the (platform) footplate or the individual footplates. Take hold of the sling from underneath in the middle and pull upwards. The wheelchair will then fold up. To fold your wheelchair until it is as small as possible, e.g., to put into a car, you can remove the footrests (depending on the model). For this purpose, open the latch from the out-side and swing the footrest to the side. Then the footrest only has to be drawn out of the frame tubing.

3. Folding back-rests (models 315 and 310) (fig. 1) To fold down the back-rest, press the black buttons on the components on the back-rest tubing and pull back the push handles.

4. Reclining back (4° - 55°) (Fig. 2): By pulling the 2 levers at the same time the backrest is released and can be set to the required position. Releasing the 2 levers will lock the backrest automatically in the new position. To bring the backrest back to the upright position, just move it (no release necessary) to the required position.

Caution: The reclining back is only to be used in combination with the active anti-tips setting

Caution: The reclining back is only to be used in combination with a stabilizer bar. Please always make sure that the stabilizer bar is locked in place on both sides.

Caution: Before folding the wheelchair, make sure that the stabilizer bar is removed

Caution: In order to avoid any reduction in the cable function, you must make sure that the cables are not completely taut.

Caution: When adjusting the back angle, be careful not to get your fingers caught.

Caution: Before adjusting the back angle, put the anti tip tubes in to the active position

5. elevating legrest (Fig. 3) To elevate: Lift your lower leg and press the release button. Pull the legrest upwards and adjust it to the desired height. As soon as you release the button, the legrest will click back into place. To lower: Lift your lower leg and press the release button to unlock the legrest to lower it. As soon as you release the button, the legrest will click back into place. CAUTION: The elevating legrest is not allowed to be adjusted into the lowest position to avoid any interference with the casters. CAUTION: Keep hands clear of the adjustment mechanism between the frame and the movable parts of the legrest while elevating or lowering the legrest.

6. Hemiplegic propulsion wheels: To move forwards, push both rims forwards at the same time. To move backwards, push both rims backwards at the same time. To move to the right, push the smaller, external rim. To move to the left, push the larger, internal rim.

## D. ASSEMBLY

1. Footrests (Fig. 4): The footrests can be rotated inwards or outwards by pressing the release latch and turning the footrest in the desired position. To remove the footrest, press the release latch, rotate the footrest approximately 90° and lift.

2. Armrests (Fig. 5): The armrest can either be flipped back or can be detached from the chair. To flip arms back, press the release latch and pull arm toward back of the wheelchair. To detach the armrest press the button located on the back of the armrest and lift it. CAUTION: Make sure that the armrest is properly locked into place. Never lift the wheelchair by the armrests as these parts are detachable and lifting the wheelchair by them may cause damage to the chair and/or injury to the user.

3. Head-rest (Fig. 6): To install the head-rest, push the head-rest tubes into the tubes at the top of the back-rest until the locking pins click into position in the holes. To remove it, press the locking pins and pull the head-rest up until it comes off.

4. Reclining back-rest stabilising bar (Fig. 6): Installation: Insert the bar into the tubes on the back of the back-rest. Press the locking pins and hold them down while inserting the stabilising bar, until the pins lock into place in the holes in the bar. To remove the bar, reverse the process.

## 5. Lap belt

Before using your wheelchair ensure the seat belt is worn and correctly adjusted before use.

The lap belt is fitted to the wheelchair as shown in the illustrations. The lap belt is fitted between the armrest support and the backrest tube. The belt should then be run across the outside of the back upholstery.

It is important that the belt is routed around the back of the wheelchair. The buckle end can be fed between the armrest and backrest.

Adjust the belt position so that the buckles are in the centre of the seat.

Adjust the lap belt to suit the user's needs as follows:

To increase the belt length:

Feed the belt through slide adjusters and male buckle to provide more belt length.

To reduce the belt length:


Feed the belt back through male buckle and slide adjusters.

Ensure the belt is not looped at the male buckle

When fastened check space between the lap belt and the user, when correctly adjusted it should be possible to insert the flat of the hand between the lap belt and the user.


Generally the Lap Belt should be fixed so that the straps sit at an angle of approximately 45°, and when correctly adjusted should not allow user to slip down in the seat.

To fasten buckle:

Firmly push the male buckle into female buckle.

To release belt:

Press the exposed sides of the male buckle and push towards the centre whilst gently pulling apart.

Do not rely on the lap belt only when the wheelchair and occupant are transported in a vehicle, use the separate occupant lap and diagonal restraints provided in the vehicle.

### Advice to Client

The lap belt must be checked on a daily basis to ensure it is adjusted correctly and it is free from any obstruction or adverse wear.

Failure to make sure that the lap belt is secure and adjusted prior to use could cause serious injury to the user. E.g. too loose a strap may allow the user to slip down in the chair and risk suffocation.

### Maintenance

Check lap belt, and securing components; at regular intervals for any sign of fraying or damage. Replace if necessary. Clean the lap belt with warm soapy water and allow to dry.

### Note


The lap belt should be adjusted to suit the end user as detailed above

Sunrise Medical also recommend that the length and fit of the belt is checked on a regular basis to reduce the risk of the end user inadvertently re-adjusting the belt to an excessive length

If in doubt about the use and operation of the seat belt then ask your healthcare professional, wheelchair dealer, carer or attendant for assistance.

### 6. Transport Kit

## E. ADJUSTMENTS

1. Footrest heights (Fig. 7) To adjust the height of the footplates, loosen the retaining bolt in the frame tube, reposition footplate and tighten bolt. For these adjustment use a 10 mm. wrench. A minimum interval of 2.5 centimetres from the ground should always be maintained.

2. Wheel locks (Fig. 8) For a slight adjustment, loosen the nut of the screw on the clamp, using a 10 mm wrench. Slide the whole wheel lock forwards or backwards, through the groove of the clamp. For a larger adjustment, remove the screw and reset the wheel lock in the appropriate hole in the frame. Fix the new position with the screw and nut.

3. Height-adjustable push handles (option only available on 305 and 312 models) (Fig. 9): To adjust the height of the push handles, loosen the knob that holds the tube of the push handle, turning it to the left. Move the tube up or down. To fix the position, turn the handle clockwise.

4. Backrest height (Fig. 9) (Only with the adjustable push handles option): Turn the handles to the left, until they are fully released. Remove the push handles, and remove the back upholstery from the tubes of the backrest. These tubes have 4 holes, which allow the height adjustment. To slide the tube upwards or downwards, remove the nut and the screw, using a 10 mm wrench and a 4 mm Allen key. Once they are removed, reposition the tube and fix it with the screw through the holes. The wheelchair leaves the factory with a pre-set back height of 40 cm.

## F. CONTROLS AND MAINTENANCE

**MAINTENANCE:** In addition to the maintenance and checks, general cleaning of both frame and upholstery can be undertaken. To maintain its appearance the upholstery should be cleaned regularly. Light soiling can be removed by sponging with a mild detergent and warm water, wiping in the direction of the texture. Rinse off with clean water before allowing to dry. Do not put upholstery in a dryer. If dirt is ingrained, the surface should be scrubbed with a soft brush. The frame should be cleaned regularly with a non-abrasive cleaner followed by polishing with a non-abrasive automotive wax.

**CONTROLS:** A periodic check of the following procedures will help sustaining the life and efficiency of your wheelchair.

**Weekly:** Check de pressure of the wheels following the instructions indicated on them. Check that the cross brace's central screw is tight and in good condition / Clean the upholstery with soapy water or adequate products / Clean the metallic parts of the chair with a dry cloth or with special products for chromium plated or painted surfaces / Check the wheel locks and adjust if necessary / Check wear of front and rear tyres.

**Quarterly:** Check and adjust if necessary: Castor play / Play in the castor forks / Push wheels' play / Offcentering of the push wheels / Lubricate with industrial oil or grease all parts subject to friction or wear / Retighten screws and nuts / A general revision of the chair's condition, repairing whatever is necessary.

**Yearly :** We recommend a complete revision once a year in a specialized establishment

## G. SPECIFICATIONS

Footrests: 70° Swing away, detachable and height adjusting. Elevating leg-rests optional

Armrests: Detachable and flip back

Standard Castors: 200 mm solid

Standard Rear Wheels: 600mm pneumatic in models: 305, 315 and 321. 300 mm. pneumatic in models 312, 310 and 341

**WHEEL LOCKS:** HIGH PUSH-TO-LOCK.

Spare Parts: Each wheelchair component is available as an spare part.

Accessories : Foot-rests: elevating. Wheels: Solid tyres. Wheel locks: extension handles and drum brake. Others: belt, anti-tip tubes, table, crutch and dripping device supports.

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## GUARANTEE

SUNRISE MEDICAL guarantees this product for two years from the date indicated on the invoice, receipt or delivery note. For products containing batteries, the battery warranty is 6 months. This warranty covers parts in need of replacement, labour and transport.

In warranty claims, the user must submit the product and the purchase invoice as proof of delivery date.

This warranty does not cover:

- NORMAL WEAR AND TEAR, e.g. on upholstery, tyres, inner tubes, rims, wheels, fabric, etc.
- ACCIDENTS e.g. bent structures, scratches, flat tyres, etc.
- IMPROPER USE and damage caused by users who weigh more than the weight limit indicated in the user's manual.
- LACK OF MAINTENANCE and failure to follow the storage instructions in the user's manual.
- MODIFICATIONS or REPAIRS carried out by persons not authorised by Sunrise Medical.

PRODUCTS WITH NO INVOICE OR SERIAL NUMBER (when appropriate).

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Sunrise Medical AB
Sunrise Medical B.V.
G2 A/S