
Leica HI 1220 – Instrucciones de uso

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RUN/STOP
(CON/DES)

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
RUN/STOP

SET

Arrow buttons for
temperature
selection

1. **Connect** the power cable.
2. **Turn on** the main switch on the back of the instrument. The green **POWER** indication on the control panel at the front will illuminate.
3. Push the **RUN/STOP** button on the control panel to heat the Leica HI1220 up to the previously entered set temperature. The display reads the actual value.

All surface materials are resistant to detergents commonly used in laboratories. The instrument must be switched off and cool prior to cleaning.

Do not use organic solvents such as xylene for cleaning!

4. The required temperature is selected using the **ARROW** buttons. At this stage, the value indicated in the display is the set temperature. Once the appropriate temperature is indicated, the **ARROW** button should be released. The indicated set value will be memorised automatically. After about 2 seconds, the display returns to the actual temperature.

To select the required temperature hold the corresponding ARROW button depressed. For the first 8 digital steps the indication will change slowly. After that, the speed of the step indication will increase.

5. The set temperature can be verified pushing the **SET** button.

The set temperature is **maintained** even if the Leica HI1220 is switched off (with the RUN/STOP button or main switch), disconnected from the power supply or in the event of a power failure.

6. **To switch off** push the **RUN/STOP** button. In daily routine operation, it is not necessary, however, to switch the instrument off with the mains switch.

7. Cleaning

Lateral ventilation slits are avoided on the Leica HI1220. The control panel is protected with a foil. This ensures smooth, easy-to-clean surfaces.

8. Caution

A change in the set temperature while the Leica HI1220 is in operation can cause a slight temperature increase over the the selected value for a moment. Therefore, specimens that are very susceptible to temperature should not be placed on the hot plate before the appropriate temperature is reached.

Safety precautions

- The instructions for use must be read carefully prior to operating the Leica HI1220.
- It must be ensured that the power supply conditions in the laboratory correspond to the data indicated on the label of your Leica HI1220 prior to connecting it to mains.
The Leica HI1220 should be connected to shockproof sockets only.
- The surface of the instrument can be very hot during operation.

Risk of burning!
Combustible and inflammable substances must not be stored near the instrument.

- Prior to repair or maintenance, the Leica HI1220 must be switched off and disconnected from mains.

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Leica Biosystems Nussloch GmbH