
Installation and Operation Manual | Installations- und Bedienungshandbuch |
Mode d'emploi | Manuale per l'installazione e l'uso | Manual de instrucciones

arium[®] 61316

Version 1.0 | Version 1.0 | Version 1.0 | Versione 1.0 | Versión 1

About This Manual

English	Page 2
Deutsch	Seite 29
Français	page 56
Italiano	pagina 83
Español	página 110

This manual instructs you on how to install, maintain and safely operate the arium® 61316 reverse osmosis system.

Sartorius Stedim Biotech has carefully designed your arium® 61316 water purification system for reliability, economy and safe operation. To ensure this, you must read these operating instructions carefully before operating the system. The manual contains important instructions that will help you ensure the reliability of your arium® 61316 and avoid potential hazards.

Please read the Safety Information on page 3.

If you have any questions about the correct use of arium® 61316, please contact us at the address below or call your local Sartorius Stedim Biotech dealer:

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
www.sartorius-stedim.com

Contents

About This Manual	2
Contents	2
Safety Information	3
Intended Usage	4
Product Description	5
arium® 61316 System	5
Control Panel Display	7
Power Connection Main Switch	7
Unpacking and Installation	8
Unpacking	8
Display Installation	9
Bench Mounting	10
Wall Mounting	10
Water Inlet Connection	11
Concentrate Outlet Tube	11
Permeate Product Water Tubing	11
Pretreatment Cartridge Installation	12
RO Module Installation	13
Pressure Tank Installation	14
Open Tank Installation	16
Initial Operation	17
System Settings	17
Language	17
Date and Time	17
Measurement Units	18
Set Point	18
Tank System	18
Audio Alarm	19
Cleaning	19
Further Options	20
System Information	20
Connecting a Printer	20
Transferring Data to PC	20
Maintenance and Servicing	21
System Cleaning	21
Tank Cleaning	22
Pretreatment Cartridge Replacement	23
RO Module Replacement	24
Fuse Replacement	24
Appendix	25
Information and Instructions on	
Disposal and Repairs	25
Troubleshooting	26
System Specifications	27
System, Consumables and Accessories	28

Safety Information

Please read through the following safety information carefully and follow the instructions exactly. This information is designed to ensure your own safety and will help prevent damage to the arium® 61316 unit.

The following symbols are used in this manual:

Warning!

Warnings alert you to a possibility of personal injury or property damage. Safety data sheets on specific products are available from Sartorius Stedim Biotech on request.

Caution!

Caution signs alert you to a possibility of damage to the equipment.

Note!

Notes alert you to pertinent facts and conditions that are important for economic operation of the equipment.

Warning!

Danger of personal injury! Servicing and repairs may only be performed by trained and qualified personnel.

Warning!

Severe electrical shock hazard or danger of electrocution!

- Use a properly grounded electrical outlet of correct voltage (100–240 VAC, 50–60 Hz) and current handling capacity to plug in the arium® 61316 unit.
- Do not place the arium® 61316 on top of electrical equipment. Routine maintenance of this unit may involve water spillage and subsequent electrical shock hazard, if it is improperly located.
- The units may only be mounted in rooms equipped with a floor drain to prevent any water damage.

Warning!

Danger of fire or explosion!

- Do not use the arium® 61316 in the presence of flammable or combustible materials; fire or explosion may result. This device contains components which may ignite such materials.
- The arium® 61316 is to be only fed with potable water of a suitable grade. Contact Sartorius Stedim Biotech whenever there is doubt.

Warning!

Danger of injury to eyes and skin!

- Avoid splashing sanitizing | cleaning solutions on clothing, eyes or skin.
- Wear protective eyeglasses and rubber | latex gloves when performing the sanitization | cleaning procedure. Laboratory chemical safety procedures should be followed when handling any laboratory chemicals.
- Ensure that all tubing connections are tight to avoid leakage.
- Carefully follow all the manufacturer's safety instructions on labels of the sanitizing | cleaning solutions.

Caution!

Danger of irreversibly damaging the arium® 61316 components!

- Be sure to replace defective fuses with those of the same type and rating.

Intended Usage

The arium® 61316 system is exclusively designed to provide reverse osmosis water for laboratory use. To ensure that the unit works properly, use only the filter media and other auxiliary media that are listed in this manual. Using this system for any other purpose shall be considered improper usage.

- The arium® 61316 may only be operated by trained personnel.
- Do not open the rear panel of the system. If the seal is broken, the warranty claim is null and void. If you encounter any problems with the system, please contact your local Sartorius Stedim Biotech Service Center.
- Operate the arium® 61316 only using original accessories or replacement parts. If you modify this water purification system on your own without consulting Sartorius Stedim Biotech, the performance and operating safety of the system are no longer guaranteed, thus constituting a safety hazard for the operator.
- Please take all pertinent precautions to prevent accidents and always observe the generally valid technical and occupational safety rules and regulations of your country. Only connect to grounded electrical outlets.
- Use only Sartorius Stedim Biotech-listed materials (such as connectors, gaskets, tools, sanitizing | cleaning agents, prefilter and cartridges).

Product Description

arium® 61316 system

The arium® 61316 is a high performance reverse osmosis (RO) system comprised of four principle parts: housing (incorporating electronics, pump and tubings), RO modules, pretreatment cartridge and pressure tank, illustrated in the flow chart.

A safety solenoid valve is upstream from the actual water purification. This valve is only open when the system is in the operating mode. After the feed water has passed the built-in pressure regulator, it is pre-filtered with help of the pretreatment cartridge. A pump maintains the proper operating pressure and forces the water through the parallel connected two RO modules. Product water can both be stored in a pressure tank or in an open tank.* The concentrate is fed directly into the waster water drain. For quality control the system performs a temperature compensated conductivity measurement in both the feed water and the product water.

If a pressure tank (optional) for storing water has been installed, it will shut off the RO system when the system pressure reaches 3 bar ("preoperate" mode) and will automatically initiate a backflushing cycle of the RO-modules with permeate.

The automatic backflush program dissolves residues that built up on the untreated water side and flushes them out to the drain. This technology guarantees an exceptionally long life of the RO modules. The system is automatically flushed whenever the system switches from the "operate" (tank <100%) mode to the "preoperate" (tank level = 100%) mode or by pushing the "standby" button.

The reverse osmosis system should be cleaned regularly. Sartorius Stedim Biotech supplies suitable sanitizing | cleaning solutions. Cleaning takes place automatically.

Feed water quality and permeate quality are continuously monitored: product water quality and rejection rate of the RO modules are continuously indicated in the display. When water quality is out of tolerance with the manually entered <Set point>, the appropriate visual alarm appears. An audio signal can be selected optional for all alarm messages.

If the system has a depressurized tank (optional) connected downstream, this tank will be controlled via a float switch (level regulator).*


The two RO modules and the pretreatment cartridge are consumables which require replacement from time to time. The interval between replacements will depend upon water consumption and feed water quality. A note for replace any consumables is indicated in the display.

Automatic sanitization programs are available for cleaning both the reverse osmosis system and the tank.

The various individual parameters can be set in the <Maintenance> mode.


An RS232 serial interface allows the user to document the product water quality in accordance with GLP guidelines (Good Laboratory Practice).

* on request


Flow chart

Item No.	Description
1	Left and right door of the arium® system
2	Display and control panel
3	Permeate outlet
4	15-pin D-sub port for connecting a remote display unit (only in conjunction with a below-bench unit)
5	Reject water outlet
6	Fuse drawer
7	Main power switch
8	Power cord receptacle
9	Serial interface RS232 with 9-pin D-sub port for connecting a printer PC
10	PLC interface for external communication
11	Tank port
12	Feed water inlet

Front view

Side view, left

Control Panel | Display

Operate the arium® system using the control panel that features four function keys and two control keys for the cursor.


Standby Key

Press the **Standby** key to switch the unit into “Standby” mode. The production of reverse osmosis water will be interrupted and a 4 minute automatic backflush cycle will rinse the RO modules. By flushing the modules, any residues that build up on the untreated water side during operation will be passed into the drain. The life time of the cartridges will be expanded.

Operate Key


With the **Operate** key you can interrupt the “standby” mode and activate the unit. Afterwards the system will start an internal initialisation and check the feed pressure of the system and the tank. The arium® 61316 will switch in the “operate” mode to produce RO water, if the tank pressure is < 2.5 bar.

When the system is in the “operate” mode, the display will show:

Menu Key

This key is used to view the System Menu and all its options.

Cursor Keys

Use these two arrow keys to move within the menu and highlight individual menu items. With the  key, move the cursor up a line, with the  key, move the cursor down a line.

Enter Key

Use the Enter key to confirm a menu item you have selected with the cursor.

Power Connection | Main Switch

Note


The arium® 61316 can be connected to single phase power sources ranging from 100-240 VAC with 50-60 Hz.

1. Use the supplied power cord and plug it into the outlet on the lower left side of the system.
2. The system power is located on lower left side. Press (“I”) on to activate the system and (“O”) to deactivate the system.

Note

All installations have to be finished before activating the system.

Control Panel | Display

AC power connector port, left-hand side of housing

Unpacking and Installation

Unpacking

1. To reduce accidental damage occurring during shipment some components have been packed separately.
2. Inspect all packages for damage. Open the boxes carefully. Remove the system and place it on a level surface. Remove the plastic bag from the unit. Inspect the system for any damage caused during transport.
3. Remove the accessories and take inventory of the equipment supplied. The accessories are shipped inside the unit and in the top cover of the box. Refer to the Equipment Supplied Section.

A wall mounting bracket is optionally available.

Note!

If you have only ordered the system (order no.: 61316), then you have to order the prefilter cartridges, RO modules and tank separately. When you order a complete system (for example, order no.: 61316030F05M1A) all components are included.

Equipment Supplied

Box 1

(order no. 61316)

- arium® system 61316
- Inlet water tubing with quick-fit adapter, PE, 3/8" outer diameter, 2.40 m length (7.9 ft.)
- Concentrate waste tube, PE, 1/4" outer diameter, 3.10 m length (10.2 ft.) (with red band at one end)
- Permeate tube, PE, 1/4" outer diameter, 3.10 m length (10.2 ft.) (with blue band at one end)
- Cleaning | Sanitizing syringe for storage tank (Order no. 611CDS)
- Cleaning | Sanitizing syringe for RO unit (Order no. 612CDS)
- Wall mounting bracket (optional)
- Power cord
- Operation manual
- QA Certificate
- Kit with accessories for installation

Box 2

(Order no. 613CPM4-----V delivered separately)

- 2 × RO modules

Box 3

(Order no. 613CPF05-----V delivered separately)

- 2 × Pretreatment cartridges

Note!

One pretreatment cartridge is required to complete initial installation of the assembly and one spare is intended for future replacement.

Box 4

(Order no. 613APV31 delivered separately)

- 30 liter storage tank* or

(Order no. 613APV50, delivered separately)

- 50 liter storage tank* or

(Order no. 613APV70, delivered separately)

- 70 liter storage tank* or

(Order no. 613APV100, delivered separately)

- 100 liter storage tank*

* at 3 bar

Note!

All installations have to be finished before "Initial Operation".

Display Installation

The display unit has been designed to be mounted in 3 ways:

At the top of the right door for bench systems, at the bottom of the right door for wall-mounted systems or in a remote location up to 2.8 m from the system (additional kit required).


To install:

- Route the supplied cables through the large opening of the right door (top or bottom).
- Use the 4 screws and washers provided to mount the display unit on the door by routing the screws from inside the door through the holes in the door to the outside and screw to the assembly.
- Connect the data and power cord of the display to the plug-in, inside the system.


Follow the course of dotted lines in the figure:

Position 1 for top mounted display
Position 2 for bottom mounted display

- Place rubber cover over the remaining recess (top or bottom) in the doors.

Mounting the display on the top

Mounting the display on the bottom

**Warning!**

Severe electrical shock hazard or danger of electrocution!

- Do not place arium® 61316 on top of electrical equipment. Routine maintenance of this unit may involve water spillage and subsequent electrical shock hazard, if located improperly.

**Warning!**

Danger of fire or explosion!

- Do not use in the presence of flammable or combustible materials; fire or explosion may result. This device contains components which may ignite such materials.
- Place the arium® system on a flat surface.
- Ensure that a feed water supply (pre-pressure minimum 1.2 bar and maximum 6.8 bar), an atmospheric drain and a suitable power supply (100–240 VAC, 50–60 Hz) are available.

Bench Mounting

Place the system on a suitable workbench. Make sure that the specified bench space is available and that the surface supports the system's operating weight 20 kg (44 lb).

Wall Mounting

A wall bracket (2) allow you to attach the system securely to the wall with a minimum amount of space. A clear wall area of 63 × 63 cm is required to mount the arium® system. Two sturdy keyholes (1) are located at the top of the back of the housing to enable you to hang the unit on the wall bracket using two screws (5). Use the predrilled holes in the bracket to fasten the bracket to the wall. The metal profile (6) located on the lower arium® housing functions as a spacer.

**Note!**


The wall mounting bracket and the required screws are not included in the standard equipment supplied with the unit. These have to be ordered separately.

1. Fasten the bracket to the wall with suitable screws and dowels.
2. Loosen the metal support (4) using the screw (3) and slide the support to the left.
3. Push the arium® system onto the pins and lower the system to engage the keyholes (1).
4. Slide the metal support (4) to the right until the metal support rests against the arium® housing and lock it using the screw (3).


**Warning!**

The reverse osmosis system operating weight is more than 20 kg (44 lb). Make sure that the wall is stable enough for mounting the system.

- Injury or equipment damage could occur if the system is not mounted properly.

arium® backside with keyholes and spacer

Wall bracket

Water Inlet Connection

1. Connect the tubing for the feed water connection, left on the unit, to the quick-connect adapter (1). Connect the other end to a potable water supply. A shut-off valve should also be installed.
2. Make sure that the pre-pressure is sufficient, minimum 1.2 bar up to maximum 6.8 bar.

Note!

All installations have to be finished before opening the feed water line.

Concentrate Outlet Tube

1. The concentrate water is drained through the outlet (3). The outlet tube (red coloured label at the end) must be connected to the red outlet connector at the rear on the lower right side of the system. Push the red end of the tube into the connector until the beginning of the label. Ensure that the tube is fastened and a free flow to the drain is possible.
2. Connect the other end of the tube to an atmospheric drain.


Permeate | Product Water Tubing

1. Connect the permeate water tube (blue coloured label at the end) to the blue fitting (2) located on the middle right side of the system. Push the blue end of the tube into the connector until the beginning of the label. Make sure that the tubing is fastened tightly.
2. Connect the other end of the tubing to a storage tank or reservoir, as described in the following part, tank installation.


Note!

Please install an enclosed shut-off valve in front of the permeate inlet connection to the tank, this allows the water supply to be interrupted, e.g. for a cleaning cycle.

Connector for feed water tubing

Connectors for reject and permeate | product water tubing

Pretreatment Cartridge Installation

The pretreatment cartridge pack used in arium® 61316 is designed to protect the two high-grade RO modules.

Constructed of virgin polypropylene, the twin chambered casing incorporates an activated carbon and catalyst followed by a 5 µm (nominal retention rating) polypropylene depth filter. This combination reduces the levels of particulate matter, free chlorine and other oxidizing agents from the feed water.

Note!

Under normal operating conditions, pretreatment cartridges will require replacement on a regular basis. We recommend that you replace the pretreatment cartridge every 3 months. A corresponding message appears in the display. If a softener has been installed because the water is too hard, the pretreatment cartridge should be replaced every 6 months.


1. Remove one pretreatment cartridge from the box and remove packaging.
2. Remove the shipping plug from the connections on the cover of the pretreatment cartridge.
3. Check to see whether the cartridge is damaged.
4. Check that no water feed is connected.
5. Place the cartridge holder for connection on the cartridge (see picture).
6. Wet the gaskets.
7. The angle on the quick-connector prevents improper connection of the cartridge. Be aware that the safety bars are open.
8. Press the water connections firmly into the cartridge ports and push the lock lever firmly into its secure position.
9. Place the cartridge in the arium® unit and make sure that it is installed without kinking.

The installation of the prefilter cartridge has been completed.


Note!

Prior to new installation, the pretreatment cartridge has to be flushed. For further steps please follow the subsequent chapter "Initial Operation", under "Flushing".

Press the adapter onto the cartridge

RO Module Installation

RO modules are designed to remove up to 99% of salt and organic contaminants found in potable water. Configured in parallel, the two RO modules of the arium® 61316 are capable of producing up to 16 liters per hour with a 55% recovery rate. When the system shuts down, automatic membrane flushing is carried out and thereby prevents potential contamination and prolongs the life of the RO modules.

The RO modules are constructed out of thin-film composites, spiral-wound membranes encapsulated in (white) polypropylene housing.

Note!

The RO modules do not contain a preservative solution for shipment. No special sanitization is necessary.

To install them, please proceed as follows:

1. Remove both RO modules from the packaging.
2. Remove shipping plugs from inlet and outlet ports on the RO modules.
3. Check that no water feed is connected.
4. Note the color-coded connectors at the top and the bottom of the modules.
Grey = feed water inlet
Red = concentrate drain
Blue = permeate drain
5. Place the first RO module, in front of the safety ring, labeled with the corresponding blue number. Push the module in the installed safety ring. Repeat this procedure with the second RO module.

RO modules

Pretreatment cartridge

6. Install both RO modules, by connecting the coloured labeled tubes in the system with the same labeled connections on the RO modules.
7. Check to make sure that all connectors are securely connected. Pull on all the connections to make sure that they are all firmly locked into the safety ring.

The installation of the RO modules have been completed.

Note!

Prior to new installation, the RO modules have to be flushed. For further steps please follow the subsequent chapter "Initial Operation", under "Flushing".

Pressure Tank Installation

The arium® 613APV pressurized storage tank, is optional supplied together with the arium® 61316 system. It is a self-contained downstream vessel. All materials used in construction are FDA approved. The water cell contained within the tank is a flexible butyl rubber membrane which expands and collapses as it is filled or emptied. As the bag fills, the pressure will increase in the gas-filled interface between the outside of the bag and the internal surface of the tank wall. This will be indicated on the built-in pressure gauge on the tank. No external venting is required.

1. Make sure system power and feed water are switched off.
2. Remove storage tank from packaging and inspect for damage.
3. Install tank at any location up to approx. 3 meters (8") from the system.
4. Push the other end of the blue labeled permeate tube, which is connected with the system, into the product inlet on the tank until the tube bottoms out in the fitting.

Note!

Bear in mind that the supplied shut-off valve has to be installed between pressure tank and RO system. Please be sure that the shut-off valve is open.

Note!

Be aware that all previous described installations have been prepared and finished.

5. Open the feed water line to the RO system and switch on the power.

Note!

By the use of RO water or a tank pressure <2.5 bar the system starts automatically to produce reverse osmosis water. If the tank is full the system stops automatically.

6. For further settings please follow the chapter "Initial Operations", "Tank System".

After successful installation the bar graph in the display indicates the reservoir level. When the tank is filled (100 %), it has a system pressure of 3 bar.

Note!

Prior to new installation, the tank has to be flushed | cleaned. For further steps please follow the subsequent chapter "Initial Operation", under "Flushing".

7. Connect ultrapure water system or other RO water equipment to the tank.

To utilize the advantages of the system's pressure, a pressure tank should always be installed. If this is not possible, a depressurized tank may also be installed. In such cases, the system can be controlled by a float switch. A socket is located on the left side of the unit.

Note!

As supplied from the manufacturer, the arium® 613APV pressure tank is configured for upright (vertical) floor mounting. However, the tanks may be oriented in any direction to accommodate installation.

Note!

The arium® pressurized tanks have two connections for consumers: e.g. arium® 611 ultra pure water system, dish washer etc.

Note!

Refer to chapter "Service and Maintenance" of this manual for cleaning instructions if required.

arium® 613APV Pressure Tank

Usable volume of the tanks at 3 bar:

Order no.: 613APV100
100-liter volume


Order no.: 613APV70
70-liter volume

Order no.: 613APV50
50-liter volume

Order no.: 613APV31
30-liter volume

Permeate feed inlet

Quick-connect adapter, permeate outlets

Permeate feed inlet

Quick-connect adapter, permeate outlets

Pressure gauge indicator for tank

Pressure gauge indicator for tank

Open Tank Installation*

Equally, with the help of a multi-function box an open storage tank can be connected to the arium® 61316. The required accessory (610AMFB1) can be ordered via Sartorius Stedim Biotech.

Note!

Please ensure that the tank meets all the installation requirements. Further details on request.

This kind of tank is controlled via a float switch that starts the system when the tank is empty and switches off when the system's tank is full. In this case, the reservoir level display is run by the three reservoir level detectors on the tank.

1. Make sure that the system power and feed water are off.
2. Install tank at suitable location.
3. Connect the multi-function box, including level sensors, as described in the attached installation guide.
4. Press the other end of the blue labeled product water tube, which is connected with the system, into the product inlet connection of the tank.

If a Sartorius Stedim Biotech product is not available, use a cable that is closed when the reservoir level is low, in order to maintain contact between the system and production.

5. Open the feed water line to the RO system and switch on the Power.

Note!

With a low tank level the arium® 61316 will change, after system check in to the "operate" mode and starts the production of RO water.

6. For further settings please follow the chapter "Initial Operations", "Tank System".

After successful installation the bar graph in the display indicates the reservoir level.

Note!

After new installation, the tank must be flushed | cleaned.

7. Connect RO water equipment to the tank.

To utilize the advantages of the system's pressure, a pressure tank should always be installed.

Note!

Sartorius Stedim Biotech can give no guarantees for product water quality.

* on request

Initial Operation

Before the initial operation of the arium® 61316, the previously installed consumables must be put through a 60-minute rinsing process.

The arium system recognizes automatically whether or not it is an initial operation when the device is switched on.

The rinsing process is automatically activated when the device is switched on ("I") in the case of initial operation.

Please follow the notes that appear on the display.

Note!

All installation steps have to be completed before you switch on the system power and pump in feed water.

1. The main power switch on the left-hand side of the housing must be on ("I"). The arium® system will then perform a system check.

The arium® system display will show the following information (in English upon initial operation):

```

SYSTEM ARIUM RO
S/N xxxxxxxxxxxx
System Check OK
U1.XX
    
```

The display will show the model type of the system, the serial number and the software version.

Following this, there is a further system test for initial operation. In the case of initial operation, the system automatically changes to the rinsing process.

Note!

If initial operation has already been carried out, the system misses out the rinsing process.

Following this, the system carries out internal initialization as well as an inspection of the inlet pressure (minimum 1 bar, maximum 6.8 bar).

```

Initialization
---- Tank level
Time         Date
    
```

When there is a tank pressure of < 2.5 bar (in the case of depressurized tank volume < 50%), the system then changes to the "operate" mode or, in the case of a tank pressure > 2.5 bar (in the case of depressurized tank volume > 50%), to the "ready" mode.

The following display appears in the "operate" mode:

```

XX  %  Rejection Rate
XX  µS/cm comp.
----- Tank Level
Time                               Date
    
```

Basic settings

After the arium® 61316 rinsing process, set the language for the display text, date and time, the desired measurement units, the limit value and the tank system.

Note!

In some menus, not all menu items are visible on the display. Repeatedly press the arrow key (⏴) to display the required language.

Display Language

This option allows you to set the language for the menu texts in the display. The system is factory set to English.

Proceed as follows:

1. Press **Menu**: The display will show the following information:

```

MENU           Maintenance
                Setup
                Info
    
```

2. Select <Setup> and press **Enter**. The following display will appear:

```

Setup          Time/Date
                Print Mode
                Set Point
                Meas.Unit
                Tank
                Language
                Audio Alarm
    
```

3. Use the arrow key to scroll down, select <Language> and press **Enter**. The following display will appear:

```

Language       English
                German
    
```

Use the arrow keys to select the desired language and press **Enter** to confirm.

The display will then switch to "operate" mode.

Measurement Units

This feature allows you to select how measurement units are displayed. You can choose between $M\Omega \times cm$ or $\mu S/cm$, to display data, compensated to 25°C, and temperature.

To set the measurement units:

Proceed as follows:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Setup> and press **Enter** to confirm.

```
Setup     Time/Date
           Print Mode
           Set Point
           Meas.Unit
           Tank
           Language
           Audio Alarm
```

3. Select <Measurement Units> and press **Enter** to confirm.

The display will show the following:

```
Meas. Unit
           MΩ x cm comp.
           μS/cm comp.
```

4. Use the arrow keys to select the desired option and press **Enter** to confirm.

The system will then switch to the previously activated operation mode.

Set Point

The option <Set point> allows the user to select an individually defined product water quality in $\mu S/cm$ or $M\Omega \times cm$ downstream of the RO modules.

To change the set point:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Setup> and press **Enter** to confirm.

```
Setup     Time/Date
           Print Mode
           Set Point
           Meas.Unit
           Tank
           Language
           Audio Alarm
```

3. Select <Set Point> and press **Enter** to confirm. The display will show the following:

```
Set Point  XX MΩxcm comp.
           or
           XX μS/cm comp.
```

4. Select current set point (between 10–80 $\mu S/cm$ or 0.013–0.1 $M\Omega \times cm$) and use the arrow keys to change the desired set point.

5. Press **Enter** to confirm the set point.

The system will then switch to the previously activated Operation Mode

Note!

When the water quality drops up ($\mu S/cm$) or down ($M\Omega \times cm$), the defined set point in the display will show the following alarm message:

```
XX % Rejection Rate
Set Point Error RO
_____ Tank Level
Time _____ Date
```

In this case please contact your Sartorius Stedim Biotech Service.

Tank System

With this option you have to set the installed tank system (open gravity or pressurized tank) to the arium® 61316.

The factory setting is always one pressure tank.

Proceed as follows:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Setup> and press **Enter** to confirm.

```
Setup     Time/Date
           Print Mode
           Set Point
           Meas.Unit
           Tank
           Language
           Audio Alarm
```

3. Use the arrow key to scroll down, select <Tank> and press **Enter** to confirm. The display will show the following:

```
Tank      Pressure Tank
           Open Tank*
```

* on request

4. Use the arrow keys to select the desired tank system and press **Enter** to confirm.

The system will then switch to the previously activated operation mode.

Audio Alarm

The option <Audio Alarm> allows the user to switch on an additional audio alarm signal.

To choose audio alarm:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Setup> and press **Enter** to confirm.

```
Setup     Time/Date
           Print Mode
           Set Point
           Meas.Unit
           Tank
           Language
           Audio Alarm
```

3. Select <Audio Alarm> and press **Enter** to confirm. The display will show the following:

```
Audio Alarm
           On
           Off
```

4. Choose <on> to get an audio alarm signal and press **Enter** to confirm.

The system will then switch to the previously activated operation.

Cleaning

All previously installed components should be cleaned during the initial operation.

Please proceed by following the maintenance steps in consecutive order:

1. "System Cleaning" (Page 21)
2. "Tank Cleaning" (Page 22)

Upon completion of these steps, the Initial Operation is concluded.

Further Options

System Information

With this option you will get additional information such as on tank level in "bar", feed water quality, and product water quality measured after the RO modules.

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
          Setup
          Info
```

2. Select <Info> and press **Enter** to confirm. Following information are displayed numerically:

```
INFO
Tank pressure (bar)
Conductivity feed (µS/cm)
Conductivity after RO
(µS/cm)
Temp (°C/°F)
```

3. Press **Menu** to go back to the Main Menu.

Printer Connection

The arium® 61316 reverse osmosis system can be connected to an optional printer by a serial interface RS232. If you have connected an arium® 611APR1 printer on the arium® system, following information can be printed:

- Model number
- Serial number
- Date and time
- Conductivity or resistivity

Note!

The connection to the arium® system will be done by an RS232 serial interface with a 9-pin sub-D port for a serial printer or a PC. The ported printer must use a separate power supply.

Similarly, every commercial printer can be used that accepts an RS232 interface. Please always follow the operating instructions and the manufacturer's instructions.

The arium® 611APR1 printer is preconfigured as follows:

Baud rate	19200cd
Parity	Space
Handshake	DTR CTS
P1:Send CR LF	ON
P2:Print format	Right
P3:Decimal-char	Point
P4:num.function	ON
P5.Inverse from	0
P6: print mode	standard
P7:Date Time	On
P8:Print text1	Off
P9:Print text 2	Off
C1:Time format	24 h
C2:Display	Time

With other printers, the settings may vary.

Connecting the Printer

1. Check whether the printer is suitable for an interface RS232 connection.
2. Connect the printer to a suitable power supply.
3. Connect the printer to the left side of the RO-unit.

4. To select a print option press the **Menu** button on the display, to activate the Main Menu.

```
MENU      Maintenance
          Setup
          Info
```

5. Select <Setup> and press **Enter** to confirm.

```
Setup     Time/Date
          Print Mode
          Set Point
          Meas.Unit
          Tank
          Language
          Audio Alarm
```

6. Choose <Print Mode> and press **Enter** to confirm. Following display is shown:

```
Print Mode
          Interval
          Display
```

With <Interval> you can make regular print outs at intervals from 0-59 minutes from the actual display screen of the "operate" mode. The print out will be done every 10 minutes.

By choosing the <Display> mode the printer will make only one hard copy from the actual display screen of the "operate" mode.

7. Set the desired form and confirm with **Enter**.

Afterwards the system will switch to the previously activated operation.

Note!

The standard arium® 61316 RO system is not supplied together with a printer. This must be purchased as an option.

Transferring Data to a PC

Using a standard interface cable you can connect a PC to the RS232 interface and transfer the data. The data can be displayed on your PC with the Hyperterminal software included in your Windows operating system or with any equivalent program.

Maintenance and Servicing

Every reverse osmosis system is supplied with potable water. In addition to various salts the potable water also contains different microorganisms present in the potable water. It is therefore only a matter of time before sanitization | cleaning is required. The unit is factory-set with the various maintenance cycles. Prompts will be given with regard to the work to be performed and consumables to be replaced.

Factory-set Service Intervals:

Repl. Prefilter Cartridge:	3 month
Repl. RO Module	24 month
Sanitization System	3 month
Sanitization Tank	6 month

After these times have elapsed, the corresponding alarm message will appear in the second display line alternating with the current display.

Note!

The alarm message is only activated in the “operate” display.

The pretreatment cartridge and RO module installation is described with the Initial Operation work. Disassembly and re-installation are performed accordingly. Please follow the instructions that appear in the display when prompted.

Sanitization is necessary to reduce bacterial growth and to minimize and | or inhibit the build-up of a biofilm within the cartridges or piping.

The frequency with which you will need to clean your unit depends on your feed water’s characteristics, your purity requirements and your water consumption.

Warning

Severe electrical shock hazard or danger of electrocution! Only trained qualified personnel may carry out service and repair work.

System Cleaning

Note!

The arium® system should be cleaned prior to new installation and at regular intervals.

To remind you to perform cleaning, the message “Timeout System Sanitization” will appear on the display according to the factory setting.

This message will blink on your display until a complete cleaning cycle has been performed.

The audio alarm can be stopped by switching the **Enter** button.

At the end of the sanitization cycle, this reminder function will reset automatically.

Warning!

The reverse osmosis system must be cleaned with a special cleaning | sanitization solution.

Wear protective eyeglasses and rubber | latex gloves when performing the sanitization procedure. Laboratory chemical safety procedures should be followed when handling any laboratory chemicals.

Have the following ready prior to cleaning | sanitization:

- Cleaning | Sanitization kit with Order no. 612CDS
- Container for 1–2 liter to prepare the cleaning solution
- 1 liter water (if possible use RO water)

To clean the system:

1. Prepare cleaning solution.

Therefore bring the content of one sanitization syringe 612CDS in the 1–2 liter container and fill it up with 1 liter water (if possible use RO water).

2. Press **Menu** to activate the Main Menu.

MENU	Maintenance Setup Info
------	------------------------------

3. Select <Maintenance> and press **Enter** to confirm. The display will show the following:

Maintenance	Repl. RO modules Repl. Prefilter Sanitization
-------------	---

System cleaning

4. If require, use the arrow key to scroll down and select <Sanitization> in the menu. Press **Enter** to confirm.

Sanitization	System Tank
--------------	----------------

5. Select <RO-Unit> and press **Enter** to confirm. The display will show the following:

Sanitization	System
Continue?	Yes No

Note!

The program can no longer be interrupted, once the safety prompt has been confirmed.

6. Choose <Yes> to proceed system cleaning and press **Enter** to confirm. Follow the instructions that appear in the display and push **Enter** after each step is completed.

Note!

After sanitization check to make sure that all tubes are securely installed and leak tight.

After all work is completed, the unit is operational again.

Tank Cleaning

Note!

The arium® tank should be cleaned prior to new installation and at regular intervals.

To remind you to perform cleaning, the message Timeout “Tank Sanitization” will appear on the display according to the factory setting.

This message will blink on your display until a complete cleaning cycle has been performed.

The audio alarm can be stopped by switching the **Enter** button.

At the end of the sanitization cycle, this reminder function will reset automatically.

Warning!

The tank system must be cleaned with a special cleaning | sanitization solution.

Wear protective eyeglasses and rubber | latex gloves when performing the sanitization procedure. Laboratory chemical safety procedures should be followed when handling any laboratory chemicals.

Have the following ready prior to cleaning:

- Cleaning | Sanitization kit with Order no. 611CDS

To clean the tank:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Maintenance> and press **Enter** to confirm. The display will show the following:

```
Maintenance
Repl. RO modules
Repl. Prefilter
Sanitization
```

3. If require, use the arrow key to scroll down and select <Sanitization> in the menu. Press **Enter** to confirm.

```
Sanitization  System
              Tank
```

4. Select <Tank> and press **Enter** to confirm. The display will show the following:

```
Sanitization  Tank
Continue?      Yes
              No
```

Note!

The program can no longer be interrupted, once the safety prompt has been confirmed.

5. Choose <Yes> to proceed tank cleaning and press **Enter** to confirm. Follow the instructions that appear in the display and push **Enter** after each step is completed.

Note!

For the injection of the sanitization solution unscrew the luer end cap from the connection of the adapter. Inject the cleaning | sanitization liquid through the port into the system (see picture) and then remove the syringe. Afterwards reattach the luer end cap.

Note!

After sanitization check to make sure that all tubes are securely installed and leaktight.

After all work is completed, the unit is operational again.

Tank cleaning

Pretreatment Cartridge Replacement

The service life time of the cartridges depends on the quality and volume of the inlet water treated. The cartridges should be changed at regular intervals.

To remind you to perform cleaning, the message "Repl. Prefilter" will appear on the display according to the factory setting.

This message will blink on your display until the prefilter cartridge has been changed.

The audio alarm can be stopped by switching the **Enter** button.

After change of prefilter cartridge, this reminder function will reset automatically.

Note!

The pretreatment cartridge installation is described with the "Initial Operation" chapter.

To change the prefilter cartridge:

1. Press **Menu** to activate the Main Menu.

```
MENU      Maintenance
           Setup
           Info
```

2. Select <Maintenance> and press **Enter** to confirm. The display will show the following:

```
Maintenance
           Repl. RO modules
           Repl. Prefilter
           Sanitization
```

3. Select <Repl. Prefilter> and press **Enter** to confirm.

```
Repl. Prefilter.
Continue?      Yes
               No
```

Note!

The program can no longer be interrupted, once the safety prompt has been confirmed.

4. Choose <Yes> to change prefilter cartridge and press **Enter** to confirm. Follow the instructions that appear in the display and push **Enter** after each step is completed.

Note!

After replacement check to make sure that all tubes are securely installed and leak tight.

After all work is completed, the unit is operational again.

Prefilter cartridge

RO Module Replacement

Every reverse osmosis system is supplied with potable water. In addition to various salts that can block the RO modules, potable water also contains different microorganisms. Therefore the RO modules should be changed at regular intervals.

To remind you to perform cleaning, the message "Repl. RO Mod." will appear on the display according to the factory setting.

This message will blink on your display until the RO modules have been changed.

The audio alarm can be stopped by switching the **Enter** button.

After change of RO modules, this reminder function will reset automatically.


Note!


The RO module installation is described with the "Initial Operation" chapter.

To change the RO modules:


1. Press **Menu** to activate the Main Menu.

2. Select <Maintenance> and press **Enter** to confirm. The display will show the following:

3. Select <Repl. RO Modules> and press **Enter** to confirm.

RO modules

Note!

The program can no longer be interrupted, once the safety prompt has been confirmed.

4. Choose <Yes> to change RO modules and press **Enter** to confirm. Follow the instructions that appear in the display and push **Enter** after each step is completed.

Note!

After replacement check to make sure that all tubes are securely installed and leak tight.

After all work is completed, the unit is operational again.

Fuse Replacement

arium® 61316 has two fuses for the unit. These are located next to the central On | Off switch on the left side of the unit in a fuse slot. This can be pulled out of the holder. The fuse is then visible.


Warning!

Severe electrical shock hazard or danger of electrocution!

- Only qualified personnel may service and repair arium® 61316.
- Disconnect arium® 61316 from the AC power outlet before fuse replacement.
- Make sure that all tubing connections are tight to avoid leakage.


1. Turn off arium® 61316 using the main power switch ("O") and disconnect it from the AC power outlet. The fuse drawer with two safety fuses is located on the left side of the housing above the main power switch.
2. Gently press the locking device of the fuse drawer upwards and pull out the fuse drawer (see figure below).

3. Remove the old fuses and replace with new fuses of the same type and rating.
4. Place fuse drawer back into the housing until it clicks into place.
5. Reconnect the unit to the AC power outlet and turn the unit back on using the main power switch ("I").
6. The unit is now ready to operate again.

Appendix

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials

that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The European legislation requires that electrical and electronic equipment be collected and disposed of separately from other communal waste with the aim of recycling it.

In Germany and many other countries, Sartorius Stedim Biotech takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our Service technicians on location or our Service Center in Goettingen, Germany:

Sartorius AG
Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and | or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet | manual or visit our Internet website (www.sartorius.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

Troubleshooting

Immediately after you discover any malfunction, turn off your arium® system at the main power switch. Then use the troubleshooting chart below to try to find out what the problem is and take the suggested action.

Problem	Causes	Solution
Pump is not running	Tank is full No power Feed water pressure too low AC inlet fuse blown Pretreatment cartridge is blocked	No action Connect power cord Check feed pressure, minimum 1 bar Replace fuses Change pretreatment cartridge
Retention rate is poor e.g. 85 %	RO modules are old, blocked	Replace modules
System leaks	Check tubing connections	Change or firmly tighten fitting
Low product flow rate	Module blocked Pretreatment cartridge blocked Feed water temperature low	Change RO modules Change pretreatment cartridge Contact Sartorius Stedim Biotech for assistance
Pretreatment cartridge life is shorter than 3 months	Inlet water is highly contaminated	Check inlet water characteristics May need external filter Contact Sartorius Stedim Biotech for assistance
Error message "Feed pressure low"	Too low inlet pressure Pretreatment cartridge blocked	Check feed pressure, minimum 1 bar Replace pretreatment cartridge
Poor water quality in the permeate	RO Module is blocked Inlet water quality has changed	Replace RO module

System Specifications

Unit Measurements:	43 cm (17") width, 48 cm (19") height and 33.4 cm (13.5") depth
--------------------	---

Note

38 cm (15") in front of the unit must be available to open the doors and 10 cm (4") on each side for connections.

System Weight:	14 kg (31 lb)
Operating Weight:	20 kg (44 lb)
Inlet Water Requirement:	Only use potable water as feed water according to national laws and regulations (Max. total hardness: 360 ppm as CaCO ₃ , max conductivity 1,500 µS/cm) In the case of hard water, it makes sense to soften the feed water Iron maximum: 0.1 ppm Colloid index or SDI: <3 Temperature range: 2°C - 35°C Pressure: min. 1 bar (15 psi), max 6.8 bar (100 psi)

Technical Data:

Tolerance limit: product water after RO modules	Adjustable from 10–80 µS/cm or 0.013–0.1 MΩ × cm
Programmable times printer in the print mode	From 1...59 min
Rejection rate	From 85 to 99%
Display conductivity resistivity value of product water	µS/cm or corresponding MΩ × cm compensated to 25°C

Operating Specifications:	Monovalent ion rejection rate: up to 98% Polyvalent ions rejection rate: up to 99% Microorganism rejection rate: >99% Particle rejection rate: >99% TOC: <100 ppb or >99% for MW >300
---------------------------	---

Flow Rate Product Permeate:	up to 16 l/h @ 25°C
-------------------------------	---------------------

RO Reject Concentrate:	up to 12 l/h
--------------------------	--------------

Recovery RO, approximately:	55% (Depends on feed water quality and temperature)
-----------------------------	---

Note

Permeate flow will change with water temperature. Lower temperature = less permeate.

Electrical:

Voltage	100–240 Vac
---------	-------------

Frequency	47–63 Hz single phase
-----------	-----------------------

Power Consumption	74 VA max
-------------------	-----------

Fuses	2 (neutral conductor phase), 4 AT, 5×20 mm
-------	--

Ambient Operating Temperature:	5–30°C, 80% relative humidity, non-condensing
--------------------------------	---

Ambient Storage:	Temperature 5–45°C, 80% relative humidity, non-condensing
------------------	---

Customer Service

- Sartorius Stedim Biotech
– Germany +49.551.3080
– USA +1.800.3687178
– UK +44.1372.737100
– France +33.1.69192100
– Italy +39.055.505671
– Spain +34.91.3588566

or your local Sartorius Stedim
Biotech dealer.

System, Consumables and Accessories

Part No.	Description
61316	RO housing incorporating electronics, pump and tubing
613CPF05-----V	Two pretreatment cartridges
613CPM4-----V	Two RO modules
613APV31	30-liter pressure vessel storage tank
613APV50	50-liter pressure vessel storage tank
613APV70	70-liter pressure vessel storage tank
613APV100	100-liter pressure vessel storage tank
611CDS2	Storage tank cleaning sanitization kit (two 30 ml syringes)
612CDS2	RO unit cleaning sanitization kit (two 30 ml syringes)
613AMDG1	Dispense gun to deliver water directly from the RO pressure tank
611APR1	Printer for data documentation
613AKD1	Connecting cable for external display (2.8 m long)
610AMFBI	Multi-function box, including level sensors, for the connection of an open tank system

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Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission
of Sartorius Stedim Biotech
GmbH.

The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make chan-
ges to the technology, features,
specifications, and design of the
equipment without notice.

Status:
December 2008,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany