Manual de instrucciones / Operating instructions Refractómetros de mano universales / Universal hand refractometers

01 Descripción / Description

- | | |
|--|---------------------------------------|
| 1. Entrada luz | 1. Entrance light prism |
| 2. Selector escala | 2. Scale selector |
| 3. Anillo polarizador | 3. Polarizer ring |
| 4. Cubre prisma | 4. Cover plate |
| 5. Anillo de ajuste dióptrico / Ocular | 5. Diopter adjustment ring / Eyepiece |

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1. Aim the front end of the refractometer to the direction of bright light, and adjust the diopter adjustment ring (5) until the reticle can be seen clearly.
2. Polarizer ring. While measuring a light concentration solution, you will find chromatism in view, you could reduce it by rotating (3).
3. Light Prism. If the light is not enough, you could rotate and adjust (1).
4. Adjust to zero. It is not necessary, this instrument has been adjusted before in factory.
5. Measurement. Open the cover plate (4), clean the surface of prism by soft cotton cloth, and drop 1-2 drops of solution to be measured. Close the cover plate, press it lightly, the read the corresponding scale of light/blue boundary, the reading is refractive index of measured solution.
6. Rotate. There are three scales in the instrument; you could get them by rotating (2).
7. Store. After measurement, clean away the adhered on the surface of prism and cover plate by moist gauze.

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1. After usage, don't use water to wash the instrument, to avoid damages in the instrument.
2. As it a precision optical instrument, you should hand it gently and take good care of it. Don't touch and scratch the optical surfaces. It should be kept in a dry clean and non-corrosiveness environment, to prevent the surface of the humidity. Please avoid strong shock during transportation.
3. The reference of temperature in 20°C. In operation, the temperature compensation should be made according to the table.

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