Samarus™

PEEK ANATOMICAL
CERVICAL CAGE

1988-2008

20

YEARS OF INNOVATION

SCIENT'X

YOUR EXPERTISE

OUR SUPPORT

The gu of an immediate

Features and Benefits of PEEK-OPTIMA®

PEEK-OPTIMA® by Invibio™ is a pure PEEK material designed to meet spine, hip and dental implants requirements.

➤ Safety

- Elasticity prevents sinking:
less risk of endplate penetration
- Rupture modulus is lower than metal:
less risk of bone damage
- Withstands repeated autoclave sterilizations


➤ Medical imaging (X-ray, CT, MRI)

- Radiolucent material to allow accurate fusion follow-up
- Embedded gold markers to facilitate implant placement verification

➤ Biocompatibility

- No carbon fibres: no risk of inflammatory reaction
- CE-marked and cleared by FDA for long-term implantation in the human body

Young's Modulus (GPa)

➤ Enhanced bone fusion

- No stress shielding
- Optimum load repartition
- Bone growth enhancement thanks to micromovements

Guarantee stability

Interbody arthodesis is a secure and efficient technique to treat cervical spine pathologies.

The Samarys™ anterior cage answers all the requirements of this technique. Its unique design allows a maximum interbody fusion.

The Samarys™ cage anatomical design offers an optimum stability and its wide range of sizes allows a perfect adaptation to each patient's anatomy. The particular shape and the pure PEEK composition of this implant ensure an excellent bone fusion.

Increased stability and fusion

- ➔ **Original design** for an optimum primary stability to restore the interbody space and preserve cervical lordosis
- ➔ **Large graft surface** to enhance bone fusion

Anatomical shape

- ➔ **Designed** to maximise contact with the intervertebral endplates

Reinforced security

- ➔ Inferior and superior teeth on surface to **avoid risks of migration**

1 Surgical protocol / Patient positioning

The procedure is performed under general anesthesia with endotracheal intubation. The patient is positioned in supine position on a standard table.


The head is placed in a neutral position or rotated in the opposite direction to the approach. Rotation should be approximately 30° to limit stress on the sterno-cleido-mastoid muscle.

The shoulders are lowered and fixed with self-adhesive strips in order to release the inferior cervical spine.


The skin incision targets the level of the lesion to be treated. It is marked on the skin after locating the lesion with an image intensifier. It may be horizontal, along a skin fold, extended to the opposite side of the median line in the case of bilateral approach, or beyond according to the number of levels to be treated.

The platysma muscle is sectioned on the same plane as the skin, then two upper and lower detachments are created to avoid muscle tension. The aponeurotic planes are sectioned with the finger up to the anterior surface of the spine.


The muscles along the neck are separated and rugined to thoroughly release the operating area and allow proper retraction using an autostatic distractor.

Under scopic control the proper level is located, then the Caspar distractor pins are positioned in the middle of the vertebral bodies over- and under-lying the discal space to be treated and parallel to the endplates. Then the anterior vertebral ligament is incised, down to the lateral unco-vertebral articular processes, before distracting the space. Discectomy is then carefully performed, if necessary with the help of a microscope.


2 Discectomy and preparation of the surgical site

The discectomy starts with the rectangular excision of the anterior part of the annulus fibrosus, using a scalpel. The disc material is resected using curettes and rongeurs. It may prove necessary to use a microscope for resecting the posterior disc material.

After full discectomy, resecting the posterior osteophytes and the caudal part of the uncus allows foramen release. For simple disc herniations, it may prove sufficient to open the ligament toward the root and remove the herniation.

This essential surgical stage gives direct access to the anterior side of the medulla and nerve roots.

Cervical distractor
21DST02


Distraction pin Ø 3.5
21FCD35-xx

Screwdriver for distraction pin
21TRV03

3 Preparing the vertebral plates

The endplate preparation is performed with a curette. This step must be sufficiently thorough in order to extract all discal debris but not so pronounced as to collapse the cancellous bone, which would lead to secondary graft collapse. The intersomatic height may be restored in different ways : autologous, heterologous graft or interbody cages.


Plate positioning is performed after resecting the anterior osteophytes. This step is necessary to obtain a close contact between the plate and the anterior vertebral body surface. The choice of the plate and its correct positioning are essential steps of the osteosynthesis.

4 Placing the trials implant

Trial implant
21CC2Axx-xx

A set of implant trials is used to determine the ideal implant size. The width of the trials is 17 or 20 mm and their height is from 5 to 7 mm. These different sizes allow to find the implant that will fit to the anatomies and the needs of the patient. The trial is secured to the holder and is impacted into the intervertebral space. Having no retaining edges, the trial can easily be removed.

5 Graft harvesting and cage filling

Bone collection device
20TRE01



Socket for Samarys cage
21SOC02

Graft pusher
21COM03

The cage can be filled with a bone substitute or an autologous graft. In cases of autogeneous bone grafting, the graft harvesting is performed by way of a 15 mm incision on the anterior iliac crest. A specific trefine can be adapted by impact on the iliac crest and allows to obtain a transcortical harvesting on different angles. Generally, one or two samples are enough to fill the cage. The closing of the iliac bone flap is made after haemostasis with wax.



The cage is positioned on its socket during the filling phase. The upper and lower faces of the cage are widely opened to offer optimum graft volume and contact area. A graft pusher matching the cage shape is used to compact the graft.

6 Installing the cage

After filling, the cage is secured to the holder tool. The holder is then positioned onto the cage.

The end piece of the holder is inserted on the front side of the implant and ensures full locking. The holder has small dimensions and locks in continuation of the cage. The cage can be impacted or installed with slight distraction on the distractor. The stop secured at the top end of the holder allows impacting the implant. The holder is removed after implanting the cage.

*** Note:** The cervical cage ensures immediate post-operative stability of the spine but, in some cases, according to the physician's advice, internal fixation with an anterior cervical osteosynthesis plate, or external collar or brace, can be added according to the surgeon's habits.

7 Closing the approach wound

The approach wound is closed after rinsing. The haemostasis is checked and a suction drain is placed on the anterior face of the spine. The platysma muscle must be carefully restored, and the skin closed by intradermic stitches with resorbable thread.

8 Removing the implant

If the instrumentation must be removed, before fusion, the cervical approach is used, down to the instrumented area. Bone bridges between the implant and the vertebral bodies must be sectioned. The holder is placed on the cage to facilitate implant removal.

Cage holder
21CC2Axx

Instructions for use

The Samarys anterior cervical cage is an implant intended for the surgical treatment of the cervical vertebrae.

GENERAL DESCRIPTION:

The Samarys anterior cervical cage has been designed for stabilisation and arthrodesis between the cervical vertebrae. The Samarys anterior cervical cage is convex in shape on its upper face and on the openings on the upper and anterior faces. It comes in several sizes to adapt to the variable morphology of the inter-vertebral spaces. The Samarys anterior cervical cage is made from PEEK. The Samarys anterior cervical cage must not be used with components coming from other manufacturers.

INDICATIONS:

The Samarys anterior cervical cage is a surgical implant for cervical arthrodesis by anterior route, designed to optimize bone fusion. It is an intersomatic cage whose main function is radicular decompression by intersomatic distraction and a design favouring osteogenesis. The indications include:

- stabilisation after treatment of a herniated cervical intervertebral disc or osteophytosis compressing the nerve roots or spinal cord,
- degenerative intervertebral instability.

CONTRAINDICATIONS :

The contra-indications for the Samarys anterior cervical cage include:

- local infection or inflammation,
 - vertebral osteoporosis,
 - malignant vertebral disease,
 - allergy or intolerance to PEEK,
 - incompatible patient age and physical condition,
 - all cases not included in the indications.
- The Samarys anterior cervical cage has not been designed, intended or sold for uses other than those indicated.

POSSIBLE ADVERSE EFFECTS:

- infection,
- post-operative migration of the implant ,
- damage to the vertebrae adjacent to the arthrodesis,
- intolerance to the material.

> **Note:**

An additional surgical operation may be needed to correct any adverse event.

> **Warnings:**

An entirely satisfactory result is not always obtained at each and every operation. This is particularly true in spinal cord surgery where many external factors can compromise the results.

OPERATING PRECAUTIONS:

The surgeon must be thoroughly familiar with the Samarys anterior cervical cage, the method of application, the instruments and the operating technique (similar to that of the CBK cage). The size of the Samarys anterior cervical cage must be chosen in relation to the individual clinical case and the desired correction. It is recommended to use a Caspar type distractor to open the intervertebral space intended for implanting the cage. Before implanting the Samarys anterior cervical cage, a careful curetage of vertebral faces must be conducted to expose healthy tissue but without making them fragile to prevent the cage from slipping inwards. The Samarys anterior cervical cage must be filled with autologous or allogenic bone or with a bone substitute to obtain bone fusion. The correct positioning of the cage with respect to the vertebrae is confirmed by radiology. After implantation, the **batch number** and the **reference of the implanted Samarys anterior cervical cage** must always be recorded in the patient's surgical records. This can be done using one of the labels attached to the sterile packaging.

The implants are single-use. A removed implant must never be implanted again.

> **Warnings:**

No silicone, formaldehyde or latex based products should be implanted with a SCIENT'X implant.

PACKAGING:

Implants delivered sterilized are conditioned in individual sterile packages. Implants delivered non-sterile and their instruments are delivered in one or more containers, one (or more) "kits", except in the case of non-assorted kits for stocks. Each kit must be closed and sealed. The packaging must be intact upon reception. If the stocked kit system is used, complete kit composition must carefully be checked. The proper state of all implants and instruments must be checked prior to any use. Storage conditions must enable implant, accessories and packaging integrity to be maintained. Damaged products must never be used and must be returned to SCIENT'X. products must not be used but must be returned to SCIENT'X.

DECONTAMINATION/CLEANING/STERILISATION:

Products delivered sachets are not sterile.

> **For the implants delivered sterile**

The implants are sterilized by Gamma radiation at doses of 25 to 40 kGy. The expiry date is 5 years. The expiry date of sterile parts is indicated on the packaging.

> **For implants and instruments not sterile on delivery:**

All implants and instruments delivered non-sterile must be **decontaminated, cleaned and sterilised** before and after use. Implants and instruments in sachets must be removed from the original packaging for the following operations:

Recommended method:

- **Decontamination:** Plunge the implants and instruments into a bactericidal and fungicidal solution of the didecyldimethylammonium chloride type diluted to 0.5 % (5mL to 1 litre water). Length of soaking: 20 min. Rinse with demineralized water.
- **Cleaning:** Wash the implants and instruments in a LANCER type machine with suitable cleaning products, rinse and dry. Any product that might damage the equipment is forbidden (such as bleach, formol, etc.). It is essential to sterilise the kit by steam using the following conditions:
- **Sterilization:** It is essential to sterilise the kit by steam using the following conditions:
 - pre-heating 25 min at 110 °C (1 bar)
 - vacuum 5 min (0.8 bar under atmospheric pressure)
 - heating 5 min at 120 °C (1 bar)
 - vacuum 5 min (0.8 bar)
 - sterilization 18 min at 134 °C (2 bar)
 - drying 20 min return to room temperature

COMPLAINTS:

Any complaints, together with the reference and lot number of the incriminated product, should be sent to Scient'x.

Date of revision of this package insert: 2006/07/24

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EVALUATION OF AN INTERVERTEBRAL CERVICAL CAGE FILLED WITH INJECTABLE BONE SUBSTITUTE

Authors: H. Pascal-Moussellard⁽¹⁾, E. Aguado⁽²⁾, Y. Catonné⁽³⁾, J.L. Rouvillain⁽¹⁾, G. Daculsi⁽⁴⁾

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Introduction

Anterior cervical discectomy and arthrodesis with use of autogenous iliac-crest bone graft (Robinson and Cloward procedures) is the gold standard procedure to treat degenerative changes of the cervical spine resulting in neck pain and/or cervical myelopathy and radiculopathy. The main problem of this procedure remains the morbidity associated to the iliac graft harvest.

The ideal system should respond to the following specifications:

- good primary stability
- disc height conservation
- biocompatible material able to provide bony fusion with respect to both conditions: - fusion rate = fusion obtained with autologous bone graft
- morbidity rate < morbidity rate of the gold standard procedure

We intend to evaluate an arthrodesis system combining an impacted intervertebral PEEK cage and an injectable bone substitute. Thanks to its design and its mechanical strength, the cage is a self-stable implant (primary stability) and restores the local lordosis as well as the discal height. The use of a bone substitute avoids the iliac bone harvest and then the morbidity associated to this procedure. The parameter we will focus on, is therefore the fusion rate.

Objective

Our objective is to evaluate the fusion rate of the arthrodesis obtained with a PEEK intervertebral cage associated to autologous bone graft or injectable bone substitute.

Material & Methods

Implant description

• **Samaris cervical cage**
(Scient'x - Guyancourt - France):

- made of PEEK (Poly Ether Ether Keton)
- radiolucent
- lordotic
- dimensions: L13 x w17 x h5-6 mm

• **Injectable Bone Substitute IBS, MBCP gel**
(Biomatlante - Vigneux de Bretagne - France)

The MBCP gel is made of 2 phases:

- a mineral phase = granules of biphasic calcium phosphate ceramic [60% hydroxyapatite - 40% β -tricalcium phosphate]

- a polymeric phase = cellulose derivative polymer HPMC (hydroxy-propyl-methyl-cellulose), vehicle for the calcium-phosphate mineral phase

Samaris cervical cage

Animal model

The goat has been chosen because of its validation as a model for the cervical spine.

Fourteen animals underwent discectomy C3-C4 + cage insertion + osteosynthesis.

The 14 animals were divided into 3 groups:

- EMPTY CAGE n = 3
- CAGE + SPONGEOUS BONE HARVESTING FROM THE ILIAC CREST n = 6
- CAGE + BONE SUBSTITUTE IBS n = 5

Goats were killed at 6 months postoperative and the C3-C4 levels were explanted.

Fusion assessment methods:

The instrumented levels were fixed with methyl metacrylate and prepared for undecalcified sectioning:

- Histological examination: polarized light microscopy + Movat's staining
- Scanning Electronic microscopy (SEM)

We considered interbody fusion complete only if there was bone bridging the disc space continuously from endplate to endplate.

Results

				SCANNING ELECTRONIC MICROSCOPY (SEM) JEOL JSM 6300 (Tokyo, Japon)					
		Polarized Light	Stained sections (Movat Pentachrome)	Qualitative analysis	Quantitative analysis				
EMPTY CAGE	FUSION 0/3 PSEUDARTHROSIS 3/3	 Fibrous tissue	 (X 10) Pseudarthrosis	 Non mineralized cartilage zone between the endplates	26,06 %	26,06 %	0		
CAGE + BONE	FUSION 4/5 PSEUDARTHROSIS 1/5	 Trabecular bone	 (X 20) Arthrodesis	 Trabecular bone oriented perpendicular to the endplates	33,8 %	33,8 %	0		
CAGE + Bone substitute	FUSION 3/6 PSEUDARTHROSIS 0/6	 Nearly fused bone. Decrease of the IBS quantity in fibrosis	 (X 10) Arthrodesis	 (X 15) Fusion in progress	 Arthrodesis: 3/6	 Mineralization in progress: 3/6	32,6 %	33,22 %	0,2
					Newly formed bone	Mineralized tissue	Cartilage / Mineralized tissue		

Discussion

Many human clinical studies showed that fusion can be achieved in the cervical spine without any graft material. We tested thus the implantation of an empty cage. The systematic absence of fusion emphasizes the **need for a graft material into the cage**.

Average fusion rate of cervical spine arthrodesis using autologous bone graft with animal models is 46,8% (range: 38% - 62%) [1] [2] [5] [6] [7] [8] [10] [11] [12]. In our study, we observed a fusion rate of 80% for cages associated to autologous bone. **This very high rate of fusion validates the use of the cage with autologous bone for cervical arthrodesis.**

The results of the quantitative SEM analysis for the bone substitute IBS are very closed to those obtained with autologous graft. The association of BCP granules to the HPMC vehicle phase nearly gives IBS the structure of a totally interconnected macroporous ceramic [3] [9]. This kind of structure helps the dynamic of the ceramic resorption - bone substitution process and allows an early newly formed bone appearance [4]. The use of this biomaterial as a bone substitute is very promising.

Conclusion

This histological study validated the use of the impacted Scient'x cage for cervical arthrodesis. The fusion rate obtained when cage-autologous bone are associated is the highest fusion rate reported in the literature.

The use of the IBS bone substitute is very promising. Research is still in progress to improve the conditions of use of this material, particularly to control the filling of the cage volume.

Recent advances in bone substitute and intervertebral cages give the opportunity to **considerably reduce the morbidity rate at the donor site.**

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

➔ Please read carefully the instructions for use bulletin.
Devices may be subject to modification. Patented

88BTN11CC2A-A / Version 2.0 July 2006