OSTEOPROTEGERIN

ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF OSTEOPROTEGERIN IN
EDTA PLASMA, HEPARIN PLASMA, SERUM OR CELL CULTURE SUPERNATANTS
CAT. NO. BI-20402 . 12 X 8 TESTS

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www.bmgrp.com

1) INTRODUCTION

Osteoprotegerin (OPG) or Osteoclast inhibitory factor (OCIF) is a dimeric glycoprotein of the TNF receptor family with a molecular weight of 60 kD and 120 kD respectively. OPG acts as a soluble secreted receptor for RANKL that prevents it from binding to and activating RANK (osteoclast differentiation and activation receptor, ODAR) on the osteoclast surface, thus inhibiting the osteoclast development.

POSSIBLE INDICATIONS

- Postmenopausal and senile Osteoporosis
- Glucocorticoid-induced Osteoporosis
- Diseases with locally increased resorption activity
- Therapy monitoring after treatment with OPG
- Arthritis
- Oncology

2.) CONTENT OF THE KIT

CONT	KIT COMPONENTS	QUANTITY
PLATE	Monoclonal Anti-OPG Antibody, pre-coated microtiter strips in a strip holder	12 x 8 tests
WASHBUF	Wash buffer, 20x concentrated	1 x 50 ml
AB	Polyclonal biotinylated anti OPG antibody (yellow cap), ready to use	1 x 7 ml
STD	Standards, (0; 0.37; 1.1; 3.3; 10; 30 pmol/l), ready to use	6 x 500 µl
CTRL	Control (yellow cap), ready to use	1 x 500 µl
ASYBUF	Assay Buffer, ready to use	1 x 25 ml
CONJ	Conjugate, (streptavidin- HRPO), ready to use	1 x 22 ml
SUB	Substrate (TMB solution), ready to use	1 x 22 ml
STOP	Stop solution, ready to use	1 x 7 ml
OPG-STOCK	Osteoprotegerin stock solution, red cap, (500 pmol/l)	1 x 500 µl

3) ADDITIONAL MATERIAL ADDED TO THE KIT

- 2 self-adhesive plastic films
- QC protocol
- Protocol sheet
- Instruction manual for use

4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- 1,5 ml reaction vials
- Precision pipettes calibrated to deliver 10-1000 µl and disposable tips
- Elisa reader for absorbance at 450 nm (or from 450 nm to 620 nm)
- Plate washer is recommended for washing
- Distilled or deionised water

5) REAGENTS AND SAMPLE PREPARATION

Samples should be stored at -20°C, for long-term storage store at -70°C. All samples should undergo only 4 freeze-thaw cycles. Lipemic or haemolysed samples may give erroneous results. Samples should be mixed well before assaying.

Reconstitute as follows:

WASHBUF (Wash buffer): Dilute the concentrate 1:20: e.g. 50 ml WASHBUF + 950 ml distilled water. Crystals in the buffer concentrate will dissolve at room temperature. Buffer is stable at 2-8°C until expiry date stated on label.

Use only diluted WASHBUF (Wash buffer) for optimum assay performance:

Cell Culture:

Dilute the OPG-STOCK (Osteoprotegerin stock solution) with Cell Culture Medium 1:20 (1+19) for STD1 (50 µl OPG-STOCK + 950 µl Medium).

Dilute the Standard for the Standard curve 1:2 (1+1) as it is shown in the scheme:

STD1 (25 pmol/l)


250 µl STD1 + 250 µl dilution medium = STD2 (12.5 pmol/l)

250 µl STD2 + 250 µl dilution medium = STD3 (6.25 pmol/l)

250 µl STD3 + 250 µl dilution medium = STD4 (3.13 pmol/l)

Use Cell Culture Medium as STD0 (Standard 0)

6) PRINCIPLE OF THE ASSAY

7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-26°C) before use in the assay
Mark position for BLANK/STD/SAMPLE/CTRL (Blank/Standard/Sample/Control) on the protocol sheet
Take microtiter strips out of the alu bag, take a minimum of one well as Blank. Store unused strips with desiccant at 2-8° C in the alu bag. Strips are stable until expiry date stated on the label.
Add 100 µl ASYBUF (Assay Buffer) into each well. Pipette additional 100 µl into well marked as blank.
Add 50 µl STD/SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective well, except blank
Add 50 µl AB (biotinylated anti OPG antibody) into each well, except blank, swirl gently
Cover tightly and incubate over night (18-24h) at 4°C
Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting plate against paper towel after the latest wash
Add 200 µl CONJ (Conjugate) into each well
Cover tightly and incubate at room temperature (18-26°C) for 1 hour
Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting plate against paper towel after the latest wash
Add 200 µl SUB (Substrate) into each well
Incubate for 20 min at room temperature (18-26°C) in the dark
Add 50 µl STOP (Stop solution) into each well
Measure absorbance immediately at 450 nm with reference 620 nm, if available

8) CALCULATION OF RESULTS

Subtract the blank extinction from all other values. Construct the standard curve from the standard values. Use commercially available software or graph paper. Obtain sample concentration from this standard curve. The assay was evaluated with 4PL algorithm. Different curve fitting methods need to be evaluated by the user. Respective dilution factors have to be considered.

9) ASSAY CHARACTERISTICS

Normal range:	Median = 1.8 pmol/l (n = 1134; age 19-96) It is recommended to establish the normal range for each laboratory.
Standard range:	0 to 30 pmol/l
Conversion factor pg/ml to pmol/l:	1 pg/ml = 0.05 pmol/l (MW: 19.9 kD)
Sample volume:	50 µl human EDTA plasma, Heparin plasma, serum, cell culture supernatant
Detection Limit:	(0 pmol/l + 3 SD): 0.14 pmol/l
Incubation time:	18 – 24 h / 1h / 20 min

10) PRECISION

Intra-Assay: 2 samples of known concentrations were tested 16 times to assess intra-assay precision

Inter-Assay: 2 samples of known concentrations were tested in 10 Assays to assess inter-assay precision

Intra-Assay (n=16)		
Mean (pmol/l)	4.59	10.76
SD	0.46	0.41
CV%	10%	4%

Inter-Assay (n=10)		
Mean (pmol/l)	5.53	10.1
SD	0.38	0.76
CV%	7%	8%

11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date. Protect reagents from direct sunlight.
- Substrate solution should remain colourless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

12) PRECAUTIONS

All test components of human source were tested with 3rd generation tests against HIV-Ab and HBsAG; and were found negative. Nevertheless, they should be handled and disposed as if they were infectious.

All liquid reagents contain 0.01% Thimerosal or 0.01%Proclin 300 as preservative.

Thimerosal is toxic! Avoid contact with skin and mucous membrane. Proclin 300 is not toxic in concentrations used in this kit. It may cause allergic skin reactions-avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Avoid all contact with the reagents by using gloves.
- Sulfuric acid is irritating to eyes and skin. Flush with water if contact occurs. Avoid contact with skin and mucous. Irritations are possible - Flush with water after contact!!

13) LITERATURE

- "Clinical Implications of the Osteoprotegerin/RANKL/RANK System for Bone and Vascular Diseases"
Lorenz C. Hofbauer, MD; Michael Schoppet, MD; JAMA, July 28, 2004—Vol 292, No. 4
- "Comparison of 10 serum bone turnover markers in prostate carcinoma patients with bone metastatic spread: diagnostic and prognostic implications".
Jung K; Lein M; Stephan C; Von Hosslin K; Semjonow A; Sinha P; Loening SA; Schnorr D Department of Urology, University Hospital Charite, Humboldt University of Berlin, Berlin, Germany. Int J Cancer 2004 Sep 20;111(5):783-91
- "Osteoprotegerin is a Risk Factor for Progressive Atherosclerosis and Cardiovascular Disease"
Stefan Kiechl, MD; Georg Schett, MD; Gregor Wenning, MD; Kurt Redlich, MD; Martin Oberhollenzer, MD; Agnes Mayr, MD; Peter Santer, MD; Josef Smolen, MD; Werner Poewe, MD; Johann Willeit, MD
Circulation. 2004;109:2175-2180.
- "Endurance running acutely raises plasma osteoprotegerin and lowers plasma receptor activator of nuclear factor kappa B ligand" [In Process Citation];
Ziegler S; Niessner A; Richter B; Wirth S; Billensteiner E; Woloszczuk W; Slany J; Geyer G
Metabolism 2005 Jul;54(7):935-8

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BI-20402 OSTEOPROTEGERIN

ASSAY PROTOCOL AND CHECKLIST

- Bring all reagents to room temperature (18 – 26°C)
- Prepare reagents and samples as instructed
- Bring unused and prepared components to the storage temperature mentioned in the package insert
- Take microtiter strips out of the alu bag and mark positions on the protocol sheet
- Add 100 µl ASYBUF (Assay buffer) into each well. Pipette additional 100 µl into well marked as blank.
- Add 50 µl STD/SAMPLE/CTRL (standard/sample/control) into all wells except blank,
- Add 50 µl AB (biotinylated anti OPG ab) into each well, except blank, swirl gently
- Cover tightly and incubate at 4°C over night (18-24 hours)**
- Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- Add 200 µl CONJ (Conjugate) into each well
- Cover tightly and incubate at room temperature (18-26°C) for 1 hour**
- Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- Add 200 µl SUB (Substrate) into each well
- Incubate at room temperature (18-26°C) for 20 minutes, in the dark**
- Add 50 µl STOP (Stop solution) into each well
- Read Optical Density at 450 nm with reference 620 nm, if available