# Urisys 1100

cobas®

Gebrauchsanweisung  
Operator's Manual  
Manual de instrucciones

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The packaging material, the identification  
plate on the instrument and the manual  
may contain the following symbols or  
abbreviations:

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Please consult instructions for use

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Caution (refer to accompanying docu-  
ments). Please refer to safety-related  
notes in the manual accompanying this  
instrument.

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This product fulfils the requirements of  
Directive 98/79/EC on in vitro diagnostic  
medical devices.

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# Table of Contents

	Page
<b>1. Introduction</b>	30
<b>2. System Description</b>	31
2.1 Measuring Principle	31
2.2 Components and Functions	32
<b>3. Software</b>	32
3.1 Overview	32
3.2 Menu Structure (Flowchart)	33
3.3 Menu Functions	35
3.4 Results Table	37
3.5 Changing the Range Limits	38
<b>4. Installation</b>	38
4.1 Modification for Combur <sup>z</sup> Test and Combur <sup>s</sup> Test Test Strips	40
4.2 Operator ID	40
4.3 Authentication	41
<b>5. Calibration</b>	42
<b>6. Reading Test Strips</b>	43
6.1 Overview	43
6.2 Normal Mode (for Single Readings)	43
6.3 Patient Report	44
6.4 Fast Mode (for Serial Readings)	44
6.5 Strip Measurement Error	45
6.6 Entering Patient ID, Operator ID and Authentication Password	45
6.7 Data Transmission to a PC or Host Computer	46
<b>7. Cleaning and Maintenance</b>	46
7.1 Cleaning the Test Strip Tray	46
<b>8. Error Messages and Troubleshooting</b>	47
<b>9. Connecting to Other Devices</b>	50
9.1 Serial Interface	50
9.2 Barcode Reader, AT/PC Keyboard	50
<b>10. Technical Information and Notices</b>	51
10.1 Technical Data	51
10.2 Safety Notices	52
10.3 Guarantee	52
<b>11. Ordering Information</b>	53
<b>12. Alphabetical Index</b>	54
<b>13. Warranty Card</b>	83
<b>14. International Addresses</b>	85

# 1. Introduction

## Revision History

Manual version	Software version	Revision date	Amendments
04	5.x	June 2008	Operator ID, limited lock-out function, Device ID, compatibility of barcode reader, ASTM protocol

The Urisys 1100 system is a reflectance photometer designed to read and evaluate the urine test strips Combur<sup>10</sup>Test UX, Combur<sup>7</sup>Test and Combur<sup>5</sup>Test\* from Roche Diagnostics. It reads the strips under standardized conditions, saves the results to memory and outputs them via its own inbuilt printer and/or serial interface.

The Urisys 1100 system is designed for In Vitro Diagnostic (IVD) use by qualified physicians and laboratory staff.

**BIOHAZARD: Treat all samples of human origin as being potentially infectious. Always observe good laboratory practice.**

Using the Urisys 1100 system eliminates factors known to affect visual evaluation of urine test strips, such as:

- Variable lighting conditions at the workplace
- People's individual skill at matching colors and their limits of concentration
- Different reaction times for the test strips
- Clerical errors
- Strong color of the urine sample

The following symbols are used throughout this document.

Symbols/Signs	Explanation
	<b>WARNING/CAUTION:</b> Indicates a potentially hazardous situation that if not avoided could result in personal injury or damage to the instrument. This symbol is also used to highlight situations that can compromise results.
	<b>BIOHAZARD:</b> Indicates a potentially dangerous situation involving the presence of biohazardous material. All safety precautions must be taken to prevent personal injury or damage to the equipment.
	<b>ATTENTION:</b> Indicates special problems or important information. Read the accompanying text carefully.


\* Combur<sup>7</sup>Test and Combur<sup>5</sup>Test urine test strips are not available in all countries.

## 2. System Description

### 2.1 Measuring Principle

The test strip is placed on a sliding tray, and a stepping motor moves it under the reading head, which remains stationary. The analyzer reads the reference pad, followed by each of the test pads on the strip.

The reading head contains LEDs that emit light at various wavelengths. Reading is done electro-optically, as follows:


The LED (1) emits light of a defined wavelength on to the surface of the test pad (2) at an optimum angle. The light hitting the test zone is reflected more or less intensely depending on the color produced on the test pad, and is picked up by the detector, a phototransistor (3) positioned directly above the test zone. The phototransistor sends an analogue electrical signal to an A/D converter (4), which changes it to digital form. The microprocessor (5) then converts this digital reading to a relative reflectance value by referring it to a calibration standard.

Finally, the system compares the reflectance value with the defined range limits (reflectance values that are programmed into the analyzer for each parameter) and outputs a semi-quantitative result (6).

Each test pad is read photometrically after a lead (incubation) time of about 55–65 seconds. In strongly alkaline urine samples, the Urisys 1100 system automatically corrects the result of the Specific Gravity test.

## 2.2 Components and Functions

### Component

1. Printer cover
2. Display/keypad
3. Test strip tray
4. START button
  
5. On/Off switch
6. Serial interface
7. Power socket
8. 5-pin DIN socket

### Function

- Flips up for insertion of printer paper
- LCD display and three function keys for menu-driven operation and interfacing with the user
- Holds and anchors the strip
- a) Starts the reading process  
b) Closes submenus and returns to the starting menu (Ready-to-Measure status)
- Powers the unit on and off
- For connection to a personal or host computer
- Socket used to connect the analyzer to the mains adapter
- For connecting a barcode reader or AT/PC keyboard

## 3. Software

### 3.1 Overview

The Urisys 1100 system's software provides a user interface that enables all laboratory-specific settings and recurrent functions to be selected via the liquid crystal display and function keys (see Sections 3.2 and 3.3).

The three function keys assume the particular function displayed on the second line of the liquid crystal display. The first line of the display is used for system status and user information.


The user interface is self-explanatory, so that only details of the major functions are presented here.

**Pressing the START button within any submenu accepts the chosen instrument state or function, and returns the system from menu display to Ready-to-Measure status.**

The instrument switches from the Ready-to-Measure status or the status displayed, respectively, to the Standby mode after 5 minutes of function key inactivity. During Standby, the date and time are displayed. Ready-to-Measure status can be resumed by pressing the START key, except when certain error messages are displayed (see Section 8).


### 3.2 Menu Structure (Flowchart)

Pressing the **START** button confirms the setting, closes the submenu, and resumes Ready-to-Measure status.

# Menu Structure (continued)

Pressing the **START** button confirms the setting, closes the submenu, and resumes Ready-to-Measure status.

### 3.3 Menu Functions

**Self Check:** During Self Check at power-on, the analyzer automatically checks that the program chip, the tray transport mechanism, the printer connection and optical system are all OK. The tray type is checked (see Sections 4 and 7.1) to ensure that it correctly matches the strip type selected in the menu.

**Calibration:** For requesting calibration with Control-Test M (see Section 5).

**Linefeed:** Causes paper to advance. The linefeed is stopped by pressing the left function key ("Stop") (see Section 4).

**New Series:** For starting a new series of measurements at sequence number 1. It is also possible to have an Automatic Patient ID (see Section 6.1).

**Memory:** The analyzer can store up to 100 results together with date and time of measurement, sequence number and patient ID (if entered). Memory is automatically cleared every time the date changes.

After memory has been erased, "NO RESULTS STORED" is displayed when "Memory" is pressed. Pressing the START button resumes Ready-to-Measure status. "MEMORY FULL" and the options "Print/Send/Clear" are displayed when the memory is full. Memory must be cleared before the analyzer can resume Ready-to-Measure status.

**Print Results:** For generating a printout of stored results. The options are as follows:

- **"All":** All results in memory
- **"Last Series":** The most recent series of readings
- **"Last One":** The most recent reading

Printing can be repeated as often as desired. Printing can be interrupted by pressing the left key ("Stop"), for example to allow a new roll of printer paper to be inserted ("Linefeed"), and subsequently resumed ("Continue"). The analyzer resumes Ready-to-Measure status when the "Home" key is pressed or when printing has finished.

**Send Results:** For sending stored results to the serial interface. Options are the same as for "Print Results". Results can be sent as often as required (see also Sections 6.7 and 9.1).

**Clear Memory:** Erases results from memory.

**Mode:** Choice of **Print**, **Measure** and **Interface** modes.

**Printer:** Printer options are:

- **"On":** The printer is switched on. Each result is printed once.
- **"2 Copies":** In Normal mode, each result is printed twice. Note: after serial readings (Fast mode) and when the printout is a repeat printout (activated by the "Print Results" function), each result is printed only once.
- **"Off":** For switching off the printer when the printout of results is only required at the end of a series of readings (activated by the "Print Results" function), or when results are to be sent via the interface to a personal computer or host only.

**Measurement:** Choice of measuring mode:

- **"NormalMode":** Read cycle 70 seconds, incubation within the analyzer (see Section 6.2)
- **"FastMode":** Read cycle 30 seconds, incubation outside the analyzer (see Section 6.4)

**Interface:** Choice of "unidirectional", "bidirectional" or "ASTM" data transfer. For further details, see Sections 6.7 and 9.

**Logdata:** Choice of sending the log files to the host PC or printing out the last 10 logged data. (Only by supervisors in Authentication mode)

**Strip Settings:** For selecting test strip type, units and range limit settings.

**Strip Type:** Choice of test strip type:

- **Combur<sup>10</sup>Test UX** ("C-10")
- **Combur<sup>7</sup>Test** ("C-7")
- **Combur<sup>5</sup>Test** ("C-5")

The Urisys 1100 system leaves our factory configured for Combur<sup>10</sup>Test UX test strips. If you intend to use Combur<sup>7</sup>Test or Combur<sup>5</sup>Test test strips, you will require the appropriate test strip tray (see Sections 4.1 and 11).

**Units:** Options are:

- **Conventional units** (mg/dL)
- **SI units** (mmol/L)
- **Arbitrary units** (1+, 2+, 3+, 4+)

Either single units or combined units (conventional/arbitrary or SI/arbitrary) can be selected.

The operator selects a unit in which the results are to be stored, printed and/or transferred to a computer. After a new unit has been selected, the repeat printout (activated by "Print Results") and all following printouts and/or data transfer (activated by "Send Results") will be in the newly chosen unit.

**Range Limits:** For changing the reflectance range limits. Changing the range limits for the result intervals alters evaluation sensitivity (see Sections 3.4 and 3.5).

**Language:** This enables the display language to be set. The "Other" option allows the operator to choose between English, German, Italian, Spanish and French.

**Device ID:** Displays the 5 digits device ID, which is part of the factory settings and cannot be changed. The device ID will also be sent to the host.

**Operator ID:** Choice between **Normal** and **Authentication** modes for the operator identification. If activated the input of an operator ID will be required by the start of the instrument.

The Operator ID will appear in the result printout and will also be sent to the host PC. Authentication mode offers a lock-out function and requires the ASTM protocol. (Refer Section 4.2 and 6.6 for details.)

**Parameter (Par.):** For selecting the order and number of parameters in the patient report and for "unidirectional" data transfer.

To change the order or number of parameters in the patient report, press the function key on the left ("New") with the display showing "SEQUENCE PARAMETERS". This calls up the first item in the patient report ("1. PARAM.:") followed by the first parameter ("SG" in the case of Combur<sup>10</sup>Test UX test strips). Press "Yes" to confirm the parameter in the position shown. Pressing "No" causes the next test strip parameter to be offered for the first position (display = "1. PARAM.: PH"; etc.).

As soon as "Yes" has been pressed to confirm a parameter as the first item, the second item is shown with "2. PARAM.:" followed by the next, not yet confirmed, test strip parameter. Press "Yes" or "No" to select the parameter as the second item in the patient report, and so on.

If there are any test strip parameters that are to be suppressed in the patient report, terminate selection following the most recently confirmed parameter by pressing "START".

For readings with Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips, the display only offers the relevant parameters for number or sequence selection.

You can undo the selection by pressing the "Default" function key when the display shows "SEQUENCE PARAMETERS". The Urisys 1100 system always measures and memorizes all of the parameters on whichever test strip is used. If, after reading, you wish to reset the choice of parameters to the standard setting, you can reinstate the complete report through the "MEMORY/Print" or "MEMORY/Send" menu.

**Date/time:** For setting the date and time.

The factory default is the date in Day-Month-Year order and the time in hours (24-hour system) and minutes. If required, the time can be displayed in the "12-hour plus a.m./p.m." system. Pressing the "Sequence" key allows the date format to be changed to Month-Day-Year or Year-Month-Day. "Set" causes the date and time to be displayed and set. Pressing the left function key (<<<<) moves the blinking cursor to the left. The time or date unit highlighted by the cursor can then be increased or decreased by pressing the keys (+ / -).

**Pressing the START button confirms the setting, closes the submenu, and resumes Ready-to-Measure status.**

### 3.4 Results Table

The Urisys 1100 system prints the results in the following gradation of concentration:

**ATTENTION: For full details on how to interpret results, consult the following Roche Diagnostics publication, *Compendium Urinalysis with Test Strips*.**

Parameter	Conventional Units (Conv.)	SI Units (SI)	Arbitrary Units (Arbitrary)
SG (Specific Gravity)	1.000	1.000	1.000
	1.005	1.005	1.005
	1.010	1.010	1.010
	1.015	1.015	1.015
	1.020	1.020	1.020
	1.025	1.025	1.025
	1.030	1.030	1.030
pH	5	5	5
	6	6	6
	6.5	6.5	6.5
	7	7	7
	8	8	8
	9	9	9
LEU (Leukocytes)	neg	neg	neg
	25 Leu/ $\mu$ L	25 Leu/ $\mu$ L	1+
	100 Leu/ $\mu$ L	100 Leu/ $\mu$ L	2+
	500 Leu/ $\mu$ L	500 Leu/ $\mu$ L	3+
NIT (Nitrite)	neg	neg	neg
	pos	pos	pos
PRO (Protein)	neg	neg	neg
	25 mg/dL	0.25 g/L	1+
	75 mg/dL	0.75 g/L	2+
	150 mg/dL	1.5 g/L	3+
	500 mg/dL	5.0 g/L	4+
GLU (Glucose)	norm	norm	neg
	50 mg/dL	3 mmol/L	1+
	100 mg/dL	6 mmol/L	2+
	300 mg/dL	17 mmol/L	3+
	1000 mg/dL	56 mmol/L	4+
KET (Ketone)	neg	neg	neg
	5 mg/dL	0.5 mmol/L	(+)
	15 mg/dL	1.5 mmol/L	1+
	50 mg/dL	5 mmol/L	2+
	150 mg/dL	15 mmol/L	3+
UBG (Uribilinogen)	norm	norm	neg
	1 mg/dL	17 $\mu$ mol/L	1+
	4 mg/dL	70 $\mu$ mol/L	2+
	8 mg/dL	140 $\mu$ mol/L	3+
	12 mg/dL	200 $\mu$ mol/L	4+
BIL (Bilirubin)	neg	neg	neg
	1 mg/dL	17 $\mu$ mol/L	1+
	3 mg/dL	50 $\mu$ mol/L	2+
	6 mg/dL	100 $\mu$ mol/L	3+
ERY (Erythrocytes)	neg	neg	neg
	10 Ery/ $\mu$ L	10 Ery/ $\mu$ L	1+
	25 Ery/ $\mu$ L	25 Ery/ $\mu$ L	2+
	50 Ery/ $\mu$ L	50 Ery/ $\mu$ L	3+
	250 Ery/ $\mu$ L	250 Ery/ $\mu$ L	4+

### 3.5 Changing the Range Limits

The reflectance range limits for the lower concentration intervals can be altered to a limited extent for all test parameters except SG (specific gravity) and pH.

To change the range limits, choose “RANGE LIMITS” and press the left function key (“New”). The lowest interval for the first parameter is displayed next to the current range limit setting. Through the keys (+ / -), the reflectance value can be raised or lowered and evaluation sensitivity thus increased or decreased.

If the first positive concentration level is to be eliminated, the reflectance range limits of the negative/normal range can be lowered manually, and the reflectance range limit of the first positive concentration interval raised to this reflectance value.

Pressing “Next” stores that setting and moves to the next concentration interval or the next parameter.

The sequence can be terminated at any point by pressing the START button. The new range limits are printed out together with a date and timestamp and the analyzer resumes Ready-to-Measure status. All changes made to the limits are flagged with an asterisk on the printout after the concentration interval.

Pressing the “Default” key in the “RANGE LIMITS” display window reinstates the factory defaults and causes them to be printed together with a date and timestamp.

If range limits have been altered manually, all results are printed with an asterisk in the line below the Device ID.

Changed range limits apply only to the test strip type for which the changes are made.

**CAUTION: The correctness of results obtained after the user has altered the ranges or reflectance values is not warranted by Roche Diagnostics. The user is responsible for validating the consistency of results after changes have been made.**

## 4. Installation

**WARNING: Read the Urisys 1100 Operator’s Manual carefully before installation, so as to ensure proper operation of the analyzer from the outset.**

**ATTENTION: If your analyzer has been exposed to marked changes in temperature and/or humidity, wait for it to acclimatize sufficiently before operating it.**

**ATTENTION: Upon receipt, check the contents of the box for completeness. If the contents have suffered any damage in transport, contact your local Roche representative.**

Content:

- Urisys 1100 analyzer
- Adapter 100 V - 240 V, 50/60 Hz
- Power Cable
- Roll of printer paper
- 1 test strip tray, Type “C”, for reading Combur<sup>10</sup>Test UX test strips
- 1 test strip tray, Type “N”, for reading Combur<sup>5</sup>Test and Combur<sup>2</sup>Test test strips
- Operator’s Manual
- Interactive Training CD

**CAUTION:** To ensure your readings are always accurate, do not set up the Urisys 1100 system in close proximity to devices that create high-frequency fields, as they may interfere and produce false results. Such devices include, for instance, walkie-talkies, mobile telephones, microwave ovens and diathermic equipment.

## Instrument connection

### Power-on

- (1) Unpack the Urisys 1100 system and place it on a stable, level surface. Do not expose the analyzer to direct sunlight or other source of direct light (for example a spot lamp). For measurements of Combur<sup>10</sup>Test UX test strips, the test strip tray with the type mark C on the reverse side is required.
- (2) Connect the AC adapter to the power socket at the rear of the instrument and to a readily accessible AC wall socket.

### Insertion of the test strip tray

- (3) Pick up the test strip tray, Type "C", and with the retaining bar nearest to you pushed down (see Figure 1), slide it into the slot below the function keys, so that the near edge of the tray is flush with the near edge of the analyzer (see Figure 2).

1

2

**CAUTION:** Be careful not to touch the grey reference pad. Contamination of the reference pad may impair the quality of the results obtained.

Use the test strip tray with the mark N on its reverse side if you measure Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips. Additionally, select the desired test strip type in the menu Setup 1. **The Urisys 1100 system checks whether the defined test strip type matches the selected test strip tray. It may be necessary to configure the system for the desired test strip type.**

### Paper Insertion

To release the printer paper cover, press the structured area immediately below the cover (see Figure 4). The cover can then be lifted back. Place the paper roll in the compartment and pull out the first few centimetres of paper to just beyond the edge of the compartment. The thermosensitive side of the paper (the outer surface of the paper roll) should be facing downwards (see Figure 5). Close the cover again by pressing until it locks audibly into position.

4

5

### ATTENTION:

**To remove the printed test report, tear off the paper by pulling it in front across the edge.**

- (4) Switch on the Urisys 1100 system using the on/off switch at the rear of the instrument.

### Self Check

- (5) The analyzer automatically carries out a self check. The default language setting for the Urisys 1100 system is English. The instrument is factory set to read Combur<sup>10</sup>Test UX test strips. It checks that the correct tray has been inserted. On completion of the self check, the tray is transported to the start position and the retaining bar opens (see Figure 3). "SELF CHECK OK" is printed together with a date and timestamp.

3

**ATTENTION:** If the message E9 Wrong Tray! is displayed after the Urisys 1100 system is switched on, press the START button. The selection menu STRIP TYPE is displayed. If you use the Combur<sup>10</sup>Test UX test strips, select C-10 by pressing the left function key. If you use Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips, select C-7 or C-5, respectively, by pressing the corresponding function keys.

When the self test is finished a sound is heard from the printer. The display will show Insert Strip. If the paper has already been inserted, the printout "Self Check OK" is executed with date and time; this happens whenever the instrument is switched on. If an error occurs, consult the Operator's Manual or contact your local Roche Diagnostics service representative.

- (6) Upon first use of the analyzer, the message “REPEAT CALIBRATION” appears after a successful self check. Press the “No” key to attain **Ready-to-Measure** status. **This also happens if the instrument is not used for more than seven days. To bring the instrument to Ready-to-Measure mode, insert the paper and execute the calibration.**

**! ATTENTION: If Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips are used, a regular calibration is not necessary with the Control-Test M, because in this case, the instrument calibration is performed via the grey reference field in the test strip tray.**

- (7) To set the language of your choice, press the key sequence “Menu/Setup 1/Setup 2/Setup 3/Language/Other” (see Sections 3.2 and 3.3).
- (8) The Urisys 1100 system leaves the factory with the default settings for Printer, Measure and Interface modes, Test Strip Type, Units, Language, Date/time format (see Section 3.2) and Operator ID. If your laboratory works differently, individual preferences can be entered via the menus.
- (9) Calibrate the analyzer with Control-Test M (see Section 5) if you intend to use Combur<sup>10</sup>Test UX test strips for urinalysis.

#### **Power-off**

We recommend that the analyzer be switched off at the end of each working day and that the mains adapter be unplugged from the AC wall socket (see Section 7.1, Cleaning).

### **4.1 Modification for Combur<sup>2</sup>Test and Combur<sup>5</sup>Test Test Strips**

For readings with Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips a Type “N” test strip tray is required (for an order number, see Section 11), which exactly matches the geometry of the shorter test strips. Set the desired strip type from the menu. After pressing START to confirm, insert the Type “N” tray in the reader. The Urisys 1100 system is then ready (see Section 6).

**! ATTENTION: Combur<sup>2</sup>Test and Combur<sup>5</sup>Test test strips do not have a compensation pad. Consequently, the Urisys 1100 system does not then compensate for strong intrinsic urine coloration during reading, so that occasional false-positive readings may be obtained.**

**! ATTENTION: For Combur<sup>2</sup>Test or Combur<sup>5</sup>Test test strips you do not need to calibrate regularly with Control-Test M. With these strips the reader is calibrated entirely through the reference pad in the tray. Always ensure that the reference pad is clean and undamaged (see Section 7).**

### **4.2 Operator ID**

The Urisys 1100 software allows the activation/deactivation of the operator identification code, containing up to 12 alphanumeric characters. The Operator ID and Authentication mode can be activated in Setup 3.

**! ATTENTION: Please ensure that you have a barcode reader and/or an AT/PC keyboard prior to activation, as you will need one for this function.**

#### **Normal**

If activated, the operator ID is asked upon every restart of the system and coming out from sleep mode. This operator ID can contain up to 12 alphanumeric characters and will be printed out together with the test results and will be sent to the host.

**ATTENTION: If the maximum length of identification is exceeded the input cursor will skip to the first character and the identification will be overwritten.**

### **4.3 Authentication**

It is possible to download up to 300 operator IDs with corresponding passwords (up to 12 alphanumeric characters) from the host PC via the ASTM protocol.

Device can be used only by an operator with ID and password in the downloaded list. Entry of incorrect operator ID and passwords causes a lock out. This prevents access of the instrument by unauthorised users.

Operators having supervisor rights (maximum 2) have access to all results, may send the instrument log file to the host or print the last 10 actions of the log file and may deactivate the authentication mode.

## 5. Calibration

**ATTENTION: The following statements only need to be observed when the Urisys 1100 system is used in conjunction with Combur<sup>10</sup>Test UX test strips.**

The Urisys 1100 system is calibrated before leaving the factory. When installed, it must be recalibrated with Control-Test M calibration strips before the first samples are read, and thereafter every seven days. Control-Test M calibration strips consist of a grey plastic material that is standardized to give constant, defined reflectance readings. The purpose of calibrating the analyzer is to compensate for aging effects that influence the optical system and the grey reference pad in the strip tray. If the compensation needed is excessive, for example because the reference pad is badly soiled, or an LED is defective and cannot emit the required amount of light, an error message is displayed (see below).

When the Urisys 1100 system is set to Combur<sup>10</sup>Test UX test strips, it automatically requests a new calibration every week. This is why, when the instrument is installed, the "REPEAT CALIBRATION!" message is displayed following the successful self check and also one week after the latest valid calibration.

Press the left function key ("Yes") to read the "START CALIBRATION!" message.

### Procedure

1. Remove a calibration strip from the Control-Test M container. Be careful not to touch the pads and do not allow them to come into contact with urine.
2. Place the calibration strip, with the test pads facing upwards, on the tray so that its leading edge is held by the clip at the front end of the insertion slot. **The retaining bar must be open** (see Figures 6 and 7). Before calibrating, ensure that the tray is clean and dry.

6

7

**ATTENTION: It is very important that the calibration test strip locks into the instrument correctly in order to ensure the quality of the calibration.**

3. Press the START button. An acknowledging beep sounds. The tray advances slightly, the retaining bar closes, and the grey reference pad on the tray and the calibration pads are read.
4. The tray is then transported back to the start position, and the retaining bar opens. Remove and dispose of the calibration strip. Use each calibration strip once only.
5. If the calibration is valid, the result is stored together with the date and time, and printed (see figure 8). Remove the test strip and discard properly (follow test strip insert instructions).

8

**CAUTION: Regular calibration is necessary to ensure the quality of the results obtained. Roche Diagnostics cannot warrant the correctness of results if the system is not calibrated regularly.**

**ATTENTION: You can recalibrate the system at any time, even if a week has not elapsed since the last calibration, for example if a Combur<sup>10</sup>Test UX test strip has returned an implausible result. Start the calibration procedure by pressing the left function key ("Calibr.") with the analyzer in Ready-to-Measure mode. The message "START CALIBRATION !" appears. Follow the calibration procedure described page 36.**

**ATTENTION: If you reply with "No" to the weekly "REPEAT CALIBRATION !" reminder, for example because you have used up your supply of Control-Test M calibration strips, the message "Repeat calibration" is printed together with each patient report from subsequent readings. Obtain a new pack of Control-Test M and recalibrate immediately.**

### Calibration printout

If the new calibration results are within the permitted range, the message “CALIBRATION O.K.” is automatically printed with a date and timestamp and also a list of reflectance values for measuring positions 1–11 for the orange LED (middle column) and green LED (right column) (see Figure 8).

### Calibration errors

If the results obtained for the reference pad or the calibration strip are outside the programmed tolerances, the following messages appear: “REFERENCE PAD ERROR !”, “CALIBRATION INVALID !” or “CALIBRATION ERROR !”.

In the event of a calibration error or if the calibration was invalid, repeat the calibration procedure with a fresh Control-Test M calibration strip. Press the START button to return to the “START CALIBRATION” menu. Follow the same calibration procedure as above. When the message “CALIBRATION O.K.” has been printed, proceed with the reading of test strips. If you continue to receive an error message, see Section 8.

If the analyzer is switched off when one of the above-mentioned calibration status messages is displayed and then switched on again, it reverts to Ready-to-Measure status following a successful self check, provided that the most recent calibration was carried out not more than a week ago. Otherwise, the message “REPEAT CALIBRATION !” appears (see above).

## 6. Reading Test Strips

### 6.1 Overview

The Urisys 1100 system is very easy to use. Simply insert the test strip when the sample number is displayed (refer to Section 6.2), then press the START button to commence reading. The analyzer can be operated in two different modes:

1. In **Normal Mode**, the Urisys 1100 system automatically waits for the strip to incubate for 55 seconds before it reads the first test pad. Seventy seconds after the START button is pressed, the measurement is completed and the test strip tray moves to the start position. Throughput in this mode is approximately 50 strips per hour.
2. In **Fast Mode**, which can be selected from the menu, the test strip is measured directly after START is pressed. In this case, it is up to the user to time the incubation period outside the analyzer (see Section 6.4). Using **Fast Mode** enables a 30-second cycle to be maintained.

Each time a strip is read, the grey reference pad in the tray is evaluated to compensate for temperature and aging effects that may influence the optical system. If the compensation needed is excessive, for example because the reference pad is badly soiled or an LED is defective and cannot emit the required amount of light, an error message is displayed (see Sections 7.1 and 8).

The Urisys 1100 system assigns to each reading a consecutive sequence number (sample number) having a maximum of three digits. The sequence start number automatically reverts to 1 each time the date is incremented. You may, if you wish, reset the sequence number to 1 via the “New Series” function, for example when one series of measurements has been completed and another is due to begin.

In **Automatic Patient ID** mode the instrument will assign automatically unique serial numbers to the tests results which have no Patient IDs. These unique numbers are ascending serial numbers based on the total number of tests performed on the instrument and cannot be altered nor cleared.

Whenever there has been five minutes of inactivity, the analyzer automatically switches to **Standby** mode. The tray advances slightly so as to close the retaining bar, and the display shows the date and time. The analyzer resumes Ready-to-Measure status when the START button is pressed.

### 6.2 Normal Mode (for Single Readings)

**BIOHAZARD: Always wear protective gloves when handling and disposing of samples of human origin.**

The Urisys 1100 system is ready to read when the display shows a sample number and “INSERT STRIP !”.

**ATTENTION: To ensure that urinalysis is carried out correctly, read the package insert included with the test strips.**

1. Dip the test strip in the urine sample and carefully wipe off any excess urine. Briefly (approximately 1 second) dab the side edge of the strip on a clean, absorbent surface (for example kitchen roll) (see Figure 9). 9
2. Place the test strip, with the test pads facing upward, on the tray so that its leading edge is held by the clip at the insertion slot. **The retaining bar must be open** (see Figure 10). About 2 mm of strip must be held under the clip (see Figure 11). 10  
11
3. Press the START button (see Figure 12). An acknowledging beep sounds. The tray advances slightly, the retaining bar closes, and the grey reference pad on the tray is read (see Figure 13). 12  
13

**ATTENTION: Ensure that the retaining bar is locked into place and that the test strip is in the correct position. If the test strip is not correctly located in the middle of the tray, move it gently to the side until it is properly aligned (see Figure 14). Be careful not to move the tray.**

4. Fifty-five seconds after the START button is pressed, the first test pad is measured, followed by the others. After that, the tray returns to the start position and the retaining bar opens. 14
5. Remove and dispose of the test strip. Wipe any urine residues from the tray with a lint-free cloth (see Figure 16). 16
6. The result will be printed out and the next sample number will be displayed. The next test strip can be dipped, wiped off, placed in the tray and read by pressing the START button. Refer to Section 6.6 for details of what must be observed when working with patient identification numbers.

### 6.3 Patient Report

The patient report is printed out together with the sequence number, device ID, operator ID, date and time. The patient's name will also appear on the print out if entered prior to measurement (see Section 6.6). Test results which diverge from negative or normal values are flagged with an asterisk before the parameter concerned. For selection of concentration units and the order of parameters in the patient report, see Section 3.3. An asterisk printed in the line below the heading but before the first parameter of the test results list indicates that the reading was taken with changed range limits (see Section 3.5) 15

Tear off the printout, if desired, by pulling it horizontally over the edge.

**WARNING: Thermal printing paper is sensitive to light and liable to age on prolonged exposure to bright light. Patient reports are thus best kept in a place away from light.**

### 6.4 Fast Mode (for Serial Readings)

You can also perform serial readings provided that the strips are dipped and incubated for about 45 seconds outside the Urisys 1100 system. For this type of reading, "Fast Mode" should be selected from the menu. Place the test strips on the Urisys 1100 system tray after allowing them to incubate for approximately 45 seconds, pressing START to begin each reading.

The reference pad is read immediately, followed by the test pads on the test strip. Reading is complete after 30 seconds and the next sample number appears in the display. When working in Fast Mode, ensure that you have a foolproof system for matching sequence numbers to samples.

#### Remarks:

1. Results obtained in Fast Mode are saved in memory immediately; however, printing is deferred. Printout of the whole series begins automatically following 60 seconds of inactivity. Check carefully whether all sequence numbers have been printed out. If any results are missing, press "Memory/Print" to printout stored results.

2. When Fast Mode is selected, the analyzer continues to operate in Fast Mode for the rest of the calendar day, even if the analyzer is switched off and then on again.
3. Fast Mode is cancelled as soon as the date is automatically incremented.

**WARNING: When performing serial measurements in Fast Mode, allow the strips to react for approximately 45 seconds before inserting them in the Urisys 1100 system and pressing START. False-low or false-negative results may be obtained for some parameters if the reaction time is too short. Likewise, false-high results may be obtained for some parameters if the incubation time outside the Urisys 1100 system is too long.**

## 6.5 Strip Measurement Error

If "STRIP MEASUREMENT ERROR !" appears, the test strip and/or analyzer have probably been incorrectly used. Refer to Section 8 for details.

## 6.6 Entering Patient ID, Operator ID and Authentication Password

### Patient ID

When the analyzer is ready to measure ("INSERT STRIP !" displayed), you may enter a Patient ID not exceeding 13 characters in length against the currently shown sequence number by means of a barcode reader or AT/PC keyboard (see Section 9.2). The Patient ID can be verified in the display window and entered again if necessary. The last Patient ID entered is stored when the START button is pressed, i.e. when reading commences, and is printed out and/or sent to the serial interface together with the test result.

### Operator ID

If the operator ID is activated, the instrument request the input of the operator ID immediately after the self check is performed when the instrument is turned on, or when it leaves the stand-by mode. You may enter an operator ID up to 12 alphanumeric characters by means of a barcode reader or an AT/PC keyboard.

### Authentication

In the authentication mode the user will be required to input his/her apart from the operator ID list, also the corresponding password in order to have access to the instrument and the test results he/she had performed.

Operators with supervisor rights have access to all results, may send the instrument log file to the host or print the last 10 actions of the log file and may deactivate the authentication mode.

The list of operator IDs with corresponding passwords may be updated from the host PC using the "Download List" function key.

**ATTENTION: A new Patient ID should only be entered from the keyboard after the preceding measurement has been completed and results printed out (Normal Mode). This allows the operator to check correct entry directly in the display.**

**ATTENTION: If a Patient ID is entered via the barcode reader while a test strip is being read, the Urisys 1100 system assigns that ID to the next sample number in the sequence. The Patient ID can only be erased by switching the Urisys 1100 system off and then on again before starting the next reading.**

**ATTENTION: If an identification will be entered and does exceed the maximum lengths of 13 characters for the patient Id or the 12 characters for the operator ID the cursor will move to the beginning (left side) of the display and delete the previously written characters.**

When data are entered from the keyboard, each character appears immediately in the display. It can be erased by backspacing, and corrected as necessary. **Press the keyboard ENTER key to terminate input**, otherwise the Urisys 1100 system cannot start reading. Press the keyboard Escape key to delete the entire entry or turn the Urisys 1100 system off and then on again.

## 6.7 Data Transmission to a PC or Host Computer

In **unidirectional** mode, the results are transmitted immediately together with the sequence number, Patient ID (if entered), date and time.

In **bidirectional** mode, transmission can only be accomplished by using the "Send" function when "MEMORY" is displayed. If a bidirectional PC/host communication link cannot be established, the Urisys 1100 system aborts transmission after several attempts and reports an "INTERFACE ERROR !" (see Section 8).

In **ASTM** mode the results, sequence number, Operator ID, Device ID, Patient ID (if entered), date and time of the measurement, and of the last calibration will be sent to the host.

For further information on the serial interface, refer to Section 9.1.

**ATTENTION:**  
Ensure that all required data is backed up on a regular basis.

## 7. Cleaning and Maintenance

The Urisys 1100 system is designed for maintenance-free operation. Protect the instrument from extremes of temperature and high atmospheric humidity (see Section 10), and keep it out of bright light (direct sunlight, spot lamps, etc.).

Maintain hygiene by keeping the exterior parts and surfaces of the instrument clean. For cleaning we recommend applying a commercial cleaning agent or disinfectant (preferably 70 % alcohol) with a moist cloth. Take care that no liquid enters the instrument.

### 7.1 Cleaning the Test Strip Tray

**BIOHAZARD: Liquid waste and strip waste are potentially biologically hazardous. Always wear gloves if handling those materials. Dispose of the used test strips according to the regulations for handling potentially infectious material.**

When inserting and removing test strips, take care that no urine residues come into contact with the retaining bar mechanism.

To prevent carry-over and an accumulation of urinary deposit, wipe off residues of urine from the tray with a lint-free cloth after each reading (see Figure 16).

16

At the end of each working day, clean the tray with water and, if necessary, with a commercial cleaning agent and disinfectant. The procedure is as follows:

1. Switch off the instrument.
2. Pull the tray out of the instrument.
3. Rinse the contaminated parts of the tray under running water, then clean with 70 % alcohol or another suitable disinfectant. Slight crystalline deposits, especially those contaminating the retaining bar mechanism or the cogs on the underside of the test strip tray, may be removed with a soft brush. Dry the parts with a dry, lint-free cloth.

**CAUTION: Take care not to damage the grey reference pad during cleaning, and ensure that it is completely clean and dry before proceeding to read.**

**CAUTION:** Ensure that the positioning hole on the side of the tray is absolutely dry (see Figure 17). This hole is used to ensure that the test strip tray is automatically positioned correctly in the instrument.

17

- Pick up the cleaned test strip tray and, with the retaining bar nearest to you and closed (down), slide it into the slot below the function keys, so that the near edge of the tray is flush with the near edge of the analyzer.

**CAUTION:** Be careful not to touch the grey reference pad. Contamination of the reference pad may impair the quality of the results obtained.

- If you wish to proceed with the next readings directly after cleaning the test strip tray, switch the Urisys 1100 system on again. During the self check the system will verify that the reference pad is in good condition for reading and that the positioning hole in the test strip tray (see Figure 17) is free. If not, an error message will be displayed (see Section 8).

17

## 8. Error Messages and Troubleshooting

Error messages are shown in the display but are not printed out. Following 5 minutes of inactivity, the analyzer switches to Standby mode. The error message is displayed again when the START button is pressed. In the event that the instrument is defective, contact your local Roche Diagnostics representative.

### **E1 REFERENCE PAD ERROR MIDDLE!**

**Cause:** The middle portion of the reference pad on the tray is soiled or damaged.  
**Action:** Switch off the instrument. Carefully clean and dry the pad. Check if it is damaged (e.g. scratched, etc.). Insert the tray again, then wait for the self check to finish. If the error message is repeated, replace the reference pad or use the spare tray. Recalibrate with Control-Test M (when using Combur<sup>10</sup>Test UX test strips).

### **E15 REFERENCE PAD ERROR BOTTOM!**

**Cause:** The bottom portion of the reference pad on the tray is soiled or damaged.  
**Action:** see E1

### **E16 REFERENCE PAD ERROR TOP!**

**Cause:** The top portion of the reference pad on the tray is soiled or damaged.  
**Action:** see E1

### **E2 WRONG STRIP !**

**Cause:** The test strip used is different from the one for which the analyzer has been programmed (e.g. Combur<sup>10</sup>Test UX test strips).  
**Action:** Press the START button. Repeat the measurement with the type of strip for which the analyzer has been programmed.

### **E3 STRIP MEASUREMENT ERROR !**

**Causes:** No test strip is present on the tray, or the strip is incorrectly positioned on the tray, the urine on the test strip has dried, the test strip has not been dipped in urine.  
**Action:** Press the START button. Repeat the measurement with a new test strip. Ensure that all test pads are dipped in the urine sample. Insert the strip correctly and ensure that the retaining bar is closed properly after START is pressed.

### **E4 CALIBRATION ERROR !**

**Cause:** Calibration values differ from those obtained in the last valid calibration.  
**Action:** Press the START button. Repeat the calibration with a new calibration strip taken from the Control-Test M container. Ensure that the strip is properly positioned under the clip on the test strip tray (see Section 5).

### **E5 CALIBRATION INVALID !**

**Cause:** Calibration values are out of tolerance.  
**Action:** Check the reference pad for soiling or damage. Clean if necessary (see Section 7.1) or use the spare tray. Repeat the calibration with a new Control-Test M calibration strip (see Section 5). If the error message appears again, the instrument is defective.

**E6 CHIP ERROR !**

**Cause:** The program chip on the right of the analyzer underneath the printer cover (see Figure 18) is missing, is not making contact, is defective or contains an old software version.

**18**

**Action:** Switch off the Urisys 1100 system. Insert the program chip and switch the instrument on again. If "CHIP ERROR" appears again, the instrument is defective.

**E7 MISSING TRAY !**

**Cause:** No tray inserted or tray not inserted far enough to be engaged by the motor.

**Action:** Insert the tray correctly (see Section 4). Press the START button.

**E8 TRAY POSITION ERROR !**

**Cause:** The positioning hole in the tray (see Figure 17) is soiled or still wet after cleaning; the retaining bar is open while the tray is advancing or the retaining bar mechanism is fouled with urinary deposit and blocked. (See Section 7.1)

**17**

**Action:** Clean, blow through or dry the positioning hole (using a lint-free cloth) to ensure that it is completely patent. Remove urinary deposits, if there are any, including those on the underside of the tray. Insert the tray again and press the START button. Ensure that the retaining bar is down and locked into place while the reading is taking place. If the error message appears again, use the spare tray.

**E9 WRONG TRAY !**

**Cause:** The test strip tray used is not the one for the programmed test strip type, or the grey reference pad is missing from the tray, or there is an instrument error.

**Action:** Press the START button. The strip type menu is displayed. The strip type must match the tray type (see Section 11). Use the correct tray, or order a new one. If the error message appears again, the instrument is defective.

**E10 LIGHT BARRIER ERROR !**

**Cause:** The light barrier used to control the position of the test strip tray is defective or the tray transport is blocked.

**Action:** Pull out the tray and return it to the start position. Press the START button. If the error message appears again, the instrument is defective.

**E11 MOTOR STEP ERROR !**

**Cause:** The stepping of the motor is out of tolerance or the advance of the tray is blocked. This may be due to

- soiling on or between the cogs
- worn or broken cogs
- defective motor.

**Action:** Carefully clean the tray. Remove any urinary deposits, including those on the underside of the tray and the cogs. If the tray is damaged insert the spare tray. Press the START button. If the error message appears again, the instrument is defective.

**E12 OPTICS ERROR !**

**Cause:** The reference pad is missing from the tray, or an LED or the phototransistor is defective.

**Action:** Attach the reference pad or use the spare tray. Press the START button. If the error message appears again, the instrument is defective.

**CLOSE PRINTER COVER**

**Cause:** The printer cover is open.

**Action:** Close printer cover.

**NO PAPER IN PRINTER**

**Cause:** No paper has been inserted or roll is finished.

**Action:** Insert new roll of paper and close printer cover.

After elimination of printer errors the results can be printed from the instrument's memory using the "Print" function.

**E14 INTERFACE ERROR !**

**Cause:** Fault in data transfer to PC or host in **bidirectional** or **ASTM** mode.

**Action:** Check the data cable. Verify that the PC or host is ready to receive data. Use the "Send" function to transfer data or press "Home" to resume Ready-to-Measure status.

**E17 INVALID PASSWORD!**

**Cause:** The entered password doesn't match.

**Action:** Enter a correct password.

<b>E18 INVALID OP.ID!</b>	<p><b>Cause:</b> The entered operator ID is not valid.</p> <p><b>Action:</b> Enter a valid operator ID.</p>
<b>E19 LIST DOWNLOAD FAILED</b>	<p><b>Cause:</b> The new operator ID list download failed.</p> <p><b>Action:</b> No action. After 2 seconds next state starts with old list if there was.</p>
<b>E20 NO VALID LIST!</b>	<p><b>Cause:</b> There is no valid list at all in the device.</p> <p><b>Action:</b> Try to download an operator ID list from the host or continue without authenticated operator.</p>
<b>CHECK MEASUREMENT</b>	<p>Analyzer prints out software and chip version number and 3-digit numbers without naming the parameters.</p> <p><b>Cause:</b> Service function is activated.</p> <p><b>Action:</b> Press the "Back" function key to return to the main menu.</p>
<b>VALUES OBTAINED ARE IMPLAUSIBLE WHEN COMPARED WITH THOSE FROM VISUAL EVALUATION</b>	<p><b>Causes:</b> Test strip incorrectly positioned, uncharacteristic test pad colors, proper incubation intervals not kept to during serial measurements. The wrong test strip, such as Combur<sup>10</sup>Test test strips, may have been used. Electromagnetic interference from other devices (see Section 4).</p> <p><b>Action:</b> Repeat the measurement with a new instrument-compatible test strip. Follow the directions carefully and ensure the test strip is correctly inserted (see Section 6.2). Repeat calibration if necessary. Remove external sources of interference, if there are any.</p>
<b>NO PRINTOUT</b>	<p><b>Causes:</b> "Printer: Off" has been selected, or the analyzer is operating in "Fast Mode" (see Section 6.4), or the printer/software is defective, or the printer is out of paper.</p> <p><b>Action:</b> Insert paper if needed. Choose "Printer: On" to re-activate the printer. Request a patient report via the "Print" function. If this fails, activate the "Line-feed" function. If there is still no response, the instrument is defective.</p>
<b>THE ANALYZER WILL NOT READ EVEN THOUGH THE SEQUENCE NUMBER IS DISPLAYED</b>	<p><b>Action:</b> If an AT/PC keyboard is connected, press the Escape key, or switch the Urisys 1100 system off and back on again.</p>

## 9. Connecting to Other Devices

### 9.1 Serial Interface

At the rear, the Urisys 1100 system has a serial interface through which it can be connected to a PC or central host computer. This is not an RS 232 type interface.

On request, Roche Diagnostics can supply a suitable standard data cable (see also Section 11). The connected PC must satisfy the requirements with regard to electrical safety laid down in EN 60950.

The interface can be used for **unidirectional**, **bidirectional** or **ASTM** communication, selectable via the menu.

When the interface is set for unidirectional communication, the data are sent as an ASCII file and can be received via a terminal program.

The entered Patient ID appears in the Urisys 1100 system display window and is also printed out and/or sent to the PC/host along with the test results (see Section 6.6). If bidirectional communication with a PC or host computer has been selected, the maximum length of the Patient ID used (either 10 or 13 characters) must be preprogrammed via the display message "INTERFACE: BIDIR." and the function "10/13" to ensure that the correct data is sent.

For further information and specifications for operation in **bidirectional** or **ASTM** mode, e.g. for connection to a host computer, contact Roche Diagnostics.

**Interface specification:** 9600 baud, 8 bits, 1 stop bit, no parity (for unidirectional and bidirectional modes).

Selectable baud rates in ASTM mode: 1200, 2400, 4800, 9600, 19200 and 38400.

**Data cable:** D-sub, 9-pin, male on instrument side, female on PC side.

		<b>Connections:</b>	
Urisys 1100			Host (PC pinout 9-pin)
2	—————	RxD	————— 2
3	—————	TxD	————— 3
4	—————	DTR	————— 4
5	—————	GND	————— 5
6			6
7			7
8			8
9			9

**CAUTION:** The use of a data cable not meeting the Roche Diagnostics specification can cause data to be lost or corrupted.

### 9.2 Barcode Reader, AT/PC Keyboard

Sample or Patient IDs, Operator IDs and corresponding passwords can be entered against each sample sequence number displayed on the LCD, either via a barcode reader (see recommended reader) or via an AT/PC keyboard, e.g. the Reflotron keyboard from Roche Diagnostics (see ordering Information chapter 11). Power is supplied by the barcode reader interface.

**Interface specification:** 5-pin DIN socket, female

<b>Pinouts:</b>		
1	clock	
2	data	
3	n/c	
4	GND	
5	+ 5 V	

## Barcode Reader

Barcode readers suitable for use with Urisys 1100 with SW Version 5.0 and above must meet the following specifications:

- Radio frequency interference class B according to EN 61326-1
- Electromagnetic interference immunity requirements for industrial locations according to EN 61326-1

There is a recommended barcode reader to read commonly used barcodes such as Codabar, Code 39, Code 128 and Interleaved 2 of 5.

If a barcode reader and an external keyboard have to be connected a data cable CAB 322 IBM AT/XT DIN is needed.

For questions regarding the operation of the Urisys 1100 system with barcode readers, please contact your local Roche Diagnostics representative.

# 10. Technical Information and Notices

## 10.1 Technical Data

<b>Dimensions:</b>	Width:	approx. 150 mm
	Depth:	approx. 290 mm
	Height:	approx. 95 mm
<b>Weight:</b>		≤ 0,8 kg
<b>Power supply:</b>	External mains adapter, Model SA 125A-0735U-S (Sino-American)	
	Input: 100 - 240 V AC, 50/60 Hz, 800 mA	
	Output: 7.5 V DC, 3000 mA	
	Polarity: - ———⚡———— +	
<b>Consumption:</b>	Operating:	max. 15 W
	Standby:	1.3 W
<b>System description:</b>	Type:	reflectance photometer
	Light source:	6 LEDs (light emitting diodes)
	Wavelengths:	565 nm (green) 3x 610 nm (orange) 3x
	Reader head:	1 head with 6 LEDs
	Measuring cycle:	Normal Mode: approx. 70 sec Fast Mode: approx. 30 sec
	Incubation time:	55-65 sec
	Printer:	thermal printer
	Display:	liquid crystal display, 2 lines of 24 characters
	Memory:	100 samples
	Date, time:	integrated clock
<b>Operating conditions:</b>	<u>Operating</u>	<u>In storage</u>
Temperature:	+15 ° to + 32 °C	- 20 ° to + 70 °C
Relative humidity:	20 % to 80 %	20 % to 85 %
Optimum operating conditions:	Temperature:	+ 20 ° to + 26 °C
	Relative humidity:	30 % to 60 %
<b>Interfaces:</b>		
PC/HOST:	serial, D-Sub socket, 9-pin, female, <b>unidirectional</b> , <b>bidirectional</b> or <b>ASTM</b> protocol (selectable)	
AT/PC keyboard:	Barcode reader 5-pin DIN socket, female	
<b>Certification marks:</b>	UL, cUL	

## 10.2 Safety Notices

This analyzer was designed and manufactured to comply with the following international regulations “Safety requirements for electrical equipment for measurement, control and laboratory use” and left the factory in a safe condition. In order to keep the instrument in a perfect and safe condition, it is up to the user to observe all instructions and warnings included in this manual.

This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.

Issued by Underwriters Laboratories Inc. (UL) for the USA and Canada.

The instrument must only be operated with the prescribed power supply unit (Class II protection).

The instrument is classified as Category II for overvoltage and Degree 2 for pollution according to IEC 664.

Opening covers or removing parts of the instrument, except where this can be achieved manually without the use of any tools, may expose voltage-carrying components. Connectors can be live, too. Never try to maintain or repair an open instrument which is carrying voltage.

If you suspect that the instrument can no longer be operated safely, turn it off and take steps to ensure that no one will subsequently attempt to use it. Make sure that only trained members of staff operate the Urisys 1100 analyzer.

Any personal computer to which the analyzer is connected must meet the EN 60950, UL 60950/CSA C22.2 No. 60950 requirements for data processing equipment.

If the instrument is to be taken out of operation entirely and disposed of, it must be disposed of in conformity with the relevant legal regulations and in co-ordination with your local authority, if appropriate.

Please note that the instrument may potentially be infectious. It should therefore be decontaminated before disposal, e.g. by cleaning the housing and the test strip tray with 70 % alcohol.

**ATTENTION: The data and information contained in this manual are accurate at the time of going to press. Any substantial changes will be incorporated in the next edition. In case of conflict between this manual and information given in package inserts, the package inserts shall take precedence.**

## 10.3 Guarantee

The statutory provisions on rights in consumer goods sales in the country of purchase shall apply.

## 11. Ordering Information

The analyzer, consumables and accessories are:

Catalogue No.  
(REF)

03 617 548	Content: Urisys 1100 analyzer, Adapter Model SA 125A-0735U-S, 100 V - 240 V, 50/60 Hz, Program Chip, Test Strip Trays C and N, Printer Paper, Operator's Manual, Quick Reference Guide, power cable, Interactive Training CD.
03 666 735	Test strip tray Type C for Combur <sup>10</sup> Test UX test strips
03 666 913	Test strip tray Type N for Combur <sup>2</sup> Test/Combur <sup>9</sup> Test test strips
03 617 572	Urisys 1100 Operator's Manual
11 544 373	Combur <sup>10</sup> Test UX (100 test strips)
11 008 552	Combur <sup>2</sup> Test (100 test strips)
11 893 467	Combur <sup>9</sup> Test (100 test strips)
11 379 194	Control-Test M (50 calibration strips)
03 666 751	Thermal paper (5 rolls)
11 906 186	Interface connection cable
11 248 685	Reflotron keyboard for patient ID entry (German)
11 248 723	Reflotron keyboard for patient ID entry (English)
11 248 707	Reflotron keyboard for patient ID entry (French)
11 248 715	Reflotron keyboard for patient ID entry (Italian)
11 248 995	Reflotron keyboard for patient ID entry (Spanish)
11 428 667	Reflotron keyboard for patient ID entry (Schwedish)
11 428 675	Reflotron keyboard for patient ID entry (Norwegian)

The following replacement parts are available:

11 907 131	Spare reference pads (5 pieces)
05345618001	Program Chip

## 12. Alphabetical Index

	<b>Section</b>		<b>Section</b>
<b>A</b>		<b>O</b>	
Acoustic signal (beep)	6.2	Operating parts	2.2
<b>B</b>		<b>P</b>	
Barcode reader	9.2	Parameters in patient report, setup	3.3
Basic settings	3.3	Pat. Id. entry	6.6
Beep (acoustic signal)	6.2	Patient identification	6.6
<b>C</b>		Positioning hole	7.1
Cable configuration	9.1	Positioning of test strip	6.2
Calibration	4., 4.1, 5.	Power supply	4.,10.1
Calibration printout	5.	Printer paper	11.
Cleaning	7.	Printer paper, insertion of	4.
Components and Functions	2.2	Printer, setting	3.3
Computer interface	3.3, 6.7, 9.1	Printing results from memory	3.3
Consumables and accessories, ordering	11.	Printout with copy	3.3
Control-Test M	3.3, 5., 11.	<b>R</b>	
<b>D</b>		Range limits	3.3, 3.4, 3.5
Data transfer	3.3, 6.7, 9.1	Range limits, change of	3.5
Date setting	3.3	Reference pad	4., 4.1, 7.1, 11.
Deleting data	3.3, 6.6	Results, deletion of	3.3
Deletion of results	3.3	Results, retrieval of	3.3
Design of result report form	3.3	Results, storage of	3.3
Display messages	3.2, 8.	<b>S</b>	
<b>E</b>		Sample identification	6.6
Error messages	8.	Self Check	3.3
<b>F</b>		Sensitivity, adaptation of	3.3, 3.5
Flow Chart Software (Menu Structure)	3.2	Sequence number	6.1
<b>G</b>		Sequence number reset	3.3, 6.1
Gradation of concentration in the result printout	3.4	Sequence of parameters in report form, setup	3.3
<b>I</b>		Siting, proper siting of the instrument	4.
Installation	4.	Software	3.2, 3.3
Interface connection cable	9.1, 11.	Standby mode	3.1
Interface, configuration of	3.3, 6.7, 9.	<b>T</b>	
Interface specification	9.1	Taking a reading	6.2, 6.4
<b>K</b>		Technical data	10.1
Keyboard	6.6, 9.2	Test strip measurement	6.
<b>L</b>		Test strip tray	2.2, 3.3, 4., 4.1, 11.
Language selection	3.3	Test strips	1, 3.3, 4., 4.1, 11.
Linefeed	3.3, 4.	Time, setting the	3.3
<b>M</b>		Troubleshooting	8.
Mains voltage	10.1	<b>U</b>	
Measuring principle	2.1	Units, setting of	3.3, 3.4, 3.5
Memory	3.3		
Menu structure	3.2, 3.3		

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Phone: +54 11 5129-8000  
Fax: +54 11 5129-8105

### **Australia**

Roche Diagnostics Australia Pty. Ltd.  
31 Victoria Avenue  
Castle Hill, NSW 2154

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Phone: +56 2 4413200  
Fax: +56 2 4626687

### **China**

Greater China (China, Hong Kong) Headquarter  
Roche Diagnostics (Shanghai) Ltd  
12F, Huaihai Plaza  
No. 1045 Central Huaihai Road  
Shanghai 200031

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Phone: +91 22 56974900  
Fax: +91 22 56974909

**Indonesia**

PT Roche Indonesia Diagnostics Division  
Artha Graha Building 21st Floor,  
Sudirman Central Business District - Lot 25,  
Jl. Jend. Sudirman Kav. 52 - 53  
12190 Jakarta

**Indonesia**

Phone: +62 21 514 000 91/81  
Fax: +62 21 514 000 90

**Italy**

Roche Diagnostics S.p.A.  
Viale G.B. Stucchi 110  
20052 Monza (MI)

**Italia**

Phone: +39 039 2817 1  
Fax: +39 039 2817 654

**Japan**

Roche Diagnostics K.K.  
DC Product Department  
5F Nippon Roche Building  
6-1, Shiba 2-chome, Minato-ku,  
Tokyo 105-0041

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Phone: +421 7 54 78 85 079  
Fax: +421 7 54 78 85 12

**South Africa**

Roche Products (Pty) Ltd - South Africa  
Diagnostics Division  
9 Will Scarlet Road/Ferndale  
PO Box 1927  
Randburg 2125

**South Africa**

Phone: +27-11-504 4600  
Fax: +27-11-886 0942

**South Korea**

Roche Diagnostics Korea Co., Ltd  
15th Floor, Sam-Wha Building, 144/17, Samsung-dong,  
Kangnam-Ku  
135-092 Seoul

**Republic of Korea**

Phone: +82 2 550 3300  
Fax: +82 2 566 1878

**Spain**

Roche Diagnostics S.L.  
Avda. de la Generalitat, s/n  
08174 Sant Cugat del Vallés (Barcelona)

**España**

Phone: +34 935 834 000  
Fax: +34 935 834 340

**Sweden**

Roche Diagnostics Scandinavia AB  
Karsbodavägen 30, Box 147  
16126 Bromma

**Sverige**

Phone: +46 (0)8 404 88 00  
Fax: +46 (0)8 98 44 42

**Switzerland**

Roche Diagnostics (Schweiz) AG  
Industriestrasse 7  
6343 Rotkreuz

**Schweiz/Suisse/Svizzera**

Phone: +41 41 799 61 61  
Fax: +41 41 799 65 52

**Taiwan**

Roche Diagnostics Ltd.  
11F, No. 35 Sec.  
3 Min Quan Road, 104 Taipei

**Taiwan**

Phone: +886-2-2505 2533  
Fax: +886-2-2502 1441

**Thailand**

Roche Diagnostics (Thailand) Ltd  
18th Floor, Rasa Tower  
555 Phaholyothin Road  
10900 Ladyao,  
Chatuchak, Bangkok

**Thailand**

Phone: +66-2-937 0840  
Fax: +66-2-937 0850

**Turkey**

Roche Diagnostik Sistemleri Ticaret A.Ş.  
Gazeteciler Sitesi – Matbuat Sokak No 3  
34394 Esentepe – İstanbul

**Türkiye**

Phone: +90 212 213 32 80  
Fax: +90 212 216 73 51

**United Kingdom**

Roche Diagnostics Ltd  
Charles Avenue  
Burgess Hill, RH15 9RY

**United Kingdom**

Phone: +44 1444 256000

**Uruguay**

Roche Diagnostics  
Solferino 4096  
P.O. Box 6468  
11400 Montevideo

**Uruguay**

Phone: +598 2 613 7888  
Fax: +598 2 613 7879

**USA**

Roche Diagnostics Corporation  
9115 Hague Road  
Indianapolis, IN 46256

**USA**

Phone: +1 317 521 2000  
Fax: +1 317 521 4295

**Vietnam**

Roche Chemical Establishments Ltd  
Viet Nam Representative Office  
180C Hai Ba Trung Street, Da Kao Ward  
District 1, Ho Chi Minh City

**Viet Nam**

Phone: +84 8 822 3955  
Fax: +84 8 822 3953

 Roche Diagnostics GmbH  
D-68298 Mannheim  
Germany

[www.diavant.com](http://www.diavant.com)  
[www.roche-diagnostics.com](http://www.roche-diagnostics.com)