DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:
AUXILAB, S.L.

Declara que el equipo | Declare that the equipment | Declare que l'appareil:
MICROSCOPIOS | MICROSCOPES | MICROSCOPES

Código | Code | Code:
50135006, 50135007, 50148007 (digital)
50158007, 50158008
50122147, 50122148, 50122107, 50122108
50100106, 50100107, 50100127, 50100117, 50100204
50128000, 50128001, 50128002, 50128022
50116000

Series | Series | Séries:
135, 148, 158, 122, 100, 128, 116

Cumple las siguientes directivas | Meet the following directives | Accomplit les directives suivantes:

2006/95/CE Directiva sobre material eléctrico destinado a utilizarse con determinados límites de tensión (Directiva sobre Baja Tensión)
Directive relating to electrical equipment designed for use within certain voltage limits (Low Voltage Directive)
Directive relative aux matériels électriques à employer dans certaines limites de tension (Directive de basse tension)

Cumple las siguientes Normas: | Meet the following Standards | Accomplit les normes suivantes:

EN 61010-1:2010 Requisitos de seguridad de equipos eléctricos de medida, control y uso en laboratorio. Parte 1: Requisitos generales
Safety requirements for electrical equipment for measurement, control and laboratory use -- Part 1: General requirements
Règles de sécurité pour appareils électriques de mesurage, de régulation et de laboratoire -- Partie 1: Exigences générales

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ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

5. MAINTENANCE AND CLEANING

All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

Please follow the processes of use of this manual.

- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the microscope (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.
- ◆ Bulb and fuse replacement: The power wire must be disconnected
 - **Bulb replacement:** Loosen the knurled screw on the underside of microscope and open the panel to expose the bulb. Remove the old bulb after it becomes cool. (The bulb will become very hot when using or after using.) Don't touch the new bulb with finger, if there is a fingerprint and dirt, that will decrease the brightness and shorten the life of the bulb, wipe it with clean and soft cloth. Hold the new bulb with the same specification with clean gloves or gauze and vertically insert the pins to the jack. Close the panel and tighten the knurled screw with finger.
 - **Fuse replacement:** Open the fuse holder with a "↔" screwdriver in the back of the microscope. Remove the old fuse and install a new fuse with the same specification. Replace fuse holder in place.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminum, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ If there is dust on the lens, blow it by rubber ball blower, after that clean the lens gently with a soft brush pen, carefully wipe off oil or fingerprints on the lens surface with lens tissue or absorbent cotton moistened with a few organic solvent (ether and alcohol 7:3)
- ◆ Don't wipe the lens surface regularly, or else the lens will be scraped, reduce the quality of the transmission and imaging. Please keep the instrument clean.

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- ◆ Turn coarse & fine focusing knob to focus the specimen till you see clear image of specimen when observing the fixed eyepiece with eye. Then rotate the diopter adjusting ring, if the image is unclear when observing the another eyepiece with another eye, also still you see clear image of specimen (Remember your eye's diopter, so that you could use next time). When using two eyes to observe, hold the base of the prism and rotate them around the axis until there is only one field of view.

Security

- ◆ The microscope must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the microscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Follow all required **applicable laws about safety and health** (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ Working environment: Indoor temperature: 0 - 40°C and maximum relative humidity: 85 %.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the microscope (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 100-240V 50-60Hz $\pm 10\%$

◆ If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Setting up

- ◆ Turn on the power switch, adjust the brightness adjusting knob to make the brightness 70% of the full load.
- ◆ Place the slide to be viewed smoothly onto the stage, cover slip to face to the objective. Clamp specimen slide carefully with the movable spring clip.
- ◆ The magnitude of incident beam of light can be changed when adjusting the aperture diaphragm. The highest resolution of the objectives can reach when the fitted aperture diaphragm is adjusted. When the objectives are changed, in order to get the best resolution of the objective, please take off the eyepiece to observe the size of aperture diaphragm in the eyepiece tube. It is better to adjust aperture diaphragm till it is a little smaller than the aperture of the objective.

Note; aperture diaphragm is not for adjusting the brightness, the brightness is adjusted through brightness adjusting knob.

- ◆ Swing out the filter holder, according to user needs put filter in the filter holder and then backtrack.
- ◆ Turn the nosepiece when changing the objective 4x or 10x, and make sure the objective is shift in the light path until hear a "click".
- ◆ When adjusting the focus⁺ in order to prevent objective touch the slide, turn the coarse focusing knob until the slide is approximately 3 mm from the objective.
- ◆ Slowly turn the coarse focusing knob until a clear image is obtained, then use the fine focusing knob to enhance the observation of the specimen to its clearest image. If the magnification is increased, here you can obtain clear image under other higher magnification objectives with a little fine adjustment.
- ◆ When using objective 100x to observe, lift the condenser to the highest position, then drop a little oil on surface of objective 100x and cover slip. If there's air bulb in oil, it will influence observation. The 100x oil immersion objective and specimen should be wiped off with a piece of soft clean cloth or lens tissue to remove the cedar oil with xylene immediately after using.
- ◆ If you find to lift the mechanical stage too tension or loosen in use. Turn the tension adjusting ring. Coarse focusing knob would be tightening if it turns in the clockwise direction, on the other hand it would be loosen.
- ◆ Turn transversal and longitudinal direction adjustable knobs located just below the stage, the specimen may be moved to the center of the eyepiece's viewing field for observation.

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4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the microscope, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or provider immediately, so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied

- ◆ Please check that all the accessories are enclosed with the equipment:

- | | |
|-------------------------------|----------------------------------|
| - Cleaning kit | - Fuse 1A |
| - CD with software and driver | - Blue, green and yellow filters |
| - Cable Schuko | - Cover |
| - Immersion oil | - USB wire |
| - Scale microscope slide | |


Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT

- ◆ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Put the binocular head into the adapter, tighten the knurled screw with finger.
- ◆ Familiarize yourself the mechanical parts of your microscope. Gently operate each part by hand to see how it behaves and what result it produces.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The microscope is supplied with a Schuko standard wire
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire and to the other end to the microscope connector.

2. DESCRIPTION

1. Eyepiece
2. Camera
3. Binocular head
4. Tornillo de ajuste
5. Stand
6. Tension adjusting ring
7. Coarse focusing knob
8. Fine focusing knob
9. Transversal adjusting knob of stage
10. Longitudinal adjusting knob of stage
11. Collector
12. Condenser
13. The adjusting rod of aperture
14. Mechanical stage
15. Objective
16. Nosepiece

3. TECHNICAL SPECIFICATIONS

Code	50148007
Viewing head	Compensation free binocular head, inclined at 30° (55 – 75 mm)
Camera	3 MP
Eyepiece	WF10x
Nosepiece	Quadruple
Stage	Double layer mechanical stage (150 x 180 mm) with moving clamp
Condenser	N.A. 1,25 Abbe condenser with iris diaphragm & filter
Focusing	Coaxial coarse & fine focusing adjustment with rack and pinion mechanism. Fine focusing scale value 0.002 mm
Light source	White LED with intensity regulation
Collector	High brightness Kohler illumination
Objectives	Achromatic

Magnification	Numerical aperture	Working distance (mm)	Cover slip thickness (mm)
4 x	0.1	37.5	0.17
10 x	0.25	6.54	0.17
40 x (S)	0.65	0.63	0.17
100 x (S) Oil	1.25	0.195	0.17

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Thank you for choosing this equipment. We sincerely wish that you enjoy your microscope Zuzi 148. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the microscope Zuzi 148, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

INDEX OF LANGUAGES

Spanish 2 - 9
 English 10 - 17

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT..... 11
 2. DESCRIPTION..... 12
 3. TECHNICAL SPECIFICATIONS..... 12
 4. INSTALLATION / SETTING UP 13
 5. MAINTENANCE AND CLEANING 16
 ANNEX I: CE CERTIFICATE 20

1. USES OF THE INSTRUMENT

Microscope Zuzi model 148 is specially designed for clinical examination and teaching demonstrations in laboratories, colleges and medical fields. It is also ideal instrument for biological, bacteriological, pathological and pharmaceutical research. Because of higher capability with lower price and humanized configuration design, it also becomes the first choice of individual, schools and scientific research organization.