MS 10310030104

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ENGLISH

SUMARY

The aim of blood coagulation is to form a clot in the vascular lesion region. Besides the formation of this clot, which acts as a blocker to avoid hemorrhage, this pro-coagulant action should be restricted to the lesion region, that is, it cannot be widespread in the circulatory system. Thus, the moment that occurs the vascular lesion, it also starts the fibrinolytic system action and coagulation physiological inhibitors. The "cascade" model, proposed for Macfarlane and David & Ratnoff, in 1964, in order to explain the blood coagulation physiology, still used today for didactic reasons, divide the coagulation in extrinsic and intrinsic pathways which converge on the activation point of X factor (start of the final pathway or common). APTT is sensitive to deficiencies or inhibitors of extrinsic coagulation pathway, therefore,

evaluates the VIII, IX, XI and XII factors (extrinsic pathway) and X, V, II, and I (common pathway). APTT is the method of choice for monitoring intravenous heparin anticoagulant therapy and for screening of coagulationlupus inhibitor. It is also sensitive to fibrin/fibrinogen degradation. It is also used, together with PT, in the pre-surgery routines and hemorrhagic conditions.

PRINCIPLE OF METHOD
The capacity of blood to form a fibrin clot by way of the coagulation intrinsic pathway requires coagulation factors I, II, V, VII, VIII, IX, X, XI and XII, platelet lipids and calcium. **COAGULATION APTT** offers phospholipid (cephalin), since the plasma used is deficient in platelets; consequently it does not have platelet lipids, a contact activator (ellagic acid) and calcium which is removed by chelation through the use of sodium citrate anticoagulant. Therefore, APTT reagent is added to a normal citrated plasma to determine FXII activation which by cleavage will activate FXI.

Starting from this point, it is necessary the presence of calcium for the activation of other factors. The time for fibrin clot formation, after the adding of calcium chloride, is measured in seconds and represents the activated partial thromboplastin time (APTT). The APTT prolonged level is proportional to the level of heparin in the plasma.

APTT REAGENT
The APTT Reagent from **WAMA Diagnóstica** is a preparation of rabbit brain cephalin and ellagic acid activator with buffer, stabilizers and preservatives. The reagent is provided ready to use. Calcium chloride is provided separately and ready for use.

KIT PRESENTATION
REF 731040-C – 40 Determinations
REF 731080-C – 80 Determinations
REF 731200-C – 200 Determinations
1. Reagent/APTT: 1 x 4 ml
2. Calcium chloride (0.02M): 1 x 4 ml
3. Instructions for use
REF 731040-C – 80 Determinations
1. Reagent/APTT: 2 x 4 ml
2. Calcium chloride (0.02M): 2 x 4 ml
3. Instructions for use
REF 731200-C – 200 Determinations
1. Reagent/APTT: 5 x 4 ml
2. Calcium chloride (0.02M): 5 x 4 ml
3. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED
- Tube
- Automatic pipette
- Double boiler– 37°C
- Timer
- Container for material disposal

PREPARATION AND STABILITY
APTT Reagent (1): Ready for use. The reagent is stable to the expiry date shown on the label when stored in the original vial at 2 to 8°C.
Calcium Chloride (0.02M) (2): Ready for use. It is stable to the expiry date shown on the label when stored in the original vial at 2 to 8°C.
Obs.: Kit presents results after its first handling and it is stable up to expiration date if stored at 2-8°C.

SPECIMEN COLLECTION
1. **Blood collection using syringe method:** Draw venous blood into a plastic or siliconized syringe. Immediately transfer 9 ml of blood into a tube containing 1 ml of 3.2% or 3.8% sodium citrate solution alternatively (4.5 of blood to 0.5 of sodium citrate, keep the ratio of 9 parts blood to 1 part sodium citrate).
2. **Blood collection using an evacuated blood collection tube:** Draw venous blood into a vacuum tube containing 3.2% or 3.8% sodium citrate solution (blue cap). Insure that a full draw has been obtained since the ratio of 9 parts blood to 1 part citrate is critical. It is recommended that the collection be performed by using the first tube without anticoagulant (dry tube). In case of collection only for coagulation tests, collect two tubes with sodium citrate and use the second tube for coagulation tests. If multiple collections are needed, the sample for coagulation test should be obtained after the collection in tube without anticoagulant (dry tube), prior to collection in tube containing EDTA.
3. **Plasma preparation:** Mix well by inversion and centrifuge at 2.500xg (approximated 3000rpm) for 15 minutes soon after the blood collection. Unless the samples are to be processed immediately, transfer the plasma into a plastic tube. Plasma that is clearly hemolyzed or contains >10.000 platelets per cubic meter or red cells is not suitable for coagulation testing.
4. **Plasma storage:** Plasma samples may be stored at room temperature (18 a 26°C) for up to 2 hours; refrigerated (2 a 8°C) for up to 4 hours; frozen (- 20°C) for up to 2 months or at -70°C for up to 6 months. Plasma may be re-centrifuged prior to freezing to assure that all cells are removed. When testing frozen samples, they should be quickly thawed and tested immediately. The samples must not have any contact with glass. Do not incubate at 37°C for longer than 5 minutes to avoid the loss of factors V and VII. Samples stored for more than 4 hours at room temperature can prolong the APTT, since FVIII is labile.

PROCEDURE
This procedure is intended to use with manual testing. For semi-automated and automated systems, please refer to your instrument manual.
1. Pre-incubate Calcium chloride (0,02M) (2) at 37°C for at least 10 minutes.
2. Pipette 100 µl of test or control plasma into a test cuvette. Incubate at 37°C for 1 to 2 minutes.
3. Pipette 100 µl of the APTT reagent (1) to the cuvette containing plasma after mixing the suspension of the reagent by inversion.
4. Incubate the mixture at 37°C for 3 minutes.
5. Rapidly add 100 µl of the pre-incubated Calcium chloride (0,02M) (2) and simultaneously start the timer.
6. Slightly mixing by inversion the tube contents, record the clotting time in seconds.

CALCULATION OF RESULTS
It is recommended to calculate the mean clotting time of duplicate samples and controls. Differences between duplicate results should be less than 5%. Repeat the test if necessary. The APTT result may be reported in seconds from the adding of the Calcium chloride up to form a fibrin clot.

APTT and ratio of patient clotting time to mean normal clotting time (obtained from APTT from plasma or pooled normal plasma) may also be reported in the test result.

EXPECTED VALUES
APTT results are influenced by the method of clot detection and can vary from laboratory to laboratory. In general, APTT tests performed on a photo-optical coagulometer will give clotting times for normal plasma in the range of 24 to 38 seconds. However, each laboratory should establish a normal range using individuals representative of its patient population. A new range should be established with any change of instrumentation, blood collection techniques, or anticoagulant.

Therapeutic ranges for monitoring oral anticoagulation therapy will vary from laboratory to laboratory. Therefore, it is essential that each laboratory establish relevant APTT ranges for its

respective patient population.
The ratio patient APTT/Mean Normal APTT should be up to 1.3. Higher or equal values mean an APTT prolonged.

Abnormal results obtained with a plasma from a patient not on anticoagulant therapy may indicate a factor deficiency or the presence of an inhibitor. The result may also be due to the effects of certain drugs and medications. Additional procedures such as the TP tests and mixing studies using factor deficient plasma are usually required.

QUALITY CONTROL

The test performance should be monitored by using normal and abnormal coagulation control plasmas such as **Control Plasma Level 1 (COD 733001-C)**, **Control Plasma Level 2 (COD 733002-C)** and **Control Plasma Level 3 (COD 733003-C)** from **WAMA Diagnóstica**. Each laboratory should establish control ranges in order to determine the allowed variation on control plasma performance.

PERFORMANCE CHARACTERISTICS

Precision: Within-run precision was assessed using **Control Plasma Levels 1, 2 and 3** from **WAMA Diagnóstica**, on manual, semi-automated and fully automated instruments. The results are shown in the following table.

Within-run Precision Results

Sample	Manual	MLA Electra 900C (semi-automated)	Fully automated
Level 1	3,3%	2,0%	1,2%
Level 2	1,3%	2,2%	0,6%
Level 3	4,2%	2,3%	0,8%

Correlation: Correlation studies were performed against the APTT reagent of a competitor using a semi-automated and a fully automated coagulometer. The results are shown in the following table.

Correlation results

Instruments	Regression coefficient	Slope	Intercept
Semi-Automated	0,938	0,727	6,96
Automated	0,930	0,952	3,76

Factor Sensitivity: APTT times obtained with **APTT Coagulation** were evaluated on factor-deficient plasmas using a semi-automated coagulometer. The results are shown the following table.

Percent Factor	APTT-Clotting Times (in seconds)			
	Factor VIII	Factor IX	Factor XI	Factor XII
100%	28,4	25,7	25,8	25,9
60%	28,2	27,8	28,6	29,6
50%	28,7	28,3	30,2	31,5
40%	29,6	30,7	32,0	33,2
30%	33,5	31,1	34,3	35,1
20%	35,0	32,9	36,5	38,7
10%	41,4	36,7	42,4	45,7

These values should be used as guidelines. Each laboratory should establish factor sensitivity using their instruments, reagents and techniques.

Heparin Sensitivity: The heparin sensitivity of **APTT COAGULATION** from **Wama Diagnóstica** can be determined using normal plasma with known levels of heparin and by measuring the APTT. The following results were obtained using the **APTT COAGULATION** and a semi-automated coagulometer.

Heparin Concentration (units/ml)	APTT - Clotting Time (in seconds)
0,0	23,0
0,1	28,0
0,2	38,1
0,3	59,0
0,4	76,1
0,5	99,4

LIMITATIONS

- To prevent discrepant results ensure the blood to anticoagulant ratio is 9:1.
- Grossly lipemic, jaundiced and hemolyzed samples may produce erroneous APTT values.
- Delay in testing, difficulty in specimen collection and venipuncture above the site of a thrombus may result in falsely prolonged APTT results.
- APTT results can vary with anticoagulation therapy depending upon the type and dosage of anticoagulant, the route for administration and the time of administration of the last dose.

PRECAUTIONS AND WARNINGS

- Reagents for in vitro diagnostic use only.
- Do not use the reagents after the expiration date.
- Follow the instructions for use prior to testing.
- The reagents should be well mixed by inversion and sealed immediately after use.
- The samples and controls to be tested should be used simultaneously in order to maintain the same test conditions.
- Do not ingest. Avoid skin and eyes contact.
- Handle all samples as potentially infectious.
- Dispose in accordance with local regulations.
- Use Good Laboratory Practice at storing, handling and material disposal.

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds evidence that running, handling and storage were not properly follow ed.

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- Arkin, CF; Adcock, DM; Ernst, DJ et al. National Committee for Clinical Laboratory Standards. *Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays: Approved Guideline – Fourth Edition*. Wayne, PA: NCCLS; 2003, NCCLS document H21-A4.
- Human Blood Coagulation, Hemostasis and Thrombosis*, 3rd ed. R Biggs, CR Rizza, Editors. Blackwell Scientific Publications, London, 1984.
- Cole E, Hall ER, Wu KK, *Principles of Antithrombotic Therapy*. In *Wu KK, Thromboembolic Disorders*, PSG Publishing Co. Inc. Itteton, 1984, p 91.
- Harker LA, *Hemostasis Manual*, FA Davis Co, Philadelphia, 1974, p 62.
- Hougie C, *The Biochemistry of Blood Coagulation, In Laboratory Evaluation of Coagulation*, American Society of Clinical Pathologists Press, Chicago, 1982, p 2.
- Karnal, AH; Telfer, A.; Pruthi, RK. How to Interpret and Pursue an Abnormal Prothrombin Time, Activated Partial Thromboplastin Time, and Bleeding Time in Adults. *Mayo Clin Proc.*, 2007, 82(7): 864-873.
- Owen CA, Bowie EJM, Thompson JH. *The Diagnosis of Bleeding Disorders*, Little Brown and Company, Boston, 1975 p 110.
- Triplett DA, Harms CS, *Procedures for the Coagulation Laboratory*, American Society of Clinical Pathologists Press, Chicago, 1981, p 7.
- Triplett DA, *Heparin: Clinical use and Laboratory Monitoring*. In Triplett DA, *Laboratory Evaluation of Coagulation*, American Society of Clinical Pathologists Press, Chicago, 1982, p 272.
- Young DS, Pestaner LC, Gibberman V, *Effects of Drugs on Clinical Laboratory Tests*, Clin Chem, 1975, 21: 1D.

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