WIELAND
D E N T A L

Wieland Dental + Technik
GmbH & Co. KG
Schwenninger Straße 13
75179 Pforzheim, Germany
Fon +49 72 31/37 05 - 0
Fax +49 72 31/35 79 59
www.wieland-dental.de
info@wieland-dental.de

WIELAND
D E N T A L

WIELAND
D E N T A L

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso

Zenotec PMMA Cast

- Fräsröhling zur Herstellung rückstandsfrei
verbreimbarer Kunststoffmodelle
- Milling blank for the fabrication of
plastic models which burn without trace
- Bloc fraisable pour la confection de maquettes
en résine calcinant sans laisser de résidus
- Grezzo per la realizzazione di modelli in
materia sintetica totalmente calcinabile
- Pieza bruta a fresar para la confección de
modelos de resina acrílica totalmente quemables

EXPECT THE DIFFERENCE! BY WIELAND.

- Größte Auswahl an Materialien in höchster Qualität
- Widest range of materials in highest quality
- Grand choix de matériaux. Tous d'excellente qualité
- Vasta scelta di materiali di primissima qualità
- Gran variedad de materiales de máxima calidad

Deutsch

Gebruiksaanweisung Zenotec PMMA Cast

1 – Eigenschaften

Chem. Charakterisierung		100 % PMMA (Polymethylmethacrylat)
Dichte	g/cm ³	1,17
E-Modul	MPa	2400
Vicat Punkt	°C	102
Kugeldruckhärte	MPa	145

2 – Allgemeines

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English

Instructions for use Zenotec PMMA Cast

1 – Properties

Chemical description		100 % PMMA (polymethyl methacrylate)
Density	g/cm ³	1.17
Modulus of elasticity	MPa	2400
Vicat softening temperature	°C	102
Ball indentation hardness	MPa	145

2 – General information

Zenotec PMMA Cast can be machined with the Zenotec milling system. For this purpose the accompanying coding must be read into the CAM module. Only recommended milling tools for the milling system should be used for processing and milling. Please always make sure customers follow the instructions for use applicable to the product.

3 – Storage

Protect from direct sunlight. Protect from moisture.

4 – Areas of application

Zenotec PMMA Cast is not intended for intraoral use. The material can be used to make frameworks for crown and bridge models. These models can then be used in the Casting process, e.g. for the fabrication of dental restorations.

5 – Finishing

To separate the milled frameworks from the blank, use crosscut carbide burs suitable for use with acrylic materials, or suitable cutting discs.

5a – Spruing

The rules for designing and making sprues depend on the alloy or the press ceramic used. You can find detailed information on the alloys in the Casting Handbook available from Wieland Dental + Technik (in German). As a rule, acrylic models must be coated in wax.

5b – Investing

Use a ceramic ring liner. Mix the investment material under vacuum and invest in the usual way. For this, Wieland recommends the use of the investment materials WILAVEST quick, WILAVEST universal, IMAGINE PressX or PressX Zr. Speed heating to 580 °C is recommended, together with a holding time of 30 min. Then the mould should be heated to the end temperature appropriate for the alloy or press ceramic. Always follow the instructions for use pertaining to the investment material and the alloy or press ceramic being used.

6 – Disposal

Zenotec PMMA Cast is made from pure PMMA. Waste can be recycled along with other acrylics. It can also be disposed of through the domestic refuse disposal system. PMMA is waterinsoluble, inert and presents no hazard to the ground water.

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