# STORZ

**KARL STORZ — ENDOSKOPE**

20 1326 20-1

XENON 100 SCB

20 1326 20-1

XENON 100 SCB

20 1326 20-1

XENON 100 SCB

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Instruction manual

Manual de instrucciones

20 1326 20-1

XENON 100 SCB

20 1326 20-1

XENON 100 SCB

20 1326 20-1

XENON 100 SCB

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## **Important information for users of KARL STORZ devices**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the XENON 100 SCB cold light source. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully. Insert this manual in its protective wallet and keep it available in a convenient and conspicuous location near the equipment.

## **General warning**

Users of KARL STORZ devices should be thoroughly trained in their use and the applicable medical procedures. These devices and accessories should not be used until the user has fully familiarized himself/herself with the instructions for use, assembly and care. Instruction manuals should be carefully studied and be available to the surgical team during the procedure; it is essential to follow the instructions contained in the instruction manual, paying particular attention to the cautions and warnings.

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
**Controls, displays,  
connectors  
and their functions**

- ① Power switch (lit up)
- ② Lamp/unit warning and pilot light (lights up red when the lamp/electronics are defective and/or from 450 h...499 h lamp operating time - it flashes as of 500 h)
- ③ Warning light (lights up when light source is configured for general purpose applications (not to be used for videoscopes for gastroenterology))
- ④ Lamp On/Off button, LED lights up when 'On' in the standard configuration / video gastroenterology
- ⑤ Up-key for light intensity
- ⑥ Down-key for light intensity
- ⑦ Display for light intensity setting  
green = standard mode / standard range  
yellow = standard mode / boost range  
red = intensive range: applicable only in general purpose configuration
- ⑧ Insufflation pump On/Off button, LED lights up when 'On'
- ⑨ Display for pump speed setting (three stages)
- ⑩ Down-key for pump speed
- ⑪ Up-key for pump speed
- ⑫ Receptacle for the fiberoptic light cable adaptor and insufflation (videoscope for gastroenterology)
- ⑬ Quick-release coupling mechanisms
- ⑭ SCB input/output connector
- ⑮ Potential equalization connector
- ⑯ Power connector and fuse holder
- ⑰ Holder for water bottle (on left of housing)

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General warning	.II
Images of the equipment	.III
controls, displays, connectors and their functions	.V
Symbols employed	.VI
<b>General information</b>	<b>.2</b>
Description of the unit	.2
<b>Safety instructions</b>	<b>.4</b>
Warnings and cautions	.4
Normal use	.8
User qualification	.9
Safety precautions at the site of installation	.10
Safety precautions when operating the unit.	.10
<b>Installation and operating instructions</b>	<b>.11</b>
Unpacking the equipment	.11
Basic equipment	.11
Connecting the unit	.12
Operating the cold light source	.15
Light	.16
Pump	.18
Increased light output	.19
General purpose application	.19
<b>Maintenance</b>	<b>.22</b>
Lamp replacement	.22
Fuse replacement	.26
General purpose configuration	.27
Cleaning, disinfection and care of the unit	.30
Maintenance	.31
Disposal	.31
Repair program	.32
Limitation of liability	.33
Manufacturer's warranty	.33
<b>Technical description</b>	<b>.34</b>
Troubleshooting	.34
Technical data	.36
Technical documentation	.38
<b>Spare parts, recomm. accessories</b>	<b>39</b>
List of spare parts/accessories	.39
<b>Appendix</b>	<b>.41</b>
Electromagnetic Compatibility (EMC) Information	.41

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## Description of the unit

The cold light source XENON 100 SCB is a powerful light system for endoscopic use, particularly in combination with videoscopes for gastroenterology. In a special configuration and with the aid of an adaptor (Art. no.: 495KS), the light source is capable of satisfying the requirements of even demanding general purpose applications. With a color temperature of 6000 K, the 100 watt xenon lamp has the same color temperature as sunlight and therefore provides an outstanding bright illumination.

The SCB feature provides remote control of the light source, system status information and easy configuration upon initialization. Modern design guarantees easy operation and maintenance.

After the unit has been installed properly and the unit and the lamp have been switched on (via the power switch and the 'lamp on/off-key'), it takes a few seconds until the lamp unfolds its full potential. By inserting the plug of the video-scope (for gastroenterology) into the light source, the antiglare device gives way and light is coupled into the light guide. The light intensity can be adjusted via the up/down keys of the intensity control either while the lamp is on or off (preselection). For standard operation with a KARL STORZ video-scope for gastroenterology, the light intensity can be adjusted from 5% to 100% without any harm to the instrument. If required, but only for short periods of time and under the conditions of increased attention to overheating of instrument or harm to the patient due to high intensity levels, the intensity may be boosted beyond 100%. A locking scheme prevents unintended activation of this mode.

The expiry of the normal lamp service interval, which amounts to 500 hours, is signalled via a flashing indicator ② on the front panel. This indicator is lit permanently 50 hours before that. A pump system for insufflation and/or water irrigation\* is integrated into the system. The pump is started and stopped via the on/off button ⑧. The pump speed is adjustable over three levels via the +/- ⑩ / ⑪ :

Level 1+2: Normal range

Level 3: For 'difficult' conditions and severe lens soiling (see also p. 18) - short periods of time only!

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**Function of the warning/pilot light on the front panel**

The pilot light is highlighted permanently when the light source is turned on and

- the lamp operating hours are between 450 h and 500 h as described above
- the ignition fails due to a defective or improperly installed lamp
- the light source is overheated due to obstructed ventilation in-/outlets or a defective fan.

The pilot light flashes if the lamp has been in service for 500 h and beyond that. Reset the signal by

- either replacing the lamp and carrying out the appropriate measures (please see the section 'Lamp replacement' page 22, for resetting please see page 25)
- or eliminating the reason(s) for ignition failure or overheating of the device.

**Property rights**

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

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### Warnings and cautions

Please read this instruction manual carefully and pay attention to it. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, read the text that follows carefully to ensure the safe and effective operation of the device.

#### Definitions

**Warning:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**Caution:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

**Note:** A Note indicates additional information about the safe operation of the unit.

**Warning:** Read this instruction manual thoroughly and be familiar with its contents prior to using the unit. Read the section on safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Before using the unit on the patient, it is imperative that you be acquainted with how the unit operates and how it is controlled.

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**WARNINGS:**

**Warning:** The instructions and interface specifications for medical devices used in combination must be observed precisely. Additional equipment which is connected to one of the interfaces of the unit must be proven to comply with the relevant IEC specifications (e.g., IEC 60601-1 for medical electrical equipment and IEC 60601-2-18 for endoscopic units). Furthermore, all the configurations must satisfy the applicable version of the system standard IEC 60601-1-1. Anyone who connects additional units becomes the system configurator and is thus responsible for ensuring that the applicable version of the system standard IEC 60601-1-1 is observed.

**Warning:** Combinations of medical devices are only assured to be safe if the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

**Warning:** Test this unit prior to each surgical procedure to ensure that it works properly.

**Warning:** Do not use in the presence of flammable anesthetics. There is a risk of explosion.

**Warning:** Danger of electric shock! Do not open the unit! Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

**Warning:** Danger of burns! The lamp and the housing may be very hot.

**Warning:** The insufflation pump serves to insufflate air and/or to supply water (e.g., for lens cleaning) in combination with the appropriate videoscopes.

**Warning:** In order to avoid embolisms, only carry out insufflations within physiological limits. **Important:** When using level 3 for a long period of time, it is essential to remember to extract the air once the procedure has been completed.

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**Warning:** Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**Warning:** The electrical installations in the operating room or examination room in which the unit is connected and operated must comply with the applicable IEC standards.

**Warning:** Keep unit out of reach of patients.

**Warning:** Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

**Warning:** Danger of glare! Never look into the light outlet of the light adapter or the air vent of the unit when the light system is switched on.

**Warning:** Do not handle the XENON lamp without wearing protective gear - at the very least gloves and safety goggles. Lamp explosion is an inherent danger and can take place anytime without a noticeable reason.

**Caution:** Always operate the cold light source with the minimum brightness setting necessary to achieve optimum illumination.

**Note:** If the cold light source is switched on and off frequently, this increases wear to the xenon lamp. Therefore, during short interruptions, the brightness should be reduced.

**Caution:** Only operate the unit at the line voltage stated on the manufacturer's identification plate.

**Caution:** When replacing fuses, only use fuses of the specified type and interrupt rating.

**Caution:** Do not store liquids on or above the unit. Avoid allowing fluids to enter the housing.

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**Videoscope**  
(for gastroenterology)

**Warning:** Danger of glare! Never look into the distal end of a videoscope.

**Warning:** High-intensity light radiated through the endoscope can give rise to high temperatures at the tip of the endoscope. For standard operation, do not exceed intensity beyond the green section of the intensity display. A higher intensity must only be used for a short time and only under direct visual control.

**Warning:** To avoid burns, do not place the videoscope onto or near patients.

**Warning:** When configured for general purpose applications (warning signal is on and the light source beeps when the power is switched on) do not use the light source in combination with videoscopes for gastroenterology and in corresponding applications.

**Caution:** Do not place a connected videoscope on combustible materials such as textiles (surgical drapes) or in the direct vicinity of swabs soaked with combustible fluids, such as disinfectant. These could catch fire from heat development at the distal end of the videoscope.

**Note:** Any damage to the unit resulting from incorrect operation or misuse is not covered by the manufacturer's warranty.

**Note:** Lamps do not have an indefinite service life, so the lamp may fail. Although it is highly unlikely that a lamp will fail during operation, but it cannot be ruled out. As a result the following applies: a) For purely diagnostic use, we recommend that you have a replacement lamp ready. b) For therapeutic endoscopic interventions we strongly recommend always having a suitable, ready-for-use replacement light source available.

**Caution:** Do not touch the lamp with bare hands. Use a protective cover or gloves! Touching the lamp with your bare hands can lead to problems during operation.

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Before using the unit on patients for the first time, it is imperative that you be acquainted with how the unit operates and how it is controlled.

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## Normal use

The XENON 100 SCB light source is to be used in combination with **KARL STORZ videoscopes for gastroenterology**. A particular feature of these endoscopes is an integrated connector for light insufflation and/or irrigation. The XENON 100 SCB supports this via an appropriate socket and an integrated air pump. In addition to this, the light source offers the following features:

- SCB connectivity/functionality for application inside a SCB installation;
- an optional configuration as a standard light source for general purpose application (alongside gastroenterology with the respective videoscopes).

The XENON 100 SCB light source shall only be operated and used in compliance with the following terms and conditions:

- use restricted to human medical applications (especially gastroenterology);
- use only in combination with KARL STORZ videoscopes for gastroenterology (gastrosopes);
- or with the aid of an adaptor (Art. no. 495KS) in combination with other KARL STORZ endoscopes (flexible and rigid);
- use in rooms for medical applications (operations);
- use only under standard environmental conditions;
- use only while connected to public mains, 100 - 240 VAC, 50/60 Hz, via hospital grade power sockets installed acc. to the applicable national standards and guidelines;
- operation only in upright position, free convection provided;

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- use strictly prohibited in the presence of flammable fluids or gases or in the presence of increased oxygen concentration;
- use strictly prohibited in the presence of increased dust concentration;
- use prohibited under the conditions of increased air pressure (diving stations) or low air pressure i.e., at high altitudes (> 3000 m above sea level);
- contact prohibited with aggressive liquids (acids, leaches, solvents) or aggressive gases;
- use prohibited outdoors or under the conditions of a condensing atmosphere (damp rooms);
- installation and use prohibited in vehicles (land, water and air) or in mobile compartments (containers). Mobile videocarts are excluded from this;
- installation and use prohibited close to heaters or radiators or other locations where the device is exposed to intensive heat or light;
- installation prohibited in the presence of intense electromagnetic fields (electrical engines, HF generators, X-ray devices, MRI scanners, etc.);
- installation and use prohibited in the presence of nuclear radiation (devices for nuclear diagnosis and treatment, etc.).

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
**User qualification**

The cold light source XENON 100 SCB may only be used by physicians and medical assistants who have a corresponding specialized qualification, proven by an officially recognized certificate, and who have been instructed thoroughly in the use of the unit.

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### Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with the applicable national regulations.

The unit is equipped with a connector for attaching a ground line. It should be connected in compliance with the applicable national regulations.

It is not intended for use in areas subject to explosion hazards. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

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### Safety precautions when operating the unit

It is the user's responsibility to make sure that the unit is safe and operating properly before using it.

During treatment using the cold light source XENON 100 SCB, the patient must be treated and kept under observation with the usual medical care. This includes, in particular, monitoring sterile application conditions where required by the type of intervention.

The XENON 100 SCB light source is fail-safe up to a very high degree. Nevertheless, it is mandatory to have a spare light source immediately available as a backup system at the site of operation in case the type of operation requires a 100% fail-safe situation as regards the light source. If so, the spare light source has to satisfy all requirements to handle the exception (startup time, connectivity, optical power, etc.). These requirements depend on the application. Therefore, it is the responsibility of the user to define these requirements and make sure that they are met by the installed backup system.

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### Unpacking the unit

Carefully unpack the cold light source and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

### Basic equipment

- 1 Cold light source XENON 100 SCB (201326 20 -1)
- 1 Power cord 400 A
- 1 Water bottle 13992 BS (incl. cap)
- 1 Irrigation bottle stand 201000 31
- 1 Instruction manual 96206137D
- 1 Wallet for instruction manual
- 1 Operator's quick reference guide (self-adhesive) 96206137X

For general use

- adaptor 495KS (optional)
- light cable 495NA (optional)

Initial configuration: When shipped, the XENON 100 SCB light source is configured for standard applications (video gastroenterology). Until further notice, the following descriptions and instructions relate to this configuration.

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
### Connecting the unit

**Note:** The cold light source XENON 100 SCB may only be used in medical rooms with electrical installations conforming to the applicable national regulations

**Warning:** This unit is not intended for use in areas subject to explosion hazards. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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### Fitting/replacing a label strip on the front panel if necessary

**Note:** The appropriate controls/indicators for the unit are designated on the label strip. The label strips are currently available in 6 languages. The 'English' label strip is pre-fitted.

Undo the quick-release coupling mechanisms on the back of the unit (turn 90° counterclockwise).

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Push the unit cover backward (see also the section 'lamp replacement' page 22) and remove the front side cover which is locked into place on the left-hand side.

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Remove the pre-fitted 'English' label strips and slide the selected label strips between the plastic and the metal front.

Fold the strips along the marked line so that the side cover can be refitted without any obstruction.

Reattach the side cover.

Close the unit cover again and secure it with the quick-release coupling mechanisms.

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Place the unit on a flat surface.

**! Caution:** Do not block any air vents. Danger of overheating.  
Do not place the unit directly against the wall since the air vents are on the rear and the right-hand side of the unit (viewed from front) (minimum clearance of 10 cm). When overheated, the unit switches off automatically.

**! Warning:** Insufficient ventilation of the unit can cause an internal build-up of heat, which leads to safety shut-off of the lamp (the pilot light ② lights up; for troubleshooting, see p. 34/35).

**Note:** Cool-down time: approx. 10 to 15 minutes.

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The unit is equipped with a connector for attaching an equipotential wire (redundant grounding wire) line ⑮. The unit's equipotential line should be installed by a qualified electrician.

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Connect the power cord, insert the power plug into the power socket ⑯, up to the stop.

**Warning:** Only insert the power cord into and remove it from electrical outlets located outside areas subject to explosion hazards.

**Caution (US / CAN only):** If the power distribution provides an output of 240 V~, a center-tapped power installation is required. Installation shall follow the National Electric Code (NEC).

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Connect the water bottle for irrigation/insufflation to the XENON 100 SCB (see also the instruction manual for the videoscope).

**Warning:** Danger of contamination! To avoid reverse flow of water into the unit, the water bottle must be fixed in the stand provided. The filled water bottle must not be placed upside down or positioned at a higher level than the unit.

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IMPORTANT: To fill/empty the water bottle, remove the bottle from the light source and turn the cap. For further instructions, see the instruction manual for the videoscope in use.

**Note:** Only use water bottles that are 2/3 full (1/3 of the volume is required as reservoir for the compressed air). For further instructions regarding irrigation/insufflation applications, see the instruction manual for the videoscope in use.

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Push in the connector (light and insufflation) of the videoscope in the corresponding socket ⑫ up to the stop. Only pick up the videoscope by its handle. Never pull on the cord.

**Warning:** *Danger of glare. Never look into the light outlet of the light adaptor, the air vent of the XENON 100 or the distal end of the videoscope when the light system is switched on.*

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Connect the tube from the water bottle (integrated into the cap) to the endoscope supply connector (irrigation bottle connection).

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**Warning:** *Danger of burns. Contact with the distal end of a videoscope connected to a cold light source may cause burns.*

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### Operating the cold light source

Switch on the unit (power switch lights up).

**Note:** *The lamp and the insufflation pump are switched off when the XENON 100 SCB is switched on.*

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### Light

#### Switch on light (lamp)

To switch the lamp on or off, press button ④. The status of the lamp is indicated by the LED button. In addition, the plastic ring of the receptacle of the videoscope is illuminated when the lamp is on.

**Note:** After the lamp has been switched off, the integrated fan runs for approx. 30 s at high speed, then the fan is automatically switched to a lower speed.

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### Changing the luminous intensity

The light intensity is set by the up/down buttons ⑤ and ⑥. The intensity setting is displayed on the strip display ⑦. Tuning the intensity setting is independent of lamp on/off status. When the lamp is off, the light intensity can be preselected.

**Note:** The intensity can be adjusted from 5% to 100% either in steps by pressing and releasing the + and - button for each step or by holding the button depressed to employ the repeat function. The increase of the intensity stops at 100% in any case.

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### Standard operation of the videoscope:

To prevent damage to the videoscope do not set the intensity to higher levels than required. On no account should 100% be exceeded for any long period of time.

**Warning:** Danger of burns! Regardless of the light intensity, always keep the tip of endoscope clear from tissue. Burns may be induced either by the high light intensity due to direct contact or by the overheated instrument tip. Never leave your instrument unattended when the light is on.

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**Boosted illumination (with videoscope):**

If required, the luminous intensity may be boosted beyond 100% for short periods of time. To enter this mode, press the up button ⑤ for two seconds while the current setting is already 100% and the lamp is switched on. The setting then skips one increment (first yellow LED in the bargraph is on) into the boost range. At the same time a beep is audible for 1/2 second. To increase the intensity further, release the up button ⑤ and press it again, either repeatedly for step by step adjustment or sustained for continued increase. In each case, this will stop at the upper end of the yellow indicated boost range at 130% luminous intensity.

In the reverse direction, there is no automatic stop at the 100% level: repeated decrease of intensity will not be interrupted at this stage.

**Warning:** Danger of burns and/or damage to the videoscope

*It is highly recommended to use boosted luminous intensity only for short periods of time. Be aware of the increased risk of harm to the patient due to the high intensity of the illumination. In addition, the instrument may become overheated (peripheral parts) beyond permissible values or even damaged by overstressed internal parts. Use boosted illumination only while under direct visual control.*

When the intensity is set above 100% and the lamp is switched off via button ④, the intensity is automatically reduced to 100%. This avoids the lamp's being restarted at unintentionally high temperatures. Intensities higher than 100% are thus only possible when the lamp is switched on and cannot be preselected.

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**Pump:**

**Turning the pump on/off:**

The internal pump for insufflation and/or irrigation is turned on/off by button (8). The pump's operating status is indicated by the button's LED status.

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**Adjusting the pump speed (flow):**

The pump speed corresponds first and foremost to the flow. Its correlation with the pressure depends on the system-resistances (instrument, hoses and connectors) and the involved volumes (water-bottle, stomach) and losses via leaks. Flow of irrigation in turn depends on the pressure. Three levels of pump speed can be selected. Stages 1 and 2 are for normal use, stage 3 is for 'difficult' conditions and severe lens soiling.

Adjustment is performed using the + button (10) and - button (11). In contrast to the luminous intensity adjustment, an autorepeat function is not implemented by pressing the button. The pump speed can also be preselected when the pump is off. The current speed setting is displayed on the strip display (9). Switching the power on resets the speed to level 2.

**Warning:** In order to avoid embolisms, pump level 3 should only be used for short periods of time (i.e., in case of difficult lens soiling).

**Note:** The insufflation pump can be used for air insufflation or water irrigation. Details depend on the videoscope in use. Please read the instruction manual for the videoscope carefully.

The following text has been taken as an example from the KARL STORZ instruction manual for the gastroscope:

“Water supply and/or lens cleaning are activated by pressing the air/water valve on the videoscope. The various valve positions mean there are different possible combinations.”

IMPORTANT: Only use water bottles that are 2/3 full (1/3 of the volume is required as a reservoir for the compressed air).

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**Note:** The insufflation pump is only operational in the standard configuration / video gastroenterology. In general purpose configuration, the insufflation pump is shut off permanently because no videoscopes for gastroenterology should be used then. However, the air outlet of the light source is not accessible without the integrated light plug.

### Increased light output

#### General purpose mode

For general purpose applications with instruments other than videoscopes for gastroenterology the light source can be configured accordingly. This does not just refer to other videoscopes, but also to all standard endoscopes. In this case, the light intensity can be adjusted to levels far beyond what is suitable for the videoscopes for gastroenterology. This mode is therefore indicated by a permanently illuminated warning symbol (when on). When switching the power-on, this symbol flashes for two seconds and a beep is audible.

**Warning:** Danger of burns and damage to devices. Never use the light source in combination with videoscopes for gastroenterology and in related applications when in this general purpose configuration. The patient may be harmed and / or the instrument may be damaged in a matter of seconds.

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**Changing the light intensity**

The light intensity in the general purpose configuration is set as described below as of page 27. According to the higher intensity levels, the bargraph display comprises a separate red section.

**Note:** The intensity can be adjusted over the full range. There is no stop at the 100% level as in the standard mode.

The intensity display is non-linear. With respect to an equal display in different modes (i.e., 100% is the same light intensity in both configurations) and the limited space available for the display on the front panel, the precision of the display decreases at higher settings. Therefore the scale ends at 260%.

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**Connecting standard KARL STORZ endoscopes (general purpose):**

Standard KARL STORZ endoscopes and XENON 100 SCB light cables can be connected to the light source with the aid of the adaptor 495KS. Instead of the plug of a videoscope for gastroenterology, this adaptor is inserted into the socket of the light source. All light cables fit into this adaptor, except the 495NCS.

**Note:** The use of the adaptor 495KS will not automatically switch the light source into the general purpose configuration. If this is not done according to instruction at a subsequent point (as of page 27), no higher levels of intensity can be achieved other than those for the standard operation with videoscopes for gastroenterology. The opposite is true depending on the losses cause by the adaptor. For demanding general purpose applications, configure the light source accordingly as described below. This will not only compensate the losses of the adaptor, but this increase of intensity will contribute to meeting the application's requirements.

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**General purpose applications:**

The XENON 100 SCB light source is capable of satisfying the requirements of a wide range of applications. Which luminous intensity setting is appropriate depends on the individual application, instrument and light cable used. Always keep the intensity as low as possible.

**Warning:** Danger of burns and/or defects of the instrument/endoscope

*It is highly recommended to use the light source in general purpose configuration with the same respect and attention to safety precautions as any other high power XENON light source. Damage to instruments or injury to the patient due to increased light intensity is always a risk of such light sources.*

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**Lamp replacement**

Switch off the unit, unplug the power cable and, if present, the SCB cables and the potential equalization cable. If hot, please let the unit cool down before you continue.

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Open the housing cover. For this purpose, first undo the 2 quick-release coupling mechanisms on the back of the unit by turning them 90° counterclockwise.

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**Note:** Slide back the housing cover just far enough so that the lamp compartment is easily accessible.

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**Opening the lamp compartment lid**

**!** **Warning:** Danger of explosion. The lamp is filled with xenon gas under high pressure (approx. 60 bar). Therefore the lamp must be handled with care and should only be replaced by electrotechnically qualified service personnel. High mechanical loads can cause the lamp to break. Protective eyewear and clothing, at least protective goggles and gloves, must be worn while handling the lamp. Explosion is a **permanent** hazard and can happen anytime without noticeable reason. Pay attention to the the accompanying instruction leaflet for the spare lamp. Turn the quick-release coupling mechanisms 90° counterclockwise.

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Lift up the lid of the lamp compartment.

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**Warning:** *Danger of burns. The lamp and its surroundings may be hot. Please wait until the lamp has cooled down.*

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

Release the locking mechanism of the lamp connector by pressing on both sides and removing the connector evenly.

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Using both index fingers, push the lamp holder sharply approximately 2 mm forward (see direction of arrow in photo) until it is released from the lock . Immediately after unlocking, use the index fingers to lift the lamp holder slightly so that the spring force does not cause it to return to its original position automatically.

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Lift up the lamp holder carefully and hold the lamp in place with your second hand.

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
Remove the old lamp.

Dispose of the old lamp in accordance with the applicable regulations and pay attention to the instructions of the leaflet inside the package of the spare lamp.

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
**!** **Warning:** Do not touch the new lamp (Art. no. 20 1325 26) with bare hands. Use protective gloves! This protects against fingerprints, which can cause problems during operation and from injury in case of an accidental lamp explosion.

Insert the new lamp. Important: Pay attention to the position of the pin and the slot. The lamp can only be installed in this position!

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Place the holder carefully over the lamp.

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Close the lamp holder by pushing it down until it locks. It slides forward and back as it is a latch. Support this movement if necessary.

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
Insert the lamp connector into the lamp socket in the unit (only possible in one position) until the retaining hooks click into place.

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Move the slide switch into the alternate position as shown. Do not reverse this until the next lamp replacement.

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Close the lid of the lamp compartment and secure with the quick-release coupling mechanism (turn 90° clockwise).

**Note:** Make sure that the clip of the quick-release coupling mechanism is not sticking up. If necessary, push it down.

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Close the housing again and secure with quick-release coupling mechanism (turn 90° clockwise).

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**Fuse replacement**

Switch off the unit and remove the power plug.

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
Remove the fuse holder ⑩ with a screwdriver or other tool.

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Insert new fuses.

**⚠ Caution:** Only use fuses of the specified type and interrupt rating.

	110...120 V~ / 220...240 V~
Line fuse	2 x T 3.15 AL 250V

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Refit the fuse holder ⑩.  
Reconnect the power plug.  
Test the instrument for proper operation.

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**General purpose configuration**

If you want to use the XENON 100 SCB in applications apart from gastroenterology and using other instruments than the corresponding KARL STORZ videoscopes dedicated for this purpose and the light-source does not satisfy the requirements in terms of brightness, you may enhance the intensity by configuring the light source for general purpose application (see below).

**Warning:** It is highly recommended to use the general purpose configuration only when other options fail. Never use the general purpose configuration in combination with videoscopes for gastroenterology.

**Warning:** If your application demands enhanced brightness and your equipment is capable of withstanding the additional stress and you configure the light source for general purpose applications, undo the modification as soon as you have finished. If the modification is intended to be permanent or it is likely that other users are involved while the unit is configured for general purpose, make sure that everybody is notified about the modification. Attach a warning sticker.

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Switch off the unit and unplug all cables. If hot, let the unit cool down before you continue. Then open the unit as shown for lamp replacement. Do not open the lamp compartment lid.

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The lamp compartment is accessible from the right-hand side of the device (viewed from the front). Pull out the attenuator slide as shown until it stops. Do not touch the lamp or other optical parts.

Then mount the cover again and lock it with quick-release coupling mechanisms on the back of the unit.

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Turn on the light source.

It will not start as usual. Instead, the light source requires the user to confirm the changed configuration. This is signalled by a flashing warning symbol ③.

- Press button ⑤ (light intensity +) twice for confirmation.

The light source then moves into normal operation (general purpose mode)

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### General purpose configuration

When in general purpose configuration, several alarms are active.

When switching the power on:

- the warning symbol ③ flashes for two seconds;
- a beep is audible at the same time.

After this, the warning symbol ③ remains highlighted while the Xenon 100 SCB is in operation.

The first time button ⑤ (intensity +) is pressed after the unit is turned on, a second audible beep is heard

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In general purpose mode, the lamp is on all the time (when the device is switched on).

The system comes up with the minimum intensity setting. Light intensity can be adjusted by buttons ⑤ and ⑥ up to the absolute maximum of the device (red LEDs).

In general purpose configuration, button ④ acts as standby button which is different from its lamp on/off function in standard configuration / video gastroenterology. Due to the fact that the lamp is on all the time, the dimmer occludes the light channel completely (no light output) when standby is on. According to this, the status LED of button ④ is on when standby is active and vice versa. When standby is off, light is emitted acc. to the intensity setting. Another difference is the intensity display ⑦: When standby is active, all the LEDs of the intensity display are off.

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
**Warning:** In general purpose configuration, the achievable light intensity is far beyond the maximum (even in boost range) in standard mode. Therefore, heat generated by the emitted light, either in endoscopic equipment or the illuminated tissue might cause severe damage.

In contrast to standard configuration / video gastroenterology, no notch function prevents from adjusting the light intensity to such high levels. Be careful.

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**Return into standard configuration / video gastroenterology:**

If you no longer need general purpose configuration, reset the light source to standard configuration / video gastroenterology

- Disconnect the light source from all lines and open it as described above.
- Replace the fixed attenuator (slide) in the active position (redo the actions for general purpose configuration; see above).
- Close the light source, reconnect all lines and turn on.

The light source then requires the user for confirmation of the new configuration by flashing the warning symbol c.

- To confirm, press the - button f for the light intensity twice.

The light source will then return to normal operation (standard mode / video gastroenterology).

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### Cleaning, disinfection and care of the unit

 **Warning:** Always pull out the power plug before cleaning! Make sure the lamp system has cooled down.

 **Warning:** Avoid allowing fluids to enter the housing of the XENON 100 SCB. When preparing the videoscope or its channel, it must always be disconnected from the XENON 100 SCB.

The exterior of the unit can be surface disinfected with a disposable cloth soaked in disinfectant (preferably alcohol-free).

#### Preparation of the water bottle 13992 BS

Empty the water bottle and then carefully clean and disinfect the container, the lid and the tube, before rinsing thoroughly and drying. We recommend sterilization using the fractionated prevacuum procedure at 134°C or EO/FO gas.

Further preparation instructions can be found in the instruction manual for the videoscope in use.

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
## Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. The warnings, cautions and notes contained in the safety instructions must be observed. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3583 (English version).

**Note:** To ensure safe operation of the system, you should perform the following procedures every 12 months:

1. Check leakage current from chassis < 100  $\mu$ A.
2. Check leakage current from videoscope < 100  $\mu$ A.
3. Check ground impedance < 0.1 ohms.
4. Check power consumption  $\leq$  rated power.

## Disposal

To dispose of this medical product, no special measures are necessary. National laws and regulations must be observed. At the end of its useful operating life, dispose of the unit as electronic scrap.

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**Disposal**

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

**Repair program**

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs directly to

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

**Important information**

In order to prevent infection, it is strictly **prohibited** to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operation reliability, and performance of this equipment only subject to the following conditions:

- assembly, system expansions, readjustments, modifications, or repairs have been performed by persons authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its instruction manual at all times.

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### Manufacturer's warranty

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen, Germany

**Note:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its reliable performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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### General

In the event of the unit's malfunctioning, the troubleshooting list provided in this chapter will assist the user in finding the possible cause. The guide explains how to rectify the various malfunctions.

 **Warning:** Always unplug the unit before all maintenance work!

 **Warning:** Maintenance work that is not mentioned in this chapter must only be carried out by KARL STORZ authorized technical personnel.

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Complete failure of the unit.

#### Possible causes:

- 1 No power from the power line.
- 2 Defective line fuse.

#### Remedy:

- 1 Check that there is electricity to the wall outlet.
- 2 Change fuses as described in the instruction manual. Make sure you use the correct type of fuse.

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#### Symptom:

No light is emitted, red pilot light is on.

#### Possible causes:

- 1 Lamp is defective.
- 2 Electronics are defective.
- 3 Overheating of the XENON 100 SCB (lamp has switched off due to overheating of the unit caused by covered air vent).

#### Remedy:

- 1 Have the lamp replaced by qualified service personnel.
- 2 Notify Service.
- 3 Remove objects from in front of the air vents. Switch off the Xenon 100 SCB, let it cool down (approx. 10–15 minutes) and switch it back on. Ensure that it is sufficiently ventilated.

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*Symptom:*

Too little light – red pilot light ② is on/flashing.

*Possible causes:*

- 1 The pilot light ② lights up red when the lamp has been operated for between 450 h and 499 h (450 h corresponds to approx. 90% of the guaranteed service life of the lamp).
- 2 The lamp is worn out (the pilot light ② flashes red when the lamp has been operated for more than 500 h).

*Remedy:*

- 1 Increase the luminous intensity with the regulator ⑤ (see warning on p. 7 'higher luminous intensity').
- 2 Have the lamp replaced by qualified service personnel.

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*Symptom:*

Too little light – red pilot light ② is not on / not flashing.

*Possible causes:*

- 1 The lamp is worn out (prove warranty claim using 'lamp log book').
- 2 Luminous intensity is set too low.
- 3 End faces of the integrated light cable of the videoscope are soiled.

*Remedy:*

- 1 Have the lamp replaced by qualified service personnel.
- 2 Increase the luminous intensity with button + ⑤.
- 3 Clean end faces of the integrated light cable and light outlet ends of the endoscope.

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
*Symptom:*

Pump on – no pressure/pressure too low.

*Possible causes:*

- 1 Insufflation pump worn out.
- 2 Kink in internal/external cable.
- 3 Air filter blocked/soiled.

*Remedy:*

- 1 Notify Service.
- 2 Correct the 'kink' (if it is internal, then call in qualified service personnel).
- 3 Unscrew the air filter (inside lamp compartment) and replace it (Art. no. 20132580).

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**Technical description**

**Technical data**

XENON 100 SCB	<b>20132620-1</b>
Line voltage	110...120/220...240 V~ (US / CAN only: if connected to 240 V~ center-tapped installation required)
Line frequency	50/60 Hz
Lamp voltage	12...14 VDC
Lamp power	100 W
Pump delivery rate	max. 5 l/min *
Pump pressure	max. 0.51 bar *
Power consumption	250 VA
Line fuses	2 x T 3.15 AL250V
Dimensions	320 mm x 110 mm x 380 mm (w x h x d)
Weight	5.4 kg
Operating temperature	10°C...40°C
Ambient pressure	+500 hPa...+1080 hPa
Storage/transport conditions:	
Air humidity (RH, non-condensing)	5%...95%
Temperature	0°C...60°C

\* Only valid for standard conditions. Output  
falls as ambient pressure decreases (= with  
increasing altitude of the place of use).

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**Standard compliance** (for 20132620-1)

**IEC 60601-1, IEC 60601-2-18, UL 2601.1 and CSA 22.2 No. 601.1-M90:**

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF

**IEC 60601-1-2 :2001:**

Please read the Electromagnetic Compatibility Information in the appendix (p. 41-53).

**Directive compliance** (for 20132620-1)

This medical device belongs to Class II a

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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### Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized spare parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

**Note:** We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Suitable videoscopes are listed in our GASTROENTEROLOGY catalog.

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More accessories are listed in our TELEPRESENCE catalogue.

For detailed information please contact:  
KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen, Germany

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**Appendix****Electromagnetic Compatibility (EMC)  
Information**

**WARNING:** Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The XENON 100 SCB model 201326 20-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC/EMI requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electro-magnetic interference to be expected in a medical environment. The XENON 100 SCB model 201326 20 is a Group 1 unit (as per CISPR 11).

The classification depends on the way in which the equipment or system uses RF energy. Group 1 includes all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device;
- increase the separation between the equipment;
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**WARNING:** The XENON 100 SCB model 201326 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

**WARNING:** The use of cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the XENON 100 SCB model 201326 20-1.

The cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Cables which have been shown to comply with EN/IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
Power cord	No	3	No	Power supply

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<b>Table 201</b>		
<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The XENON 100 SCB model <b>20 1326 20-1</b> is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 100 SCB model <b>20 1326 20-1</b> should ensure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The XENON 100 SCB model <b>20 1326 20-1</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The XENON 100 SCB model <b>20 1326 20-1</b> is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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<b>Table 202</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The XENON 100 SCB model <b>201326 20-1</b> is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 100 SCB model <b>201326 20-1</b> should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>EN/IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ * (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Complies <5% $U_T$ * (>95% dip in $U_T$ ) for 0.5 cycle  Complies 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  Complies 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  See guidelines	Mains power quality should be that of a typical commercial or hospital environment.       If the user of the equipment or system requires continued operation during power mains interruptions, or operation at <5% $U_T$ (>95% dip in $U_T$ ) for 5 seconds, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
The XENON 100 SCB, model <b>20</b> 1326 20-1, is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 100 SCB model <b>20</b> 1326 20-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V<sub>rms</sub> 150 kHz to &lt; 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V<sub>rms</sub></p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the XENON 100 SCB model <b>20</b> 1326 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = [3.5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
<p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than <b>3 V/m</b>.</p>			

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<b>Table 206</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the XENON 100 SCB model 20132620-1</b>			
The XENON 100 SCB model 20132620-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance <i>d</i> [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	$d = \left[ \frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <b>P</b> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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# WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:  
Company stamp/signature:

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Tipo del aparato: Número de serie:

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Firma/Fecha:

To be filled out by instrument owner:  
Return address/company stamp:

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KARL STORZ GmbH & Co. KG  
Resident Representative Office  
80/33 (44/19) Dang Van Ngu,  
F10-Q, Phu Nhuan  
Ho Chi Minh City, Vietnam  
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.  
#05 – 08 San Centre, 171 Chin Swee Road  
Singapore 0316, Singapore  
Phone: +65 6376-1066, Fax: +65 6376-1068  
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.  
174 Parramatta Rd.  
Camperdown, N.S.W. 2050, Australia  
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199  
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.  
Hong Kong Representative Office  
Unit 1601, Chinachem Exchange Square  
1 Hoi Wan Street, Quarry Bay,  
Hong Kong, People's Republic of China  
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Beijing Representative Office  
Room 610, China Life Tower  
No. 6, Chaowai Street  
Beijing, 100020, People's Republic of China  
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Unit 3901-3904, Tower 1 Grand Gateway  
No. 1 Hong Qiao Road  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Chengdu Representative Office  
F-5, 24/F., Chuanxing Mansion,  
No. 18 Renming Road South  
Chengdu, Sichuan, 610016, People's Rep. of China  
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Guangzhou Representative Office  
Room 1119-20, Dongshan Plaza  
69 Xianlie Road Middle, Dongshan District,  
Guangzhou, Guangdong, 510095,  
People's Rep. of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd  
3791 Jalan Bukit Merah  
06-11 e-Centre @ Redhill  
Singapore 159471, Singapore  
Tel. No. +65 63761066, Fax. No. +65 63761068  
Email : infoasia@karlstorz.com.sg  
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
3791 Jalan Bukit Merah  
10-18 e-Centre @ Redhill  
Singapore 159471  
Tel. No. +65 65325548, Fax No. +65 65323832  
Email : infosp@karlstorz.com.sg  
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.  
Bois Hongo Building 6Fl, 3-42-5 Hongo  
Bunkyo-ku, Tokyo 113-0033, Japan  
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)