Serie / Series 120

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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ANNEX I: CE CERTIFICATE

AUXILAB S.L.

CE DECLARATION OF CONFORMITY
ZUZI MICROSCOPES SERIES 120 of Auxilab, S.L for the Directive of Machines (89/392/CEE modified) and the regulations adopted for their transposition

NAME OF THE MANUFACTURER / IMPORTER:

AUXILAB, S.L.

ADDRESS:

**Polígono Morea Norte, 8
31191 Beriáin (Navarra)**

WE STATE THAT:

**ZUZI MICROSCOPES SERIES 120
CODES 50120006, 50120007 and 50120008**

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
- ◆ Directive 73/23/CEE modified over the electric security
- ◆ Directive 89/336/CEE modified over the electromagnetic compatibility

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BERIAIN 11th May 2006

Signed by: ALFONSO AINCIBURU SANZ
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
 Internet: www.auxilab.es · Email: correo@auxilab.es

10. ACCESSORIES

This professional series can be completed with several Zuzi accessories that increase the range of use:

Phase contrast

- ◆ We have available two models of phase contrast kits:
 - Simple phase contrast (code 90100022).
 - Turret (code 90137022): If we change the condenser, we can choose the magnifications we need using the most suitable objectives.

Thermostatic stage (code 90100597)

- ◆ It allows to keep the sample in a suitable temperature, it is very useful for the observation of living samples (e.g. cells, microorganisms, semen).

Fluorescent attachment (code 90137035)

Eyepieces

- ◆ We have available a wide range of eyepieces of different magnifications, wide field with pointer, with micrometer, etc. Ask your distributor/dealer or enter the web www.auxilab.es.

Objectives

- ◆ There is also a wide range of objectives of different magnifications: Achromatic, semi-plan, plan, etc. Ask your distributor/dealer or enter the web www.auxilab.es.

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


- a.1) Focus the preparation with a low magnifying objective (4x).
- a.2) Put a drop of immersion oil in the cover glass.
- a.3) Rotate the nosepiece and place the immersion objective, then adjust focusing with the micrometric control.
- b) You should avoid the appearance of bubbles in the oil film, since it spoils the resulting image.
- c) It is advisable to dry the frontal lens after each observation. **Never use a yellow duster!** It is better to dry it with a special paper for the cleaning of glasses or with a piece of linen cloth. The objective must never be submerged in any kind of liquid, since it can provoke the detachment of the lens.

Picture 10: Tension adjustment ring

- ♦ Depending on the observant the tension of the macrometric focal controls can be modified, since they are provided with a tension adjustment ring that is placed next to the right macrometric control of the microscope (pic. 10). However, you must avoid to loosen it excessively since this could produce that the stage falls and therefore focusing accuracy problems.

Picture 11: Focusing limit control

- ♦ There is also a device for the control of the focusing limit. All models of this series are provided with a ring that avoids that the stage raises beyond a certain point, that can be determined once the sample has been focused, and this fact can avoid possible accidental collisions while observing (pic. 11). This mechanism also allows the possibility of getting a quick focusing, accelerating the observation process when you are working with different samples.

8. TROUBLESHOOTING

Before sending the microscope to our Technical Assistance Department you should check the following points:

If the lamp does not light:

- ◆ Check that the plug is in good conditions and that it is on a suitable earth connection.
- ◆ Check the fuse and the lamp are in good conditions.

If the field of view appears to be cut out:

- ◆ Check that the nosepiece is correctly fitted in its position. To check it, rotate the nose-piece slightly clockwise and counterclockwise until it fits properly.
- ◆ Check the condenser is well centered. You should do this by adjusting three small screws placed around the condenser mount.

If there is visible dust or dirtiness on the field of view:

- ◆ Check if there is dust in the precondenser lens, upper lens of the condenser and eye-piece, as well as the cleanliness of the sample specimen. Once you have found the area causing the problem, clean it as previously instructed.

If the preparation takes out of focus while we carry out the observation:

- ◆ Check the tension adjustment ring, since it may be too weak and it could provoke an involuntary fall of the stage.

If a section of the field of view is taken out of focus:

- ◆ Check that the objective is properly placed in the light beam path. If not, rotate it until it is properly fitted.
- ◆ Check if the sample is properly placed on the stage and make sure that it is firmly held by the clip.

If the frontal lens of the objective touches the preparation while focusing or while changing from a low magnifying objective to a higher one:

- ◆ Check if the sample specimen is upside down (the microscope cover glass on top of the microscope slide) and if no proceed to place it correctly. Check also if the cover glass has the suitable thickness (0.17mm).

9. USEFUL ADVICE

- ◆ As all professional microscopes series 120 have a rotary head, it is advisable to proceed to observation by the frontal part of the stage, as it enables accessing to the mechanical controls.
- ◆ When you use the immersion objective you have to take into account the following:
 - a) To use all the numerical aperture of the immersion objective, both the objective and the sample specimen should be immersed in oil. To do so:

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6. CHOOSING OBJECTIVES AND EYEPIECES

The observed image loses surface and definition as the magnifications are higher. This increase of magnifications must be obtained by means of more powerful objectives and not from higher magnifying eyepieces, since the eyepiece only increases the image produced by the objective, therefore the more magnifications the eyepiece has, the more the loss of definition, brightness and surface of the image it provokes.


For routine observations you must use the lower magnifying eyepiece with the most powerful objectives. The higher magnifying eyepieces will be reserved for especial observations, taking into account that they reduce the definition and do not increase the resolution.

7. CAUSES OF A DEFECTIVE IMAGE

In case of a defective image you must check that:


- ♦ The lighting has been properly done, the lighting intensity cannot be excessive nor too weak. Never adjust it with the condenser's diaphragm. Both the condenser and the lamp must be well centered. Check that there is no diffusor filter between the field's diaphragm and the aperture one.
- ♦ Eyepieces fit well and objectives are properly tightened.
- ♦ The cleanliness of the whole optical system. The process should be as follows:
 - a) Make the eyepieces rotate and if the specks move, clean them.
 - b) Make the head rotate. You must never disassemble the head. However, you can clean carefully those accesible surfaces of the prisms blowing them with a rubber pear.
 - c) Rotate the objective. If the parasit images also rotate, clean it with a small dry brush trying to remove the dust.
 - d) Observe the frontal surface with a magnifying glass or an inverted eyepiece.
- ♦ There should be enough immersion oil and it must not contain bubbles or impurities.
- ♦ The thickness of microscope slides, microscope cover glasses and assembly medium so as to stop the focusing at medium or high magnifications. There are standard dimensions for both microscope slides and cover glasses. Both microscope slides and cover glasses should be clean. You should check that the slides are properly placed and that there are no superimposed cover glasses.

Picture 8: Changing the lamp

- ◆ If it is necessary to change the fuse, loosen the protective lid that is on the base of the microscope (Fuse, pic. 9), change the fuse 0.5A and place it in the same position, putting again the protective lid.
- ◆ Please use the plastic cover provided to keep the microscope away from dust lying on the optical parts when it is not used for a long period of time.
- ◆ Please keep the original packaging to transport the equipment, either when it is not being used for a long time or in case you send it for an overhaul.

Picture 9: Changing the fuse

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, you should clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xilol or toluene.
- ◆ To remove the dust settled on the lens blow with a rubber pear or clean it with a small soft brush.
- ◆ You should use non-corrosive lubricants in metallic parts, being careful of not touching optical parts.

ATTENTION!

IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.

◆ When you are not using the microscope for a long period of time please make sure it is unplugged from the net in order to avoid possible accidents.

◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).

◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).

◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

◆ Please follow the processes of use of this manual.

◆ This manual should be available for all users of this equipment.

◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.

◆ The microscope is supplied with a Schuko standard wire and it should be connected to a current wire provided with an earth wire, it should be handy to be disconnected in case of emergency.

◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the microscope (burns, hurts...) or damages in nearby equipments.

◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department. Make sure you do not use lamps with higher power, as you may provoke an overheating or any other malfunction.

◆ If the lamp burns out, change it with an Zuzi halogen lamp 6V 20W, being careful not to touch the bulb with your hands.

IMPORTANT: Before changing the lamp bulb or the fuse, make sure that the microscope is disconnected from the net.

◆ In order to enable the changing of the lamps these microscopes are provided with a device at the bottom of its base (pic. 8). Loosen the screw and take it off, then loosen the screws of the lamp holder and replace it, holding it firmly. Center the lamp holder again, if necessary.


IMPORTANT: You should *always* start operating by using the objective with lesser magnification, as it facilitates focusing and makes impossible to ruin the slides or get the objectives dirty. Furthermore, it is essential to adjust the previous lighting to using the higher magnification, and it helps you to see the whole topographic structure of the sample so as you can concentrate on the more appealing parts to observe them at higher magnifications.

- ◆ By turning the nosepiece you will pass to an objective with higher magnification. As they are paracentric you will only have to adjust focusing with the micrometric control.
- ◆ You should adjust the condenser's height to obtain a homogeneous lighting of the field of view (the higher magnifying objective is used, the smaller will be the field of view). Thus, the condenser's precision will be higher as we use high magnification values. At 40x (objective 4x and eyepiece 10x) you should descend the condenser to obtain a wider field of view and avoid having excessive lighting. The condenser will have to be raised progressively using higher objectives to concentrate the light beam and get a better illumination and contrast. With the immersion objective the condenser will have to be raised to its maximum height.
- ◆ Iris diaphragm enables decreasing the numerical aperture to a value similar to objective's. Thus you avoid marginal illumination that decreases contrast and therefore it will be more closed when you use a lesser magnification objective and with a small numerical aperture. You will only use its total aperture when you use the oil immersion objective.
- ◆ In order to properly carry out the adjustment you should start with the maximum aperture. As contrast increases because the diaphragm is closing, diffraction phenomena appear, thus limiting resolution. The correct adjustment will be obtained when we have the maximum contrast without altering resolution.
- ◆ A useful way of making the adjustment is the following one:

a) Take off the eyepiece and observe the objective's exit pupil through the empty tube, adjusting the diaphragm aperture. It is normally convenient to keep the diaphragm aperture about approx. 70-80% of the objective's numerical aperture (pic. 7).

b) The closing of the diaphragm, reducing the numerical aperture of the optical system, increases the focusing depth or penetration power. The diffraction limit must never be surpassed under pretext of increasing the focusing depth.

c) The diaphragm must not be used to reduce the lighting intensity.

Picture 7: Diaphragm's aperture

Security

- ◆ The microscope must be used by previously qualified staff that knows how to operate it thanks to the user manual.
- ◆ You should put the microscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.

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
- ◆ Assemble the head by adjusting it with the built-in screw (pic. 4). Insert the eyepieces (WF10x). The binocular tube is normally placed in the same direction as the microscope front, but if necessary it can be placed in any other direction.


Picture 4: Placing the eyepieces

- ◆ Always use the microscope dust cover when the microscope is not being used in order to avoid dust settling on the optical parts, especially for long periods of time (weekend, holidays, ...).

- ◆ Models provided with binocular observation system must have the eyepiece tubes properly adjusted. This will depend on:

a) Interpupillary distance (pic. 5). Adjust the interpupillary distance by separating or joining the eyepieces until there is a total fusion of the two images.


Picture 5: Interpupillary distance

b) Dioptic compensation of the tubes: You will appreciate a difference in the focusing more or less marked, closing alternatively one eye, then the other. To correct this defect you should adjust the dioptic ring to zero. Looking through the right eyepiece you have to focus the sample specimen using the macro and micrometric controls. Once you have focused the sample, look through the left eyepiece and then adjust the dioptic adjusting ring until you can obtain a clear, sharp image. In case the observer wears glasses he must do the adjustment with them on, making sure he/she maintains a certain distance from the eyepiece lens so as to avoid grazing the glasses.

Once the values are found it will be very useful to learn them by heart, mostly if the microscope is shared with other users. This is highly advisable to avoid repeating the adjustment each time the microscope is used.

- ◆ Before placing the sample specimen on the stage, this one should be descended to a certain distance, higher than the working distance of the objective with the lower magnifying values. Then, you should place the specimen on the mechanical stage fitting it into the set square with the help of the moveable clips. Finally you should centre the specimen sample regarding the stage aperture by looking from the outside.

Picture 6: Focusing

- ◆ Looking by the eyepiece and with the lower magnifying objective, raise the stage with the macrometric control until the image appears. Then, move the micrometric control until obtaining a correct focusing (pic.6).

Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The microscope is supplied with a Schuko standard wire
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire and the other extreme to the microscope connector.


Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

- ◆ If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Setting up

- ◆ The observer should be in a comfortable position, properly sit down and with the back in a straight position. It is advisable to work on a dark table in order to avoid the light that dazzles and reduces the good definition of the images, thus avoiding ludicrous fatigue. Pay attention to background lighting and install the microscope away from windows.

- ◆ Take down the mechanical stage by means of the macrometric (coarse) controls. Screw the objectives into the nosepiece in an ascendant order (4x-10x-40x-100x) in clockwise direction (pic. 3).

Picture 3: Placing the objectives

VERY IMPORTANT: You must never hold the microscope by the mechanical stage or the tube, thus all weight would be supported by the micrometric screw, whose notches would be worn out slowly. You should always hold it by the base.

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- R: Retractable.

- I: Immersion objective.

- N.A.: Numerical aperture, it determines the objective's properties. The numerical aperture of a microscope objective is a measure of its ability to gather light and resolve fine sample detail at a fixed object distance. Thus, the bigger it is, the brighter and better resolved the image will be.

- Working distance: It is the distance (mm) between the sample specimen and the frontal lens of the objective when the microscope is focused.

- Focal distance: It is the distance (mm) between the main plan image of the system and its image focus.

- Resolution: It is the reciprocal value of the separating power, that is to say, the smallest distance between two points on a sample that can still be distinguished as two separate entities. It is calculated in lines/mm.

- Field's number: It expresses the diameter (mm) of the field diaphragm composed by the eyepiece.

- Field of view: Dimensions (mm) of the real field we are observing.

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the microscope, take off the involving plastic and take off the poliespan protection in which it comes fitted. Take off all the protective items and, without connecting the microscope to the net, make sure that it does not present any damage because of the shipment. In case the microscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ◆ Please check that all the accessories are enclosed with the equipment:

- Base	- Spare halogen lamp of 6V 20W
- Head (depending on the model supplied)	- Immersion oil
- Objectives 4X,10X, 40X (R), 100X (R) (I)	- Fuse
- Eyepieces (WF10X)	- Plastic cover
- Standard Shuko cable	- User manual
(10A, 250V, UNE 20356, EN 60 320).	- Warranty certificate
- Set of filters (green and blue)	

We will only accept the equipment's return within 15 days after delivery provided it comes in its original wrapping.

Stage

- ◆ Coaxial control, vertical controls.
- ◆ Dimensions: 130 x 140 mm
- ◆ Movement range: 75(X) x 50(Y) mm.
- ◆ Graduated scales at 0,1 mm.

Iris diaphragm

- ◆ It is placed under the condenser and the mechanical stage. It controls the light going through the condenser.

Lighting system

- ◆ Low tension halogen lamp 6V 20W, built-in power input (220-230V, 50-60Hz \pm 10%), with intensity control potentiometer.

Optical specificationsAchromatic objectives

They correct optical aberrations. This optical system is able to transmit white light without separating constituent colours.

The characteristics of each objective are codified with some marks as shown above:

Picture 2

- 40:** Objective magnification
0.65: Numerical aperture
160: Tube's length
0.17: Microscope slide's thickness

Achromatic objectives

Magnifications	4X	10X	20X	40X	60X	100X
Working distance (mm)	17.50	7.32	1.26	0.63	0.41	0.19
Numerical aperture	0.10	0.25	0.40	0.65	0.85	1.25

Total magnifications

It is the result of multiplying the eyepiece magnification and the objective magnifications.

Objectives	4X	10X	20X	40X (R)	60X (R)	100X(R)(I)
Eyepieces WF10X						
Total magnific.	40	100	200	400	600	1000
Field of view (mm)	4.50	1.80	0.90	0.45	0.30	0.18

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3. TECHNICAL SPECIFICATIONS

Technical specifications

		120/6	120/7	120/8
Head	Monocular	♦		
	Binocular		♦	
	Trinocular			♦
Eyepieces	WF 10x	♦	♦	♦
Nosepiece	Quadruple	♦	♦	♦
	4X N.A: 0.10			
	10X N.A: 0.25			
	40X (R) N.A: 0.65			
	100X (R) I N.A:1.25			
Objectives	Achromatic (anti-fungus)	♦	♦	♦
Condenser	Abbe N.A 1.25, diaphragm & filter holder	♦	♦	♦
Stage	Fixed, mechanical with mov. 130x140 mm	♦	♦	♦
	Movement range 75 x 50 mm			
	Clips for samples			
Lighting	Transmitted, halogen lamp 6V 20W	♦	♦	♦
Dimensions	180 mm x 380 mm x 230 mm	♦	♦	♦

Base

It is metallic, very stable, provided with coaxial controls at both sides to enable macro-metric and micrometric focusing. The right control has a friction control ring and the left one has a stage-blocking handle. The mentioned handle stops the stage at the higher point in order to avoid that the objective may crash the slide and so damage either the sample specimen or the objective itself.

Head

We can supply you with different heads, depending on the model:

- ◆ Monocular, inclined 45° and rotary 360°. It can be fixed in any position.
- ◆ Binocular, inclined 45° and rotary 360°. Symmetrical interpupillary graduation with scale between 55 and 75mm and dioptic compensation ring ± 5 dioptres.
- ◆ Trinocular, it is provided with a binocular head, inclined 45° and rotary 360°. Symmetrical interpupillary graduation with scale between 55 and 70mm, dioptic compensation ring ± 5 dioptres and a vertical tube to be adapted to a photographic camera or a video.

Condenser

- Abbe, double lens.
- ◆ Numerical aperture: 1.25.
- ◆ Moveable filter holder. It includes blue and green filters (32 mm diameter).
- ◆ Vertically variable by pinion and rack.

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	18
2. DESCRIPTION	18
3. TECHNICAL SPECIFICATIONS	19
4. INSTALLATION / SETTING UP	21
5. MAINTENANCE AND CLEANING	25
6. CHOOSING OBJECTIVES AND EYEPIECES	27
7. CAUSES OF A DEFECTIVE IMAGE	27
8. TROUBLESHOOTING	28
9. USEFUL ADVICE	28
10. ACCESSORIES	30
ANNEX I: CE CERTIFICATE	31


1. USES OF THE INSTRUMENT

Zuzi professional microscopes series 120 are very hard equipments designed to provide the highest satisfaction in a wide range of demands. They are equipped with high quality optics and a wide range of accessories to complete the microscope.

Zuzi series 120 are suitable for education as well as for veterinary, clinical fields, etc.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your professional microscope Zuzi series 120. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user. The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of *Maintenance*, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the professional microscope Zuzi series 120, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.com or by telephone: +34 807 117 040 (0,30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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