Rose Bengal Brucellosis Antigen

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For veterinary use only.

### Name and Intended Use

Pourquier\* Rose Bengal Ag is IDEXX's Rose Bengal assay for serological detection of antibodies against Brucellosis (*B. melitensis*, *B. abortus* and *B. suis*) by rapid slide agglutination.

### Descriptions and Principles

Pourquier\* Rose Bengal Ag uses the principle of agglutination on plates using an antigen stained with Rose Bengal in a buffered acidified medium (pH  $3.65 \pm 0.05$ ). The acidified medium avoids non specific agglutinations. In presence of antibodies against Brucellosis in the sample, the Brucellosis antigen will agglutinate to the sample. Samples showing presence of agglutination (even slight) are considered positive samples (presence of anti-Brucellosis antibodies). Samples showing no agglutination are considered negative samples (absence of anti-Brucellosis antibodies).

The suspension of Brucella antigen used in this kit is a concentrated suspension of *Brucella abortus* antigen (Weybridge strain 99) inactivated by heat and phenol (0,5%), dispersed in an acid buffer and stained by Rose Bengal. It is calibrated to give a positive reaction at a dilution 1:45 and a negative reaction at a dilution 1:55 with the OIE International Standard Serum (OIEISS) or National secondary standard, in accordance to OIE requirements.

### Reagents

Store all reagents at 2–8°C, protected from light. **Do not freeze.**

Reagents		Volume	
1	Rose Bengal Brucellosis Antigen	5 x 10 mL	100 mL

NOTE: see table on page 17 for the description of international symbols used on the kit labels.

### Materials Required but Not Provided

- Plates to dispense reagents and samples
- Precision Micropipettes and Multi-dispensing micropipettes (reagents volumes listed in the “Test Protocol” require pipette precision of  $\pm 10\%$  for volumes lower or equivalent to 10  $\mu\text{l}$  and  $\pm 5\%$  for the other indicated volumes)
- Disposable pipette tips
- Plate shaker

## Precautions and Warnings for Users

- Handle all biological material as potentially infectious.
- Do not pipette by mouth.
- Do not eat, drink or smoke where specimens or kit reagents are being handled.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations.
- Care should be taken to prevent the contamination of kit components.
- Do not use components past their expiration dates and do not mix components from different serials.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing and washing throughout this procedure are necessary to maintain precision and accuracy.

## Test Procedure

The Rose Bengal Brucellosis Antigen and the serum samples to be tested must be allowed to come to 18–26°C for 30–60 minutes before use, then mixed by inverting and swirling.

Only sufficient antigen for the day's tests should be removed from the refrigerator.

Use a separate pipette tip for each sample.

1. Obtain plates to dispense reagents and samples and record the sample position on a worksheet.
2. Dispense 25–30  $\mu\text{l}$  of each serum to be tested on the plate.
3. Gently shake the Rose Bengal Brucellosis Antigen vial and dispense the same volume (25–30  $\mu\text{l}$ ) besides each sample.
4. Carefully mix the Rose Bengal Brucellosis Antigen and the serum to produce a circular or oval zone approximately 2 cm in diameter.
5. Shake the plate gently for exactly 4 minutes and read immediately.

Note: In order to standardize the reading, it is better to introduce as control for each series of analysis a positive control serum for Brucellosis (available from IDEXX, product code P00130) and a known-negative serum sample (available from IDEXX, product code P00131).

## Results

Check that the results of controls are in accordance with the expected results.

## Interpretation of Results

- Samples showing presence of agglutination (even slight) are considered positive samples (presence of anti-Brucellosis antibodies).
- Samples showing no agglutination are considered negative samples (absence of anti-Brucellosis antibodies).

## Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action
1. <b>Sample distribution</b>	<ul style="list-style-type: none"><li>– Obtain plates to dispense reagents and samples and record the sample position on a worksheet.</li><li>– Dispense 25–30 <math>\mu</math>l of each serum sample on the plate.</li><li>– Gently shake the Rose Bengal Brucellosis Antigen vial and dispense the same volume (25–30 <math>\mu</math>l) besides each serum.</li><li>– Carefully mix the Rose Bengal Brucellosis Antigen and the serum to produce a circular or oval zone approximately 2 cm in diameter.</li></ul>
2. <b>Sample incubation and reading</b>	Shake the plate for exactly 4 minutes and read immediately.
3. <b>Interpretation</b>	No agglutination: absence of antibodies (negative). Agglutination (even slight): presence of antibodies (positive).

**Manufactured by:**

IDEXX Montpellier SAS  
326 rue de la Galéra  
34090 Montpellier - France

**For technical assistance:**

Contact your IDEXX area manager or distributor  
or visit: [www.idexx.com/production/contact](http://www.idexx.com/production/contact)  
IDEXX Technical Support: 00-800-727-43399

\*Pourquier is a trademark or a registered trademark of IDEXX Laboratories, Inc. in the United States and/or other countries.

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IDEXX Europe B.V.