AUXILAB S.L.

Material de laboratorio
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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POWER SUPPLIES FOR ELECTROPHORESIS

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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
- ◆ When you are not using the power supply for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the power supply by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the equipment (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the power supply it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The power supply is supplied with a Schuko wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the power supply by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the power supply (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the power supply and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

Time

Nahita power supplies present a time selecting range from 1 min to 99 h 59 min.

- ◆ Press the selection key (2.5); a cursor (←) will be displayed. Move the cursor by pressing the selection key until selecting the parameter to be changed (the operating time in this case, T_S).
- ◆ Press the increase (2.6) and decrease (2.7) keys until adjusting the selected parameter to the desired value.
- 5. Press the selection key (2.5) to move the cursor to "T=0:00 Start".
- 6. Press Start/Pause (2.4); 4 beeps will sound and then the power supply will start running.

Note: If pressing Start/Pause (2.4) without previously selecting "T=0:00 Start", 5 beeps will sound and the cursor will automatically move to "T=0:00 Start". Then, press Start/Pause (2.4) to start the power supply.

Note: During operation the value of the different parameters can be changed by selecting the desired parameter with the selection key (2.5) and then increasing (2.6) or decreasing (2.7) its value. E.g., if you need to work at a constant voltage of 200 V and the power supply only reaches to 150 V and keeps a constant current of 100 mA, the current value can be increased during operation to make the power supply reach 200 V and keep this value constant.

Note: By pressing the selection key during operation the selected values of voltage (U_S), current (I_S) and time (T_S) will be displayed as well as their real values (U , I y T).

Note: To stop temporarily operation, press Start/Pause (2.4); the power supply will stop immediately and "STOP" will be displayed. To restart operation, press the selection key (2.5) to move the cursor to "GO ON" and then press Start/Pause (2.4).

If you are not using the power supply for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Security

- ◆ The power supply must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the power supply on top of a horizontal, plane and stable table making a free space of at least 30 cm per side.
- ◆ Do not put the power supply near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations, dust and dry environments.
- ◆ During operation, dangerous material such as flammable or infectious liquids must be out of the working area..

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Red lead	Red terminal, positive
Black lead	Black terminal, negative

Note: In case of connecting two or more cells at the same time, it will be necessary to work at a constant voltage, since the current displayed is the sum of the currents of the different electrophoresis cells while the voltage is the same for all the cells.

2. Press power on/off (2.2). The display will light up and four beeps will sound; finally the last parameters selected will be shown.

3. To work with the same parameters, press Start/Pause (2.4) and the equipment will start running.

4. To change the value of voltage, current and running time proceed as follows:

Voltage

Nahita power supplies present a voltage selecting range of 6-600 V (Model 301/6) and 2-300 V (Model 301/7) with increments of 1 V.

◆ Press the selection key (2.5); a cursor (←) will be displayed. Move the cursor by pressing the selection key until selecting the parameter to be changed (the output voltage in this case, U_S).

◆ Press the increase (2.6) and decrease (2.7) keys until adjusting the selected parameter to the desired value.

Note: To work at a constant voltage, adjust the voltage to the working value and then adjust the current to the maximum value allowed in the range of the power supply. Thus, when the power supply reaches the selected voltage value will keep it constant though the current value is lower than the selected.

Current

Nahita power supplies present a current selecting range of 4-400 mA with increments of 1 mA (Model 301/6) and 5-2000 mA with increments of 2 mA (Model 301/7).

◆ Press the selection key (2.5); a cursor (←) will be displayed. Move the cursor by pressing the selection key until selecting the parameter to be changed (the output current in this case, I_S).

◆ Press the increase (2.6) and decrease (2.7) keys until adjusting the selected parameter to the desired value.

Note: To work at a constant current, adjust the current to the working value and then adjust the voltage to the maximum value allowed in the range of the power supply. Thus, when the power supply reaches the selected current value will keep it constant though the voltage value is lower than the selected.

4. INSTALLATION /SETTING UP

Preliminary inspection

Unwrap the power supply, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all protections and without connecting the power supply to the electrical current, make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Power supply
- Schuko wire
- Warranty
- User's manual

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using the power supply, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

Please put the power supply on top of a horizontal, plane and stable table making a free space of at least 30 cm per side. Do not put the power supply near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. During operation, dangerous material such as flammable or infectious liquids must be out of the working area.

The power supply is supplied together with a Schuko wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz \pm 10% provided with earth wire and to the other end to the power supply connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz \pm 10%.

1. Plug the leads of the electrophoresis cell lid to the power supply. Make sure that the leads are connected to the proper output terminal:

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
- 2.1 Output terminal (positive, red; negative, black)
- 2.2 Power on/off
- 2.3 Operation light: it switches on when the power supply starts operating.
- 2.4 Start/Pause
- 2.5 Selection key: it moves the cursor to the different parameters that can be adjusted.
- 2.6 Increase: it increases the value of the parameter selected with the selection key.
- 2.7 Decrease: it decreases the value of the parameter selected with the selection key.
- 2.8 LCD display: it shows the value of voltage, output current, selected time, etc.

3. TECHNICAL SPECIFICATIONS

Code	53010006	53010007
Output voltage	6-600 V	2-300 V
Output current	4-400 mA	5-2000 mA
Output power	240 W	300 W
Resolution	1 V, 1 mA	1 V, 2 mA
Timer	1 min - 99 h 59 min	
Stability	≤ 1% (constant voltage); ≤ 2% (constant current)	
Precision	≤ 2% (constant voltage); ≤ 3% (constant current)	
Ripple factor	< 2%	
Working conditions	5-40 °C / =80% RH	
Power	AC, 220 V ± 10% (50 Hz ± 2%)	
Dimensions	315x290x128 mm	
Weight	5 Kg	
Applications	DNA, RNA, proteins, 2D electrophoresis	Western blotting semi-dry blotting

INDEX OF LANGUAGES

Spanish 3-10
 English..... 11-18

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT 12
 2. DESCRIPTION 12
 3. TECHNICAL SPECIFICATIONS 13
 4. INSTALLATION / SETTING UP..... 14
 5. MAINTENANCE AND CLEANING 17
 ANNEX I: CE CERTIFICATE 20


1. USES OF THE INSTRUMENT

Nahita power supplies are specially designed to operate together with the horizontal and vertical electrophoresis cells. Controlled by microprocessor, power supplies can work at constant voltage or constant current output and they also present timer control to work in concrete periods of time with the capability to pause and resume a run.

Nahita power supplies present 4 parallel output jacks to be connected to up to 4 electrophoretic cells at the same time being ideal for clinical diagnostic and research laboratories.

2. DESCRIPTION

- 1.1 Fuse
- 1.2 AC plug
- 1.3 Fan: it starts operating when the voltage or current reaches the selected value.
- 1.4 Vent

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita electrophoresis power supply. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita electrophoresis power supply, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the power supply.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0,30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.