# STORZ

## KARL STORZ — ENDOSKOPE

201123 20

Kaltlicht-Fontäne HALOGEN 150

201123 20

Cold light fountain HALOGEN 150

201123 20

Fuente de luz fría HALOGEN 150

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**Instruction Manual**

**Manual de instrucciones**

	201123 20	Kaltlicht-Fontäne HALOGEN 150
	201123 20	Cold light fountain HALOGEN 150
	201123 20	Fuente de luz fría HALOGEN 150

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HALOGEN 150 cold light fountain. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG  
Utilization and duplication by third parties have to be expressly authorized by KARL STORZ GmbH & Co. KG  
All rights reserved.

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**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

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



**Controls, displays,  
connectors and  
their uses**

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Photographs of the unit ..... A 1  
Controls, displays, connectors,  
and their uses ..... A 2  
Symbols employed ..... A 3

**General information**

Device description ..... 2

**Safety instructions**

Warnings and cautions ..... 3  
Intended use ..... 6  
User qualification ..... 7  
Safety precautions at the site of installation ..... 7  
Safety precautions when operating the  
unit ..... 7

**Installation and operating instructions**

Unpacking the equipment ..... 8  
Basic equipment ..... 8  
Installation ..... 8  
Operating the cold light fountain ..... 9  
Adaptor for light cables of other manufacturers ..... 12

**Maintenance**

Lamp replacement ..... 13  
Fuse replacement ..... 15  
Cleaning and care of the unit ..... 16  
Maintenance ..... 16  
Servicing and repair ..... 16  
Repair program ..... 17  
Limitation of liability ..... 18  
Manufacturer's warranty ..... 18

**Technical description**

Troubleshooting ..... 19  
Technical data ..... 21  
Technical documentation ..... 22  
General circuit diagram ..... 23

**Spare parts,  
recommended accessories**

List of spare parts ..... 24  
Accessories ..... 24

**Appendix**

Cleaning agents and disinfectants ..... 30  
Electromagnetic Compatibility  
(EMC) Information ..... 32

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**Device description**

The cold light fountain HALOGEN 150 is suitable for all types of endoscopic examination. Color temperature is approx. 3.400 K.

A high-quality condenser system ensures high luminous efficiency while modern design guarantees easy operation and unproblematic servicing.

Adaptors for light cables of other manufacturers are available as accessories.

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**Warnings and cautions**

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

**Warning:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**Caution:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**Note:** A Note indicates special information about operating the product, or clarifies important information.

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**Warning:** Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**Note:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

**Warning:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

**Warning:** Risk of explosion if used in the presence of flammable anesthetics.

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**Warning:** The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

**Warning:** Keep out of reach of patients.

**Warning:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**Warning:** Grounding reliability can only be achieved when the equipment is connected to a properly installed "Hospital Only" or "Hospital Grade" receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**Warning:** Always unplug the unit before doing any maintenance work on it (e. g. cleaning).

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**Warning:** Danger of electric shock. Do not open the unit. Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

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**Caution:** Always adjust the cold light fountain to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, either by direct vision or when coupled to a video camera.

**Caution:** Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

**Caution:** Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

**Caution:** Use only fuses of the correct rating.

**Caution:** Do not touch glass bulb with your bare fingers, use the protective cover. Touching the lamp with bare fingers may result in starting difficulty.

**Caution:** Federal law restricts this device to sale by or on the order of a physician (USA only).

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**Warning:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and at the tip of the endoscope.

**Warning:** Never leave an illuminated light cable with or without an endoscope connected, on or near a patient; the risk of patient burns can result.

**Caution:** Do not place connected light cables on combustible materials such as textiles (drapes) or in the direct vicinity of swabs, etc. soaked with combustible fluids, such as disinfectant. These could catch fire from heat development at the open end of the light cable.

Therefore always attach the open end of a connected light cable to the provided retainer ④.

**Note:** Lamps do not have an indefinite service life, so it is possible that one may fail during use. For this reason a spare lamp system or cold light source should always be kept available during therapeutic endoscopic interventions.

**Note:** This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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### **Intended use**

The cold light fountain HALOGEN 150 Model 20 1123 20 is used as a light source in conjunction with rigid and flexible endoscopes and other instruments that are designed for use with a cold light fountain. Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The cold light fountain HALOGEN 150 may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

**Note:** KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
**User qualification**

The cold light fountain HALOGEN 150 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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**Safety precautions at the site of installation**

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. It is not intended for use in hazardous zones.

This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected before power is applied to the unit in compliance with current national regulations.

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**Safety precautions when operating the unit**

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment using the cold light fountain HALOGEN 150 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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### Unpacking the equipment

Carefully unpack the cold light fountain and remove the unit and its accessories from their packaging.

Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
### Basic equipment

- 1 Cold light fountain HALOGEN 150
- 1 Power cord 400 A
- 1 Instruction manual
- 1 Wallet for instruction manual

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### Installation

**Note:** The cold light fountain HALOGEN 150 may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

**Warning:** This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Set the unit on a flat surface.

**Caution:** Do not cover the air vent. Danger of overheating. Make certain there is free air circulation towards the top and to the side (distance of rear panel to wall: 10 cm). If the unit overheats it will shut down automatically. Turn off the main power switch ① and let the light source cool down (see error list).

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
## Installation and operating instructions

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The unit is equipped with a connector ⑥ for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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
### Operating the cold light fountain Setting of the system voltage

Remove the power fuse holder ⑧ (see section "Fuse replacement", page 15) and take out the voltage selection module.

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Insert the voltage selection module into the power fuse holder in such a way that the system voltage displayed in the viewing window on the back of the power fuse holder corresponds to the local power supply.

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Check the inserted power fuses using the following table:

100 VAC	2 x T 3.15 AL250V
120 VAC	2 x T 3.15 AL250V
230 VAC	2 x T 1.6 AL250V
240 VAC	2 x T 1.6 AL250V

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Connect power cord. Insert power cord into power cord receptacle ⑨ as far as it will go.

**Warning:** The unit should only be connected to or disconnected from the power supply outside potentially explosive locations.

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Plug the light cable into socket ⑤ as far as it will go. To disconnect the light cable pull on the connector, never on the cable itself.

**Warning:** Danger of glare. Never look into the open end of a light cable.

**Warning:** Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

**Warning:** Light cable socket may be very hot; particularly if the temperature switch has been triggered and the unit has switched off.

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Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base.

**Note:** KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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If you do not plan to conduct the light cable to the endoscope right away, plug the loose end into the dummy socket ④.

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**Note:** The entrance and exit faces of the prepared light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Remove heavy soiling with a lint-free disposable cloth, disinfect the entrance and exit faces with a cotton tip applicator soaked in alcohol and dry.

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Switch on the unit. The power switch ① and the operating pilot lamp ③ are lit.

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**Note:** The light exit ⑤ is provided with an anti-glare flap which does not allow any direct outlet of light. Light only appears when the light cable is attached.

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Adjust the required luminous intensity with the intensity switch ② (3 levels).

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**Note:** Frequently switching the cold light fountain on and off causes the halogen lamp to wear out more quickly. For relatively brief interruptions, set the luminous intensity switch (2) to Level 1.

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### Adaptor for light cables of other manufacturers

If necessary light exit (5) can be unscrewed and replaced by adaptors for light cables of other manufacturers (see section “Spare parts, recommended accessories”, page 29).

To do so, unscrew light exit (5) by rotating counterclockwise.

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Place the adaptor on the screw thread and tighten up by rotating clockwise.

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**Lamp replacement**

Switch off the unit. Disconnect the power plug.

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Open the housing cover. For this purpose, first undo the 2 thumbscrews at the bottom of the unit.

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Turn the unit back up the right way and slide the housing cover towards the back.

**Note:** Slide back the housing cover just far enough that the lamp compartment is easily accessible.

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Remove the lid from the lamp compartment. To do this, loosen the thumbscrew.

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Pull out the defective halogen lamp.


**Warning:** Danger of burns. Lamp may be hot.

**Caution:** Do not touch the glass of the new bulb with your fingers.

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**Note:** Do not touch glass bulb with your bare fingers. The glass bulb should be free of grease and other residue build-up. (If necessary, wipe the bulb with a cotton swab moistened with alcohol before operation).

Please use the protective cloth provided with the new lamp when inserting the lamp. Insert new halogen lamp (15 V/150 W) into the contact socket until firmly in place. The lamp is automatically adjusted.

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Reinstall the lamp compartment lid. Close the housing and tighten the thumbscrews. Reconnect the power cord. Carry out function test.

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**Fuse replacement**

Switch off the unit and remove the power cord.

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
Remove the power fuse holder ⑧ with a suitable tool.

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**⚠ Caution:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

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Refit the power fuse holder ⑧.  
Reconnect the power cord.  
Test the instrument for proper operation.

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### Cleaning and care of the unit

**Warning:** Always pull out power plug before cleaning! Make sure the lamp has cooled down.

**Caution:** Avoid allowing fluids to enter the unit.

The exterior of the unit can be surface disinfected with a disposable cloth soaked in disinfectant (preferably alcohol-free).

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
### Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3515 (English version).

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### Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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**Disposal**

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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**Repair program**

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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**Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations in the room where the unit is connected and operated comply with the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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### Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply. Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen  
Germany

**Note:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen  
Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**Troubleshooting**

**Warning:** Always unplug the unit before all maintenance work.

*Symptom:*  
Complete failure of the unit.

*Possible causes:*  
– No power from the power line.  
– Defective power fuse.

*Remedy:*  
– Check that there is electricity to the wall outlet.  
– Change fuses as described in the instruction manual. Make sure you use the correct type of fuse.

*Symptom:*  
No light output, operating pilot lamp ③ is off. Power switch is lit (on).

*Possible causes:*  
– Lamp defect.  
– Temperature safety device triggered.

*Remedy:*  
– Insert new lamp.  
– Disconnect unit from power and allow to cool down. Switch on again after approx. 5 minutes. Make certain there is sufficient ventilation.

*Symptom:*  
Too little light.  
*Possible causes:*  
– Luminous intensity switch is set too low (Stage 1 or 2).  
– Ends of light cable and/or endoscope are soiled.  
– Light cable is defective.  
– Dimmer defective.

*Remedy:*  
– Set the luminous intensity switch to a higher stage (2 or 3)  
– Clean ends of light cables, and light entrances and exits on the endoscope.  
– Replace light cable.  
– Call service.

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*Symptom:*  
No light output in spite of changing the bulb.

*Possible causes:*

- Bulb inserted is defective.
- Automatic cutout ⑩ has operated.
- Bulb base is defective.

*Remedy:*

- Fit a new bulb.
- Activate the automatic cutout ⑩ by pressing.
- Have the unit repaired.

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**Standard compliance** (for 201123 20)

**According to IEC 60601-1:**

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF 

**According to IEC 60601-1-2:**

Please read the Electromagnetic Compatibility Information in the appendix (p. 32-43).

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**Directive compliance** (for 201123 20)

**According to Medical Device Directive (MDD):**

This medical product belongs to Class I

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

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### Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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
**Accessories**

**Fiberoptic light cable**

Fiberoptic light cables contain a large number of single glass fibers. The spectral transmission of a fiberoptic light cable is approximately the same for all of the wavelengths of visible light. It decreases, however, towards the blue end of the spectrum. This causes the light to take on a warmer shade – an effect which becomes more intensive the longer the light cable. Fiberoptic light cables are opaque to ultraviolet light.

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**Cleaning, disinfection and sterilization of fiberoptic light cables**

**Note:** Comprehensive information on cleaning light cables can be found in the 'Light cable' instructions.

1. **Manual (preliminary) cleaning:** The light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Wipe the outside of the light cables with a sponge or cloth soaked in cleaning solution.\*
2. **Disinfection:** immerse in a disinfection solution (see Appendix, page 30, 31).\*

\* re 1./2. The light cable must always be rinsed until it is free of all residue after it has been cleaned or disinfected manually. Microbiologically pure/sterile water must be used for the final rinsing.

3. **Thermal decontamination by machine** (up to 93°C)

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4. **Manual final cleaning:** Clean the fiber surfaces and optical end faces with a soft cloth, sponge or cotton swab, soaked with 70% alcohol. Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).
5. **Sterilization:** Steam sterilization at 134°C (+3°C) STERIS SYSTEM1® processor or gas or plasma sterilization.

**Note:** The instructions for 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

Use a cleaning agent/disinfectant approved by KARL STORZ for cleaning/disinfection (see Appendix). You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

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### Fiberoptic light cable





#### Combination of light cable with endoscope

For safety reasons, the diameters of fiberoptic/ fluid\* light cables and endoscopes should be used in the combination shown in the table opposite. \* see page 27

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**Fluid light cables**

In fluid light cables the light is not transmitted through glass fibers, but through a special liquid contained in the light cable. Fluid light cables are stiffer than glass fiber light cables and bending them too much will cause a severe loss of light. However, they do transmit a brighter light than fiberoptic light cables with the same diameter. For this reason, fluid light cables are ideally suited for endoscopic photography. Fluid light cables are permeable to ultraviolet light.


**Caution:** Damage to the elastomer coat will impair the light cable.

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**Cleaning, disinfection and sterilization of fluid light cables**

1. Manual (preliminary) cleaning: The fluid light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Remove heavy soiling with a moist sponge or lint-free disposable cloth.\*  
2. Disinfection: immerse in a disinfection solution (see Appendix, page 30, 31).\*

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3. Machine decontamination

**Caution:** Carry out chemical disinfection by machine at max. 65°C.

**Caution:** When preparing light cables in a cleaning and disinfection machine, particular care is necessary in the choice of chemical agents and machine programs.

4. Manual final cleaning: Clean the fiber surfaces and optical end faces with a soft cloth, sponge or cotton swab, soaked with 70% alcohol.

Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).

5. Sterilization: Gas sterilization or STERIS SYSTEM1® processor.

**Caution:** Fluid light cables must not be steam or plasma sterilized under any circumstances.

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More accessories are listed in our TELEPRESENCE catalogue.

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**Typ**

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
- (A)= Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets**  
**Manual preparation**

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Hersteller/Manufacturer/Fabricante	Handelsname/Commercial name/Nombre comercial	Typ/Type/Typo
Ecolab GmbH & Co. OHG	Sekusept Plus Sekusept Pulver Classic Ultrasemit AF	1, 2, 3 1, 2, 3 1, 2, 3
Estee Pharma GmbH	Ultrasreptin Aktiv Ultrasreptin Classic	1, 2, 3 1, 2, 3
Fresenius AG	Afid Afid plus Sporcid FF	1, 2 2 1, 2
Holifa Polska Sp. zo.o	Poisept Holifa	1, 2
Indeba Industria E Comercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido Darodor 3000 Darodor Sinaldehyd 2000	1 2 2
Lysoform Dr. Rosemann GmbH	Aldasan 2000 Almyrol Desoform Lysoformin 3000	1, 2 1, 2 1, 2, 3 1, 2
Medichem International	MedDis MediZyme Edisonite Super	2 1 1, 3
Merz Hygiene GmbH	Mucadont-IS Mucadont-Zymaktiv Mucocot-1 neu	2, 3 1, 3 1, 2, 3
METASYS Medizintechnik GmbH (formerly LONZA)	ID 50	2
Orochemie	A 20	1, 3, A
Promagent AB	Wavacide	2
Schumacher, Dr. GmbH	Descoton Extra Descoton forte Perfektan Endo Perfektan Neu	1, 2 2 2, 3 1, 2
Schülke & Mayr GmbH	Perfektan TB Gigasept Gigasept AF forte Gigasept FF Gigasept Instru AF Gigasept Med forte Gigazyme Lysetol V	1, 2, 3, A 2, 3 1, 2 1, 2 1, 2, 3, A 1, 2 1, 3 2
Steris	EnzyCare 2	1, 3
Weigert, Dr. GmbH & Co.	neodisher LM 2 neodisher mediclean neodisher mediclean forte neodisher medizym neodisher SeptoClean	1, 3 1 1, 3 1, 3 1, A
Whiteley Industries PTY.Ltd	neodisher SeptoPreClean Aidol Plus Matrix	1, 2 1 2

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**\* Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination sheath, **not** the housing/handle part.

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**Type**

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Machine preparation**

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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**Appendix****Electromagnetic Compatibility (EMC) Information**

**WARNING:** Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The HALOGEN 150 model 201123 20 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HALOGEN 150 model 201123 20 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**WARNING:** The HALOGEN 150 model 201123 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

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<b>Table 201</b>		
<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
HALOGEN 150 model <b>20</b> 1123 20 is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 150 model <b>20</b> 1123 20 should ensure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	HALOGEN 150 model <b>20</b> 1123 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The HALOGEN 150 model <b>20</b> 1123 20 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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<b>Table 202</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
<p>HALOGEN 150 model <b>20 1123 20</b> is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 150 model <b>20 1123 20</b> should ensure that it is used in such an environment.</p>			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ * (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Complies <5% $U_T$ * (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 0.5 cycle See guidelines 70% $U_T$ (30% dip in $U_T$ ) for 1 cycle See guidelines Complies <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the equipment or system requires continued operation during power cuts, short term power interruptions or if power supply is erratic, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

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<b>Table 204</b> <b>Guidance and manufacturer's declaration – electromagnetic immunity –</b> <b>for equipment and systems that are not life-supporting</b>			
HALOGEN 150 model <b>20 1123 20</b> is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 150 model <b>20 1123 20</b> should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6  Radiated RF IEC 61000-4-3	$3 V_{rms}$ 150 kHz to < 80 MHz  $3 V/m$ 80 MHz to 2.5 GHz	$3 V_{rms}$  $3 V/m$	Portable and mobile RF communications equipment should be used no closer to any part of the HALOGEN 150 model <b>20 1123 20</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz $d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz  Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
<p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than <b>3 V/m</b>.</p>			

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<b>Table 206</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the HALOGEN 150 model 201123 20</b>			
The HALOGEN 150 model 201123 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	$d = \left[ \frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance <b>d</b> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <b>P</b> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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# WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen/Germany

To be completed by the supplier/importer:  
Company stamp/signature:

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\_\_\_\_\_  
Tipo del aparato: Número de serie:

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\_\_\_\_\_  
Firma/Fecha:

To be filled out by instrument owner:  
Return address/company stamp:

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KARL STORZ GmbH & Co. KG  
Resident Representative Office  
80/33 (44/19) Dang Van Ngu,  
F:10-Q. Phu Nhuan  
Ho Chi Minh City, Vietnam  
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.  
#05 – 08 San Centre, 171 Chin Swee Road  
Singapore 0316, Singapore  
Phone: +65 6376-1066, Fax: +65 6376-1068  
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.  
174 Parramatta Rd.  
Camperdown, N.S.W. 2050, Australia  
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199  
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.  
Hong Kong Representative Office  
Unit 1601, Chinachem Exchange Square  
1 Hoi Wan Street, Quarry Bay,  
Hong Kong, People's Republic of China  
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Beijing Representative Office  
Room 610, China Life Tower  
No. 6, Chaowai Street  
Beijing, 100020, People's Republic of China  
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Unit 3901-3904, Tower 1 Grand Gateway  
No. 1 Hong Qiao Road  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Chengdu Representative Office  
F-5, 24/F., Chuanxing Mansion,  
No. 18 Renming Road South  
Chengdu, Sichuan, 610016, People's Rep. of China  
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Guangzhou Representative Office  
Room 1119-20, Dongshan Plaza  
69 Xianlie Road Middle, Dongshan District,  
Guangzhou, Guangdong, 510095,  
People's Rep. of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd  
3791 Jalan Bukit Merah  
06-11 e-Centre @ Redhill  
Singapore 159471, Singapore  
Tel. No. +65 63761066, Fax. No. +65 63761068  
Email : infoasia@karlstorz.com.sg  
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
3791 Jalan Bukit Merah  
10-18 e-Centre @ Redhill  
Singapore 159471  
Tel. No. +65 65325548, Fax No. +65 65323832  
Email : infosgpp@karlstorz.com.sg  
Email : servicesgpp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.  
Bois Hongo Building 6F1, 3-42-5 Hongo  
Bunkyo-ku, Tokyo 113-0033, Japan  
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)