



Yes, you can.®

At'm QT



Take Along Chair

Anything else is a compromise.™



Yes, you can.®

Wherever you go, you're up and At'm in 60 seconds.

Wherever you're going, At'm QT goes, too. It's the power chair lightweight and compact enough to fit in the trunk of a car. Simple enough to assemble in 60 seconds. Powerful enough for a full day of adventure.

How easy is At'm QT? If you can set up a lawn chair, you're ready. No mobility lifts, no tools, just three easy pieces to assemble. Open the seat, snap it on the lightweight base, add the quick disconnect battery and go. No wonder we call At'm QT the Take Along Chair!



The At'm QT goes wherever you go.

The three pieces of At'm QT fit in the smallest car trunk. And it's so lightweight, you don't need a mobility lift. That means you can travel in anyone's car – no van or special transport required. At'm QT can even be checked as airline baggage.



The At'm QT assembles in seconds.

How easy is it to get your At'm QT up and at'm? First, place the base on the ground – the wheel locks will hold it steady. Next, fold open the seat like a lawn chair and lock it onto the base. Finally, put the battery in place and plug it in. It's that simple!



The At'm QT was made easy for you to drive. Just flip the switch and go.

A touch of the joystick has you moving up to 4 mph. And At'm QT is spunky and maneuverable enough to take you wherever you want to go. Around the kitchen, across thresholds, over the curb or through the mall. So what are you waiting for? With At'm QT, nothing's stopping you from going, doing, and enjoying life.

Dimensions At'm QT

Length	
w/footrests folded up	28.5"
w/footrests down	35.5"
Width	
base w/o joystick	25.5"
over arms and joystick	25"
Seat-to-floor	20"
Weight capacity	250 lb.
Speed range	0 - 4 mph
Incline capability	6'
Weight	
Base only	34 lb.
Seat only	14 lb.
Battery pack	19 lb.
Ground clearance	2" minimum
Operating range*	Up to 8 miles
Warranties	
Frame	3 years
Electronics	1 year

* Range varies with battery condition, surface, terrain and operator weight

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

The Invacare® At'm QT is not intended for street or highway use.

Invacare, The Blue Medallion Design and Yes, you can. are registered trademarks. Specifications are subject to change without notification.

© 2005 Invacare Corporation
Form No. 05-128



"Invacare says 'Yes, you can.' and I couldn't have said it better myself."